

Endotoxin Measurement System

Toxinometer[®] ET-7000

- ▶ Supports "the Bacterial Endotoxin Test" in the United States / European / Japanese Pharmacopoeias
- ▶ Compliant with FDA 21 CFR, Part 11
- ▶ Conforms to the international certification standards for C-UL (CSA) and CE



Since 1985

State-of-the-art analysis system configured to ensure data integrity*

* Data integrity refers to the completeness, accuracy, and consistency of data.




Endotoxin

What is Endotoxin?

Endotoxin is a lipopolysaccharide (LPS) that comprises the cell wall of Gram-negative bacteria. Endotoxin induces various biological reactions, such as fever, when even a small amount (i.e. pg) enters the bloodstream. Endotoxins exist in environments where gram-negative bacteria reside and remain even after the bacteria die. They cannot be deactivated completely by autoclaving because of their resistance to heat. According to the Pharmacopoeia, **dry heat sterilization for at least 30 minutes at a temperature equal to 250°C or higher is required for complete deactivation of endotoxins.**

Heat stabilization of endotoxin (LPS: 500 ng/vial)





* The make and model of the PC pictured may differ from those actually supplied.

High Reliability

High quality specialized glass tubes make it possible to maintain...

- Complete Sterilization
- Endotoxin Free Environment
- No Cross Contamination

High Sensitivity

When combined with our Chromogenic reagent can offer...

- A maximum sensitivity of 0.0002 EU/mL.
- More than 3 digits dynamic range of 0.0002 - 0.5 EU/mL.

(when using Limulus Color KY single test.)

Flexibility

ALL 3 methods available in just one system.

- Gel-clot technique
- Turbidimetric technique
- Chromogenic technique

Allows for continuous sample measurement
Measures Max. 128 samples simultaneously

Expandability to Support Future Growth

- Controlled remotely from a Windows® PC
- Up to eight measurement modules can be additionally installed

Wako's Toxinometer[®], Highly Advanced Technology for Bacterial Endotoxin Testing

Gel-clot technique

Chromogenic technique


Turbidimetric technique




Offering reliable systems that satisfy global needs

- Meets Global Standards (C-UL (CSA) / CE).
- Pre-installed BET compliant software

All 3 techniques for BET are available on one system



Recovery of Spiked-Endotoxin in Products




Single Test Reagent for More Specific Measurements

The easiest and most beneficial way to use the Toxinometer[®] is by utilizing our Single Test Reagent vials, which allow you to just spike your sample and go!



Principal of Measurement

Light from an LED goes through the reaction tube filled with reaction mixture via aperture diaphragm 1. The light passes through the reaction mixture and then, while coming through aperture diaphragm 2, is detected by the silicon photocell.




Determination of Kinetic Measurement



The Toxinometer® reports your Reaction Time as either Gelation Time (Tg for Turbidimetric technique) or Activation Time (Ta for Chromogenic technique) based on the methodology used. The higher the concentration of endotoxin present, the shorter the reaction time.

The Toxinometer® measures the Transmittance Ratio of each well independently and simultaneously.

The Reaction Time is determined when the Transmittance Ratio falls below the Threshold value.



Time-Course-Graph & Standard Curve Example: Using Limulus Color KY Single Test Wako



Reagent Examples for Toxinometer

| Technique | Code No. | Product Name | Quantitative Range (EU/mL) | Gel-clot Sensitivity (EU/mL) | Quantity | Kit Contents |
|--------------------------|-------------|---|----------------------------|------------------------------|-----------|-------------------------------|
| Turbidimetric & Gel-clot | WPESK-0015 | PYROSTAR™ ES-F Single Test, 0.015 EU/mL | 0.001-10 | 0.015 | 25 tests | 25 vials + CSE 1 vial |
| | WPEK4-50015 | PYROSTAR™ ES-F Multi Kit, 0.015 EU/mL | 0.001-10 | 0.015 | 200 tests | 4 vials × 5.2 mL + CSE 1 vial |
| Chromogenic | 291-53601 | Limulus Color KY Single Test Wako | 0.0002 - 5 | — | 25 tests | 25 vials + CSE 1 vial |
| | 291-53101 | Limulus Color KY Test Wako | 0.0005 - 5 | — | 60 tests | 3 vials × 2 mL + CSE 1 vial |

* CSE : Control Standard Endotoxin


* A wide variety of other reagents are available. Please contact us for the other reagents.

Application

Wide range of Applications


Bacteria are broadly divided into two categories: Gram-negative and Gram-positive bacteria. Every bacterium has peptidoglycan as a component of the cell wall. A Gram-negative bacterium contains endotoxin in the outer membrane of the cell wall. The cell wall of a fungus contains β -glucan.

In combination with dedicated reagents, the Toxinometer® can be used for a wide range of applications such as research and monitoring of microbial contamination.



SLP Reagents

SLP reagent is a freeze-dried product prepared from silkworm larvae plasma. The reaction mechanism is shown in the following figure. When the reagent reacts with peptidoglycan and β -glucan, it eventually forms melanin, resulting in a black coloration of the sample. As when utilizing lysate reagent, a highly-sensitive measurement of this coloration is possible with the Toxinometer®.




SLP-HS Single Reagent Set II

| Code No. | Product name | Test# |
|-----------|------------------------------|-------|
| 296-81001 | SLP-HS Single Reagent Set II | 20 |

SLP Reagent reaction mechanism

LAL Reagent Bacterial Endotoxin Test (BET)

A lysate reagent prepared from the amebocytes of the Atlantic horseshoe crab (*Limulus polyphemus*) is used to detect bacterial endotoxins. As shown in the Figure below, the cascade reactions begin due to the presence of endotoxin, whereby Factor C, a serine protease precursor, is initially activated. There follows the sequential activation of Factor B, also a serine protease precursor and a pro-clotting enzyme, which hydrolyzes coagulogen into coagulin, forming an insoluble gel. In the Bacterial Endotoxin Test, endotoxin can be quantified in three ways: measurement of gel formation, increased turbidity, or release of a synthetic substrate. Endotoxin-specific LAL reagents are not activated by (1 \rightarrow 3)- β -D-glucan, as opposed to other BET compliance tests.



 : activated factors


LAL Reagent reaction mechanism

Basic configuration

Toxinometer® ET-7000/E Part 11 Set [Code No. 299-35821]
Toxinometer® ET-7000/U Part 11 Set [Code No. 294-35871]

[Contents]

- Toxinometer® ET-7000 (1 unit)
- Toximaster® QC8 Part 11 (5 licenses)
- Personal Computer (1 unit)
- Toxinometer® ET-7000 Part 11 System Validation Documents




The system complies with FDA 21 CFR, Part 11.

Toxinometer® ET-7000/E Non Part 11 Set [Code No. 296-35831]
Toxinometer® ET-7000/U Non Part 11 Set [Code No. 291-35881]

[Contents]

- Toxinometer® ET-7000 (1 unit)
- Toximaster® QC7
- Personal Computer (1 unit)



The system doesn't comply with FDA 21 CFR, Part 11.

Options

Toxinometer® ET-7000/E [Code No. 292-35811]
Toxinometer® ET-7000/U [Code No. 297-35861]

※ Additional module for Part 11 Set or Non-Part 11 Set.

Specification (*1)

| Item | Explanation |
|---------------------|---|
| Functions | Transmitted light quantity measuring function (capable of measuring 16 samples simultaneously) Temperature control function Automatic light quantity check function |
| Light source | High intensity blue LED Central wavelength: 430 nm |
| Detector | Silicon photocell |
| Temperature control | Dry bath: 30±1.0 degrees C/37±1.0 degrees C (can be changed by software for ET-7000) Warmup time: 20 minutes (when preset temperature is 37 degrees C and surrounding temperature is 25 degrees C) |
| Display | The LED indicates measurement The LED indicates errors and information during checking |
| Weight | 6.3 kg (±10%) |
| Size | W 190 mm × D 420 mm × H 130 mm (protrusions not included) |

| Item | Explanation |
|-------------------|---|
| Power source | 220-240 (±10%) VAC (Model: ET-7000/E) 100-120 (±10%) VAC (Model: ET-7000/U) |
| Frequency | 50/60 Hz |
| Power consumption | Max 120 W |
| Environment | During operation |
| | When temperature is set at 37 degrees C Ambient temperature: 15 to 30 degrees C Humidity: 30 to 85%, non-condensing When temperature is controlled at 30 degrees C Ambient temperature: 15 to 25 degrees C Humidity: 30 to 85%, non-condensing |
| | During stored |
| | Ambient temperature: -20 to 60 degrees C Humidity: 30 to 85%, non-condensing |
| Location | Indoor |
| Altitude | 2000 m or lower |

(*1) These specifications are common to the Toxinometer® ET-7000 Non Part 11 Set and Part 11 Set. Up to 8 measurement modules can be connected to the unit to enable simultaneous measurement of 128 samples.


Related products

Endotoxin-free tip BioCleanTip Wako®

| Code No. | Product Name | Package |
|-----------|--------------------------------------|---------|
| 291-35021 | BioCleanTip Wako® 200 II 200 µL | 100 pcs |
| 298-35031 | BioCleanTip Wako® 1000 II 1000 µL | 100 pcs |
| 294-35011 | BioCleanTip Wako® Extend S II 200 µL | 100 pcs |

Endotoxin-free, Test Tube for Endotoxin Test and Aluminum Cap

| Code No. | Product Name | Size | Quantity |
|-----------|---------------------------------------|--------------|-------------|
| 292-32751 | Limulus Test Tube-S with Aluminum Cap | φ 12 × 75 mm | 10 pcs × 8 |
| 293-26551 | Limulus Test Tube-S | φ 12 × 75 mm | 10 pcs × 10 |
| 293-28251 | Aluminum Cap-S | φ 15 × 18 mm | 10 pcs × 10 |



"Microsoft® and Windows® are registered trademarks of Microsoft Corporation in the United States and other countries.
All other marks and names mentioned herein may be trademarks of their respective companies."

Contents of this catalog are subject to change without notice.

FUJIFILM Wako Pure Chemical Corporation

<http://fwk.fujifilm.co.jp/en/index.html>
1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
fwk-cservice@fujifilm.com Online Catalog: www.e-reagent.com

FUJIFILM Wako Chemicals U.S.A. Corporation

www.wakousa.com
1600 Bellwood Road, Richmond, VA 23237, U.S.A.
Toll-Free (U.S. only): +1 877 714 1920
Tel: +1 804 714 1920 Fax: +1 804 271 7791
wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH

www.wako-chemicals.de
Fuggerstr 12, D-41468 Neuss, Germany
Tel: +49 2131 311 0 Fax: +49 2131 311 100
labchem@wako-chemicals.de