High Molecular Amyloid β Oligomer ELISA Kit Wako Ver.2

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Amyloid β ($A\beta$) is a peptide consisting of about 40 amino acids and accumulates in the brain of Alzheimer disease patients forming a structure known as senile plaques. Therefore, $A\beta$ is considered as a causal factor of Alzheimer disease. In recent years, an $A\beta$ oligomer hypothesis has been focused. The hypothesis is that soluble aggregates called $A\beta$ oligomers resulting from polymerization of $A\beta$ monomers, causes synaptic failure in nerve cells, not insolubilized and accumulated $A\beta$ fibrils. There are several types of $A\beta$ oligomers due to the difference in molecular weight, and low molecular $A\beta$ oligomers (such as low-n oligomer) and high molecular $A\beta$ oligomers (such as ADDL, $A\beta$ *56, proto-fibril) have been reported. This product is an ELISA kit that can specifically measure high molecular $A\beta$ oligomers.

【Performance】

Standard curve range	Human cerebrospinal fluid (CSF) : 0.41 ~ 100 pM Human plasma (EDTA), serum : 0.16 ~ 40 pM *Based on 16-mer MAP peptide
Reactive $A\beta$ oligomers	\geq 9-mer $A\beta$ oligomer
Sample	Human CSF Human plasma (EDTA) Human serum in vitro $A\beta$ oligomer
Sample volume	Human CSF : 25 μ L (4 times dilution) Human plasma (EDTA) or serum : 50 μ L (2 times dilution)
Assay time	4.5 hours
Detection method	Luminescence method

【Reagents supplied】

Components	Condition	Volume
(1) Antibody (BAN50)-coated Plate	Use as it is.	1 plate, 96 wells (8 × 12)
(2) Standard Solution (10 nM)	Use after preparation.	100 μ L/1 vial
(3) Buffer	Use as it is.	60 mL/1 bottle
(4) Biotin-conjugated Antibody (BAN50) Solution	Use after preparation.	100 μ L/1 vial
(5) Peroxidase-conjugated Streptavidin Solution	Use after preparation.	100 μ L/1 vial
(6) Luminescent Reagent 1	Use as it is.	6 mL/1 bottle
(7) Luminescent Reagent 2	Use as it is.	6 mL/1 bottle
(8) Wash Solution (10 ×)	Use after preparation.	100 mL/1 bottle
(9) Sample Buffer for Serum/Plasma	Use as it is.	50 mL/1 bottle
(10) Plate Seal	Use as it is.	4 sheets
Instruction manual		1 copy

2.56	6.40 pM standard solution : 200 μ L	300 μ L
1.02	2.56 pM standard solution : 200 μ L	300 μ L
0.41	1.02 pM standard solution : 200 μ L	300 μ L
0.00	–	300 μ L

For human plasma (EDTA) or serum

Concentration (pM) of standard solution	Volume of the standard solution	Buffer
100	Standard Solution (10 nM) : 10 μ L	990 μ L
40.0	100 pM standard solution : 200 μ L	300 μ L
16.0	40.0 pM standard solution : 200 μ L	300 μ L
6.40	16.0 pM standard solution : 200 μ L	300 μ L
2.56	6.40 pM standard solution : 200 μ L	300 μ L
1.02	2.56 pM standard solution : 200 μ L	300 μ L
0.41	1.02 pM standard solution : 200 μ L	300 μ L
0.16	0.41 pM standard solution : 200 μ L	300 μ L
0.00	–	300 μ L

- (4) Biotin-conjugated Antibody (BAN50) Solution
Dilute the solution 100 times with the Buffer.
- (5) Peroxidase-conjugated Streptavidin Solution
Dilute the solution 100 times with the Buffer.
- (6) Luminescent Reagent 1 and (7) Luminescent Reagent 2
Mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.) 15 to 30 minutes before use. Store the reagent solution away from light until use.
Example : Mix 6 mL of Luminescent Reagent 1 and 6 mL of Luminescent Reagent 2 (when all of 96 wells are used)
- (8) Washing Solution (10 \times)
Dilute the solution 10 times with purified water (distilled water).
Example : Dilute 100 mL of the concentrated wash solution (10 \times) with 900 mL of purified water (distilled water) (when all of 96 wells are used)
- Use other reagents as they are.

[Storage and stability of reagents]

- (1) Antibody (BAN50)-coated Plate
Return unused (under a sealed and refrigerated condition) antibody-coated strips in the zip-sealed pack in the container and store at 2 to 10°C. The stability is maintained within the expiration date.
- (2) Standard Solution (10 nM)
When the kit is divided and used, take out the standard solution from the refrigerator and dilute it immediately before use and keep the remaining stock solution with the cap tightly closed at 2 to 10°C without returning to room temperature. The stability is maintained within the expiration date. Use the respective diluted standard solutions immediately and do not store them after use.
- (3) Buffer
If a part of the solution is used, transfer a volume slightly greater than the necessary volume in another container. Promptly close the cap of the container containing the remaining solution tightly and store at 2 to 10°C

without returning to the room temperature. The stability is maintained within the expiration date.

(4) Biotin-conjugated Antibody (BAN50) Solution and (5) Peroxidase-conjugated Streptavidin Solution

When the kit is divided and used, take out the solutions from the refrigerator and dilute them immediately and keep the remaining stock solutions with the caps tightly closed at 2 to 10°C without returning to room temperature. The stability is maintained within the expiration date. Discard the diluted solutions remaining after use.

(6) Luminescent Reagent 1 and (7) Luminescent Reagent 2

When the kit is divided and used, take out the reagents from the refrigerator and prepare them. Keep the remaining luminescent reagents with the caps tightly closed at 2 to 10°C without returning to room temperature. The stability is maintained within the expiration date.

(8) Washing Solution (10 ×)

Store the Wash Solution (10 ×) at 2 to 10°C with its cap tightly closed. The stability is maintained within the expiration date. Discard the diluted wash solution remaining after use.

(9) Sample Buffer for Serum/Plasma

If a part of the solution is used, transfer a volume slightly greater than the necessary volume in another container. Promptly close the cap of the container containing the remaining solution tightly and store at 2 to 10°C without returning to the room temperature. The stability is maintained within the expiration date.

[Preparation for sample]

[CSF] dilute 2-4 times CSF with buffer at room temperature.


[Plasma (EDTA) or serum] dilute 2 times plasma (EDTA) or serum with Sample Buffer for Serum/Plasma and incubate 30 min at room temperature.

[Measurement procedure]

1. Wash the Antibody (BAN50)-coated Plate 4 times by filling the wells with the prepared wash solution. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
2. Dispense 100 μ L each of the standard solutions at the respective concentrations in the reference standard measurement wells.
3. Dispense 100 μ L each of the diluted samples in the wells.
4. Shake the plate using a microplate shaker.
5. Attach the Plate Seal and leave to stand at room temperature (20 to 25°C) for 2 hours.
6. After completion of the reaction, discard the reaction solution, fill the wells with the wash solution, and wash the plate 4 times. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
7. Dispense 100 μ L each of the biotin-conjugated antibody (BAN50) solution to the wells. Shake the plate using the microplate shaker.
8. Attach the Plate Seal and leave to stand at room temperature (20 to 25°C) for 2 hours.
9. After completion of the reaction, discard the reaction solution, fill the wells with the wash solution, and wash the plate 4 times. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
10. Dispense 100 μ L each of the peroxidase-conjugated streptavidin solution to the wells. Shake the plate using the microplate shaker.

11. Attach the Plate Seal and leave to stand at room temperature (20 to 25°C) for 30 minutes.
12. After completion of the reaction, discard the reaction solution, fill the wells with the wash solution, and wash the plate 4 times. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
13. Dispense 100 μ L each of the mixed luminescent reagents to the wells. Shake the plate for 1 minute using the microplate shaker.
14. After agitation, determine the luminescent intensity using the 96-well microplate reader (for luminescence measurement). It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.
15. Create a standard curve by plotting the concentrations (pM) of the standard solutions along the X axis and the luminescence intensity along the Y axis. Read the concentration (pM) corresponding to the luminescence intensity of the diluted sample. Multiply the concentration read with the sample dilution factor to obtain a measured value. Convert the measured value to the concentration based on the 16-mer MAP peptide.
 - * Use of a tertiary polynomial equation with 4 or 5 parameters is recommended for the arithmetic processing using computer software.

Calibration curve

The concentration is a value converted on the basis of 16-mer MAP peptide.

[Precautions for use]

1. Use polypropylene containers when handling the sample. Storage in a freezer at -35°C or lower is recommended when a sample is stored for a long period of time. Do not repeat freezing and thawing. Thaw the frozen sample and thoroughly shake just before measurement. Dilute the sample just before use.
2. Do not use a sample with hemolysis or containing high lipid.
3. Centrifuge a sample with opacity and insoluble matter before measurement.
4. With a sample suspected of the influence of an interfering substance, check the linearity of dilution of the same sample using 2 or more different dilution factors. When a sample is diluted, dilute it with the Buffer in advance using test tubes (PP, PE) and dispense in the measurement wells. The dilution factor is 2 ~ 4-time dilution in the standard operating procedure.
5. This kit should be used by a person who has completed the training for the ELISA or under his/her supervision.
6. For measurement by the manual operation, the kit should be used by a person who is well-trained for stable reproducibility of pipetting operation.
7. Wear gloves, goggles and protective gown during preparation and operation

of the kit.

8. Take care not to touch the reagents to the skin. If any of the reagents of this kit has adhered to the eye, mouth, wound or skin by error, promptly take emergent measures such as washing thoroughly with tap water and visit a physician if necessary.
9. Do not eat, drink or smoke at a place where this kit is used.
10. Handle samples with thorough care as a substance with a risk of infection. This kit contains animal-derived ingredients.
11. Immerse used samples and used consumables in 1% formalin, 2% glutaraldehyde or sodium hypochlorite solution at 0.1% or higher for at least 1 hour. Or, autoclave these before disposal. Discard the used consumables and unused drugs in accordance with the rules of the testing site and local regulations.
12. Do not mix the reagent of different lot numbers. Be sure to affix the Plate seal during standing for reaction in each step to prevent well drying, contamination of foreign materials, temperature bias and evaporation of the dispensed reagents.
13. The ELISA is affected by the measurement environment. Strictly observe measurement procedure and room temperature of 20 to 25 °C (on the experimental table or incubator) at the place for reaction on standing. Never perform the measurement in an environment with air blowing (including circulation from the air-conditioner) and low humidity.

[Example of measurement]

○ Spiked recovery test


Human plasma	Not spiked	+0.5 pM Std.		+0.1 pM Std.		+5 pM Std.		+30 pM Std.	
	Assay value	Assay value	Recovery	Assay value	Recovery	Assay value	Recovery	Assay value	Recovery
	1.85 pM	2.36 pM	102%	2.85 pM	100%	7.11 pM	105%	31.2 pM	97.8%

Human serum	Not spiked	+0.5 pM Std.		+0.1 pM Std.		+5 pM Std.		+30 pM Std.	
	Assay value	Assay value	Recovery	Assay value	Recovery	Assay value	Recovery	Assay value	Recovery
	1.58 pM	2.07 pM	98%	2.59 pM	101%	6.64 pM	101%	31.9 pM	101%

Human CSF1	Not spiked	+3 pM Std.		+6 pM Std.		+15 pM Std.	
	Assay value	Assay value	Recovery	Assay value	Recovery	Assay value	Recovery
	5.16 pM	8.01 pM	95%	10.6 pM	91%	18.1 pM	86%

Human CSF2	Not spiked	+10 pM Std.		+15 pM Std.		+35 pM Std.	
	Assay value	Assay value	Recovery	Assay value	Recovery	Assay value	Recovery
	29.1 pM	38.8 pM	97%	44.0 pM	99.3%	59.9 pM	88%


○ Dilution linearity

○ Assay of clinical samples

Cerebrospinal fluid samples from subjects without Alzheimer's disease (non-AD), patients with mild cognitive impairment (MCI), and those with Alzheimer's disease (AD) were assayed using this ELISA kit.

	Age	Sex	MMSE	CSF	
				A β oligomer (pM)	
Control	69	M	30	2.37	
	60	F	30	1.65	
	59	M	30	2.82	
	46	M	28	2.24	
	60	F	29	1.34	
MCI	69	F	24	3.6	
	76	M	23	2.62	
	82	M	26	2.33	
	69	F	23	2.67	
	73	F	25	2.34	
AD	79	F	10	3.56	
	88	M	24	3.13	
	62	M	16	4.13	
	78	M	13	3.95	
	63	M	14	3.88	
	74	F	16	3.25	
	76	M	19	3.69	
	84	M	16	4.34	

[Outline of the measurement procedure]

- Return the plate and reagents thoroughly to room temperature (20 - 25°C).
- Dilution of the concentrated wash solution : Dilute 10 times with purified water returned to room temperature.
- Dilute CSF 2-4 times with buffer at room temperature (example : sample (60 μ L) + Buffer (180 μ L)).
- Dilute serum or plasma (EDTA) 2 times with Sample Buffer for Serum/Plasma and incubate 30 min at room temperature (example : sample (120 μ L) + Sample Buffer for Serum/Plasma (120 μ L)).
- Dilution of the standard solution (example) : Dilute with the Buffer returned to room temperature.

For CSF									
Concentration (pM)	100	40.0	16.0	6.40	2.56	1.02	0.41	0	
Standard solution (μ L)	10nM: 10	200*	200*	200*	200*	200*	200*	200*	0
Buffer (μ L)	990	300	300	300	300	300	300	300	300

For plasma (EDTA) or serum									
Concentration (pM)	100	40.0	16.0	6.40	2.56	1.02	0.41	0.16	0
Standard solution (μ L)	10nM: 10	200*	200*	200*	200*	200*	200*	200*	0
Buffer (μ L)	990	300	300	300	300	300	300	300	300

* : Standard solution at a 1-step higher concentration

- Antibody (BAN50)-coated Plate**
- ↓ Washing 4 times (* (i))
- Sample or standard solution** **100 μ L/well**
- ↓ Shake(* (ii)), 2-hour reaction at room temperature (20 - 25°C), standing(* (iii))
- ↓
- *Preparation of the biotin-conjugated antibody (BAN50) solution (dilute 100 times with the Buffer returned to room temperature)
- ↓ Washing 4 times (* (i))
- Biotin-conjugated antibody (BAN50) solution** **100 μ L/well**
- ↓ Shake(* (ii)), 2-hour reaction at room temperature (20 - 25°C), standing(* (iii))
- ↓

- *Preparation of the peroxidase-conjugated streptavidin solution (dilute 100 times with the Buffer returned to room temperature)
- ↓ Washing 4 times (* (i))
- Peroxidase-conjugated streptavidin solution** **100 μ L/well**
- ↓ Shake(* (ii)), 30-minute reaction at room temperature (20 - 25°C), standing (* (iii))
- ↓
- *Preparation of the luminescent reagent (mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.))
- ↓ Washing 4 times (* (i))
- Luminescent reagent** **100 μ L/well**
- ↓ Shake for 1 minute at room temperature (20 - 25°C)
- Measurement of luminescence intensity (to be measured between 10 and 20 minutes)**

(* (i)) In each washing, dispense the wash solution in the wells, shake lightly for 10 seconds on the palm and discard. After 4 consecutive washings, invert the plate on a paper towel and pat to eliminate the wash solution completely. After elimination of the wash solution, promptly dispense the subsequent solution taking care not to dry the wells. The volume of the wash solution to be pipetted to the plate is about 300 μ L/well.

(* (ii)) Shake 3 times each at about 600 - 800 rpm for 10 seconds.

(* (iii)) After completion of agitation, attach a Plate Seal and leave to stand. Peel the protective paper from the Plate Seal and attach the adhesive side to the plate. Never re-use the Plate Seal used once.

【Storage condition】 Store at 2 to 10°C.

【Expiration date】 Indicated on the label

【Packaging】 For 96 tests

【References】

- 1) Fukumoto, H., *et al.*: *FASEB J.*, 24, 2716 (2010).
- 2) Kasai, T., *et al.*: *Biochem. Biophys. Res. Commun.*, 422, 3, 375 (2012).
- 3) Kasai, T., *et al.*: *Neurosci. Lett.*, 551, 17 (2013).
- 4) Chambers, J. K., *et al.*: *Acta. Neuropathol. Commun.*, 3, 78 (2015).
- 5) Kasai, T., *et al.*: *PLoS One*, 12, 4, e0174630 (2017).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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