LabAssay™ Ammonia (Modified Fujii-Okuda Method)

For Research Use Only Not for use in diagnostic procedure

[Introduction]

Ammonia is converted to urea in the urea cycle and excreted.

LabAssay™ Ammonia is the reagent kit for assay of ammonia based on a colorimetric method using the Berthelot reaction. This kit is used for the quantitative determination of ammonia in mouse, rat and human blood. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

[Kit contents]

	Kit contents	State	Amount
(1)	Deproteinizing Reagent Sodium Tungstate Dihydrate Phosphoric Acide	liquid	100 mL/1 bottle
(2)	Chromogen Reagent A Phenol Sodium Pentacyanonitrosylferrate (III) Dihydrate	liquid	50 mL/1 bottle
(3)	Chromogen Reagent B Potassium Hydroxide	liquid	25 mL/1 bottle
(4)	Chromogen Reagent C Potassium Carbonate Sodium Hypochlorite	liquid	50 mL/1 bottle
(5)	Standard Solution (NH ₃ -N 400 µg/dL)	liquid	15 mL/1 bottle
(6)	Dilute Solution for Standard	liquid	20 mL/1 bottle


※Assay in a microplate format; 700 tests. Assay in a test tube format; 60 tests.

[Storage and expiration]

When the complete kit is stored at 2 °C~10 °C, the kit is stable until the expiration date shown on the label on the container.

[Assay principle]

The Deproteinizing Reagent deactivates the enzymes and removes the components that inhibit color reaction in a sample. After adding Deproteinizing Reagent to a sample and centrifuging, add Phenol, Sodium Pentacyanonitrosylferrate(III) Dihydrate, and Sodium Hypochlorite to the supernatant. The indophenol blue pigment is obtained by the following schematic reaction. The amount of ammonia contained in the sample is determined by measuring the absorbance of the blue color.

※Sodium pentacyanonitrosylferrate(III) dihydrate acts as catalyst

[Materials and apparatuses to be prepared]

(For microplate method)

- 96 wells microplate (transparent type) ·Micropipette ·Incubator maintaining at 37 °C *
- Plate mixer *
- Microplate reader with a 630 nm wavelength filter (* if the microplate reader is not equipped.)

(For test tube method)

- Test tube or microtube ·Pipette ·Incubator maintaining at 37 °C
- Spectrophotometer with a 630 nm wavelength filter

[Preparation of reagents to be used]

Before the measurement, transfer the required reagents to the tubes and adjust the reagents to the following temperatures.

- Return Ammonia Standard Solution, Deproteinizing Reagent, and Dilute Solution for Standard to room temperature*.
- Incubate Chromogen Reagent A, B, and C at 37 °C for an hour.

* Room temperature: 20 °C~25 °C

(1) Standard Solution (Microplate method)

Standard Solution is prepared by dilution of the provided Standard.

No.	Standard Solution	Dilute Solution for Standard	Concentration
1	50 μL	150 μL	100 μg/dL
2	100 μL	100 μL	200 μg/dL
3	150 μL	50 μL	300 μg/dL
4	undiluted	—	400 μg/dL

(2) Standard Solution (Test tube method)

Standard Solution is prepared by dilution of the provided Standard.

No.	Standard Solution	Dilute Solution for Standard	Concentration
1	200 μL	600 μL	100 μg/dL
2	400 μL	400 μL	200 μg/dL
3	600 μL	200 μL	300 μg/dL
4	undiluted	—	400 μg/dL

[Procedure]

Before the measurement, transfer the required reagents to the tubes and adjust the reagents to the following temperatures.

- Return Standard Solution, Deproteinizing Reagent and Dilute Solution for Standard to room temperature*.

- Incubate Chromogen Reagent A, B and C at 37 °C for an hour.

* Room temperature: 20 °C~25 °C

(1) Assay in a microplate format

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
An experiment can be conducted in a microtube			
Deproteinizing Reagent	Deproteinizing reagent : Sample = 4 : 1 (v/v)	Deproteinizing reagent : Standard solution = 4 : 1 (v/v)	70 µL
Sample			—
	Mix well. Centrifuge 5000×g, 4°C, for 15 min. Subsequent operations should be performed in the wells	Mix well. Subsequent operations should be performed in the wells.	
	Supernatant 70 µL	Mixed solution 70 µL	—
Chromogen Reagent A	70 µL	70 µL	70 µL
Chromogen Reagent B	35 µL	35 µL	35 µL
Chromogen Reagent C	70 µL	70 µL	70 µL
After adding each reagent to a microplate, mix well and incubate at 37 °C, for 20 min. Mix well, then let stand at room temperature* for 30 min. Mix well, then measure the absorbance of the test sample and standard solution at 630 nm with the blank as the control.			

* Room temperature : 20 °C~25 °C

(2) Assay in a test tube

Perform the assay in the tubes according to the following table scheme.


	Test	Standard	Blank
Operations should be performed in the micro tubes.			
Deproteinizing Reagent	Deproteinizing reagent : Sample = 4 : 1 (v/v)	Deproteinizing reagent : Standard solution = 4 : 1 (v/v)	800 μL
Sample			—
	Mix well. Centrifuge 5000×g, 4 °C, for 15 min.	Mix well.	
	Transfer the supernatant (800 μL) to another tube.	Transfer mixed solution (800 μL) to another tube.	—
Chromogen Reagent A	800 μL	800 μL	800 μL
Chromogen Reagent B	400 μL	400 μL	400 μL
Chromogen Reagent C	800 μL	800 μL	800 μL
After adding each reagent to a tube, mix well and incubate at 37 °C, for 20 min. Mix well, and place in water to cool (20 °C~25 °C) for 30 min. Mix well, then measure the absorbance of the test sample and standard solution with the blank solution as the control. Spectrophotometer : 630 nm			

*Room temperature : 20 °C~25 °C


※Assay in a test tube; 60 tests.

[Standard curve]

[assay in a microplate]

[assay in a test tube]

[Performance]

(1) Sensitivity [assay in a test tube method]

- The absorbance is 0.010 to 0.070, when measuring purified water as a sample.
- The absorbance is 0.200 to 0.370, when measuring 400 µg/dL ammonia as a sample.

(2) Specificity

The ammonia concentration is less than $\pm 35\%$, when measuring the known concentration of control blood as a sample.

[Usage Notes]

(1) Sample

- Ammonia in whole blood increases at room temperature. Therefore, collect the blood and transfer to test tubes or microtubes containing Deproteinizing Reagent and centrifuge the samples to separate supernatant as quickly as possible.
- Store the supernatant on ice or at $2\text{ }^{\circ}\text{C}\sim 10\text{ }^{\circ}\text{C}$ and measure the sample within an hour.
- Do not use anticoagulants such as heparin and double oxalate.
- Amino acid preparations can cause a slightly negative effect on the assay.
- Amino acids in the blood may affect the assay.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purpose other than those described herein.
- Reagents which are kept open for a long time cause dispersion of the measured values by absorbing ammonia. The opened vials should be capped and store at $2\text{ }^{\circ}\text{C}\sim 10\text{ }^{\circ}\text{C}$.
- Avoid exposing reagents to direct sunlight during operation.
- Do not use the reagents at other reaction temperatures and reaction times other than described herein.
- Wash the Equipment with a weak alkaline solution. Wash them the with tap water. Rinse them well with deionized water. Tap water and distilled water may contain a lot of volatile ammonium salts. Pay attention to washing of experimental instruments.
- When collecting the supernatant after centrifuging, pay attention not to mix with suspended matter. Suspended matter may cause dispersion of the measured values.
- This kit is for research use only. Not for use in diagnostic procedure.

(3) Safety precautions

- If reagents come into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult with a physician if necessary.
- Deproteinizing Reagent is an acid solution which is less than pH 3.
- Chromogen Reagent B and C are alkaline solutions which are more than pH 11.
- A pipette with a safety pipette filter should be used.

LabAssay™ Ammonia (Modified Fujii-Okuda method) (#LABNH3-M1)

(4) Waste

- The waste should be processed appropriately according to the laws and regulations.
- Chromogen Reagent A contains 0.015 % Sodium Pentacyanonitrosylferrate(III) Dihydrate and 40g/L Phenol.
- All the devices including, reagents and vials, that come into contact with the specimen should be considered potentially infectious.

LabAssay™ Ammonia (Modified Fujii-Okuda Method)

[Storage condition] Store the kit at 2 °C ~ 10 °C

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