## Caffeine ELISA Kit Wako

### [1. Introduction]

Caffeine is a type of alkaloid found in coffee beans and tea leaves. It has excitatory effects on the central nervous system and can stimulate the motor center and medullary respiratory center. It also exhibits vasoconstrictive, cardiotoxic, diuretic, and gastric secretion-promoting effects. Due to these properties, caffeine has been reported to be associated with various diseases including dementia, Parkinson's disease, migraines, and angina. This product is an ELISA kit that can measure caffeine in human saliva, serum or plasma and mouse/rat serum or plasma.

### [2. Assay Principle]

Each well of the assay plate is coated with an anti-caffeine antibody.

A standard solution or a sample is added to a well along with peroxidase-conjugated caffeine, and the mixture is allowed to react.

Finally, by measuring the peroxidase activity (absorbance) in the well, the concentration of caffeine in the sample can be determined.

### [3. Kit Performance]

Calibration curve range	0.244 ~ 1,000 ng/mL
Analyte	Caffeine
Sample	Mouse serum/plasma Rat serum/plasma Human saliva/serum/plasma (EDTA and Hep)
Sample amount required	55 $\mu$ L (2-fold dilution, duplicate)
Assay time	about 2 hours and 20 minutes
Detection method	Colorimetry
Measurement wavelength	450 nm/Ref, 620 nm

### [4. Precautions for Use]

1. It is important to bring the kit reagents to room temperature (20-25°C) before use (approximately 2 hours).
2. Prepare the reagents to be used in the next step in advance, before starting each step.
3. Use this kit after completing training in ELISA methods or under the guidance of an instructor.
4. For manual measurements, individuals who have demonstrated consistent reproducibility in pipette operation should perform the assay.
5. Wear gloves, glasses, and protective clothing during preparation and when using this kit.
6. Prevent reagents coming into contact with the skin. In the event the kit's reagents accidentally come into contact with

eyes, mouth, wounds, skin, etc., immediately wash thoroughly with running water and seek medical attention from a physician if necessary.

7. Do not eat, drink, or smoke in the area where this kit is being used.
8. Handle samples with particular care as they may pose a risk of infection. This kit contains animal-derived components.
9. Used specimens and consumables should be soaked in a sodium hypochlorite solution ( $\geq 0.1\%$ ) for at least one hour, or autoclaved and disposed of. Used consumables and unused chemicals should be disposed of in accordance with the regulations of the facility and local laws.
10. Do not mix reagents from different lots.
11. It is important to cover the plate with a plate seal to prevent drying of the wells, contamination, temperature fluctuations, and evaporation of dispensed reagents during incubation at each step.
12. ELISA methods can be affected by the assay environment. Strictly adhere to the following: Room temperature should be 20-25°C (on the bench or in the incubator) during operation and incubation. Do not perform the assay in environments with wind speeds (including air conditioning) exceeding 0.4 m/sec and humidity levels below 30%.

### [5. Materials supplied]

#### 5-1. Kit Components

Item	State	Amount
Antibody-coated Plate	Use after washing	96 wells (8 × 12)/1 plate
Caffeine Standard	Use after preparation	100 $\mu$ L/1 vial
Buffer	Ready-to-use	60 mL/1 vial
Peroxidase-conjugated Caffeine Solution	Use after preparation	50 $\mu$ L/1 vial
TMB Solution	Ready-to-use	12 mL/1 vial
Stop Solution	Ready-to-use	12 mL/1 vial
Wash Solution (10 ×)	Use after preparation	100 mL/1 vial
Plate Seal	—	4 seals
Instruction Manual	—	1 copy

#### 5-2. Handling of Unused Reagents

When using the kit in divided portions, refer to the following for the stability and storage of unused reagents.

If unopened, they remain stable through the expiration date.

- (1) Antibody-coated Plate  
The plate can be used in divided portions.  
Unused antibody-coated strips should be returned into the zip-lock bag provided in the kit and stored at 2-10°C.
- (2) Caffeine Standard  
Store unused portion at 2-10°C.  
Note : Do not store diluted standard solutions.
- (3) Buffer  
When using a portion of the solutions, transfer a slightly larger amount than necessary to separate containers.

The remaining unused solutions should be tightly capped and stored at 2-10°C immediately ; do not bring them to room temperature.

(4) Peroxidase-conjugated Caffeine Solution

When using the kit in divided portions, take out the solution from the refrigerator when ready to prepare dilutions. The remaining undiluted solution should be tightly capped and stored at 2-10°C immediately ; do not bring it to room temperature.

Discard any remaining diluted solutions.

(5) TMB Solution

When using a portion of the solution, transfer a slightly larger amount than necessary to another container.

The remaining unused solutions should be tightly capped and stored at 2-10°C immediately ; do not bring them to room temperature.

(6) Stop Solution

Store remaining unused solution at 2-10°C.

(7) Wash Solution (10×)

Cap tightly and store at 2-10°C.

Discard any remaining diluted Wash Solution.

**[6. Apparatus, equipment, and materials required]**

- Purified water (distilled water)
- Tubes for dilution of standard solutions and samples
- Glassware for dilution of Wash Solution (graduated cylinders, beakers)
- Micropipettes (for 10 μL disposable tips and 200 to 500 μL disposable tips)
- Continuous dispensing pipette, capable of dispensing 100 μL continuously (preferred if available)
- Paper towels or other absorbent material (to remove liquid left on the plate after washing)
- Vortex-type mixer
- Microplate shaker (range approx. 600 to 800 rpm)
- Automatic washer for 96-well plate (if available) or washing bottle
- Microplate reader for 96-well plates (for absorbance measurement)
- Software for data processing

**[7. Sample Preparation]**

As per the standard procedure, dilute the samples with the buffer solution provided in the kit and use for the assay.

It is recommended that you optimize the dilution ratio to fall within the range of standard curve concentrations.

Sample	Dilution ratio
Mouse serum/plasma	2-fold
Rat serum/plasma	
Human serum/plasma	
Human saliva	4-fold

Recommendations and precautions before sample preparation

- Sample collection
  - Samples should be human/mouse/rat serum and plasma

collected and separated according to the standard method.

- The passive drool method is recommended for the collection of saliva samples. (the passive drool method requires subjects to salivate into a vial or tube.)

It has been confirmed that samples collected using cotton swabs or inert polymer swabs can be used similarly to samples collected by the passive drool method.

Note : caution should be taken because some swab materials may absorb target analytes.

- Sample storage
  - For long-term storage, samples should be frozen at -35°C or below. Avoid repeated freezing and thawing.
  - Frozen samples should be thawed just before the assay and thoroughly mixed.
- Sample dilution
  - Use the buffer solution provided in the kit as the dilution solvent.
  - Use test tubes (PP, PE) to dilute the samples.
  - Dilute the samples just before use.
- Other issues
  - When turbidity or insoluble materials are present in the samples, remove them by centrifugation or other methods before the assay.
  - If it is suspected that the sample contains an interfering substance, perform a dilution linearity check with two or more different dilutions of the same sample.
  - When using a serum-separating agent, check for compatibility in advance.
  - Do not use hemolyzed or high-lipid samples.

**[8. Reagent Preparation]**

Reagents in the kit must be brought to room temperature (20-25°C) before use (about 2 hours).

8-1. Preparation of Standard Solution

Before starting the procedure, Prepare a caffeine standard stock solution (15 μg/mL) with the buffer solution supplied with the kit pre-equilibrated to room temperature.

The following is an example.

Concentration (ng/mL)	Volume of Standard Solution	Buffer
1,000	Standard stock solution : 30 μL	420 μL
250	1,000 ng/mL Solution : 50 μL	150 μL
62.5	250 ng/mL Solution : 50 μL	150 μL
15.6	62.5 ng/mL Solution : 50 μL	150 μL
3.91	15.6 ng/mL Solution : 50 μL	150 μL
0.977	3.91 ng/mL Solution : 50 μL	150 μL
0.244	0.977 ng/mL Solution : 50 μL	150 μL
0.00 (B <sub>0</sub> )	—	150 μL

8-2. Peroxidase-conjugated Caffeine Solution

Dilute 10,000-fold with the buffer solution that has been equilibrated to room temperature.

Note : A two-step dilution is recommended.

### 8-3. Wash Solution (10×)

Before starting the procedure, dilute 10-fold with purified water (distilled water) that has been equilibrated to room temperature. Example : 100 mL of Wash Solution (10×) + 900 mL of purified water (distilled water) (when all 96 wells are used)

• Other reagents are ready to use.

### [9. Assay Procedure]

• Prepare the reagents to be dispensed in the next step before starting each step.

• For dispensing reagents in steps (2), (9), and (12), continuous dispensing using a multichannel pipette is recommended.

(1) Fill each well of the antibody-coated plate with the diluted Wash Solution and remove. Repeat four times.

Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.

Note : Steps (2) through (5) below should be completed within 30 minutes.

(2) Dispense 50  $\mu\text{L}$  of diluted Peroxidase-conjugated Caffeine Solution into each well.

(3) Gently agitate the plate using a microplate shaker.

Shaking conditions : 600-800 rpm for 10 seconds. Repeat 3 times.

After each cycle of shaking, stop and then restart again.

(4) Dispense 50  $\mu\text{L}$  each of diluted standard solutions into the wells designated for standards.

(5) Dispense 50  $\mu\text{L}$  each of diluted samples into the wells designated for samples.

(6) Agitate the plate using a microplate shaker.

Shaking conditions : 600-800 rpm for 10 seconds. Repeat 3 times.

After each cycle of shaking, stop and then restart again.

(7) Cover with the plate seal and incubate at room temperature (20-25°C) for 2 hours without shaking.

Note : Peel the protective paper off the plate seal and seal with the adhesive side facing the plate.

Do not reuse a plate seal that has already been used.

(8) After the reaction is complete, remove the reaction mixture and fill each well with Wash Solution and wash 4 times.

Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.

(9) Dispense 100  $\mu\text{L}$  of TMB solution pre-equilibrated to room temperature into each well.

(10) Agitate the plate using a microplate shaker.

Shaking conditions : 600-800 rpm for 10 seconds. Repeat 3 times.

After each cycle of shaking, stop and then restart again.

(11) Cover with the plate seal and incubate at room temperature (20-25°C) for 20 minutes.

Note : Peel the protective paper off the plate seal and seal with the adhesive side facing the plate.

Do not reuse a plate seal that has already been used.

(12) Dispense 100  $\mu\text{L}$  of Stop Solution pre-equilibrated to room

temperature into each well.

(13) After shaking, measure the absorbance at 450 nm (secondary wavelength 620 nm) with a plate reader.

[Supplementary Notes]

### Washing Procedure

For each wash, dispense the Wash Solution into the wells, then agitate gently by hand for 10 seconds and discard.

After 4 consecutive washes, place the plate upside down on a paper towel and tap to completely remove the Wash Solution.

After removing the Wash Solution, immediately dispense the next solution to prevent the wells from drying out. The recommended volume of Wash Solution to be dispensed is 300  $\mu\text{L}$ /well.

If foam remains during washing, when manually washing using a wash bottle, allow the Wash Solution to overflow during the final cycle of washing to prevent the formation of foam.

### [10. Calculation Method]

(1) Create a standard curve by plotting the concentration of the standard solution (ng/mL) on the X-axis against the  $B/B_0$  (%) value on the Y-axis.

(2) Read the concentration (ng/mL) corresponding to the  $B/B_0$  (%) of the diluted sample.


(3) Multiply the read concentration by the sample dilution factor to obtain the measured concentration.

[Note]

For calculations using computer software, we recommend the use of a cubic polynomial, and 4 or 5 parameters.

$B/B_0$  (%) = (Absorbance of each standard solution or sample / Absorbance of the zero standard ( $B_0$ ))  $\times$  100

### [11. Typical standard Curve]

### [12. Assay Procedure Summary]

Ensure that you read the instruction manual and check the

sample conditions, assay conditions, and assay procedure before beginning the assay.

- Plates and reagents must be brought to room temperature (20 to 25°C) before use (about 2 hours).
- Dilution of Wash Solution (10×) : Dilute 10-fold with purified water pre-equilibrated to room temperature.
- Preparation of standard solutions (example) : Prepare a caffeine standard stock solution (15 μg/mL) with the buffer solution supplied with the kit pre-equilibrated to room temperature. The following is an example.

Dilution example	Concentration(ng/mL)	1,000	250	62.5	15.6	3.91	0.977	0.244	0
Standard Solution (μL)	Stock Solution :	30	50*	50*	50*	50*	50*	50*	—
Buffer (μL)		420	150	150	150	150	150	150	150

\* : Standard solution at 1 step higher concentration

- Preparation of Peroxidase-conjugated Caffeine Solution : Dilute 10,000-fold with Buffer pre-equilibrated to room temperature.
- Prepare samples according to [7. Sample Preparation].

**Antibody-coated Plate**

- ↓ Wash 4 times (\*1).

**Diluted Peroxidase-conjugated Caffeine Solution 50 μL/well**

- ↓ Agitate using a microplate shaker (\*2).

**Standard solutions or diluted samples 50 μL/well**

- ↓ Agitate using a microplate shaker (\*2), then incubate at room temperature (20-25°C) for 2 hours (\*3).

- ↓ Check TMB Solution is equilibrated to room temperature.

- ↓ Wash 4 times (\*1).

**TMB Solution 100 μL/well**

- ↓ Agitate using a microplate shaker (\*2), then incubate at room temperature (20-25°C) for 20 minutes (\*3).

**Stop Solution 100 μL/well**

- ↓ Agitate using a microplate shaker (\*2).

**Absorbance Measurement (main wavelength, 450 nm ; sub-wavelength, 620 nm ; 600 to 650 nm)**

(\*1) For each wash, dispense the Wash Solution into the wells, then agitate gently by hand for 10 seconds and discard. After 4 consecutive washes, place the plate upside down on a paper towel and tap to completely remove the Wash Solution. After removing the Wash Solution, immediately dispense the next solution to prevent the wells from drying out. The recommended volume of Wash Solution to be dispensed is 300 μL/well. If foam remains during washing, when manually washing using a wash bottle, allow the Wash Solution to overflow during the final cycle of washing to prevent the formation of foam.

(\*2) Shaking should be done at 600-800 rpm for 10 seconds. Repeat 3 times.

(\*3) After shaking is complete, cover with a plate seal and incubate without shaking. Peel the protective paper off the plate seal and seal with the adhesive side facing the plate. Do not reuse a plate seal that has already been used.

**[Storage]**

Store at 2°C to 10°C

**[Expiration date]**

Indicated on the label

**[Package size]**

For 96 assays

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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