## LBIS™ Human IL-7 ELISA Kit

Please, read this instruction carefully before use.

### 1. Intended use

Interleukin 7 (IL-7) is a glycoprotein cytokine with a molecular weight of approximately 25000 Da consisting of 152 amino acid residues and is an essential cytokine for the adaptive immune system. It is mainly produced by non-hematopoietic cells, including thymocytes, bone marrow and lymphoid tissue stromal cells, keratinocytes, and intestinal epithelial cells.

IL-7 is considered a hematopoietic cytokine that contributes to host defense by its critical and nonredundant role in inducing differentiation and homeostasis of immune cells, including NK cells, B cells, and T cells. Therefore, it is a promising candidate for an immunotherapeutic agent against chronic infectious diseases. It is also associated with rheumatoid arthritis, chronic colitis, and asthma. In addition, the substantial immune reconstitution effect has been demonstrated. Its function in promoting antiviral activity and as a booster factor for anticancer drugs and immune checkpoint inhibitors has been revealed ; thus, the high-sensitivity IL-7 measurement is expected to be helpful in immunotherapy.

This kit is a sandwich ELISA system for quantitative measurement of human IL-7 and can precisely measure IL-7 in human serum (plasma) and with high sensitivity.

Research use only. Not for use in diagnostics procedures.

### 2. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBISTM Human IL-7 ELISA Kit, standards or samples are incubated in polyclonal antibody-coated wells to capture IL-7. After 1 hour incubation and washing, biotinylated anti-IL-7 antibody is added and incubated further for 1 hour to bind with captured IL-7. After washing, HRP (horse radish peroxidase)-conjugated streptavidin is added and incubated for 30 minutes. After washing, bound HRP-conjugated streptavidin is reacted with a TMB Solution for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to IL-7 concentration. The standard curve is prepared by plotting absorbance against standard IL-7 concentrations. IL-7 concentrations in unknown samples are determined using this standard curve.

### 3. Performance Characteristics

• Precision of assay (Within assay variation)

2 serum samples, 5 replicates assay

n/ID	Sample1 (pg/mL)	Sample2 (pg/mL)
1	2.00	30.4
2	2.13	31.6
3	2.21	31.5
4	2.14	32.1
5	2.14	31.3
mean	2.13	31.4
SD	0.075	0.64
CV (%)	3.5	2.0

• Reproducibility (Between assay variation)

3 serum samples, 4 days, 3 replicates assay

Day/ID	Sample1 (pg/mL)	Sample2 (pg/mL)	Sample3 (pg/mL)
0	124	7.72	4.08
1	126	7.93	3.99
2	129	7.86	4.18
3	121	7.64	3.88
mean	125	7.79	4.03
SD	3.3	0.13	0.13
CV (%)	2.6	1.7	3.2

- Recovery test (Standard human IL-7 was added in 4 concentrations to serum or plasma sample and assayed.) The recoveries were 92.0% - 105%.


Serum sample

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
0	8.37	–	–
12.5	19.9	11.5	92.0
25.0	31.9	23.5	94.0
50.0	57.5	49.1	98.2
100	113	105	105

Plasma sample (EDTA)

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
0	8.33	–	–
12.5	19.9	11.6	92.8
25.0	32.4	24.1	96.4
50.0	57.7	49.4	98.8
100	109	101	101

- Dilution test (Serum sample and plasma sample were serially diluted by 3 steps.) The dilution curves showed linearity with  $R^2=0.998 - 1$ .

#### 4. Technical Tips / Precautions

- ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20°C - 25°C. Avoid performing assay under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- As the Antibody-coated Plate is module type of 8 wells × 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in a place where this kit is used.
- The materials must not be pipetted by mouth.
- Do not use reagents with different lot numbers together.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.

- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

## 5. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells/1 plate
(B) Human IL-7 Standard	Freeze-dried. Use after reconstitution.	1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody	Freeze-dried. Use after reconstitution.	1 bottle
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	150 $\mu$ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 $\times$ )	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	–	4 sheets

### **[Storage and Stability]**

#### [(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip lock originally used for well-plate container and store at 2°C - 10°C.

#### [(B) Human IL-7 Standard]

The reconstituted standard solution (original standard solution) should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted standard solutions after use.

#### [(C) Buffer] & [(F) TMB Solution]

Use only volume you need for your assay. Remaining reagents should be stored at 2°C - 10°C closing the cap tightly. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

#### [(D) Biotin-conjugated Antibody]

The reconstituted (D) solution should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted working solution.

#### [(E) Peroxidase-conjugated Streptavidin Solution]

Remaining working solution (already diluted) should be disposed. The rest of the undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

#### [(H) Stop Solution]

Close the cap tightly and store at 2°C - 10°C.

#### [(I) Wash Solution (10 $\times$ )]

The rest of undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Dispose any remaining diluted buffer.

## 6. Equipments or supplies required but not supplied Use as a check box

Deionized water (or Distilled water)  Test tubes for preparation of standard solution series.  Glassware for dilution of Wash Solution (10 $\times$ ) (a graduated cylinder, a bottle)  Pipettes (disposable tip type). One should be able to deliver 10  $\mu$ L - 100  $\mu$ L precisely, and another for 100  $\mu$ L - 1000  $\mu$ L.  Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 100  $\mu$ L.  Paper towel to remove washing buffer remaining in wells.  A vortex-type mixer.  A shaker for 96 well-plate (600 rpm - 1200 rpm)  An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle.  A 96 well-plate reader (450 nm  $\pm$  10 nm, 620 nm : 600 nm - 650 nm)  Software for data analysis.

## 7. Preparation of Samples

- This kit is intended to measure IL-7 in human serum or plasma.
- We recommend using EDTA-2Na, EDTA-2K or heparin as anticoagulant.
- Samples should be immediately assayed or stored below -35°C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- Dilution of a sample should be made in a PP or PE test tube using (C) Buffer prior to adding them to wells.

## 8. Preparation of Reagents

- ◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

### [Concentrated reagents]

#### [(B) Human IL-7 Standard]

Reconstitute the (B) Human IL-7 Standard with purified water in the amount specified in "Reconstitution of standard"\* to prepare a original standard solution (2360 pg/mL).


Then prepare using the kit-attached (C) Buffer that has been allowed to warm up to room temperature.

\*Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer	Concentration
Original standard solution 50 $\mu$ L	450 $\mu$ L	236 pg/mL
236 pg/mL solution 200 $\mu$ L	350 $\mu$ L	85.8 pg/mL
85.8 pg/mL solution 200 $\mu$ L	350 $\mu$ L	31.2 pg/mL
31.2 pg/mL solution 200 $\mu$ L	350 $\mu$ L	11.3 pg/mL
11.3 pg/mL solution 200 $\mu$ L	350 $\mu$ L	4.13 pg/mL
4.13 pg/mL solution 200 $\mu$ L	350 $\mu$ L	1.50 pg/mL
Blank	350 $\mu$ L	0

} Standard Curve

#### [(D) Biotin-conjugated Antibody]

Add 150  $\mu$ L of deionized water (or distilled water) of room temperature to (D) Biotin-conjugated Antibody (Freeze dried) (Original (D) solution). Prepare working solution by dilution of original (D) solution with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

#### [(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

#### [(I) Wash Solution (10 $\times$ )]

Dilute 1 volume of the concentrated (I) Wash Solution (10 $\times$ ) to **10 volumes** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated (I) Wash Solution (10 $\times$ ) and 900 mL of deionized water (or distilled water).


## 9. Assay Procedure

Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the wells with washing buffer and discard 4 times (\*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells (\*⑤).
  - (2) Pipette **100  $\mu$ L** of standard solution to the wells designated for standards.
  - (3) Pipette **50  $\mu$ L** of (C) Buffer to the sample wells, and then add **50  $\mu$ L** of sample to the wells.
  - (4) Shake the plate gently on a plate shaker (\*②).
  - (5) Stick a Plate Seal (\*③) on the plate and incubate for 1 hour at 20°C - 25°C.
  - (6) Discard the reaction mixture and rinse wells as step (1).
  - (7) Pipette **100  $\mu$ L** of Biotin-conjugated antibody solution to all wells and shake as step (4).
  - (8) Stick a Plate Seal (\*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
  - (9) Discard the reaction mixture and rinse wells as step (1).
  - (10) Pipette **100  $\mu$ L** of Peroxidase-conjugated Streptavidin Solution to all wells and shake as step (4).
  - (11) Stick a Plate Seal (\*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
  - (12) Discard the reaction mixture and rinse wells as step (1).
  - (13) Pipette **100  $\mu$ L** of (F) TMB Solution to wells and shake as step (4).
  - (14) Stick a Plate Seal (\*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
  - (15) Add **100  $\mu$ L** of the (H) Stop Solution to all wells and shake as step (4).
  - (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (\*④)) immediately using a plate reader.
- \*Refer to 13. Summary of Assay Procedure for notes on \*①, \*②, \*③, \*④, and \*⑤.

## 10. Calculations

- (1) Create a standard curve by reducing the data using computer software capable of generating 3rd order regression curve, or 4 or 5 parameters. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IL-7 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.
- (2) Using the standard curve, read the IL-7 concentration of samples at its absorbance, and multiply the assay value by dilution factor if the sample has been diluted. The standard operation method is 5-fold dilution.

## 11. Calibration

The measured value obtained by this kit can be converted to the unit concentration based on NIBSC/WHO standard IL-7 (Code : 90/530) by multiplying by the conversion factor specified in a separate sheet according to the following calculation formula.

NIBSC/WHO (90/530) unit concentration (U/mL) = conversion factor (\* 1)  $\times$  measured IL-7 value (pg/mL) (\* 2)

(\* 1) The exchange rate is different depending on the lot. Please see the specified conversion factor in the annex document.

(\* 2) Use the assay value that multiplied by dilution rate.

## 12. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
  - 2) Reagents necessary for coloration such as Biotin-conjugated Antibody, Peroxidase-conjugated Streptavidin Solution, or TMB Solution might not be added.
  - 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody or Peroxidase-conjugated Streptavidin Solution.
  - 4) Contamination of enzyme inhibitor(s).
  - 5) Influence of the temperature under which the kits had been stored.
  - 6) Excessive hard washing of the well plate.
  - 7) Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.
- Blank OD was higher than that of the lowest standard concentration (1.50 pg/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found 96 well-plate is empty when I opened the box.

A-2 : As this kit is dried type, not preservation stabilizer is added.

## 13. Summary of Assay Procedure : Use as a check box

\*First, read this instruction manual carefully and start your assay after confirmation of details.

Bring the well-plate and all reagents back to **20°C - 25°C for 2 hours**.

( I ) Wash Solution (10×) must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C -25°C.

Standard solution dilution example

Reconstitute the (B) Human IL-7 Standard with purified water in the amount specified in "Reconstitution of standard"\* to prepare a original standard solution (2360 pg/mL).

Then prepare using the kit-attached (C) Buffer that has been allowed to warm up to room temperature.

\*Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Preparation of diluted IL-7 standard solutions :

Conc. (pg/mL)	236	85.8	31.2	11.3	4.13	1.50	0
Std. sol. (μL)	50	200*	200*	200*	200*	200*	—
Buffer (μL)	450	350	350	350	350	350	350

\*One rank higher standard.

- Antibody-coated Plate
- ↓ Washing 4 times (\*①), (\*⑤)
- Samples (diluted sample : Buffer **50 μL** + Sample **50 μL**), or Standards 100 μL
- ↓ Shaking (\*②), Incubation for 1 hour at 20°C - 25°C. (Standing (\*③))
- Add deionized water (or distilled water) of room temperature to (D) Biotin-conjugated Antibody (Original (D) solution). Prepare working solution by dilution of original (D) solution to **1 : 100** with the (C) Buffer returned to 20°C - 25°C. (This should be prepared during incubation.)
- ↓ Washing 4 times (\*①), (\*⑤)
- Biotin-conjugated Antibody solution (Diluted) 100 μL
- ↓ Shaking (\*②), Incubation for 1 hour at 20°C - 25°C. (Standing (\*③))
- Dilute Peroxidase-conjugated Streptavidin Solution to **1 : 100** with the (C) Buffer returned to 20°C - 25°C. (This should be prepared during incubation.)
- ↓ Washing 4 times (\*①), (\*⑤)
- Peroxidase-conjugated Streptavidin Solution (Diluted) 100 μL
- ↓ Shaking (\*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (\*③))
- ↓ Washing 4 times (\*①), (\*⑤)
- TMB Solution 100 μL
- After dispensing, the color turns to blue depending on the concentration.
- ↓ Shaking (\*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (\*③))
- Stop Solution 100 μL
- After dispensing, the color turns to yellow depending on the concentration.
- ↓ Shaking (\*②), Immediately shake.
- Measurement of absorbance (450 nm, Ref 620 nm (\*④)) immediately.  
Ref. wave cancels the dirt in the back of plate.

\*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 sec. and remove the buffer. Guideline of washing volume : 300 μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipettes, sometimes the background tends to be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.

Standard of plate-washing pressure : 5 mL/min. - 25 mL/min. (Adjust it depending on the nozzle's diameter.)


\*② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.

\*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

\*④ 600 nm - 650 nm can be used as reference wavelength.

\*⑤ After removal of wash buffer, immediately dispense the next reagent.

## Assay Worksheet

### 14. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

#### LBIS™ Human IL-7 ELISA Kit

[Storage] Store the kit at 2°C - 10°C (Do not freeze)  
[Expiration date] Indicated on the container  
[Package] For 96 tests  
[Cat #] 290-87501

---

### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

#### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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