LBIS™ Mouse IL-17A ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

Mouse Interleukin 17A (IL-17A) is secreted as a disulfide-linked homodimer of 30 - 35 kDa glycoproteins consisting of subunits with 133 amino acid residues and functions as a pro-inflammatory cytokine. The IL-17 family consists of six structurally related isoforms (IL-17A to F), but IL-17A is generally known as IL-17. Recently, it has been shown that IL-17 is produced by Th17 cells, a subset distinct from other CD4+ T cells such as Th1 and Th2 cells. IL-17A efficiently induces the expression of inflammatory mediators in various cells, including macrophages, stromal cells, epithelial cells, endothelial cells, and fibroblasts.

IL-17A has been reported to play essential roles in autoimmune diseases such as rheumatoid arthritis (RA), inflammatory diseases, allergic reactions, and protection against bacterial infection. It is attracting attention as a new therapeutic target.

This kit is a sandwich ELISA system for quantitative measurement of mouse IL-17A. It can precisely measure IL-17A in mouse serum (plasma) and with high sensitivity.

Research use only. Not for use in diagnostics procedures.

2. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse IL-17A ELISA Kit, standards or samples are incubated in polyclonal antibody-coated wells to capture IL-17A. After 2 hours' incubation and washing, biotinylated anti-IL-17A antibody is added and incubated further for 1 hour to bind with captured IL-17A. After washing, HRP (horse radish peroxidase)-conjugated streptavidin is added and incubated for 30 minutes. After washing, bound HRP-conjugated streptavidin is reacted with a TMB Solution for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to IL-17A concentration. The standard curve is prepared by plotting absorbance against standard IL-17A concentrations. IL-17A concentrations in unknown samples are determined using this standard curve.

3. Performance Characteristics

- Precision of assay (Within assay variation)

2 serum samples, 5 replicates assay

n/ID	Sample1 (pg/mL)	Sample2 (pg/mL)
1	261	18.6
2	248	17.8
3	243	18.4
4	264	18.2
5	274	19.2
mean	258	18.4
SD	12	0.53
CV (%)	4.8	2.9

- Reproducibility (Between assay variation)

3 serum samples, 4 days, 3 replicates assay

Day/ID	Sample1 (pg/mL)	Sample2 (pg/mL)	Sample3 (pg/mL)
0	281	40.6	5.92
1	309	39.4	5.65
2	308	41.2	6.11
3	297	41.8	5.98
mean	299	40.7	5.92
SD	13	1.0	0.20
CV (%)	4.4	2.4	3.3

- Recovery test (Standard mouse IL-17A was added in 3 concentrations to serum or plasma sample and assayed.) The recoveries were 95.4% - 105%.

Serum sample

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
–	1.41	–	–
10	11.2	9.79	97.9
50	49.1	47.7	95.4
250	264	263	105

Plasma sample (EDTA)

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
–	0.328	–	–
10	9.99	9.66	96.6
50	50.6	50.3	101
250	257	257	103

- Dilution test (Serum sample and plasma sample were serially diluted by 3 steps.) The dilution curves showed linearity with $R^2=0.9998$.

4. Technical Tips / Precautions

- ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20°C - 25°C. Avoid performing assay under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- As the Antibody-coated Plate is module type of 8 wells × 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in a place where this kit is used.
- The materials must not be pipetted by mouth.
- Do not use reagents with different lot numbers together.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

5. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells/1 plate
(B) Mouse IL-17A Standard	Freeze-dried. Use after reconstitution.	1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody	Freeze-dried. Use after reconstitution.	1 bottle
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	100 μ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 \times)	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	—	4 sheets

[Storage and Stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip lock originally used for well-plate container and store at 2°C - 10°C.

[(B) Mouse IL-17A Standard]

The reconstituted standard solution (original standard solution) should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted standard solutions after use.

[(C) Buffer] & [(F) TMB Solution]

Use only volume you need for your assay. Remaining reagents should be stored at 2°C - 10°C closing the cap tightly. Once

opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(D) Biotin-conjugated Antibody]

The reconstituted (D) solution should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted working solution.

[(E) Peroxidase-conjugated Streptavidin Solution]

Remaining working solution (already diluted) should be disposed. The rest of the undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(H) Stop Solution]

Close the cap tightly and store at 2°C - 10°C.

[(I) Wash Solution (10×)]

The rest of undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Dispose any remaining diluted buffer.

6. Equipments or supplies required but not supplied Use as a check box

Deionized water (or Distilled water) Test tubes for preparation of standard solution series. Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle) Pipettes (disposable tip type). One should be able to deliver 10 μL - 100 μL precisely, and another for 100 μL - 1000 μL. Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 100 μL. Paper towel to remove washing buffer remaining in wells. A vortex-type mixer. A shaker for 96 well-plate (600 rpm - 1200 rpm) An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. A 96 well-plate reader (450 nm ± 10 nm, 620 nm : 600 nm - 650 nm) Software for data analysis.

7. Preparation of Samples

- This kit is intended to measure IL-17A in mouse serum or plasma.
- We recommend using EDTA-2Na, EDTA-2K or heparin as anticoagulant.
- Samples should be immediately assayed or stored below -35°C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- Dilution of a sample should be made in a PP or PE test tube using buffer prior to adding them to wells.

8. Preparation of Reagents

- ◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Mouse IL-17A Standard]

Reconstitute the (B) Human IL-17A Standard with purified water in the amount specified in "Reconstitution of standard"* to prepare a original standard solution (8 ng/mL).


Then prepare using the kit-attached (C) Buffer that has been allowed to warm up to room temperature.

*Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer	Concentration
Original standard solution 50 μL	450 μL	800 pg/mL
800 pg/mL solution 200 μL	340 μL	296 pg/mL
296 pg/mL solution 200 μL	340 μL	110 pg/mL
110 pg/mL solution 200 μL	340 μL	40.6 pg/mL
40.6 pg/mL solution 200 μL	340 μL	15.1 pg/mL
15.1 pg/mL solution 200 μL	340 μL	5.58 pg/mL
5.58 pg/mL solution 200 μL	340 μL	2.06 pg/mL
Blank	340 μL	0

} Standard Curve

[(D) Biotin-conjugated Antibody]

Add 100 µL of deionized water (or distilled water) of room temperature to (D) Biotin-conjugated Antibody (Freeze-dried) (Original (D) solution). Prepare working solution by dilution of original (D) solution with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

[(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

[(I) Wash Solution (10×)]

Dilute 1 volume of the concentrated (I) Wash Solution (10×) to **10 volumes** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated (I) Wash Solution (10×) and 900 mL of deionized water (or distilled water).

9. Assay Procedure


Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the wells with washing buffer and discard 4 times (*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells (*⑤).
- (2) Pipette **50 µL** of standard solution to the wells designated for standards.
- (3) Pipette **25 µL** of (C) Buffer to the sample wells, and then add **25 µL** of sample to the wells.
- (4) Shake the plate gently on a plate shaker (*②).
- (5) Stick a Plate Seal (*③) on the plate and incubate for 2 hours at 20°C - 25°C.
- (6) Discard the reaction mixture and rinse wells as step (1).
- (7) Pipette **50 µL** of Biotin-conjugated Antibody solution to all wells and shake as step (4).
- (8) Stick a Plate Seal (*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
- (9) Discard the reaction mixture and rinse wells as step (1).
- (10) Pipette **50 µL** of Peroxidase-conjugated Streptavidin Solution to all wells and shake as step (4).
- (11) Stick a Plate Seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
- (12) Discard the reaction mixture and rinse wells as step (1).
- (13) Pipette **50 µL** of (F) TMB Solution to wells and shake as step (4).
- (14) Stick a Plate Seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
- (15) Add **50 µL** of the (H) Stop Solution to all wells and shake as step (4).
- (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (*④)) immediately using a plate reader.

*Refer to 12. Summary of Assay Procedure for notes on *①, *②, *③, *④, and *⑤.

10. Calculations

- (1) Create a standard curve by reducing the data using computer software capable of generating 3rd order regression curve, or 4 or 5 parameters. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse IL-17A concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.
- (2) Using the standard curve, read the IL-17A concentration of samples at its absorbance, and multiply the assay value by dilution factor if the sample has been diluted. The standard operation method is 5-fold dilution.

11. Trouble Shooting

- Low absorbance in all wells
Possible explanations :
 - 1) The standard or samples might not be added.
 - 2) Reagents necessary for coloration such as Biotin-conjugated Antibody, Peroxidase-conjugated Streptavidin Solution, or TMB Solution might not be added.
 - 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody or Peroxidase-conjugated Streptavidin Solution.
 - 4) Contamination of enzyme inhibitor(s).
 - 5) Influence of the temperature under which the kits had been stored.
 - 6) Excessive hard washing of the well plate.
 - 7) Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.
- Blank OD was higher than that of the lowest standard concentration (2.06 pg/mL).
Possible explanations :
Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.)
 - High coefficient of variation (CV)
Possible explanation :
 - 1) Improper or inadequate washing.
 - 2) Improper mixing of standard or samples.
 - 3) Pipetting at irregular intervals.
- Q-1 : Can I divide the plate to use it for the other testing?
A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.
- Q-2 : I found 96 well-plate is empty when I opened the box.
A-2 : As this kit is dried type, not preservation stabilizer is added.

12. Summary of Assay Procedure : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

Bring the well-plate and all reagents back to **20°C - 25°C for 2 hours**.

(I) Wash Solution (10×) must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.

Standard solution dilution example

Reconstitute the (B) Human IL-17A Standard with purified water in the amount specified in "Reconstitution of standard"* to prepare a original standard solution (8 ng/mL).

Then prepare using the kit-attached (C) Buffer that has been allowed to warm up to room temperature.

*Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.


Preparation of diluted IL-17A standard solutions :

Conc. (pg/mL)	800	296	110	40.6	15.1	5.58	2.06	0
Std. sol. (μL)	50	200*	200*	200*	200*	200*	200	—
Buffer (μL)	450	340	340	340	340	340	340	340

*One rank higher standard.

- Antibody-coated Plate
- ↓ Washing 4 times (*①), (*⑤)
- Samples (diluted sample : Buffer **25 μL** + Sample **25 μL**), or Standards 50 μL
- ↓ Shaking (*②), Incubation for 2 hours at 20°C - 25°C. (Standing (*③))
- Add deionized water (or distilled water) of room temperature to (D) Biotin-conjugated Antibody (Freeze-dried)
- (Original (D) solution). Prepare working solution by dilution of original (D) solution to **1 : 100** with the (C) Buffer returned to 20°C - 25°C. (This should be prepared during incubation.)
- ↓ Washing 4 times (*①), (*⑤)
- Biotin-conjugated Antibody solution (Diluted) 50 μL
- ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- Dilute HRP-conjugated streptavidin to **1 : 100** with the (C) Buffer returned to 20°C - 25°C.
- (This should be prepared during incubation.)
- ↓ Washing 4 times (*①), (*⑤)
- Peroxidase-conjugated Streptavidin Solution (Diluted) 50 μL
- ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- ↓ Washing 4 times (*①), (*⑤)
- TMB Solution 50 μL
- After dispensing, the color turns to blue depending on the concentration.
- ↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))
- Stop Solution 50 μL
- After dispensing, the color turns to yellow depending on the concentration.
- ↓ Shaking (*②), Immediately shake.
- Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately.
- Ref. wave cancels the dirt in the back of plate.
- *① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 sec. and remove the buffer. Guideline of washing volume : 300 μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipettes, sometimes the background tends to be high. If so, change washing frequency from 4 times to 5 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.
- Standard of plate-washing pressure : 5 mL/min. - 25 mL/min. (Adjust it depending on the nozzle's diameter.)
- *② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.
- *③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.
- *④ 600 nm - 650 nm can be used as reference wavelength.
- *⑤ After removal of wash buffer, immediately dispense the next reagent.

Assay Worksheet

13. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

LBIS™ Mouse IL-17A ELISA Kit

[Storage] Store the kit at 2°C - 10°C (Do not freeze)

[Expiration date] Indicated on the container

[Package] For 96 tests

[Cat #] 293-94701

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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