## LBIS™ Human Apo B-48 ELISA Kit

Please, read this instruction carefully before use.

### 1. Intended use

LBIS™ Human Apo B-48 ELISA Kit is a sandwich ELISA system for quantitative measurement of human Apolipoprotein B-48 (Apo B-48). This is intended for research use only.

### 2. Introduction

Plasma lipoproteins are complex substances composed of certain ratios of lipids and proteins, and are important factors for serum lipids transportation. On the surface of the lipoprotein structure, apolipoproteins are present and playing roles in stabilization of lipoprotein structure, in activation of enzymes relating to lipoprotein metabolism, and in binding with lipoprotein receptors on the cell surface. Apolipoprotein-B48 (ApoB-48) has 48% amino acid sequence of apolipoprotein B-100 which is present in lipoproteins, VLDL, LDL and HDL of liver origin. ApoB-48 is an apolipoprotein specific to a lipoprotein, chylomicron (CM) which is formed in intestine and carries exogenous lipids derived from foods to the liver and peripheral tissues. Measurement of ApoB-48 is useful in observation and pursuit of postprandial dynamics of lipoprotein, lipid-soluble nutrients and drugs. Co-estimation of LDL-, HDL-cholesterol and B-48 will give important information about exogenous and endogenous cholesterol (see reference), and also estimation of CM-remnant that is thought to be a risk factor of cardiovascular atherosclerosis.

Reference : Otokozawa, S., *et al.* : *Atherosclerosis*, **205**, 197-201 (2009).

### 3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Human Apo B-48 ELISA Kit, standards or samples are incubated in monoclonal antibody-coated wells to capture Apo B-48. After 1 hour's incubation and washing, biotinylated anti-Apo B-48 antibody is added and incubated further for 1 hour to bind with captured Apo B-48. After washing, Peroxidase-conjugated Streptavidin Solution is added, and incubated for 30 minutes. After washing, bound HRP-conjugated streptavidin is reacted with a TMB Solution for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to Apo B-48 concentration. The standard curve is prepared by plotting absorbance against standard Apo B-48 concentrations. Apo B-48 concentrations in unknown samples are determined using this standard curve.

### 4. Performance Characteristics

- Assay range : The assay range of the kit is 2.5 ng/mL - 160 ng/mL.
- Specificity : The antibodies used in this kit are specific to Apo B-48.  
Human Apo B-100 had no cross reactivity with this kit.
- Precision of assay  
Within assay variation (3 samples, 5 replicates assay), Mean CV was within 10%.
- Reproducibility  
Between assay variation (3 samples, 3 days, 3 replicates assay), Mean CV was within 10%.
- Recovery test  
Standard Apo B-48 was added in 4 concentrations to 2 serum samples and were assayed.  
The recoveries were 94.9% - 101%.
- Dilution test  
Serum sample was serially diluted by 4 steps.  
The dilution curves showed linearity with  $R^2 = 0.996 - 0.999$ .

### 5. Reference Assay Data

Human apo B-48 assay data

Mean assay value : 4.60  $\mu\text{g/mL}$ , SD : 1.54  $\mu\text{g/mL}$ , n = 18, fasting

These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for Apo B-48 levels independently.

### 6. Precautions

- For professional use only. Beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.

- Do not drink, eat or smoke in the areas where assays are carried out.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the stop solution. The stop solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Avoid contact with the acidic stop solution and TMB Solution containing hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- The materials must not be pipetted by mouth.
- Unused samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- Use clean laboratory glassware.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

## 7. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells/1 plate
(B) Human Apo B-48 Standard	Frozen-dried. Use after dilution.	1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody Solution	Concentrated. Use after dilution.	100 $\mu$ L/1 bottle
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	100 $\mu$ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 $\times$ )	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	–	4 sheets

### **[Storage and Stability]**

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Human Apo B-48 Standard]

**The reconstituted standard solution (original solution) should be used within 24 hours if stored in a refrigerator. Dispose remaining prepared solution. For longer storage, freeze it immediately after reconstitution, and store it under -35°C.**

[(C) Buffer] & [(F) TMB Solution]

Use only volume you need for your assay. Remaining reagents should be stored at 2°C - 10°C fastening the cap tightly. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(D) Biotin-conjugated Antibody Solution] & [(E) Peroxidase-conjugated Streptavidin Solution]

Unused working solution (already diluted) should be disposed. If not opened, store at 2°C - 10°C.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10 $\times$ )]

If not opened, store at 2°C - 10°C. Dispose any unused diluted buffer.

## 8. Equipments or supplies required but not supplied Use as a check box

Deionized water (or Distilled water)  Test tubes for preparation of standard solution series.  Glassware for dilution of Wash Solution (10 $\times$ ) (a graduated cylinder, a bottle)  Pipettes (disposable tip type). One should be able to deliver 50  $\mu$ L precisely, and another for 200  $\mu$ L and 400  $\mu$ L.  Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 50  $\mu$ L.  Paper towel to remove washing buffer remaining in wells.  A vortex-type mixer.  A shaker for 96 well-plate (600 - 1200 rpm)  An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle.  A 96 well-plate reader (450 nm  $\pm$  10 nm, 620 nm : 600 nm - 650 nm)  Software for data analysis.

## 9. Preparation of Samples

This kit is intended to measure Apo B-48 in human serum or plasma (**citric acid for plasma is not suitable**). We recommend to use EDTA-2Na, EDTA-2K or heparin as anticoagulant. Don't use citric acid. Samples should be immediately assayed or stored below -35°C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

If presence of interfering substance is suspected, examine by dilution test at more than 2 points.

Dilution of a sample should be made in a PP or PE test tube using buffer solution prior to adding them to wells (Be careful with pipetting since samples diluted with buffer have highly viscosity). Mix well and pipette 50  $\mu\text{L}$  of diluted sample into a well for assay.

Recommended dilution is 100 $\times$ . For high value samples, more than 200 $\times$  dilution may be needed. When diluting samples, vortex them well (1800 rpm - 2000 rpm 3 times for 10 seconds) with care not to make bubbling.

Hemolytic and hyperlipemic samples are not suitable.

\* To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 160 mg/dL with this kit.

[Storage and Stability]

Apo B-48 in samples will be inactivated if stored at 2 $^{\circ}\text{C}$  - 10 $^{\circ}\text{C}$ . If it is necessary to store samples in refrigerator (2 $^{\circ}\text{C}$  - 10 $^{\circ}\text{C}$ ), add aprotinin at final concentration of 100 - 500 KIU/mL (KIU : kallikrein inhibitor unit). If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35 $^{\circ}\text{C}$ . Defrosted samples should be mixed thoroughly for best results. Avoid repeated freeze-thaw cycles.

## 10. Preparation of Reagents

◆ Bring all reagents of the kit to room temperature (20 $^{\circ}\text{C}$  - 25 $^{\circ}\text{C}$ ) before use.

◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

### [Concentrated reagents]

[(B) Human Apo B-48 Standard]

**Add 400  $\mu\text{L}$  of deionized water (or distilled water) of room temperature to Standard human Apo B-48 (Frozen-dried) and let it stand still for 10 minutes (Original solution).** Avoid bubbling of the solution, and shake it sufficiently at 1800 rpm - 2000 rpm 3 times for 10 seconds. Dissolve the standard in vial completely, and let it be clear solution.

Caution for vortexing : Put a vial on Vortex and start mixing. Do not put vial on Vortex while moving. That caused bubbling.

**Reconstituted solution (original solution) is highly viscous. Before preparing standard solutions, get accustomed to treat highly viscous solution (for pipetting, we recommend co-washing procedure). In order to prepare the standard solution precisely, it is very important to pipette correctly the first 200  $\mu\text{L}$ .**

Sufficient shaking with Vortex (1800 rpm - 2000 rpm 3 times for 10 seconds) in every step of dilution is very important. Be careful not to make bubbles while mixing. Bubbling causes bad reproducibility.

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer	Concentration (ng/mL)
Original solution 200 $\mu\text{L}$	200 $\mu\text{L}$	160
160 ng/mL solution 200 $\mu\text{L}$	200 $\mu\text{L}$	80
80 ng/mL solution 200 $\mu\text{L}$	200 $\mu\text{L}$	40
40 ng/mL solution 200 $\mu\text{L}$	200 $\mu\text{L}$	20
20 ng/mL solution 200 $\mu\text{L}$	200 $\mu\text{L}$	10
10 ng/mL solution 200 $\mu\text{L}$	200 $\mu\text{L}$	5.0
5.0 ng/mL solution 200 $\mu\text{L}$	200 $\mu\text{L}$	2.5
0 (Blank)	200 $\mu\text{L}$	0

[(D) Biotin-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**.

10 mL of the diluted solution is enough for 96 wells.

[(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 100**.

10 mL of the diluted solution is enough for 96 wells.

[(I) Wash Solution (10 $\times$ )]

Dilute I volume of the concentrated Wash Solution (10 $\times$ ) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated Wash Solution (10 $\times$ ) and 900 mL of deionized water (or distilled water).

## 11. Assay Procedure

**Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.**

- (1) Wash the (A) Antibody-coated Plate by filling the wells with washing buffer and discard 4 times (\*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells (\*⑤).
- (2) Pipette 50  $\mu\text{L}$  of properly diluted samples to the designated sample wells.
- (3) Pipette 50  $\mu\text{L}$  of standard solution to the wells designated for standards.
- (4) Shake the plate gently on a plate shaker (\*②).

- (5) Stick a Plate Seal (\*③) on the plate and incubate for 1 hour at 20°C - 25°C.
- (6) Discard the reaction mixture and rinse wells as step (1).
- (7) Pipette 50  $\mu$ L of Biotin-conjugated Antibody Solution to all wells, and shake as step (4).
- (8) Stick a Plate Seal (\*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
- (9) Discard the reaction mixture and rinse wells as step (1).
- (10) Pipette 50  $\mu$ L of Peroxidase-conjugated Streptavidin Solution to all wells, and shake as step (4).
- (11) Stick a Plate Seal (\*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
- (12) Discard the reaction mixture and rinse wells as step (1).
- (13) Pipette 50  $\mu$ L of (F) TMB Solution to wells, and shake as step (4).
- (14) Stick a Plate Seal (\*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
- (15) Add 50  $\mu$ L of the (H) Stop Solution to all wells and shake as step (4).
- (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (\*④)) immediately using a plate reader immediately.


\* Refer to 16. Summary of Assay Procedure for notes of \*①, \*②, \*③, \*④ and \*⑤.

## 12. Technical Tips

- In manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of stop solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antibody-coated Plate is module type of 8 wells  $\times$  12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

## 13. Calculations

- (1) Prepare a standard curve by plotting standard concentration on X-axis and absorbance on Y-axis.
- (2) Using the standard curve, read the Apo B-48 concentration of a sample at its absorbance\*, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.
- (3) Higher reaction temperature than 25°C tends to produce high absorbance on the whole. Though it depends on assay instruments, if coloring proceeds to area where there's low reliability of absorbance and slope of standard curve decreases, the area should not be used for calculation. Repeat assay by setting reaction temperature between 20°C - 25°C.

- \* We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.
- \* Physiological or pathological situation of human should be judged comprehensively taking other examination results into consideration.
- \* Human Apo B-48 assay standard curve (an example). Absorbance may change due to assay environment.

## 14. Trouble Shooting

- Low absorbance in all wells  
Possible explanations :

- 1) The standard or samples might not be added.
  - 2) Reagents necessary for coloration such as Biotin-conjugated Antibody Solution, Peroxidase-conjugated Streptavidin Solution, or TMB Solution might not be added.
  - 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody Solution or Peroxidase-conjugated Streptavidin Solution.
  - 4) Contamination of enzyme inhibitor (s).
  - 5) Influence of the temperature under which the kits had been stored.
  - 6) Excessive hard washing of the well plate.
  - 7) Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.
- Blank OD was higher than that of the lowest standard concentration (2.5 ng/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
  - 2) Improper mixing of standard or samples.
  - 3) Pipetting at irregular intervals.
- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2 : When we manufacture 96 well-plate, we insert preservation stabilizer in wells.

## 15. References

1 : Determination of apolipoprotein B-48 in serum by a sandwich ELISA.

Kinoshita, M. *et al.* : *Clinica Chimica Acta*, **351**, 115-120 (2005).

2 : Effects of intensive atorvastatin and rosuvastatin treatment on apolipoprotein B-48 and remnant lipoprotein cholesterol levels.

Otokozawa, S. *et al.* : *Atherosclerosis*, **205**, 197-201 (2009).


## 16. Summary of Assay Procedure : Use as a check box

\* First, read this instruction manual carefully and start your assay after confirmation of details.

Bring the (A) Antibody-coated Plate and all reagents back to 20°C - 25°C for 2 hours.

Wash Solution (10×) must be diluted to 10 times by deionized water (or distilled water) that returned to 20°C - 25°C.

**Standard solution dilution example**

- Antibody-coated Plate
- ↓ Washing 4 times (\*①), (\*⑤)
- (Diluted) Samples, or Standards 50 μL
- ↓ Shaking (\*②), Incubation for 1 hour at 20°C - 25°C. (Standing (\*③))
- Dilute Biotin-conjugated Antibody Solution to **1 : 100** with (C) Buffer returned to 20°C - 25°C.  
(This should be prepared during incubation.)
- ↓ Washing 4 times (\*①), (\*⑤)
- Biotin-conjugated Antibody Solution (Diluted) 50 μL
- ↓ Shaking (\*②), Incubation for 1 hour at 20°C - 25°C. (Standing (\*③))
- Dilute Peroxidase-conjugated Streptavidin Solution to **1 : 100** with (C) Buffer returned to 20°C - 25°C.  
(This should be prepared during incubation.)
- ↓ Washing 4 times (\*①), (\*⑤)
- Peroxidase-conjugated Streptavidin Solution (Diluted) 50 μL
- ↓ Shaking (\*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (\*③))
- ↓ Washing 4 times (\*①), (\*⑤)
- TMB Solution (After dispense, the color turns to blue depending on the concentration.) 50 μL
- ↓ Shaking (\*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (\*③))
- Stop Solution (After dispense, the color turns to yellow depending on the concentration.) 50 μL
- ↓ Shaking (\*②) (Immediately shake.)
- Measurement of absorbance (450 nm, Ref 620 nm (\*④)) immediately.  
Ref. wave cancels the dirt in the back of plate.

\*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300 μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution. Standard of plate-washing pressure : 5 - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

\*② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.

\*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.


\*④ 600 nm - 650 nm can be used as reference wavelength.

\*⑤ After removal of wash buffer, immediately dispense the next reagent.

### Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	160 ng/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	80 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	40 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	20 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	10 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	5.0 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	2.5 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

## Assay Worksheet

### 17. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

LBIS™ Human Apo B-48 ELISA Kit

[Storage] Store at 2°C - 10°C (Do not freeze).

[Expiration date] Indicated on the label

[Package] For 96 tests

[Cat #] 298-88901

## Attention

### Very Important Points to Read before Starting Assay.

#### 1) Temperature

Bring all the reagent solutions up to room temperature (20°C - 25°C) before starting assay.

For this purpose, place and keep all the reagent solutions on a laboratory table for about two hours under the room temperature set at 22°C. Keep the reaction temperature between 20°C and 25°C.

#### 2) Reconstitution of the standard solution

The standard preparation attached to the kit is a frozen-dried material. Open the cap of the container, and add purified water as indicated (400 µL, room temperature). Fasten the cap, and let it stand for 10 minutes. Then vortex the container (1800 rpm - 2000 rpm 3 times for 10 seconds) under enough care not to form bubbles. Bubbling may cause denaturation of the protein. Confirm that the material is completely dissolved and solution is clear. Vortexing should be started from still-standing state. If you put the container on the rotating stirrer, it may cause bubbling.

#### 3) Samples

For collecting plasma samples, use heparin (final concentration : 10 µg/mL (1.2 U/mL) - 100 µg/mL), EDTA-2Na (final concentration : 1.0 mg/mL - 1.5 mg/mL), or EDTA-2K (final concentration : 1.1 mg/mL - 1.7 mg/mL). **Do not use citrate.** Use PP or PE tubes when samples are stored frozen. **(Do not use glass tubes.)**

#### 4) Preparation of standard solutions and dilution of samples

As the original reconstituted standard solution is highly viscous, start preparation of the standard solutions after getting used to measure highly viscous solution through enough training. Accurate measurement and delivery of the first 200 µL original standard solution is a very important step for preparation of the standard series. Co-washing method of pipetting written in the instruction paper of Eppendorf pipette is recommended. In every step of serial dilution, enough vortexing (1800 rpm - 2000 rpm 3 times for 10 seconds) is necessary. As the buffer used for dilution contains a surface active agent, enough attention must be paid to avoid bubble formation. Use a properly calibrated micropipette. Accurate pipetting is essential for excellent precision. Use exclusive (or matched) tips for your pipette. For the treatment of highly viscous solution, use tips with rather wide opening, and fix repeating of co-washing, pushing and releasing speed of pipetting, and angle of the pipette-holding. Refer to instruction paper of the pipette of your use for the exact operation and care. As a rule, preparation of the standard solutions, and preparation of assay samples must be done right before starting assay, and dilute the reagents (biotin-labelled anti-Apo B-48 antibody, peroxidase-conjugated avidin) during 2 hours' incubation for antigen-antibody reaction.

#### 5) Addition of standard solution, sample and reagent solution to wells

In pipetting above solutions, be careful not to cause bubble formation. Bubble formation may influence assay variation. For addition of reagent solutions, we recommend the use of a repeating dispenser like Eppendorf's Multipette Plus. In blowing out of the solution from the pipette, do not release the solution at the center of the well but at wall of the well, and then take out the pipette through "touch and go".

After addition of any reagent solution, shake the well plate on a plate shaker (800 rpm 3 times for 10 seconds). Adjust rotation speed to get the best mixing result without bubbling.

#### 6) Reaction period

The time counting should be started from the pipetting to the first well. But if the pipetting to a plate should take more than 15 minutes due to unexperienced situation, start time counting from the end of pipetting to the last well.

#### 7) Washing of the well plate

When you wash 96-well plate using an automatic plate washer, adjust the speed and pressure of pouring and sucking of washing buffer not to cause removal of solidified antibody due to excessive water stream, or increase the blank absorbance owing to incomplete washing, or causing big variation in coloring due to uneven washing. We recommend rather weak stream or pressure and more repeat. Never forget the maintenance of the apparatus with proper intervals. Be careful not to dry wells.

#### 8) Quality controls

We have LBIS™ Human Apo B-48 Control Set (Cat # 292-89009) as an option. This set contains high and low controls (250 µL each) and should be stored below -35°C. Expiration is 12 months from production and dilution is not necessary (already diluted). This will be helpful for your quality control in each assay.

---

### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : +81-6-6203-3741  
Facsimile : +81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

#### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bollwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : +1-804-271-7677  
Facsimile : +1-804-271-7791  
<http://www.wakousa.com>

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