LabAssay™ ATX

1. Intended use

This kit is a research reagent for measuring autotaxin (ATX).

ATX is a glycoprotein with a molecular weight of 125 kDa that was isolated as a cell migration promoting factor from the culture supernatant of human malignant melanoma cells. It is known that liver damage such as fibrosis causes metabolic inhibition of ATX, which leads to retention of ATX in the blood, resulting in an increase in its blood concentration.

2. Storage and expiration date

Store at 2°C - 10°C and do not freeze. The expiration date is indicated on the label on the outer box of the kit.

3. Kit Component Reagents

Components		Use Status	Amount
(1)	Reacting Solution	Solution	16 mL/1 bottle
(2)	Substrate	Freeze-dried	For 8 mL/1 bottle
(3)	Substrate Dissolving Solution	Solution	8 mL/1 bottle
(4)	ATX Standard	Freeze-dried	2 bottles
(5)	Buffer	Solution	5 mL/1 bottle
(6)	Stop Solution	Solution	16 mL/1 bottle

*For 150 tests

4. Principle of the method

Lysophosphatidylcholine is hydrolyzed by the lysophospholipase D activity of autotaxin (ATX) in the sample to form choline. The resultant choline is oxidized by the action of choline oxidase (COD) to form hydrogen peroxide. The resultant hydrogen peroxide oxidatively condenses n-ethyl-N-(3-sulfopropyl)-3-methylaniline (TOPS) and 4-aminoantipyrine quantitatively through the action of peroxidase (POD) to produce a blue-violet pigment. The reaction is stopped by addition of the Stop Solution and absorbance is measured to determine ATX activity in the sample.

5. Equipment or supplies required but not provided in the kit

- 96-well microplate (transparent type)
- Micropipette
- Microtube
- Pipette
- Incubator maintained at 37°C
- Plate mixer
- Microplate reader with 546 nm/700 nm wavelength filter

6. Preparation of reagents

- ① Reacting Solution : Ready to use. After opening the bottle, store at 2°C - 10°C and use within 4 weeks.
- ② Substrate Solution : Reconstitute 1 bottle of Substrate with 8 mL of Substrate Dissolving Solution to make the Substrate Solution. After preparation, store at 2°C - 10°C and use within 4 weeks.
- ③ Dilution series of ATX Standard : Reconstitute ATX Standard with purified water in the amount specified in "Reconstitution of standard"* to prepare the original standard solution (55 U/L). Then prepare a dilution series using the Buffer included in the kit that has been allowed to warm up to room temperature.

*Find and check "Reconstitution of standard" on this product page. As the amount of purified water to be added varies by lot, be sure to check it for every lot.

After preparing each concentration of the standard solution, use immediately and do not store.

Example of preparation of dilution series of standard solution

ATX activity (U/L)	55.0	27.5	13.8	6.88	3.44	1.72	0
ATX Standard (μL)	Original standard solution : 40	50*	50*	50*	50*	50*	-
Buffer (μL)	-	50	50	50	50	50	50

*One rank higher standard solution.

④ Stop Solution : Ready to use. After opening the bottle, store at 2°C - 10°C .

*If crystals have formed, return to room temperature, stir, and completely dissolve before use.

7. Preparation of specimen

Serum/Plasma

- Analyze specimens immediately after collection. If it is difficult to analyze immediately, store frozen until analysis.
- Do not use EDTA, an anticoagulant, because it influences the assay.
- The anticoagulant heparin does not significantly influence the assay when used in normal amounts.
- Hemolysis may influence the assay.
- Ascorbic acid does not affect the assay up to 10 mg/dL.
- Dilute specimen with saline and repeat the assay if the measured value exceeds the measurable range, and multiply the result by the dilution factor.

8. Assay procedure

Ensure that the reagents are brought up to room temperature (20°C - 25°C) before use.


- (1) Dispense $100\ \mu\text{L}$ of Reacting Solution into each well of a 96-well microplate.
 - (2) Dispense $10\ \mu\text{L}$ of each concentration of standard solution into the well in which the standard is to be measured.
 - (3) Dispense $10\ \mu\text{L}$ of sample into the well in which the sample is to be measured.
 - (4) Shake the 96-well microplate on a plate mixer (*①).
 - (5) Allow to react in a incubator at 37°C for 20 minutes.
 - (6) Remove the 96-well microplate from the incubator and dispense $50\ \mu\text{L}$ of Substrate Solution into each well.
 - (7) Shake the 96-well microplate on a plate mixer (*①).
 - (8) Allow to react in a incubator at 37°C for 20 minutes.
 - (9) Remove the 96-well microplate from the incubator and dispense $100\ \mu\text{L}$ of Stop Solution into each well.
 - (10) Immediately after shaking (*①), measure the absorbance of each well at 546 nm (reference wavelength 700 nm) using a spectrophotometer for microplates.
- (*①) Guideline for shaking : 500 rpm - 600 rpm for 10 seconds, repeated 6 times.

9. Calculation

- (1) Create a calibration curve for each measurement, with the ATX activity (U/L) on the X-axis and absorbance on the Y-axis.
 - (2) From the calibration curve, read the ATX activity (U/L) corresponding to the absorbance of the sample. If a diluted sample is used, the concentration reading is multiplied by the sample dilution ratio to obtain the measurement value.
- *If the absorbance of a sample is outside the standard curve, dilute the sample with saline to an appropriate dilution ratio and repeat the assay.
- *In computer software calculations, we recommend using a cubic polynomial, or 4 or 5 parameters.

10. Standard Curve

- Measurement range : 1.72 - 55 U/L

11. Performance

- Sensitivity
 - When Buffer is measured as a sample, the absorbance is less than 0.05.
 - When a standard solution of a specific concentration (ATX 55 U/L) is measured as a sample, the absorbance ranges from 0.39 - 0.72.
- Accuracy
 - When control serum of known concentration is measured, the measured value is within $\pm 15\%$ of the known concentration.
- Simultaneous repeatability
 - When the same sample is measured simultaneously 5 times, the CV (%) of measured values is 15% or less.

Notes

- If any reagent accidentally comes into contact with the mouth, eyes, or skin, flush immediately with a large amount of water. Consult a physician if necessary.
- Handle samples being cognizant of the risk of infection (e.g., viruses).
- Wear disposable gloves to avoid the risk of infection when performing a test.
- Store reagents under the specified conditions and do not use expired reagents.
- Do not use reagents that have been accidentally frozen because accurate results may not be obtained.
- After opening, use the reagent as soon as possible. If it is stored, close the lid and store under the specified conditions.
- Do not top up any reagent.
- Be careful regarding contamination when collecting reagents with a pipette.
- Do not mix and use reagents with different lot numbers.
- As the pressure inside the Substrate vial is reduced, open it gently so that the contents do not spill out.
- Open vials carefully to avoid injury from the aluminum cap.
- When disposing of the reagents, dispose of them according to Waste Management and Public Cleansing Law (Waste Disposal and Cleaning Act) and wastewater standards. Substrate Solution contains cyano complex (156 mg/L as cyan).
- Any reagent or reagent bottle that comes into contact with a specimen should be considered potentially infectious and handled accordingly.
- Any scattering/splashes of samples or waste fluids should be removed with a disinfectant such as sodium hypochlorite (effective chlorine concentration 1,000 ppm) or glutaraldehyde (2%).

[References]

- (1) Masuda A. et al. : Clin Chim Acta., 412 (21-22) : 1944-1950 (2011).
- (2) Tokuhara Y. et al. : PLOS ONE, 10 (6) : e0130074 (2015).
- (3) Shao Y. et al. : Medicine (Baltimore), 98 (13) : e14973 (2019).
- (4) Weng J. et al. : J Cell Mol Med., 23 (2) : 1050-1058 (2019).
- (5) Shimamoto R. et al. : TOSOH Research & Technology Review, 61 : 99-103 (2017).

LabAssay™ ATX

[Storage] Store at 2°C - 10°C (Do not freeze)
[Term of validity] Indicated on the label.
[Package] For 150 tests
[Cat #] 293-96901

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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