Vasopressin ELISA Kit *Wako*

[1. Introduction]

This ELISA kit is designed for quantitative determination of (Arg8)-vasopressin. Vasopressin is a 9-amino acid peptide hormone produced in the hypothalamus and released mainly from the posterior pituitary gland. It is called (Arg8)-vasopressin (AVP) because in humans and many mammals the eighth amino acid is arginine. Vasopressin has three types of receptors, V1a, V1b and V2. The V2 receptor is reported to be abundantly expressed in the kidney and involved in its antidiuretic action. On the other hand, V1a and V1b are expressed predominantly in the central nervous system. Recently, it was found that vasopressin is also widely distributed in the brain, and its function in the central nervous system has been attracting attention. Vasopressin has been reported to be involved in stress, social behavior, information processing, spatial learning, and aggressive behavior. Furthermore, vasopressin has been linked to various psychiatric disorders (depression, schizophrenia, autism, anxiety disorders, eating disorders, etc.), and research on vasopressin may be useful in elucidating these disorders and developing therapeutic agents.

[2. Assay Principle]

The microplate is coated with anti-vasopressin antibody. In each well, the standard solution or sample and biotin-conjugated vasopressin are incubated to proceed with the antigen-antibody reaction. Furthermore, peroxidase-conjugated streptavidin is added to proceed with the biotin-streptavidin binding. Finally, peroxidase activity in each well is measured to determine (Arg8)-vasopressin in the sample.

[3. Kit Performance]

Calibration curve range	1.00 - 2,000 pg/mL
Analyte	(Arg8)-vasopressin
Sample	Human : Serum, Plasma (EDTA), Saliva, Urine Mouse : Serum, Plasma (EDTA) Rat : Serum, Plasma (EDTA)
Sample amount required	50 μ L (Measurable at n=2) ^{*1}
Assay time	Approx. 2.5 hours
Detection method	Luminescent detection ^{*2}

*1 Sample preparation is required for measurement. Please refer to [7. Sample Preparation].
*2 Requires luminescence plate reader for measurement.

[4. Precautions for Use]

1. It is important to bring the kit reagents to room temperature (20 - 25°C) before use (approximately 2 hours).
2. Prepare the reagents to be used in the next step in advance, before starting each step.
3. Use this kit after completing training in ELISA methods or under the guidance of an instructor.
4. For manual measurements, individuals who have demonstrated consistent reproducibility in pipette operation should perform the assay.
5. Wear gloves, glasses, and protective clothing during preparation and when using this kit.
6. Prevent reagents coming into contact with the skin. In the event the kit's reagents accidentally come into contact with eyes, mouth, wounds, skin, etc., immediately wash thoroughly with running water and seek medical attention from a physician if necessary.
7. Do not eat, drink, or smoke in the area where this kit is being used.
8. Handle samples with particular care as they may pose a risk of infection. This kit contains animal-derived components.
9. Dispose of used specimens, used reagents, consumables, and unused reagents in accordance with your institution's regulations and local laws and regulations.
10. Do not mix reagents from different lots.
11. It is important to cover the plate with a plate seal to prevent drying of the wells, contamination, temperature fluctuations, and evaporation of dispensed reagents during incubation at each step.
12. ELISA methods can be affected by the assay environment. Strictly adhere to the following : Room temperature should be 20 - 25°C (on the bench or in the incubator) during operation and incubation. Do not perform the assay in environments with wind speeds (including air conditioning) exceeding 0.4 m/sec and humidity levels below 30%.

[5. Materials Supplied]

5-1. Kit Components

Components	Use Status	Amount
Antibody-coated Plate	Use after washing.	96 wells (8×12) / 1 plate
Vasopressin Standard	Freeze-dried. Use after reconstitution.	2 bottles
Buffer	Ready to use.	60 mL / 1 bottle
Biotin-conjugated Vasopressin	Freeze-dried. Use after reconstitution.	2 bottles
Peroxidase-conjugated Streptavidin Solution	Use after dilution.	100 μL / 1 bottle
Luminescent Reagent 1	Ready to use.	6 mL / 1 bottle
Luminescent Reagent 2	Ready to use.	6 mL / 1 bottle
Wash Solution (10×)	Use after dilution.	100 mL / 1 bottle
Sample Buffer 1	Ready to use.	10 mL / 1 bottle
Sample Buffer 2	Ready to use.	5 mL / 1 bottle
Sample Buffer 3	Ready to use.	15 mL / 1 bottle
Plate Seal	—	4 sheets

5-2. Handling of Unused Reagents

When using the kit in divided portions, refer to the following for the stability and storage of unused reagents.

If unopened, they remain stable through the expiration date.

- (1) Antibody-coated Plate
Antibody-coated Plate can be used in divided portions.
Unused antibody-coated strips should be returned into the zip-lock bag provided in the kit and stored at 2 - 10°C.
- (2) Vasopressin Standard
Store unused portion at 2 - 10°C. After dissolving, store at 2 - 10°C and use within 1 week.
Note : Do not store diluted standard solutions.
- (3) Buffer
When using a portion of the solutions, transfer a slightly larger amount than necessary to separate containers. The remaining unused solutions should be tightly capped and stored at 2 - 10°C immediately ; do not bring them to room temperature.
- (4) Biotin-conjugated Vasopressin
After reconstitution, store at 2 - 10°C and use within 1 weeks.
Do not store Biotin-conjugated Vasopressin solution diluted 100-fold with Buffer.
- (5) Peroxidase-conjugated Streptavidin Solution
When using the kit in divided portions, take out the solution from the refrigerator when ready to prepare dilutions. The remaining undiluted solution should be tightly capped and stored at 2 - 10°C immediately ; do not bring it to room temperature.
- (6) Luminescent Reagent 1/Luminescent Reagent 2
When using the kit in several divided portions, take out the solution from the refrigerator, prepare, immediately close the lid tightly, and store the remaining stock solution at 2 - 10°C (do not return to room temperature).
- (7) Wash Solution (10×)
When storing the Wash Solution (10×), close the lid tightly, and store at 2 - 10°C.
- (8) Sample Buffer 1/Sample Buffer 2/Sample Buffer 3
When storing the Sample Buffer 1/Sample Buffer 2/Sample Buffer 3, close the lid tightly, and store at 2 - 10°C.

[6. Apparatus, Equipment, and Materials Required]

- ☐ Purified water (distilled water)
- ☐ Tubes for dilution of standard solutions and samples
- ☐ Glassware for dilution of Wash Solution (graduated cylinders, beakers)
- ☐ Micropipettes (for 10 μL disposable tips and 100 to 500 μL disposable tips)
- ☐ Continuous dispensing pipette, capable of dispensing 100 μL continuously
- ☐ Paper towels or other absorbent material (to remove liquid left on the plate after washing)
- ☐ Vortex-type mixer
- ☐ Microplate shaker (range approx. 500 to 800 rpm)
- ☐ Automatic washer for 96-well plate (if available) or washing bottle
- ☐ 96-well microplate reader for luminescence detection
- ☐ Software for data processing

[7. Sample Preparation]

7-1. Recommendations and Precautions Before Sample Preparation

- Sample collection
Add 100 KIU*/mL aprotinin and 1/750 amount (v/v) Proclin950 to sample.
* Kallikrein Inhibitor Unit
- Sample storage
 - It is recommended to store the sample with 100 KIU/mL aprotinin and 1/750 amount (v/v). Proclin950 at -80℃ or lower for long-term storage. Avoid repeated freezing and thawing.
 - Frozen samples should be thawed just before the assay and thoroughly mixed.
If samples are to be diluted, they should be prepared on an as-needed basis.
- Sample dilution
 - Please dilute the sample before the following preparation process (7-2. Sample Preparation).
 - Use saline as the dilution solvent.
 - Use test tubes (PP, PE) to dilute the samples.
 - Dilute the samples just before use.
- Other issues
 - When turbidity or insoluble materials are present in the samples, remove them by centrifugation or other methods before the assay.
 - If it is suspected that the sample contains an interfering substance, perform a dilution linearity check with two or more different dilutions of the same sample.

7-2. Sample Preparation

- (1) Add the room temperature Sample Buffer 1 to the sample as shown in the table below and shake.
- (2) Add Sample Buffer 2 at room temperature to (1), shake, and incubate for 5 minutes.
Shake it again and incubate 5 minutes.
- (3) Add Sample Buffer 3 which gel is mixed evenly to (2).
Shake it and incubate 5 minutes. Shake it again and incubate 5 minutes.
- (4) After shake it, centrifuge (5,000 - 6,000 g, 10 minutes, 4℃) and use the supernatant as a sample for measurement. Do not suck the gel in the buffer when collecting it.

Sample Preparation		
Sample	Saliva	Serum/Plasma/Urine
Sample Volume	50 μL	50 μL
Sample Buffer 1	50 μL	50 μL
Sample Buffer 2	5 μL	10 μL
Sample Buffer 3	75 μL	150 μL
Dilution Factor	2.85	3.70

[8. Reagent Preparation]

Reagents in the kit must be brought to room temperature (20 - 25℃) before use (about 2 hours).

8-1. Preparation of Standard Solution

Reconstitute Vasopressin Standard with purified water in the amount specified in “Reconstitution of standard”* to prepare the original standard solution (200 ng/mL). Then prepare a dilution series using the Buffer included in the kit that has been allowed to warm up to room temperature.

*Find and check “Reconstitution of standard” on this product page. As the amount of purified water to be added varies by lot, be sure to check it for every lot.

Examples :

Concentration (pg/mL)	Volume of the standard solution	Volume of Buffer
2,000	Standard stock solution : 30 μL	2,970 μL
500	2,000 pg/mL solution : 50 μL	150 μL
125	500 pg/mL solution : 50 μL	150 μL
25.0	125 pg/mL solution : 50 μL	200 μL
5.00	25.0 pg/mL solution : 50 μL	200 μL
1.00	5.00 pg/mL solution : 50 μL	200 μL
0 (B0)	—	150 μL

8-2. Biotin-conjugated Vasopressin

Before starting the procedure, dissolve the lyophilized product in 100 μ L of purified water and dilute the required amount 100-fold with Buffer.

The stock solution after dissolution can be stored at 2 - 10°C for up to 1 week.

Discard any remaining diluted solutions.

8-3. Peroxidase-conjugated Streptavidin Solution

During the 2 hours competitive reaction, dilute 100-fold with room temperature Buffer.

8-4. Luminescent Reagent 1/Luminescent Reagent 2 * extemporaneously prepared

Mix Luminescent Reagent 1 and Luminescent Reagent 2 in a 1 : 1 ratio (v/v) and stored protected from light 30 minutes before use*.

(e.g.) When solution for 96-well reaction is prepared.

Mix 6 mL of Luminescent Reagent 1 and 6 mL of Luminescent Reagent 2.

* This means during the reaction of Peroxidase-conjugated Streptavidin Solution.

8-5. Wash Solution (10 $\times$)

Wash Solution (10 $\times$) is diluted 10-fold with purified water (distilled water).

(e.g.) When solution for 96-well reaction is prepared.

Add 100 mL of Wash Solution (10 $\times$) to 900 mL of purified water (distilled water) to prepare 1,000 mL of Wash Solution (1 $\times$).

- Other reagents are ready-to-use.

[9. Assay Procedure]

- Prepare the reagents to be dispensed in the next step before starting each step.
- For dispensing reagents in steps (2), (9) and (13), continuous dispensing using a multichannel pipette is recommended.
- (1) Fill each well of the Antibody-coated Plate with the diluted Wash Solution and remove. Repeat four times.
Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.
Note : Steps (2) - (5) below should be completed within 10 - 15 minutes.
- (2) Dispense 50 μ L of Biotin-conjugated Vasopressin solution to each well.
- (3) Agitate the plate on a microplate shaker.
Lightly agitation (for about 2 to 3 seconds once) until the liquid level in the well is even.
- (4) Dispense 50 μ L of standard solution into the wells designated for standards.
- (5) Dispense 50 μ L of pretreated-samples (please refer to [7. Sample Preparation]) into the wells designated for samples.
- (6) Agitate the plate using a microplate shaker.
Note : Shaking conditions : 600 - 800 rpm for 10 seconds.
Repeat 3 times.
After each cycle of shaking, stop and then restart again.
- (7) Cover with a Plate Seal.
Incubate at room temperature (20 - 25°C) for 2 hours without shaking.
Note : Peel the protective paper off the plate seal and seal with the adhesive side facing the plate.
Do not reuse a plate seal that has already been used.
- (8) After the reaction is complete, remove the reaction mixture and fill each well with Wash Solution and wash 4 times.
Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.
- (9) Dispense 100 μ L of diluted Peroxidase-conjugated Streptavidin Solution to each well.
- (10) Agitate the plate using a microplate shaker.
Note : Shaking conditions : 600 - 800 rpm for 10 seconds.
Repeat 3 times.
After each cycle of shaking, stop and then restart again.
- (11) Cover with a Plate Seal.
Incubate at room temperature (20 - 25°C) for 30 minutes without shaking.
Note : Peel the protective paper off the plate seal and seal with the adhesive side facing the plate.
Do not reuse a plate seal that has already been used.
- (12) After the reaction is complete, remove the reaction mixture and fill each well with Wash Solution and wash 4 times.
Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.
- (13) Dispense 100 μ L of Luminescent Reagent prepared at room temperature into each well.
- (14) Shake the plate for 1 minute using the microplate shaker.

(15) After agitation, determine the luminescent intensity (RLU) using the 96-well microplate reader (for luminescence measurement). It is recommended to perform the measurement 10 to 20 minutes after addition of the Luminescent Reagent.

[Supplementary Notes]

Washing Procedure

For each wash, dispense the Wash Solution into the wells, then agitate gently by hand for 10 seconds and discard. After 4 consecutive washes, place the plate upside down on a paper towel and tap to completely remove the Wash Solution. After removing the Wash Solution, immediately dispense the next solution to prevent the wells from drying out. The recommended volume of Wash Solution to be dispensed is 300 μ L/well. If foam remains during washing, when manually washing using a wash bottle, allow the Wash Solution to overflow during the final cycle of washing to prevent the formation of foam.

[10. Calculation Method]

- (1) Create a standard curve by plotting the concentration of the standard solution (pg/mL) on the X-axis against the B/B0 (%) value on the Y-axis.
- (2) Read the concentration (pg/mL) corresponding to the B/B0 (%) of the diluted sample.
- (3) Multiply the read concentration by the sample dilution factor to obtain the measured concentration.

[Note]


- Luminescence intensity varies depending on the manufacturer and model of the plate reader.
- For calculations using computer software, we recommend the use of a cubic polynomial, and 4 or 5 parameters.

Examples :

Standard solution concentration (pg/mL)	RLU ($\times 100$)	B/B ₀ (%)
2,000	3,177	28.0
500	6,111	53.9
125	8,530	75.2
25.0	9,941	87.7
5.00	10,506	92.7
1.00	10,733	94.7
0	11,336	100

$B/B_0 \text{ (\%)} = (RLU \text{ of standard or sample} / RLU \text{ of } 0 \text{ pg/mL (B}_0\text{)}) \times 100$

[11. Typical Standard Curve]

[12. Assay Procedure Summary]

Ensure that you read the instruction manual and check the sample conditions, assay conditions, and assay procedure before beginning the assay.

- ☐ Plates and reagents must be brought to room temperature (20 to 25°C) before use (about 2 hours)
- ☐ Dilution of Wash Solution (10×) : Dilute 10-fold with purified water pre-equilibrated to room temperature.
- ☐ Preparation of Standard Solution :

Reconstitute Vasopressin Standard with purified water in the amount specified in “Reconstitution of standard”* to prepare the original standard solution (200 ng/mL). Then prepare a dilution series using the Buffer included in the kit that has been allowed to warm up to room temperature.

*Find and check “Reconstitution of standard” on this product page. As the amount of purified water to be added varies by lot, be sure to check it for every lot.

The following is an example.

Dilution example	Conc.(pg/mL)	2,000	500	125	25.0	5.00	1.00	0
	Standard (μL)	stock solution :30	50*	50*	50*	50*	50*	—
	Buffer (μL)	2,970	150	150	200	200	200	150

* : Standard solution at 1 step higher concentration

- ☐ *Preparation of Biotin-conjugated Vasopressin solution: Dissolve the lyophilized product in 100 μL of purified water and dilute the required amount 100-fold with the Buffer.
- ☐ *Prepare samples according to [7. Sample Preparation].

- ☐ **Antibody-coated Plate**
- ☐ ↓ Wash 4 times (*1)
- ☐ **Diluted Biotin-conjugated Vasopressin solution** **50 μL/well**
- ☐ ↓ Agitate using a microplate shaker
*Lightly agitation (for about 2 to 3 seconds once)
- ☐ **Standard solution or sample after preparation** **50 μL/well**
- ☐ ↓ Agitate using a microplate shaker (*2), room temperature (20 - 25°C)
- ☐ ↓ Cover with a Plate Seal. (*3)
- ☐ ↓ Incubate at room temperature (20 - 25°C) for 2 hours.
- ☐ *Preparation of Peroxidase-conjugated Streptavidin Solution
(dilute 100-folds with Buffer equilibrated to room temperature)
- ☐ ↓ Washing 4 times (*1).
- ☐ **Peroxidase-conjugated Streptavidin Solution** **100 μL/well**
- ☐ ↓ Agitate using a microplate shaker (*2), room temperature (20 - 25°C)
- ☐ ↓ Cover with a Plate Seal. (*3)
- ☐ ↓ Incubate at room temperature (20 - 25°C) for 30 minutes.
- ☐ *Preparation of the Luminescent Reagent
(mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.))
- ☐ ↓ Washing 4 times (*1).
- ☐ **Luminescent Reagent** **100 μL/well**
- ☐ ↓ Agitate at room temperature (20 - 25°C) for 1 minute.
- ☐ **Measurement of luminescence intensity (to be measured between 10 - 20 minutes)**
- (*1) For each wash, dispense the Wash Solution into the wells, then agitate gently by hand for 10 seconds and discard. After 4 consecutive washes, place the plate upside down on a paper towel and tap to completely remove the Wash Solution. After removing the Wash Solution, immediately dispense the next solution to prevent the wells from drying out. The recommended volume of Wash Solution to be dispensed is 300 μL/well. If foam remains during washing, when manually washing using a wash bottle, allow the Wash Solution to overflow during the final cycle of washing to prevent the formation of foam.
- (*2) Shaking should be done at 600 - 800 rpm for 10 seconds. Repeat 3 times.
- (*3) After shaking is complete, cover with a plate seal and incubate without shaking. Peel the protective paper off the plate seal and seal with the adhesive side facing the plate. Do not reuse a plate seal that has already been used.

Vasopressin ELISA Kit *Wako*
[Storage] Store at 2°C - 10°C (Do not freeze).
[Expiration date] Indicated on the label.
[Package] For 96 tests
[Cat #] 291-97301

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