Phosphorylated Tau T181 ELISA kit Wako

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Tau is a microtubule-associated protein that is mainly expressed in neurons in the central nervous system and regulates stability of microtubules. In the brain of patients with Alzheimer's disease, the aggregates of phosphorylated Tau, and neurofibrillary tangles are formed, and the extent of formation has been reported to correlate with severity of dementia. Therefore, Tau has been studied to investigate causes of Alzheimer's disease and to develop drugs for treating the disease. As the concentration of Tau and phosphorylated Tau in cerebrospinal fluid (CSF) of Alzheimer's disease patients is higher than that of non-dementia patients, it is important to measure the concentration of Tau and phosphorylated Tau in CSF. This product is high sensitive ELISA kit to measure Tau with phosphorylated threonine 181 (T181) in CSF.

【Kit Performance】

| | |
|----------------------------|----------------------------------|
| Range of calibration curve | 4.40 - 500 pg/mL |
| Measurement subject | Tau pT181 |
| Sample | Human cerebrospinal fluid (CSF)* |
| Sample volume | 20 μ L |
| Measurement time | 15-18h + 2h |
| Detection method | Luminescent system** |

*Not recommended for plasma and serum.

**Requires luminescence plate reader for measurement

【Kit Content】

| | Component | Condition | Volume or quantity |
|----|---|--------------------------|----------------------------------|
| 1 | Antibody-coated Plate | Use as is | 1 plate 96 wells (8 \times 12) |
| 2 | Phosphorylated Tau Standard | Use after reconstitution | 1 tube |
| 3 | Buffer | Use as is | 60 mL/1 bottle |
| 4 | Biotin-conjugated Antibody Solution | Use after preparation | 100 μ L/1 tube |
| 5 | Peroxidase-conjugated Streptavidin Solution | Use after preparation | 100 μ L/1 tube |
| 6 | Luminescent Reagent 1 | Use as is | 6 mL/1 tube |
| 7 | Luminescent Reagent 2 | Use as is | 6 mL/1 tube |
| 8 | Wash Solution (20 \times) | Use after preparation | 50 mL/1 bottle |
| 9 | Standard Buffer (blue color) | Use as is | 12 mL/1 tube |
| 10 | Sample Buffer (red color) | Use as is | 50 mL/1 bottle |
| 11 | Plate Seal | | 3 sheets |
| 12 | Instruction for use | | 1 copy |

【Principle of assay】

On wells of an assay plate, anti-phosphorylated Tau T181 monoclonal antibody (epitope, Tau pT181) is immobilized. To these wells, the standard solution or sample and biotinylated anti-Tau monoclonal antibody (epitope, 211-231.a.a. or center part of Tau) are added for reaction. Then, peroxidase-conjugated streptavidin is added for reaction. Finally, by measuring the peroxidase activity in each well, a concentration of phosphorylated Tau T181 in the sample can be determined.

【Equipment or materials required but not supplied】

- Purified water (distilled water)
- Tubes for dilution of standard solution/specimen

- Glass apparatuses for dilution of washing solution (graduated cylinder and beaker)
- Pipets with a disposable tip (ones capable of accurately transferring 10 μ L and 100 to 500 μ L of liquid)
- Dispenser capable of dispensing 50 μ L or 100 μ L continuously
- Water-absorbable material such as paper towel (to remove liquid remaining on a plate after washing)
- Mixer (Vortex type)
- Microplate shaker (about 600 to 800 rpm)
- 96-well plate washer (if available) or bottle for washing
- 96-well plate reader (for luminescent assay)
- Data processing software

【Preparation methods of reagents】

Make sure that all kit reagents are brought to room temperature (20°C to 25°C) before use.

(1) Preparation of standard solutions

Reconstitute Phosphorylated Tau Standard with purified water as specified in Appendix* (to prepare the Standard Stock Solution at 10 ng/mL), and prepare the standard solutions by mixing with Standard Buffer provided in the kit as specified in the table below.

*Because the volume of purified water to be added differs from lot to lot, check the volume specified in the Appendix.

| Concentration in the standard solution (pg/mL) | Volume of the standard solution | Standard Buffer |
|--|--|-----------------|
| 500 | Stock solution (10 ng/mL) : 50 μ L | 950 μ L |
| 227 | Standard solution (500 pg/mL) : 200 μ L | 240 μ L |
| 103 | Standard solution (227 pg/mL) : 200 μ L | 240 μ L |
| 47.0 | Standard solution (103 pg/mL) : 200 μ L | 240 μ L |
| 21.3 | Standard solution (47.0 pg/mL) : 200 μ L | 240 μ L |
| 9.70 | Standard solution (21.3 pg/mL) : 200 μ L | 240 μ L |
| 4.40 | Standard solution (9.70 pg/mL) : 200 μ L | 240 μ L |
| 0.00 | — | 240 μ L |

(2) Biotin-conjugated Antibody Solution

Dilute this reagent 100 folds with Buffer.

(3) Peroxidase-conjugated Streptavidin Solution

Dilute this reagent 100 folds with Buffer.

(4) Luminescent Reagent 1 and Luminescent Reagent 2

Mix the Luminescent Reagent 1 and Luminescent Reagent 2 at the ratio of 1 : 1 (vol./vol.) 15 to 30 minutes before use.

Protect the mixture from light until use.

Example : Mix 6 mL each of the Luminescent Reagent 1 and Luminescent Reagent 2 (when all of 96 wells are used).

(5) Wash Solution (20 \times)

Dilute the Wash Solution (20 \times) 20 folds with purified water (distilled water) for use.

Example : Mix 50 mL of the Wash Solution (20 \times) with 950 mL of purified water (distilled water) (when all of 96 wells are used).

The other reagents are directly used.

【Stability of reagents and storage methods】

(1) Antibody-coated Plate

Unused antibody-immobilized strips (kept refrigerated and sealed) may be returned into the zip-seal bag provided in the kit and stored at 2°C to 10°C. They are stable within the shelf life.

(2) Phosphorylated Tau Standard

Reconstitute Phosphorylated Tau Standard with purified water as specified in Appendix* to prepare the Standard Stock Solution at 10 ng/mL. Then, prepare the standard solutions with the Standard Diluent provided in the kit, which has been brought to room temperature. For the volume of purified water to be added, see the Appendix. (*Because the volume of purified water to be added differs from lot to lot, check the volume specified in the Appendix.) Use the prepared standard solutions immediately and do not store these.

(3) Buffer

If a part of Buffer is used, transfer the volume slightly greater than needed to another container, immediately close the lid tightly without bringing the remaining Buffer to room temperature, and store at 2°C to 10°C. They are stable within the shelf life.

(4) Biotin-conjugated Antibody Solution and (5) Peroxidase-conjugated Streptavidin Solution

If the kit is divided for multiple assays, prepare these solutions by dilution just after taking out of the refrigerator, immediately close the lid tightly without bringing the remaining Stock solution to room temperature, and store at 2°C to 10°C. They are stable within the shelf life. Discard the remaining diluted solutions if any.

(6) Luminescent Reagent 1 and (7) Luminescent Reagent 2

If the kit is divided for multiple assays, use these reagents for preparation just after taking out of the refrigerator, immediately close the lid tightly without bringing the remaining Luminescent Reagents to room temperature, and store at 2°C to 10°C. They are stable within the shelf life.

(8) Wash Solution (20×)

Store the Wash Solution (20×) with the lid tightly closed at 2°C to 10°C if applicable. They are stable within the shelf life. Discard the remaining diluted washing solution if any.

(9) Standard Buffer and (10) Sample Buffer

If the kit is divided for multiple assays, use these diluents for preparation just after taking out of the refrigerator, immediately close the lid tightly without bringing the remaining diluents to room temperature, and store at 2°C to 10°C. They are stable within the shelf life.

[Preparation method of samples]


Dilute a sample 3 folds with the Sample Buffer (red) that has been brought to room temperature (example, 20 μL of sample + 40 μL of Sample Buffer).

[Assay procedure]

1. Fill each well of the Antibody-coated Plate with the prepared washing solution, and empty the wells. Repeat this washing operation 4 times. Then, reverse the plate, and tap it against paper towel lightly to remove the remaining liquid from the wells.
2. To the standard assay wells, dispense 50 μL of each standard solution as specified.
3. To the sample assay wells, dispense 50 μL of each sample diluted with the Sample Buffer (prepared sample) as specified.
4. Agitate the plate on a microplate shaker.
5. Affix a plate sealer onto the plate, and allow the plate to stand at 2°C to 10°C for 15 to 18 hours.
6. After end of the reaction, discard the reaction mixture, fill each well with the washing solution, and empty. Repeat this washing operation 4 times. Then, reverse the plate, and tap it against paper towel lightly to remove the remaining liquid from the wells.
7. To each well, dispense 50 μL of the biotin-conjugated antibody solution. Agitate the plate on a microplate shaker.
8. Affix a plate sealer onto the plate, and allow the plate to stand at room temperature (20°C to 25°C) for 90 minutes.

9. After end of the reaction, discard the reaction mixture, fill each well with the washing solution, and empty. Repeat this washing operation 4 times. Then, reverse the plate, and tap it against paper towel lightly to remove the remaining liquid from the wells.
 10. To each well, dispense 50 μ L of the peroxidase-conjugated streptavidin solution. Agitate the plate on a microplate shaker.
 11. Affix a plate sealer onto the plate, and allow the plate to stand at room temperature (20°C to 25°C) for 30 minutes.
 12. After end of the reaction, discard the reaction mixture, fill each well with the washing solution, and empty. Repeat this washing operation 4 times. Then, reverse the plate, and tap it against paper towel lightly to remove the remaining liquid from the wells.
 13. To each well, dispense 50 μ L of the prepared luminescent reagent mixture. Agitate the plate on a microplate shaker for 1 minute.
 14. After the agitation, determine luminescent intensities with a 96-well microplate reader (for luminescent assay). It is recommended to determine the intensities within 10 to 20 minutes after addition of the luminescent reagent mixture.
 15. Prepare the standard curve by plotting the concentration (pg/mL) in the standard solution on the horizontal axis and the luminescent intensity on the vertical axis. Read the concentration (pg/mL) corresponding to the luminescent intensity of the diluted sample. Multiply the read concentration by the dilution factor (3) to obtain the assay value.
- *It is recommended to use a polynomial of the 3rd degree or 4 or 5-parameter logistic fitting if the processing is performed by computer software.

Calibration curve (example)

【Precautions for use】

1. Use polypropylene containers when handling the sample. It is recommended to store the sample frozen at -35°C or lower if long-term storage is intended. Avoid repeated freezing and thawing. Thaw the frozen sample just before assay and thoroughly agitate. Prepare the sample before use.
2. Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
3. If presence of interfering substances is suspected for the sample, dilute the concerned one at multiple dilution ratios to check the dilution linearity. Dilute the prepared sample with the Sample Buffer.
4. The kit shall be used by personnel who have completed training for enzyme-linked immunosorbent assay (ELISA) or under supervision of a director.
5. If the assay is manually performed, personnel who handle a pipet in a reproducible and stable manner shall be engaged.
6. Wear gloves, glasses, and protective garment during preparation and operation with this kit.
7. Avoid skin contact with reagents. If contact of a reagent of this kit with the eye, mouth, wound, or skin occurs, provide a first-aid treatment such as immediately washing the contacted area with tap water thoroughly, and seek medical advice where necessary.


8. Do not drink, eat, or smoke in a place where this kit is used.
9. Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
10. Soak used samples and consumables in 1% formalin, 2% glutaraldehyde or 0.1% or higher sodium hypochlorite solution for at least 1 hour. Or autoclave such materials for sterilization before discarding. Discard used consumables and unused chemical agents in accordance with rules of the facility as well as applicable local laws and regulations.
11. Do not use reagents with different lot numbers together. When allowing the plate to stand in each step, always affix a plate sealer to protect the wells from drying, contamination with foreign matters, uneven temperature, and evaporation of a dispensed reagent.
12. ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20°C to 25°C (on a bench or in an incubator). Avoid performing assay under air flow (including air flow from an air-conditioning equipment) or in a low-humidity environment.

【Examples of assay operation】

○ Spiked recovery test


| | | | | | | | |
|-------------------------------|-------------|------------------------------|----------|-------------------------------|----------|-------------------------------|----------|
| Human cerebrospinal fluid [1] | None-spiked | Spiked with pTau at 20 pg/mL | | Spiked with pTau at 100 pg/mL | | Spiked with pTau at 200 pg/mL | |
| | Assay value | Assay value | Recovery | Assay value | Recovery | Assay value | Recovery |
| | 53.7 pg/mL | 73.4 pg/mL | 98.5% | 147 pg/mL | 93.3% | 271 pg/mL | 109% |
| Human cerebrospinal fluid [2] | None-spiked | Spiked with pTau at 20 pg/mL | | Spiked with pTau at 100 pg/mL | | Spiked with pTau at 200 pg/mL | |
| | Assay value | Assay value | Recovery | Assay value | Recovery | Assay value | Recovery |
| | 14.0 pg/mL | 32.8 pg/mL | 94.0% | 107 pg/mL | 93.0% | 225 pg/mL | 106% |

○ Dilution linearity

○ Assay of clinical samples

Cerebrospinal fluid samples from subjects without Alzheimer's disease (non-AD), patients with mild cognitive impairment (MCI), and those with Alzheimer's disease (AD) were assayed using this ELISA kit.

| | | | | | |
|--------------------|-------------|-----------------|-------------|----------------|------------|
| Non-AD1 | 22.8 | MCI1 | 23.5 | AD1 | 35.2 |
| Non-AD2 | 14.6 | MCI2 | 67.7 | AD2 | 36.2 |
| Non-AD3 | 9.70 | MCI3 | 38.4 | AD3 | 193 |
| Non-AD4 | 17.2 | MCI4 | ND | AD4 | 75.0 |
| Non-AD5 | 24.8 | MCI5 | 63.6 | AD5 | 204 |
| Non-AD mean | 17.8 | MCI mean | 48.3 | AD mean | 109 |

→ A significant difference was observed in assay value between non-AD subjects and AD patients.

[Outline of assay procedure]

- Bring the plate and reagents to room temperature (20°C to 25°C) adequately.
- Dilution of Wash Solution (20×) : Dilute it 20 folds with purified water that has been brought to room temperature.
- Dilution of the standard solution (example) : Reconstitute the Phosphorylated Tau Standard with purified water as specified in Appendix* (to prepare the Phosphorylated Tau Standard Stock Solution at 10 ng/mL), and prepare the standard solutions by mixing with the Standard Buffer provided in the kit as specified in the table below.

| Example of dilution | Concentration (pg/mL) | 500 | 227 | 103 | 47.0 | 21.3 | 9.70 | 4.40 | 0 |
|------------------------|-----------------------|------|------|------|------|------|------|------|-----|
| Standard solution (μL) | Stock solution: 50 | 200* | 200* | 200* | 200* | 200* | 200* | 200* | 0 |
| Standard Buffer (μL) | | 950 | 240 | 240 | 240 | 240 | 240 | 240 | 300 |

*: Standard solution at a 1-step higher concentration

- Dilute a sample 3 folds with the Sample Buffer (red) that has been brought to room temperature (example, 20 μL of sample + 40 μL of Sample Buffer).
- Antibody-coated Plate**
- ↓ Washing 4 times (*[1])
- Standard solution or prepared sample 50 μL/well**
- ↓ Agitation (*[2]) and incubation at 2°C to 10°C for 15 to 18 hours (*[3])
- ↓
- *Preparation of biotin-conjugated antibody solution (dilute 100 folds with Buffer that has been brought to room temperature)
- ↓ Washing 4 times (*[1])
- Biotin-conjugated antibody solution 50 μL/well**
- ↓ Agitation (*[2]) and incubation at room temperature (20°C to 25°C) for 90 minutes (*[3])
- ↓

- *Preparation of peroxidase-conjugated streptavidin solution (dilute 100 folds with Buffer that has been brought to room temperature)
 - ↓ Washing 4 times (*[1])
 - Peroxidase-conjugated streptavidin solution 50 μ L/well**
 - ↓ Agitation (*[2]) and incubation at room temperature (20°C to 25°C) for 30 minutes (*[3])
↓
 - *Preparation of luminescent reagent mixture (mix the Luminescent Reagent 1 and Luminescent Reagent 2 at the ratio of 1: 1 [vol./vol.])
 - ↓ Washing 4 times (*[1])
 - Luminescent reagent mixture 50 μ L/well**
 - ↓ Agitation at room temperature (20°C to 25°C) for 1 minute.
 - Measurement of luminescent intensity (measure between 10 and 20 minutes after addition)**
- (*[1]) In each washing operation, dispense the washing solution into wells, gently agitate the filled plate on the palm for about 10 seconds, and then empty the wells. After washing the wells 4 times consecutively, reverse the plate, and tap it against paper towel to remove the washing solution completely. After removal of the washing solution, immediately dispense the next solution with care not to dry the wells. Use of a pipet set at the liquid volume of 300 μ L may be appropriate for dispensing the washing solution into each well.
- (*[2]) Three repeats of agitation at 600 to 800 rpm for 10 seconds may be appropriate.
- (*[3]) After agitation, affix a plate sealer, and allow it to stand. Remove the liner from the plate sealer, and apply its adhesive side to the plate for affixation. Do not re-use any plate sealer.

[Storage]

Stored at 2- 10°C

[Expiration date]

Indicated on the label

[Package]

For 96 assays

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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