


ScreenFect™A *plus* Transfection Protocol

The optimal condition for successful transfection varies. Start any new transfection by testing the recommended two concentrations of ScreenFect™ *plus* Reagent to determine an optimum amount.

2-Step method (Forward transfection method)

1


Day 0




Pre-Cultured cells

2

Day 1




Reagent
Dilution Buffer



SFA P-reagent
Plasmid DNA
Dilution Buffer





DNA-lipid complex



DNA-lipid complex

3

Day 2~



Steps

Seed cells to be 70-90% confluent at transfection

Dilute ScreenFect™A plus Reagent※¹ in Dilution Buffer, and then mix well
※¹ Vortex the reagent before use

Dilute DNA and SFA P-reagent in Dilution Buffer, and then mix well

Add diluted DNA (+SFA P-reagent) to diluted ScreenFect™A plus Reagent, and then incubate for 5 minutes ~ at room temperature※²
※² Incubation is available until the step 4 has been completed

Add DNA-lipid complex from step 3 to well of cell culture plate from step 1

Visualize/analyze transfected cells

Procedure Details

Component	96-well	24-well	12-well	6-well
Adherent cells or suspension cells	1.0-4.0 × 10 ⁴	0.5-2.0 × 10 ⁵	1.0-4.0 × 10 ⁵	0.25-1.0 × 10 ⁶
Seed cells to be 70-90% confluent at transfection. The medium replacement may improve a transfection efficiency in a 2-Step method.				
Dilution Buffer for ScreenFect™A plus	5 μL	25 μL	50 μL	125 μL
DNA : Transfection Reagent ratio	1 : 3 1 : 4	1 : 3 1 : 4	1 : 3 1 : 4	1 : 3 1 : 4
ScreenFect™A plus Transfection Reagent	0.15 μL 0.2 μL	0.75 μL 1.0 μL	1.5 μL 2.0 μL	3.75 μL 5.0 μL
Dilution Buffer for ScreenFect™A plus	5 μL	25 μL	50 μL	125 μL
DNA (0.1-2.5 μg / μL)	50 ng	250 ng	500 ng	1250 ng
SFA P-reagent (2μL / μg DNA)	0.1 μL	0.5 μL	1.0 μL	2.5 μL
Diluted DNA (+SFA P-reagent)	5 μL	25 μL	50 μL	125 μL
Diluted ScreenFect™A plus Transfection Reagent	5 μL	25 μL	50 μL	125 μL
Final composition [per well]	96-well	24-well	12-well	6-well
DNA-lipid complex	10 μL	50 μL	100 μL	250 μL
DNA amount	50 ng	250 ng	500 ng	1250 ng
SFA P-reagent used	0.1 μL	0.5 μL	1.0 μL	2.5 μL
ScreenFect™A plus Transfection Reagent used	0.15 or 0.2 μL	0.75 or 1.0 μL	1.5 or 2.0 μL	3.75 or 5.0 μL
Medium volume	100 μL	500 μL	1000 μL	2000 μL

Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.

For support, please visit the <http://screenfect.jp>