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

Figure 1. Preparation of Pd/Fib catalyst

Scheme 1. Mechanism of Pd/Fib generation

Scheme 2. Reduction of aromatic carbonyls

Table 1. Chemoselective hydrogenation of olefin in the presence of aromatic carbonyl

Entry	Substrate	Time (h)	Product	Yield (%)
1		30		97 (99 ^a , 97 ^b)
2		3		100
3		37		99
4		46		74
5		2		quant
6		10		99
7		24		100
8		27		100

^a5 atm of pressure of hydrogen. ^bThe reaction was performed at 50°C.

Table 2. Chemoselective hydrogenation of olefin and azide functionalities using 2.5% Pd/Fib

Entry	Substrate	Time (h)	Product	Yield (%)
1		24		100
2		20		98
3		21		99
4		12		100
5		24		91
6		24		100
7		3		94

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Table 3. Chemoselective hydrogenation of olefins and azide in the presence of benzyl esters

Entry	Substrate	Solvent	Time (h)	Product	Yield (%)
1		CD ₃ OD	6		81
2		THF- <i>d</i> ₆	7		91
3		CD ₃ OD	6		69
4		THF- <i>d</i> ₆	7		93
5		MeOH	23		50
6		MeOH	18		77
7		MeOH	24		33
8		THF	24		98
9		MeOH	8		99
10		MeOH	12		100
11		MeOH	6		97
12		MeOH	17		100

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Table 4. Chemoselective hydrogenation of olefins in the presence of aromatic *N*-Cbz

Entry	Substrate	Solvent	Time (h)	Product	Yield (%)
1		MeOH	44		64
2		THF	5		92
3		MeOH	48		97
4		THF	24		48
5		THF	34		99
6		MeOH	25		50
7		MeOH	22		100
8		MeOH	32		92

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Olfactory marker protein (OMP) from PAGE band to structure and function : An overview.

Department of Anatomy and Neurobiology, Program in Neuroscience, School of Medicine, University of Maryland, Baltimore Frank L. Margolis Ph. D. and Jae Hyung Koo Ph. D.

Abstract

The peripheral vertebrate olfactory system is comprised of the olfactory sensory neurons (OSNs) that reside high in the nasal vault within the olfactory neuroepithelium. The axons of these neurons form Cranial Nerve I that innervates the olfactory bulb of the CNS. These mature chemosensory neurons can be replaced from progenitor cells throughout the life of the animal and are characterized by the expression of the olfactory marker protein (OMP) that is a unique hallmark of them. The identification and characterization of this protein will be the primary subject of this essay after we address the biology of the olfactory system. Excellent overviews of the anatomy, neurobiology, function and clinical aspects of the olfactory system are found in Doty, R. L. *et al.*¹.

Introduction

The vertebrate olfactory neuroepithelium (OE) is a generative neuroepithelium that is a site of continuing neurogenesis throughout life (Figure 1). This is an unusual property of OSNs as the vast majority of central nervous system (CNS) neurons are not capable of being replaced and, if damaged or destroyed as

a result of trauma or disease are irrevocably lost. By contrast, in all vertebrates mature OSNs that die or degenerate for whatever reason are replaced from mitotically active progenitors residing in this generative neuroepithelium. These newly formed neurons undergo maturation and can correctly reform synaptic connections with their neuronal targets in the olfactory bulb (OB) of the CNS that were lost as a result of the degeneration of the dying OSNs. This progenitor mitotic activity and subsequent neuronal maturation is associated with changes in gene expression in the OE. The accompanying process of degeneration of the synaptic terminals in the olfactory bulb is also associated with transsynaptic alterations in gene expression^{2,3}. Thus, olfactory bulb neuron phenotype is modulated by degeneration and regeneration of the OSNs and their associated synaptic activity in the OB. In addition to alterations in connectivity, the olfactory bulb neuronal phenotype can be manipulated by alterations in afferent OSN activity. This has been demonstrated by altering or restricting peripheral odor stimulation resulting in modulation of gene expression in target neurons in the olfactory bulb⁴. Current studies in my lab demonstrate that gene expression in deeper cortical areas that receive olfactory input can also be modulated by manipulation of afferent input to the OB⁵ and Kim, H. H. *et al.* in prep). Taken

together all these studies confirm that the central and peripheral portions of the olfactory system are extremely plastic and responsive to manipulations of the external odor environment. In addition, olfactory function is modulated by administration of various systemically administered drugs. These include agents used in chemotherapeutic treatments of cancer, chronic alcohol abuse, and other therapeutic and environmental agents^{1,6}. Clearly it is critical to identify unique molecular reagents that will facilitate the study of this system.

Discovery of OMP

Therefore, we began to search for neuronal cell-specific examples of gene expression in this pathway. The thought was that such proteins could serve as reagents to study the physiology and function of individual classes of neurons in the presence of their cellular neighbors. This may seem somewhat naive today but we refer to a time when many of the standard techniques of contemporary biology were still undreamed of and in the future. Thus, it is important to recall that this investigation began before the advent of monoclonal antibodies, 2D gel electrophoresis, recombinant DNA, PCR, transgenic mice, etc. Indeed we have adopted these various techniques as they arose and have applied them in order to unravel the role of this protein in the olfactory system.

We began by using the low resolving power of 1D non-denaturing gel electrophoresis to search for small acidic proteins that exhibited CNS regional specificity. Extracts of mouse brain regions exhibited quantitative and qualitative differences in staining patterns across CNS regions and in the olfactory system an apparently unique protein band was observed that was absent from other regions of the CNS. To study this in detail, and to study the function of this novel protein, it was necessary to obtain purified homogeneous protein and generate antisera to it. Olfactory tissue was collected from many rats and the protein was purified to homogeneity using the mobility in PAGE as an assay. This purified protein was then used to generate polyclonal antisera^{7,8}. These antisera to the protein were used to characterize its immunoreactivity in several species and many brain regions. It rapidly became apparent


Figure 1. Schematic of olfactory neuroepithelium.

The olfactory neuroepithelium is avascular and rests on a highly vascularized substratum studded with glands whose ducts project to the surface. Non-neuronal sustentacular (supporting) cells are interspersed among the mature olfactory sensory neurons (OSNs). Neuronal precursor cells lie at the base of the olfactory epithelium (OE). The panel illustrates the life cycle of the OSN. Precursor globose basal cells (GBC) undergo mitosis, migration, maturation and differentiation to mature OSNs that express olfactory marker protein (OMP) and odor receptors. Immature neurons express GAP43. Additional markers permit identification of various stages of OSN differentiation and of the various cell types in the OE. Axons of OSNs project to the olfactory bulb where they synapse with mitral and tufted cells and juxtglomerular interneurons. The expression of dopamine (DA), its biosynthetic enzyme tyrosine hydroxylase, and fos in tufted cells and interneurons, is regulated by OSN activity. Modified from G. Ronnett.

the protein was conserved across many vertebrate species and was restricted almost exclusively to the mature olfactory neurons^{9,10}. As a result, and in the absence of any information about its function, we called the protein OMP (Olfactory Marker Protein). Even now, over 30 years later, its function is only beginning to be unraveled. The OMP has been demonstrated to be present in all mammals tested including humans, in several species of fish, in marsupials, in amphibia, and as best as we can tell is present in mature olfactory neurons of all vertebrate species¹¹. The antisera to OMP gave no evidence of immunologically cross-reactive materials in any invertebrate species. The critical question that remained unanswered was that of function. However, even in the absence of knowledge of its function, the nearly unique specificity of the antiserum for mature olfactory neurons has proven to be an extremely valuable reagent. It has been invaluable for studies of olfactory neurogenesis in developing animals, and in response to surgical and chemical lesions of the olfactory system. Its broad phylogenetic cross reactivity has made this antiserum valuable for immunocytochemical studies of the olfactory system in species as diverse as mice (Figures 2a, b) and humans (Figure 3). In addition, the presence of the OMP throughout the cytoplasm of the OSN has also made the antiserum a valuable reagent for the study of the organization of the olfactory bulb during development^{12,13,14} and its response to disturbances of the olfactory system. Immunocytochemistry with this antiserum has facilitated these studies by both light and electron microscopy¹⁵. Curiously, the OMP was also identified in a small subpopulation of neurons in the hypothalamus^{16,17} (Figure 4).


Figure 2.
 (a) OMP immunostaining of mouse olfactory neuroepithelium visualized by confocal microscopy. The OSN cell bodies are apparent with their dendritic processes projecting to the surface of the epithelium where they terminate in dendritic knobs. The OSN axons project centrally and gather into bundles seen in round cross sections. Koo, J. H. *et al.* (2005).
 (b) OMP immunostaining of mouse olfactory bulb (OB) visualized by confocal microscopy. The OSN axons are seen penetrating the surface of the olfactory bulb where they gather into globular neuropil containing structures called glomeruli. No staining is seen deeper in the OB indicative of the specificity of the antiserum. Koo, J. H. *et al.* (2005).

Cloning and mapping of the OMP gene

In addition to the interest in understanding the function of the OMP another gnawing question was to learn why it is so highly restricted to the mature OSNs. The growing availability of recombinant DNA technology provided opportunities to learn more about the mechanisms regulating the highly selective expression of this protein. We determined the amino acid sequence of OMP from HPLC purified tryptic peptides using semi-automated peptide sequencing. Knowing the protein sequence facilitated new directions in the study of this protein and its gene. Peptide sequence information permitted prediction of degenerate oligo-nucleotides that could then be used to identify OMP cDNA clones from an olfactory neuroepithelium library. With this information it was possible to generate additional reagents to study the distribution of the OMP mRNA confirming its virtual restriction to expression in mature OSNs. Further it allowed the mapping of the OMP gene to a defined area of mouse chromosome 7 and human chromosome 11 that was near the locus for an hereditary auditory defect called shaker-1 in mouse and Usher syndrome in humans respectively. The OMP gene was isolated and characterized leading to the demonstration that the entire OMP gene was contained in a single exon and that only one copy of it existed in the genome^{18,19}.

This information and the determination of the nucleotide sequence of this genomic region enabled us to analyze the genomic elements that were responsible for the highly restricted pattern of OMP expression. Using *in vitro* DNA-binding gel-shift assays we

were able to identify and characterize the upstream regulatory elements of the OMP gene. Several of the promoter elements of the OMP gene were identified and it became clear that at least one such element was present in several genes whose expression pattern was highly selective for regulating genes preferentially expressed in the OSNs^{20,21,22}.

OMP transgenic and knock-out mice

This set the stage to take advantage of the then novel technology for generating transgenic mice to perform promoter analyses *in vivo* in OSNs^{23,24}. The OMP promoter was then used by many laboratories to drive expression of various genes into mature OSNs *in vivo* to analyze the biology of the olfactory system. One of the most powerful ways to do this was to utilize the technique of homologous recombination in mice *in vivo* by which a known gene can be inserted into the exact genomic locus normally occupied by OMP²⁵. This allowed the normal genomic regulatory elements that regulated the expression of OMP to drive the expression of an ectopic gene that could be used a probe of function. In one such example the OMP was replaced by the fluorescent protein EGFP leading to the identification of a group of OMP expressing olfactory neurons in an anatomical site that was previously unknown^{26,27,28}. This technology also offered an opportunity to overcome one glaring deficit, that of the function of OMP. After all of these studies we still did not have any insight as to the function of this protein whose expression was developmentally regulated and phylogenetically conserved. Indeed


Figure 3. Immunostaining for OMP in human olfactory tissue. Mature olfactory neurons are stained brown in the middle of the olfactory epithelium with their dendritic processes projecting to the surface of the epithelium where they end in dendritic knobs. OMP-stained axon bundles are visible deep in the epithelium. Buiakova, O. I. *et al.* (1994).

with the advent of extensive gene cloning and with the sequencing of whole genomes it became clear that the gene for OMP was present in every vertebrate species analyzed. The predicted amino acid sequences are >50% identical across all vertebrate species. The gene for this protein is absent from the genomes of invertebrate species including *Drosophila*, *C. elegans*, as well as from unicellular eucaryotes such as yeast. Nevertheless it became apparent that comparison of its amino acid sequence with all the sequences in the databases did not identify any functional domains that might provide information as to its function.

Function of OMP

Therefore, we utilized the technique of homologous recombination noted above to generate mice in which the OMP gene was deleted²⁹. This promised to provide an entrée to function as any novel phenotype observed in the absence of any OMP *in vivo* would provide clues as to function of the OMP. We were not disappointed in this hope. The OMP mice appeared superficially normal but careful electrophysiological and behavioral analyses demonstrated that the OMP-null mice were deficient in their olfactory function. The OMP-null mice were responsive to odors but required 50–100 times higher concentrations to achieve the same behavioral responses as the wild-type intact mice³⁰. In addition,


Figure 4. OMP immunostaining in mouse hypothalamus visualized by confocal microscopy.

The ventricle is apparent as a vertical dark area to the left of the figure. A subset of neurons stained with OMP are evident, as are numerous axonal processes throughout the figure. Koo, J. H. *et al.* (2005).

their electrophysiological responses exhibited a delayed recovery to baseline compared to controls indicating that they were somehow compromised in their ability to recover after odor stimulation^{29,31}. Other studies demonstrated that the OMP-null mice were generally compromised in several parameters associated with their ability to respond to odors^{32,33}. Curiously although the expression of OMP protein is stringently regulated and highly conserved phylogenetically its effect on olfactory transduction processes is quite subtle. We expected dramatic effects on various aspects of olfactory-mediated behavior such as those associated with mating, maternal care, food finding etc. but instead its role is more of a modulator. Ongoing studies to identify exactly where in the overall sensory transduction process are underway³⁴.

Where in the process was not clear as the molecular steps in the transduction cascade appear to be all well characterized without requiring the presence of an additional component. Furthermore, two olfactory subsystems, the main olfactory system and the vomeronasal system both express OMP but they utilize very different transduction cascades^{35,36}. Since OMP is abundant in both subsystems it is difficult to see how it fits in these pathways. One option is that it plays a role not in the primary response transduction pathway, but that it plays a role in regulating the mechanisms associated with return of intracellular calcium to basal levels³⁴. To date, these questions, while still unanswered, are under investigation. Nevertheless, all the evidence indicates that OMP is a participant in the overall olfactory transduction cascade.

Structure of OMP

Another complementary approach to learn something of OMP function was being pursued simultaneously. Possibly the primary amino-acid sequence of the protein is not *per se* a direct determinant of function, but is the basis of a conserved element of the three dimensional structure of OMP that is critical to its function in olfactory transduction. Therefore, we undertook to determine the 3-dimensional structure of OMP in solution^{37,38}. To address this we utilized heteronuclear NMR spectroscopy to determine the structure of OMP in solution. Another group independently determined the structure by X-ray crystallography³⁹. Both structures are essentially

identical and demonstrate that OMP exhibits the beta clam fold formed from eight beta strands and a pair of helical domains and two loops. One of the loops is in the so-called omega loop configuration that is postulated to participate in protein-protein interaction. This loop is highly flexible and could adopt a more ordered structure on interaction with a protein partner. This is consistent with the postulated role of OMP in transduction and in its interaction with the Bex protein^{15,17,37,40}.

Conclusion

These brief comments are intended to illustrate the potential importance of the OMP in olfactory function and the value of the antiserum for studies of the degeneration, regeneration and function of the olfactory system in all vertebrates including humans.

(Literature cited)

- 1) "Handbook of olfaction and gustation," ed. by Doty, R. L., Marcel Dekker, Inc. New York (2003).
- 2) Baker, H., Kawano, T., Margolis, F. L. and Joh, T. H. : "Transneuronal regulation of tyrosine hydroxylase expression in olfactory bulb of mouse and rat.", *J. of Neurosci.*, **3**, 69-78 (1983).
- 3) Margolis, F. L., Verhaagen, J., Biffo, S., Huang, F. L., and Grillo, M. : "Regulation of gene expression in the olfactory neuroepithelium: a neurogenetic matrix." *Prog Brain Res.*, **89**, 97-122 (1991).
- 4) Cho, J. Y., Min, N., Franzen, L. and Baker, H. : "Rapid down-regulation of tyrosine hydroxylase expression in the olfactory bulb of naris-occluded adult rats.", *J. Comp. Neurol.*, **369** (2), 264-76 (1996).
- 5) Kim, H. H., Puche, A. C. and Margolis, F. L. : "Odorant deprivation reversibly modulates NR2B-mediated CREB phosphorylation in mouse piriform cortex.", *A. Chem. S.*, abstract (2006).
- 6) Schiffman, S. S. and Zervakis, J. : "Taste and smell perception in the elderly : effect of medications and disease.", *Adv. Food Nutr. Res.*, **44**, 247-346 (2002).
- 7) Margolis, F. : "A brain protein unique to the olfactory bulb.", *Proc. Nat. Acad. Sci.*, **69**, 1221-1224 (1972).
- 8) Keller, A. and Margolis, F. L. : "Immunological studies of the rat olfactory marker protein.", *J. Neurochem.*, **24**, 1101-1106 (1975).
- 9) Monti-Graziadei, G. A., Margolis, F. L., Harding, J. W. and Graziadei, P. P. C. : "Immunocytochemistry of the olfactory marker protein.", *J. Histochem. Cytochem.*, **25**, 1311-1316 (1977).
- 10) Farbman, A. I. and Margolis, F. L. : "Olfactory marker protein during ontogeny : Immunohistochemical localization.", *Dev. Biol.*, **74**, 205-215 (1980).
- 11) Margolis, F. L. : "Olfactory marker protein (OMP).", *Scand. J. Immunol. Suppl.*, **9**, 181-99 (1982).
- 12) Kim, H. and Greer, C. A. : "The emergence of compartmental organization in olfactory bulb glomeruli during postnatal development.", *J. Comp. Neurol.*, **422** (2), 297-311 (2000).
- 13) Kasowski, H. J., Kim, H. and Greer, C. A. : "Compartmental organization of the olfactory bulb glomerulus.",

- J. Comp. Neurol.*, **407** (2), 261-74 (1999).
- 14) Kosaka, K., Toida, K., Margolis, F. L. and Kosaka, T. : "Chemically defined neuron groups and their subpopulations in the glomerular layer of the rat main olfactory bulb-II. Prominent differences in the intraglomerular dendritic arborization and their relationship to olfactory nerve terminals.", *Neuroscience.*, **76** (3), 775-86 (1997).
 - 15) Koo, J. H., Gill, S., Pannell, L., Menco, B. and Margolis, F. L. : "The interaction of Bex and OMP reveals a metabolically active covalent dimer of OMP.", *J. Neurochem.*, **90**, 102-116 (2004).
 - 16) Baker, H., Grillo, M. and Margolis, F. L. : "Biochemical and Immunocytochemical Characterization of Olfactory Marker Protein in Rodent Central Nervous System.", *J. Comp. Neurol.*, **285**, 246-261 (1989).
 - 17) Koo, J. H., Saraswati, M. and Margolis, F. L. : "Immunolocalization of Bex protein in the mouse brain and olfactory system.", *J. Comp. Neurol.*, **487** (1), 1-14 (2005).
 - 18) Danciger, E., Mettling, C., Vidal, M., Morris, R. and Margolis, F. L. : "The OMP gene : Its structure and olfactory neuron specific expression in transgenic mice.", *Proc. Natl. Acad. Sci. USA.*, **86**, 8565-8569 (1989).
 - 19) Buiakova, O., Rama Krishna, N.S., Getchell, T.V. and Margolis, F.L. : "Human and Rodent OMP genes : Conservation of structural and regulatory motifs and cellular localization.", *Genomics*, **20**, 452-486 (1994).
 - 20) Wang, M. M., Tsai, R. Y., Schrader, K. A. and Reed, R. R. : "Genes encoding components of the olfactory signal transduction cascade contain a DNA binding site that may direct neuronal expression.", *Mol. Cell Biol.*, **13** (9), 5805-13 (1993).
 - 21) Kudrycki, K., Stein-Izsak, C., Behn, C., Grillo, M., Akeson, R. and Margolis, F. L. : "Olf-1-binding site : characterization of an olfactory neuron-specific promoter motif.", *Mol. Cell Biol.*, **13** (5), 3002-14 (1993).
 - 22) Behrens, M., Venkatraman, G., Gronostajski, R. M., Reed, R. R. and Margolis, F. L. : "NFI in the development of the olfactory neuroepithelium and the regulation of olfactory marker protein gene expression.", *Eur. J. Neurosci.*, **12** (4), 1372-84 (2000).
 - 23) Walters, E., Grillo, M., Oestreicher, A. B. and Margolis, F. L. : "LacZ and OMP are co-expressed during ontogeny and regeneration in olfactory receptor neurons of OMP promoter-lacZ transgenic mice.", *Int. J. Dev. Neurosci.*, **14** (7-8), 813-22 (1996).
 - 24) Margolis, F. L. : "Regulation of olfactory neuron gene expression.", *Cytotechnology*, **11**, 17-22 (1993).
 - 25) Mombaerts, P., Wang, F., Dulac, C., Chao, S. K., Nemes, A., Mendelsohn, M., Edmondson, J. and Axel, R. : "Visualizing an olfactory sensory map.", *Cell*, **87** (4), 675-86 (1996).
 - 26) Storan, M. J. and Key, B. : "Septal organ of Grunberg is part of the olfactory system.", *J. Comp. Neurol.*, **494** (5), 834-44 (2006).
 - 27) Fuss, S. H., Omura, M. and Mombaerts, P. : "The Gruenberg ganglion of the mouse projects axons to glomeruli in the olfactory bulb.", *Eur. J. Neurosci.*, **22** (10), 2649-54 (2005).
 - 28) Koos, D. S. and Fraser, S. E. : "The Gruenberg ganglion projects to the olfactory bulb.", *Neuroreport*, **16** (17), 1929-32 (2005).
 - 29) Buiakova, O. I., Baker, H., Scott, J. W., Farbman, A., Kream, R., Grillo, M., Franzen, L., Richman, M., Davis, L. M., Abbondanzo, S., Stewart, C. L. and Margolis, F. L. : "Olfactory Marker Protein (OMP) Gene Deletion Alters Physiological Activity Of Olfactory Neurons.", *Proc. Nat. Acad. Sci.*, **93**, 9858-9863 (1996).
 - 30) Youngentob, S. and Margolis, F. L. : "OMP gene causes an elevation in behavioral threshold sensitivity.", *Neuroreport*, **10**, 15-19 (1999).
 - 31) Ivic, L., Pyrski, M., Margolis, J. W., Richards, L. J., Firestein, S. and Margolis, F. L. : "Adenoviral vector mediated rescue of the OMP-null mouse.", *Nature Neuroscience*, **3**, 1113-1120 (2000).
 - 32) Youngentob, S. L., Kent, P. F. and Margolis, F. L. : "OMP gene deletion results in an alteration in odorant-induced mucosal activity patterns.", *J. Neurophysiol.*, **90** (6), 3864-73 (2003). Epub 2003 Aug 13.
 - 33) Youngentob, S. L., Margolis, F. L. and Youngentob, L. M. : "OMP gene deletion results in an alteration in odorant quality perception.", *Behav. Neurosci.*, **115** (3), 626-31 (2001).
 - 34) Kwon, H. J., Leinders-Zufall, T., Koo, J. H. and Margolis, F. L. : "The absence of OMP causes compromised NCX activity in the mouse olfactory receptor neurons.", *Soc. for Neurosci.*, Abstract (2005).
 - 35) Lucas, P., Ukhanov, K., Leinders-Zufall, T. and Zufall, F. : "A diacylglycerol-gated cation channel in vomeronasal neuron dendrites is impaired in TRPC2 mutant mice : mechanism of pheromone transduction.", *Neuron*, **40** (3), 551-61 (2003).
 - 36) Berghard, A., Buck, L. B., and Liman, E. R. : "Evidence for distinct signaling mechanisms in two mammalian olfactory sense organs.", *Proc. Natl. Acad. Sci. USA.*, **93** (6), 2365-9 (1996).
 - 37) Baldisseri, D. M., Margolis, J. W., Weber, D. J., Koo, J. H. and Margolis, F. L. : "Olfactory marker protein (OMP) exhibits a beta-clam fold in solution : implications for target peptide interaction and olfactory signal transduction.", *J. Mol. Biol.*, **319** (3), 823-37 (2002).
 - 38) Wright, N.T., Margolis, J. W., Margolis, F.L. and Weber, D. J. : "Refinement of the solution structure of rat Olfactory Marker Protein (OMP).", *J. Biomol. NMR.*, **33** (1), 63-8 (2005).
 - 39) Smith, P. C., Firestein, S. and Hunt, J. F. : "The crystal structure of the olfactory marker protein at 2.3 Å resolution.", *J. Mol. Biol.* **319** (3), 807-21 (2002).
 - 40) Behrens, M., Margolis, J. W. and Margolis, F. L. : "Identification of members of the Bex gene family as olfactory marker protein (OMP) binding partners.", *J. Neurochem.*, **86**, 1289-1296 (2003).

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Wako


Figure. Immunofluorescence staining of adult mouse olfactory epithelium with goat anti-OMP (Wako Chemicals USA, Code#544-10001). Green : OMP staining was visualized with Cy2 (Jackson ImmunoResearch). Data was provided by Dr. Frank L. Margolis and Dr. Jae Hyung Koo, Department of Anatomy and Neurobiology, School of Medicine, University of Maryland.

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