

Wako Product Update

GREEN CHEMISTRY -revised version-

Please visit the Wako Online Catalog
<http://search.wako-chem.com>

Wako

Index

Organic Chemistry

1. Catalysts

A. Metal Catalysts..... 1

- a. Co [Cobalt]
- b. Cr [Chromium]
- c. Hf [Hafnium]
- d. Ni [Nickel]
- e. Os [Osmium]
- f. Pd [Palladium]
- g. Pt [Platinum]
- h. Rh [Rhodium]
- i. Ru [Ruthenium]
- j. Sc [Scandium]
- k. Ti [Titanium]
- l. W [Tungsten]
- m. Zr [Zirconium]

B. Polymer-supported Catalysts..... 3

- a. Microencapsulated Catalysts
- b. Triphenyl phosphinated polystyrene complex
- c. Palladium Catalysts
- d. HAP supported Catalyst
- e. met-DIRUX
- f. FibreCat®
- g. Custom Synthesis..... 10

C. Reaction in Water..... 11

- a. Rare-earth Triflate
- b. Pd-Nanocage
- c. Surfactant Combined Catalysts

D. Chiral Phase-Transfer Catalysts (PTC)..... 13

- a. Maruoka Catalysts
- b. TaDiAS

E. Metallocenes..... 16

F. Solid Super-Acid Catalysts..... 18

- a. Nafion®
- b. Zirconia, Sulfated and Zirconia Tungstate

G. Biocatalysts..... 20

- a. Amano Lipase/ Acylase

H. Anthracene..... 21

I. Ionic Liquid..... 21

2. Ligands for Catalytic Asymmetric Synthesis

A. Chiral Crown Ethers..... 22

B. BINAP..... 22

C. linked-BINOL..... 23

3. Others

A. Sulfur Compounds with Less Odor..... 24

B. Condensing agent..... 25

C. Tosyloxylation Reaction..... 26

D. CFC-Alternatives..... 26

E. Dehydrated Solvents..... 27

Analytical Chemistry

1. Paper Fiber-made TLC sheet..... 30

ALPHABETICAL INDEX

☞ See page #28

A. Metal Catalysts

	description		Page
 Co	Metallocene	Cobaltocene	16
	Polymer-supported Catalyst	Dichlorobis(triphenylphosphine)cobalt(II), Supported PS Resin	4
 Cr	Metallocene	Bis(pentamethylcyclopentadienyl)-chromium (II)	16
	Metallocene	Bis(butylcyclopentadienyl)hafnium (IV) Dichloride	16
 Hf	Metallocene	Bis(cyclopentadienyl)hafnium (IV) Dichloride	16
	Polymer-supported Catalyst	Dibromobis(triphenylphosphine)nickel (II), Supported PS Resin	4
 Ni	Polymer-supported Catalyst	Dipotassium Dioxotetrahydroxoosmate (VIII)-Triethylamine PE fibres	9
	Polymer-supported Catalyst	FibreCat®3003, 3004	
	Polymer-supported Catalyst	Osmium (VIII) Oxide, Microencapsulated	3
	Polymer-supported Catalyst	Osmium (VIII) Oxide, Microencapsulated PEM	
	Polymer-supported Catalyst	Osmium (VIII) Oxide-Pyridine PE fibres	9
 Os	Polymer-supported Catalyst	Acetonitrile Dichloropalladium (II)-Triphenylphosphine PE fibres	
	Polymer-supported Catalyst	Di- μ -chlorobis ([η -allyl] palladium (II)), Supported PEG-PS Resin	5
	Polymer-supported Catalyst	FibreCat® 1001, 1007, 1026	
	Polymer-supported Catalyst	Palladium (II) Acetate-Dicyclohexylphenylphosphine PE fibres	9
	Polymer-supported Catalyst	Palladium (II) Acetate-Triphenylphosphine PE fibres	
 Pd	Polymer-supported Catalyst	Palladium (II)-Hydrotalcite	5
	Polymer-supported Catalyst	Palladium (II)-Hydrotalcite (modified)	
	Stable Pd/C Catalyst	Palladium-Activated Carbon Ethylenediamine Complex (Pd: 5%)	6
	Catalyst acted in Water	Palladium-Nanocage	11
	Polymer-supported Catalyst	PEG-PS resin-supported phosphine-Palladium complex	5
	Polymer-supported Catalyst	PEP-Pd	
	Polymer-supported Catalyst	Tetrakis(triphenylphosphine)palladium(0), Supported PS Resin	4
	Polymer-supported Catalyst	Triphenylphosphine Pd (0), Microencapsulated	
 Pt	Polymer-supported Catalyst	cis-Bis(acetonitrile)dichloroplatinum (II)-Triphenylphosphine PE fibres	
	Polymer-supported Catalyst	Dihydrogen Hexachloroplatinate (IV) n-Hydrate-Pyridine PE fibres	
	Polymer-supported Catalyst	FibreCat® 4001, 4003	
 Rh	Polymer-supported Catalyst	Bis(η -norbornadiene)rhodium (I) Tetrafluoroborate-Triphenylphosphine PE fibres	9
	Polymer-supported Catalyst	Di- μ -chlorobis (η -norbornadiene)dirhodium (I)-Triphenylphosphine PE fibres	
	Polymer-supported Catalyst	FibreCat® 2003, 2006	
 Ru	Polymer-supported Catalyst	FibreCat® 3002	
	Catalyst	met-DIRUX	8
	Catalyst	Methanethiolate-Bridged Diruthenium Complex	
	Polymer-supported Catalyst	RuHAP	7
	Polymer-supported Catalyst	Ruthenium(III)-Hydroxyapatite	
 Sc	Metallocene	Ruthenocene	16
	Polymer-supported Catalyst	Sodium Ruthenate (VII)-Triethylamine PE fibres	9
	Catalyst acted in Water	Scandium (III) Trifluoromethanesulfonate	11
	Catalyst acted in Water	Scandium Tris(dodecyl sulfate)Trihydrate	12

 Sc	Polymer-supported Catalyst	Scandium Trifluoromethanesulfonate, Microencapsulated	4
	Catalyst acted in Water	STDS	12
 Ti	Metallocene	Bis(butylcyclopentadienyl)titanium (IV) Dichloride	16
	Metallocene	Bis(cyclopentadienyl) titanium Dichloride	
	Metallocene	Bis(methylcyclopentadienyl)-titanium (IV) Dichloride	
	Metallocene	Bis(t-butylcyclopentadienyl)-titanium (IV) Dichloride	
	Metallocene	Ethylenebis(1-indenyl)titanium (IV) Dichloride	
 W	Metallocene	Bis(cyclopentadienyl)tungsten (IV) Dichloride	19
	Solid Super-Acid Catalyst	Zirconia, Tungstate	
	Solid Super-Acid Catalyst	Zirconia, Tungstate, Pellet	
 Zr	Metallocene	Bis(cyclopentadienyl)dimethyl-zirconium (IV)	16
	Metallocene	Bis(cyclopentadienyl)zirconium (IV) Chloride Hydride	
	Metallocene	Bis(dodecylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(ethylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(indenyl)zirconium (IV) Dichloride	
	Metallocene	Bis(methylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(n-butylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(octylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(propylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(t-butylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Bis(trimethylsilylcyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	(Cyclopentadienyl)(indenyl)-zirconium (IV) Dichloride	
	Metallocene	Dimethylsilylenebis(cyclopentadienyl)-zirconium (IV) Dichloride	
	Metallocene	Diphenylmethylidene(cyclopenta-dienyl)(9-fluorenyl)zirconium (IV) Dichloride	
	Metallocene	meso-Dimethylsilylenebis(2-methyl-1-indenyl)-zirconium (IV) Dichloride	
	Metallocene	rac-Dimethylsilylenebis(1-indenyl)-zirconium (IV) Dichloride	17
	Metallocene	rac-Dimethylsilylenebis(2-methyl-1-indenyl)zirconium (IV) Dichloride	17
	Metallocene	rac-Dimethylsilylenebis(4,5,6,7-tetrahydro-1-indenyl)zirconium (IV) Dichloride	
	Metallocene	meso-Ethylenebis(1-indenyl)-zirconium (IV) Dichloride	16
	Metallocene	rac-1,2-Ethylenebis(2-methyl-1-indenyl)zirconium (IV) Dichloride	17
	Metallocene	rac-Ethylenebis(1-indenyl)-dimethylzirconium(IV)	
	Metallocene	rac-Ethylenebis(4,5,6,7-tetrahydro-1-indenyl)zirconium (IV) Dichloride	16
	Metallocene	Isopropylidene(3-methylcyclopenta-dienyl)(9-fluorenyl)zirconium (IV) Dichloride	
	Metallocene	Isopropylidene(cyclopentadienyl)-(9-fluorenyl)zirconium (IV) Dichloride	16
	Metallocene	Isopropylidenebis(cyclopentadienyl)-zirconium (IV) Dichloride	
	Solid Super-Acid Catalyst	Zirconia, Sulfated	19
	Solid Super-Acid Catalyst	Zirconia, Sulfated, Pellet	
	Solid Super-Acid Catalyst	Zirconia, Tungstate	
	Solid Super-Acid Catalyst	Zirconia, Tungstate, Pellet	
	Metallocene	Zirconocene Dichloride	16

B. Polymer-supported catalysts

a. Microencapsulated Catalysts


[Features]

1. Readily recoverable and reusable by filtration
2. High catalytic activity
3. Utilizing patented technology that reduces release of catalyst from resin
4. Environmentally friendly

Microencapsulated OsO₄, prepared from polystyrene based on the microencapsulated technique

Osmium (VIII) oxide, Microencapsulated is a useful, safe, non-volatile, recoverable and reusable catalyst for the asymmetric dihydroxylation of olefins on an industrial scale.

Low Toxicity due to Low Volatility



Osmium (VIII) oxide, Microencapsulated (MC OsO₄)

Cat. #153-02081 1 g

RT, Solid (black chipped mass)

MW : 254.23 (OsO₄)

CAS:20816-12-0

OsO₄ content : approx. 10 %


IATA : Not restricted

[Application]

Best applied when N-methylmorpholine-N-oxide (NMO) is used as a cooxidant.

[Reaction]

Recovery and reuse of MC OsO₄



Run	1	2	3	4	5
Yield of Product (%)	84	84	83	84	83
Recovery of Catalyst (%)	quant.	quant.	quant.	quant.	quant.

Osmium (VIII) oxide, PEM Microencapsulated

(PEM-MC OsO₄)

Cat. #158-02411 1 g

RT, Solid (black chipped mass)

MW : 254.23 (OsO₄)

CAS:20816-12-0

OsO₄ content : approx. 10 % (Fluorescent X-rays analysis)

IATA : Not restricted


[Application]

PEM-based microencapsulated osmium (VIII) oxide is hydrophilic.



Best applied when potassium ferricyanide (K₃[Fe(CN)₆]) is used as a cooxidant. When an asymmetric ligand is used, PEM- MC OsO₄ can be applied for asymmetric dihydroxylation reaction without the slow addition of olefins.

[Reaction]

PEM-MC OsO₄ (5mol %)



Run	1	2	3
Yield of Product (%)	85	66	84
ee (%)	78	78	78
Recovery of Catalyst (%)	quant.	quant.	quant.



[References]

- 1) S. Kobayashi, M. Endo and S. Nagayama: *J. Org. Chem.*, **63**, 6094 (1998).
- 2) S. Kobayashi, T. Ishida and R. Akiyama: *Org. Lett.*, **3**, 2649 (2001).

Wako will perform contract services to supply microencapsulated osmium (VIII) oxide on an industrial scale and to synthesize diols and asymmetric diols.

(a. Microencapsulated Catalysts)

Sc Scandium Trifluoromethanesulfonate, Microencapsulated

Cat. #196-12041 1 g

RT, Solid

MW : 492.15 ($C_3F_9O_9S_3Sc$)

CAS : 144026-79-9

Appearance : White mass


Sc(OTf)₃ content : approx. 10 %

[Features]

1. Readily recoverable and reusable by filtration
2. Usable in both batch and flow system
3. Hardly any elution of Sc(OTf)₃


Scandium Trifluoromethanesulfonate (Sc(OTf)₃), a Lewis acid enabling to use water as a solvent for a wide range of organic reactions. Wako offers the polystyrene resin-supported catalyst based on the microencapsulation technique.

[Reaction 1]




[Reaction 3]

Friedel-Crafts Acylation (Batch System)




[Reaction 2]



[Reaction 4]

?????? (Batch System)



[Other Reactions]

1. By using imines : AzaDiels-Alder reaction, Allylation reaction
2. Three-component condensation reaction : Mannich-type reaction, Stretcher reaction
3. By using carbonyl compounds : Aldol reaction, Michael reaction, Cyanation reaction, Diels-Alder reaction


b. Triphenylphosphinated polystyrene complex

Catalyst made by fixing metals (Pd, Ni and Co) with the triphenylphosphine structure on polystyrene. This triphenylphosphine structure is bound not by a linker but in a form including the benzene ring of polystyrene. It is therefore relatively stable in various reactions.

[Features]

1. Roughly homogeneous-sized particles with an average diameter of approximately 60 μ m, enabling efficient reaction
2. Readily recoverable and reusable by filtration

[structural formula]



Product List

Catalog No.	Product	Package Size	Content	Structure
Ni 042-28421	Dibromobis(triphenylphosphine)nickel (II), Supported PS Resin	500 mg	Ni : 0.87 mmol/g	2~10°C, Solid ①
Co 049-28431	Dichlorobis(triphenylphosphine)cobalt(II), Supported PS Resin	500 mg	Co : 0.95 mmol/g	2~10°C, Solid ②
Pd 205-15561	Tetrakis(triphenylphosphine)palladium(0), Supported PS Resin	500 mg	Pd : 0.07 mmol/g	2~10°C, Solid ③

c. Palladium Catalysts


Di- μ -chlorobis [$(\eta$ -allyl) palladium (II)], Supported PEG-PS Resin

[PEG-PS resin-supported phosphine-Pd complex; PEP-Pd]

Cat. #043-27731 500 mg




2~10°C, Solid (yellow sticky mass)

Amphiphilic Resin-Supported Pd-Phosphine Catalyst



[Features]

1. Readily recoverable and reusable by filtration for over 30 times
2. Amphiphilic resin-supported phosphine-palladium complex
3. Exhibits high catalytic activities in Tuji-Trost reaction, Suzuki-Miyaura cross-coupling reaction and Heck reaction in water by using PEP-Pd



[References]

- 1) Allylic substitution: (a) Uozumi, Y., Danjo, H., and Hayashi, T., *Tetrahedron Lett.*, **38**, 3557-3560 (1997), (b) Uozumi, Y., et al., *Tetrahedron Lett.*, **39**, 8303-8306 (1998), (c) Danjo, H., et al., *Tetrahedron*, **55**, 14341-14352 (1999).
- 2) Hydroxycarbonylation: Uozumi, Y., Watanabe, T., *J. Org. Chem.*, **64**, 6921-6923 (1999).
- 3) Cross-coupling: Uozumi, Y., et al., *J. Org. Chem.*, **64**, 3384-3388 (1999).

Palladium (II)-Hydrotalcite (Pd : 1.5 %)

Cat. #161-20543 5 g

RT, Solid (pale yellow powder)

Palladium (II)-Hydrotalcite (m*) (Pd : 0.8%)

Cat. #168-20553 10 g


RT, Solid (pale yellow powder)

* m : modified


[Features]

1. Pd(II)-hydrotalcite is a heterogeneous palladium catalyst that homogeneous palladium complex is immobilized on hydrotalcite ($Mg_6Al_2(OH)_{16}CO_3 \cdot H_2O$) which is a naturally produced basic clay mineral.
2. The noble immobilized Pd catalyst is found to be effective for the oxidation of a wide range of alcohols, using atmospheric oxygen or air.
3. In this catalytic system, various alcohols are readily converted to the corresponding aldehydes or ketones with high chemoselectivity or regioselectivity.

[Reaction 1] Oxidation of Geraniol



[Reaction 2] Pd(II)-Hydrotalcite in the Oxidation of Benzyl alcohol



a) reaction conditions: Pd(II)-hydrotalcite (900 mg, 0.15 mmol Pd), benzyl alcohol (3.0 mmol), pyridine (1.5 mmol), toluene (30 mL), 80°C, O₂

b) reaction conditions: Pd(II)-hydrotalcite (300 mg, 0.05 mmol Pd), benzyl alcohol (1.0 mmol), pyridine (0.2 mmol), toluene (10 mL), 65°C, air (1 atm)

c) reaction conditions: Pd(II)-hydrotalcite(m) (1800 mg, 0.15 mmol Pd), benzyl alcohol (3.0 mmol), pyridine (1.5 mmol), toluene (30 mL), 80°C, O₂

[References] 1) T. Nishimura, N. Kakiuchi and S. Uemura: *Chem. Commun.* **2000**, 1245-1246 (2000)

2) N. Kakiuchi, T. Nishimura, M. Inoue and S. Uemura: *Bull. Chem. Soc. Jpn.*, **74**, 165 (2001).




Palladium-Activated Carbon Ethylenediamine Complex (Pd 5%) [Pd/C (en)]

5% Pd/C (en) is a heterogeneous catalyst, whose structure is a 1:1 complex of the Pd/C's Pd and ethylenediamine. It enables a selective catalytic reduction of various functional groups under neutral conditions. It can be readily removed by filtration after the reaction.

5% Pd/C (en) is an excellent reduction catalyst that is not combustible, unlike the general Pd/C, and has good storage stability.

The development into industrial scale is anticipated.

In the catalytic reductions using 5% Pd/C (en), it is possible to easily reduce various functional groups, including olefin, azido, nitro, benzyl ester, aromatic halogen, etc., while inhibiting the reduction of protecting groups, such as benzyl ether¹⁾, Cbz (benzyloxycarbonyl) group^{1), 2)} which is an aliphatic amine, *O*-TBDMS (*t*-butyldimethylsilyl) group³⁾, epoxide⁵⁾ and benzyl alcohol⁵⁾.

[Reaction 1]**Selective reduction in the presence of benzyl ether group¹⁾****[Reaction 2]****Selective reduction in the presence of Cbz group^{3), 4)}****[Reaction 3]****Chemoselective hydrogenation of epoxides⁵⁾**

Catalog No.	Description	Grade	Package Size	Storage
163-21441	Palladium-Activated Carbon Ethylenediamine Complex (Pd 5%)	for Organic Synthesis	1 g	Protect from light
169-21443			5 g	

[References]

- 1) H. Sajiki, K. Hattori, K. Hirota : *J. Org. Chem.*, **63**, 7990 (1998).
- 2) K. Hattori, H. Sajiki, K. Hirota : *Tetrahedron*, **56**, 8433 (2000).
- 3) K. Hattori, H. Sajiki, K. Hirota : *Tetrahedron Lett.*, **41**, 5711 (2000).
- 4) H. Sajiki, K. Hattori, K. Hirota : *Chem. Eur. J.*, **6**, 2200 (2000).
- 5) H. Sajiki, K. Hattori, K. Hirota : *J. Chem Soc., Perkin Trans. 1*, 4043 (1998).

Pd-Nanocage ...See page #11

d. HAP supported catalyst

Ruthenium(III)-Hydroxyapatite (RuHAP)

for Organic Synthesis

Cat. # 182-01851 1 g

188-01853 5 g






RT, Solid

RuHAP, developed by Prof. Kiyotomi Kaneda of Osaka Univ., is a novel, bound oxidation catalyst in the presence of molecular oxygen. Hydroxyapatite is a natural, inorganic crystallized compound. Using its surface as ligand, RuHAP is obtained by treating the surface with RuCl₃ aqueous solution. Because molecular oxygen (O₂) or air is used for reaction, it is unnecessary to use toxic heavy metals such as chromic acid or manganese dioxide, and clean oxidation can be carried out.

[Features]

1. High activity catalyst, 2. can be separated easily from products, 3. recyclable, 4. clean Aerobic oxidation without using toxic reagents

[Reactions]


Oxidation of alcohols 	Intracellular competitive oxidation of primary and secondary alcohols 
Intracellular competitive oxidation of amine and alcohol 	Oxidation of amines 
Oxidation of silanes 	[Reference] <ol style="list-style-type: none"> 1) Yamaguchi, K., Mori, K., Mizugaki, T., Ebitani, K. and Kaneda, K.: <i>J. Am. Chem. Soc.</i>, 122, 7144 (2000). 2) Mori, K., Yamaguchi, K., Mizugaki, T., Ebitani, K. and Kaneda, K.: <i>Chem. Commun.</i>, 461 (2001). 3) Mori, K., Tano, M., Mizugaki, T., Ebitani, K. and Kaneda, K.: <i>New J. Chem.</i>, 1536 (2002).

e. Methanethiolate-bridged diruthenium complex

met-DIRUX

met-DIRUX is a chiral methanethiolate-bridged diruthenium catalyst in which pentamethylcyclopentadiene and methanethiol are coordinated to RuCl₃. It can catalytically react with various substrates via formation of ruthenium-allenylidene complexes with compounds such as propargyl alcohol. For example, it is applicable to cyclization of terminal diynes¹⁾, propargylic substitution reaction of propargyl alcohol^{2) 3)}, cycloaddition with phenol derivatives⁴⁾, propargylation of aromatic compounds⁵⁾ and allenylidene-ene reaction with alkenes⁶⁾.


met-DIRUX is also applicable to asymmetric synthesis and is an air-stable catalyst.




met-DIRUX :[Cp*RuCl(μ_2 -SMe)₂RuCp*Cl]
(Methanethiolate-Bridged Diruthenium Complex)

[Reaction]

Propargylic Substitution Reaction²⁾³⁾



Cycloaddition with Phenol Derivatives⁴⁾



Catalog No.	Description	Grade	Package Size	Storage
130-14581	met-DIRUX	For Organic Synthesis	200 mg	RT

[References]

- 1) Y. Nishibayashi, M. Yamanashi, I. Wakiji, M. Hidai : *Angew. Chem. Int. Ed.*, **39**, 2909 (2000).
- 2) Y. Nishibayashi, M. Yamanashi, I. Wakiji, M. Hidai : *J. Am. Chem. Soc.*, **122**, 11019 (2000).
- 3) Y. Nishibayashi, I. Wakiji, Y. Ishii, S. Uemura, M. Hidai : *J. Am. Chem. Soc.*, **123**, 3393 (2001).
- 4) Y. Nishibayashi, Y. Inada, M. Hidai, S. Uemura : *J. Am. Chem. Soc.*, **124**, 7900 (2002).
- 5) Y. Nishibayashi, M. Yoshikawa, Y. Inada, M. Hidai, S. Uemura : *J. Am. Chem. Soc.*, **124**, 11846 (2002).
- 6) Y. Nishibayashi, Y. Inada, M. Hidai, S. Uemura : *J. Am. Chem. Soc.*, **125**, 6060 (2003).



f. FibreCat®

FibreCat® catalysts, which are immobilized precious metal catalysts with ligands, use graft polymerized polyethylene resins in the form of fibres as a carrier. Their fibre formation allow functional groups to easily interact with substances.

The catalysts are readily recoverable and reusable just as with common immobilized catalysts. Other excellent features are listed below.

[Features]

1. Easy filtration → FibreCat® can be easily separated by filtration because they are insoluble in almost all solvents.
2. High safety → FibreCat® can be treated as noncombustible heterogeneous catalysts.
3. High mechanical strength → FibreCat® are easy to handle because they suffer no physical degradation on stirring.
4. Wide compatibility → FibreCat® can be used in a wide range of solvents, such as aqueous, organic, and even nonpolar hydrocarbon.
5. High heat stability → FibreCat® can be used in general at 120°C and up to 150°C depending on reaction conditions.
6. Ligand addition → FibreCat® do not require additional ligands because they already contain ligands.

[Reaction 1] 1000 Series : Pd containing catalysts (for use in coupling reactions such as Suzuki coupling reaction, Heck reaction)**[Reaction 2] 2000 Series : Rh containing catalysts (for use in hydrogenation reactions)****[Reaction 3] 3000 Series : Ru or Os containing catalysts (for use in selective oxidation reactions)**

Wako Catalog No.	description	Remark	Pkg. Size
324-56621	FibreCat® 1001 (Palladium(II) Acetate-Triphenylphosphine PE fibres)	Pd content : 2.7~4%	1g 5g
320-56623			
321-56631	FibreCat® 1007 (Palladium(II) Acetate-Dicyclohexylphenylphosphine PE fibres)	Pd content : 3~4%	1g 5g
327-56633			
328-56641	FibreCat® 1026 (Acetonitrile Dichloropalladium(II)-Triphenylphosphine PE fibres)	Pd content : 3.5~4.5%	1g 5g
324-56643			
325-56651	FibreCat® 2003 (Bis(η -norbornadiene)rhodium(I) Tetrafluoroborate-Triphenylphosphine PE fibres)	Rh content : 3~4%	1g 5g
321-56653			
322-56661	FibreCat® 2006 (Di- μ -chlorobis(η -norbornadiene)dirhodium(I)-Triphenylphosphine PE fibres)	Rh content : 2~3.5%	1g 5g
328-56663			
329-56671	FibreCat® 3002 (Sodium Ruthenate(VII)-Triethylamine PE fibres)	Ruc ontent : 3.5~4.5%	1g 5g
325-56673			
326-56681	FibreCat® 3003 (Osmium(VIII) Oxide-Pyridine PE fibres)	Osc ontent : Approx. 7.5%	1g 5g
322-56683			
323-56691	FibreCat® 3004 (Dipotassium Dioxotetrahydrooxosmate(VIII)-Triethylamine PE fibres)	Os content : Approx. 7.5%	1g 5g
329-56693			
326-56701	FibreCat® 4001 (Dihydrogen Hexachloroplatinate(IV) n-Hydrate-Pyridine PE fibres)	Pt content : 3.5~5%	1g 5g
322-56703			
323-56711	FibreCat® 4003 (<i>cis</i> -Bis(acetonitrile)dichloroplatinum(II)-Triphenylphosphine PE fibres)	Pt content : 3.5~4.5%	1g 5g
329-56713			

FibreCat® is a registered trademark of Johnson Matthey plc.


g. Custom synthesis using Microencapsulated Osmium (VIII) oxide (OsO_4)

Wako offers custom synthesis of (chiral) diol on an industrial scale using Microencapsulated Osmium (VIII) OsO_4 , which has less volatility.


[Features]

1. We offer custom synthesis on an industrial scale using Microencapsulated Osmium (VIII) OsO_4 .
2. Residual OsO_4 in synthetic products is fewer than when simple OsO_4 is used. (Actual concentration of residual OsO_4 is reported)
3. Microencapsulated Osmium (VIII) OsO_4 can be designed for various reaction systems.
4. Microencapsulated Osmium (VIII) OsO_4 can be used for asymmetric reaction with asymmetric ligands.

Low Toxicity due to Low Volatility




[Before]



[After 3 days]

Please visit
our homepages!

Wako Online Catalog



<http://search.wako-chem.com>

Wako USA homepage



<http://www.wakousa.com>

Wako GmbH homepage



<http://www.wako-chemicals.de>

C. Reaction in water

a. Rare-earth Triflate

It is important for green chemistry that auxiliary substances in organic reactions are used as little as possible and that they are innocuous when used.

Solvents are usually mostly used in organic reactions.

If the reactions can be conducted in aqueous solution instead of organic solvents which are commonly used, the environmental impacts are significantly reduced. Rare-earth triflates, as Lewis acid catalysts, allow the use aqueous solution as a solvent for a wide range of organic reactions.

Catalog No.	Description	MW	CAS	Package Size
195-11391				1 g
191-11393	Scandium (III) Trifluoromethanesulfonate	(CF ₃ SO ₃) ₃ Sc=492.16	144026-79-9	5 g

[References]

- 1) S. Kobayashi, I. Hachiya, M. Araki and H. Ishitani: *Tetrahedron Lett.*, **34**, 3755 (1993).
- 2) S. Kobayashi : *Synlett*, 689 (1994).


Palladium-Nanocage

Cat. #160-20471 1 g

RT, Solid (pale yellow powder)


MW : 2992.48 (C₈₄H₉₆F₄₈O₃₆Pd₆)

Pd-Nanocage which is self-assembled from(en)Pd(NO₃)₂ (as an adhesive) and triazine derivatives (as a molecular panel) is found to promote the aerobic, aqueous oxidation of styrene and its derivatives. It acts as a reverse phase-transfer catalyst, whereas (en)Pd²⁺ as an oxidation catalyst.




[Reaction 1]

Diels-Alder Reaction




[Reaction 2]

Wacker Oxidation Reaction



[Reaction 3]

Wacker Oxidation Reaction




** : ○ shown in the above-mentioned chemical formulae shows the space of Nanocage.

[References]


- 1) D. Oguro, M. Miyazawa, H. Oka, K. Yamaguchi, K. Ogura, M. Fujita : *Nature*, **378**, 469 (1995).
- 2) T. Kusukawa, M. Fujita : *J. Am. Chem. Soc.*, **121**, 1397 (1999).
- 3) H. Ito, T. Kusukawa, M. Fujita : *Chem. Lett.*, 598 (2000).
- 4) M. Yoshizawa, T. Kusukawa, K. Yamaguchi, M. Fujita : *J. Am. Chem. Soc.*, **122**, 6311 (2000).

b. Surfactant Combined Catalysts

Organic reactions in water without use of damaging organic solvents are of great current interest from the standpoint of environmental or ecological concerns. The noble catalysts, which enable the use of water as a solvent for a wide range of organic reactions, will contribute to progress in environmentally benign chemical processes by reducing the use of organic solvents.



STDS



DBSA (n=10-14)

[Features] Immediately forming reaction substrates and emulsion in the water, realizing high-level hydrophobic condition

Sc Scandium Tris (dodecyl sulfate) Trihydrate [STDS]

Cat. #194-12341 1 g
#190-12343 5 g


RT, Solid (white, powder)

MW : 895.17 ($\text{Sc}(\text{OSO}_3\text{C}_{12}\text{H}_{25})_3 \cdot 3\text{H}_2\text{O}$)


Assay (Titration) : 95.0+ %

STDS is a Lewis acid-surfactant combined compound. The newly designed Lewis acid has been shown to be a noble catalyst in the aldol reactions, Mannich-type reactions, allylation reactions and Michael reactions. It has very different characteristics from conventional Lewis acids. The stable dispersion systems including the catalyst and organic substrates are formed in water and highly water-sensitive ketene silyl acetics even react smoothly in the dispersion systems.

[Reaction 1] Aldol Reaction



[Reaction 2] Synthesis of α -Amino phosphonic acid



Linear-Alkylbenzenesulfonic Acid [DBSA]

Cat. #017-15065 500 mL

$\sim 25^\circ\text{C}$, protect from light, Liquid

[Class 8, UN 2584]


MW : 326.49 ($\text{C}_{18}\text{H}_{30}\text{O}_3\text{S}$)

DBSA is a Brønsted acid-surfactant combined compound and catalyzes efficiently three-components Mannich-type reaction of aldehydes, amines and silyl enolates in water. DBSA with the organic substrates forms stable colloidal particles in water and dehydrate reactions like an esterification can even be achieved in the colloidal particles.


[References]

- 1) K. Manabe, Y. Mori, T. Wakabayashi, S. Nagayama and S. Kobayashi : *J. Am. Chem. Soc.*, **122**, 7202 (2000).
- 2) T. Wakabayashi, S. Kobayashi : *Tetrahedron Lett.*, **39**, 5389 (1998).
- 3) K. Manabe, Y. Mori and S. Kobayashi : *Synlett.*, **9**, 1401 (1999).
- 4) K. Manabe, X-M. Mori, S. Kobayashi : *J. Am. Chem. Soc.*, **123**, 10101 (2001).

[Reaction 1] Mannich Reaction




[Reaction 2] Esterification



D. Chiral Phase-Transfer Catalysts(PTC)**a. Maruoka Catalysts**

Chiral phase-transfer catalysts, which are optically active spiro ammonium salts with 2 binaphthyl rings that allow for easy molecule designing, were invented by Professor Maruoka of the Kyoto University.



(R,R)-3,4,5-Trifluorophenyl-NAS Bromide shows high catalytic activity and high enantioselectivity in asymmetric alkylation of α -amino acid derivatives¹⁾. By using (R,R)-3,5-Bistrifluoromethylphenyl-NAS Bromide, β - hydroxy α -amino acid derivatives which are important chiral units of biologically active peptides are obtained in good yield by the aldol reaction of glycine derivatives with aldehydes. It was found that erythro isomers, which are major products, were obtained with high enantioselectivity²⁾.



① (R=3,5-Bistrifluoro-Ph), ② (R=3,4,5-F³-Ph)

[Features]

1. The rational molecular design of C₂-symmetric chiral quaternary ammonium salt
2. Reaction under mild organic / aqueous biphasic conditions
3. Catalyzing the asymmetric alkylations of α -amino acid derivatives and the direct asymmetric aldol reactions of glycine schiff base with aldehydes to the corresponding β -hydroxy- α -amino acid derivatives

[Reaction 1]**Asymmetric alkylation of α -amino acid derivative****[Reaction 2]****Asymmetric aldol reaction of β -hydroxy- α -amino acid derivative****Maruoka Catalysts**

Wako Catalog No.	Description	Package Size	Storage
029-14921	(R,R)-3,5-Bistrifluoromethylphenyl-NAS Bromide [Maruoka Catalyst RR-Bistrifluoromethylphenyl Br Form]	100 mg	RT
025-14923		500 mg	
201-15921	(R,R)-3,4,5-Trifluorophenyl-NAS Bromide [Maruoka Catalyst RR-Trifluorophenyl Br Form]	100 mg	RT
207-15923		500 mg	

[References]


- 1) Ooi T., Takeuchi M., Kameda M. and Maruoka K. : *J. Am. Chem. Soc.*, **122** (21), 5228-5229 (2000).
- 2) Ooi T., Taniguchi M., Kameda M., Maruoka K. : *Angew. Chem. Int.*

b. TaDiAS

TaDiAS (Tartrate-derived DiAmmonium Salt), an asymmetric phase-transfer catalyst which contains 2 quarternary ammonium salts in its molecule, was developed by the Shibasaki group of The University of Tokyo.

(It is considered that) TaDiAS is believed to produce high asymmetric yields of asymmetric compounds by fixation of the anion position in asymmetric space as a result of the collaboration of 2 quarternary ammonium salts.

TaDiAS 1a and 1b, the iodide salts, show high catalytic activity and high enantioselectivity in asymmetric alkylation of α -amino acid derivatives. TaDiAS 2a and 2b, the tetrafluoroboric acids, show high catalytic activity and high enantioselectivity in asymmetric Michael addition of α -amino acid derivatives^{1) 2)}




1a) R¹: t-Bu, R²: Me, R³: OMe, X: I

TaDiAS-[(4S,5S)-2-t-butyl-2-methyl-N,N,N',N'-tetrakis(4-methoxybenzyl)] Diiodide

2a) R¹: Pr, R²: Pr, R³: Me, X: BF₄

TaDiAS-[(4S,5S)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate)





1b) R¹: t-Bu, R²: Me, R³: OMe, X: I

TaDiAS-[(4R,5R)-2-t-butyl-2-methyl-N,N,N',N'-tetrakis(4-methoxybenzyl)] Diiodide

2b) R¹: Pr, R²: Pr, R³: Me, X: BF₄

TaDiAS-[(4R,5R)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate)

[Reaction 1]**Asymmetric alkylation of α -amino acid derivative^{1) 2)}****[Reaction 2]****Asymmetric Michael addition of α -amino acid derivatives^{1) 2)}****TaDiAS**




Catalog No.	Description	Pkg. Size
208-16151	TaDiAS-[(4S,5S)-2-t-butyl-2-methyl-N,N,N',N'-tetrakis(4-methoxybenzyl)] Diiodide	100 mg
201-16141	TaDiAS-[(4R,5R)-2-t-butyl-2-methyl-N,N,N',N'-tetrakis(4-methoxybenzyl)] Diiodide	100 mg
202-16171	TaDiAS-[(4S,5S)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate)	100 mg
205-16161	TaDiAS-[(4R,5R)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate)	100 mg

[References]

- 1) T. Shibuguchi, Y. Fukuta, Y. Akachi, A. Sekine, T. Ohshima, M. Shibasaki : *Tetrahedron Lett.*, **43**, 9539 (2002).
- 2) T. Ohshima, V. Gnanadesikan, T. Shibuguchi, Y. Fukuta, T. Nemoto, M. Shibasaki : *J. Am. Chem. Soc.*, **125**, 11206 (2003).

Various Kinds of Compounds for Organic Chemistry































These three kinds of brochures are available now !!!










F. Metallocene

Metallocene is a general term for bis(cyclopentadienyl) complex salts, $(C_5H_5)_2M$, consisting of 2 cyclopentadiene rings and various transition metals. Ferrocene is the prototypical metallocene. It has high heat stability and many derivatives are produced.

In addition to the ferrocene, many metallocenes containing various transition metals such as Hf, Ru and Ti have been developed and are on the market. Wako offers the following metallocenes and its derivatives.

Cobaltocene	Bis(pentamethylcyclopentadienyl)-chromium(II)	Bis(cyclopentadienyl)-hafnium(IV) Dichloride	Bis(butylcyclopentadienyl) hafnium(IV) Dichloride	Ruthenocene
 [1277-43-6] 320-39241 1g 320-39241 1g 320-39241 1g	 [74507-61-2] 323-27881 1g	 [12116-66-4] 029-13201 1g	 [85722-08-3] 324-39261 500mg	 [1287-13-4] 327-39251 1g
Bis(cyclopentadienyl)titanium(IV) Dichloride  [1271-19-8] 204-04281 1g 200-04283 10g 202-04282 25g	Bis(methylcyclopentadienyl)-titanium(IV) Dichloride  [1282-40-2] 324-35121 1g	Bis(butylcyclopentadienyl) titanium(IV) Dichloride  [1282-40-2] 321-27963 5g	Bis(t-butylcyclopentadienyl)-titanium(IV) Dichloride  [79269-71-9] 327-27921 1g 323-27923 5g	Ethylenebis(1-indenyl)-titanium(IV) Dichloride  [112531-76-7] 327-39371 500mg
Bis(cyclopentadienyl)-tungsten(IV) Dichloride  [12184-26-8] 322-52401 1g	Bis(cyclopentadienyl)zirconium(IV) Chloride Hydride  [37342-97-5] 327-35111 5g 321-39271 500mg 325-35112 25g	Zirconocene Dichloride  [1291-32-3] 260-00541 1g 266-00543 5g 268-00542 25g	Bis(cyclopentadienyl)-dimethyl-zirconium(IV)  [12636-72-5] 320-28011 1g 326-28013 5g	Bis(methylcyclopentadienyl)-zirconium(IV) Dichloride  [12109-71-6] 328-39301 500mg
Bis(ethylcyclopentadienyl)-zirconium(IV) Dichloride  [73364-08-6] 325-39291 500mg	Bis(propylcyclopentadienyl)zirconium(IV) Dichloride  [73364-09-7] 325-39311 500mg	Bis(n-butylcyclopentadienyl)-zirconium(IV) Dichloride  [73364-10-0] 324-27931 1g 320-27933 5g	Bis(t-butylcyclopentadienyl)-zirconium(IV) Dichloride  [32876-92-9] 328-27951 1g 324-27953 5g	Bis(octylcyclopentadienyl)-zirconium(IV) Dichloride  [191803-21-1] 322-45051 100mg
Bis(dodecylcyclopentadienyl)-zirconium(IV) Dichloride  [191803-23-3] 328-39281 500mg	Bis(trimethylsilylcyclopentadienyl)zirconium(IV) Dichloride  [60938-59-2] 322-39321 500mg	Isopropylidenebis(cyclopentadienyl)zirconium(IV) Dichloride  [138533-79-6] 321-47341 100mg	Dimethylsilylenebis(cyclopentadienyl)zirconium(IV) Dichloride  [86050-32-0] 321-39391 500mg	(Cyclopentadienyl)(indenyl)-zirconium(IV) Dichloride  [80155-52-8] 329-39331 500mg
Isopropylidene(cyclopentadienyl)-(9-fluorenyl)zirconium(IV) Dichloride  [130638-44-7] 328-47351 100mg	Isopropylidene(3-methylcyclopentadienyl)(9-fluorenyl)zirconium(IV) Dichloride  [133190-48-4] 325-47361 100mg	Diphenylmethylenecyclopentadienyl(9-fluorenyl)zirconium(IV) Dichloride  [132510-07-7] 320-39361 500mg	Bis(indenyl)zirconium(IV) Dichloride  [12148-49-1] 321-49041 1g 327-49043 5g	meso-Ethylenebis-(1-indenyl)-zirconium(IV) Dichloride  [162429-20-1] 329-45061 100mg

rac-Ethylenebis-(1-indenyl)-dimethylzirconium(IV)  [136844-77-4] 326-45071 100mg 322-45073 500mg	rac-1,2-Ethylenebis(2-methyl-1-indenyl)-zirconium(IV) Dichloride  [143278-87-9] 326-39461 100mg	rac-Dimethylsilylenebis-(1-indenyl)-zirconium(IV) Dichloride  [121009-93-6] 324-39401 500mg	meso-Dimethylsilylenebis-(2-methyl-1-indenyl)-zirconium(IV) Dichloride  [182210-68-0] 323-45081 100mg	rac-Dimethylsilylenebis-(2-methyl-1-indenyl)-zirconium(IV) Dichloride  [149342-08-5] 320-45091 100mg
rac-Ethylenebis-(4,5,6,7-tetrahydro-1-indenyl)-zirconium(IV) Dichloride  [100163-29-9] 324-39381 100mg 320-39383 500mg	rac-Dimethylsilylenebis-(4,5,6,7-tetrahydro-1-indenyl)-zirconium(IV) Dichloride  [126642-97-5] 322-47371 100mg			

CFC-Alternatives

Non-ozone-depleting fluorinated solvents were developed as Chlorofluorocarbon (CFC) -alternatives.

Global warming potential is low. We offer two mixtures of *n*- and iso-butyl isomers, 99.0+% (cGC)



Ethyl Nonafluorobutyl Ether (mixture of isomers), 99.0+% (cGC)

051-06652 25mL

055-06655 500mL

RT, Liquid

M W : 264.09 ($C_6H_5F_9O$)

[Features]

- Water-repellent
- Easily soluble in various organic solvents


Methyl Nonafluorobutyl Ether (mixture of Isomers), 99.0+% (cGC)

139-13412 25mL

133-13415 500mL

RT, Liquid

M W.: 250.06 ($C_5H_3F_9O$)



Solubility in various organic solvents at 25°C

	Solvent					
	Methanol	1-Butanol	Hexane	Dodecane	Diethylether	Acetone
Ethyl Nonafluorobutyl Ether	○	○	○	○	○	○
Methyl Nonafluorobutyl Ether	○	16.8 (w/w%)	○	5.9 (w/w%)	○	○

○ : Very soluble in the solvent.


☞ See page #26 on the related products.

G. Solid Super-Acid Catalysts

Nafion® NR-50

RT, Solid (nearly white~brown, pellet)

Strong acid polymer in a bead form (perfluorinated ion exchange polymer : PFIEP), which can be applied for a wide range of organic reactions requiring an acid catalysis.



Nafion® is a registered mark of DuPont.

[How to recycle]

- Boil for two hours in 150 mL of H₂O, followed by filtration
- Add the resin that is obtained after filtration to 200 mL of 20~25% nitric acid, and agitate it for 4~5 hours at a room temperature. Then filter it.
- Repeat nitric acid treatment several times.
- Rinse the resin that has been treated with acid with distilled water until the water becomes neutral.
- Dry the rinsed resin for 24 hours at a temperature of 100~105°C.

Catalog No.	Description	CAS	Package Size
144-05991			5g
142-05992	Nafion® NR-50	118473-68-0	25g

[References]

- S. Kanemoto, H. Saimoto, K. Oshima, H. Nozaki : *Tetrahedron Lett.*, **25**, 3317 (1984).
- M. Hino, K. Arata : *Apple. Catol.*, **18**, 401 (1985).

Nafion® Dispersion Solution

Nafion® is a Perfluorosulfonic Acid/PTFE copolymer in the (H⁺) form and has unique characters that are widely applicable. Three wt. % formulations, 5%, 10% and 20% are available.

[Physical properties]

	DE520	DE521	DE1020	DE1021	DE2020	DE2021
Polymer Content (%)	min.5.0 max.5.4	min.5.0 max.5.4	min.10.0 max.12.0	min.10.0 max.12.0	min.20.0 max.22.0	min.20.0 max.22.0
Water Content (%)	45 ± 3	45 ± 3	87~90	87~90	34 ± 2	34 ± 2
VOC Content (%)	50 ± 3	50 ± 3	< 1	< 1	46 ± 2	46 ± 2
1-Propanol	48 ± 3	48 ± 3	-	-	44 ± 2	44 ± 2
Ethanol	< 4	< 4	-	-	< 2	< 2
Mixed Ethers and Other VOCs	< 1	< 1	-	-	< 1	< 1
Specific Gravity	0.92-0.94	0.92-0.94	1.05-1.07	1.05-1.07	0.97-0.99	0.97-0.99
Available Acid Capacity (meq/g, H ⁺ polymer basis)	> 1.00	> 0.90	> 1.00	> 0.90	> 1.00	> 0.90
Viscosity (cP;at 25 °C and 40sec ⁻¹ Shear Rate)	10~20	10~20	3~6	3~6	200~400	200~400

Wako Catalog No.	Description	Package Size	Storage
326-46372 324-46373	5% Nafion® Dispersion Solution DE521	25 mL 100mL	RT
323-46382 321-46383	5% Nafion® Dispersion Solution DE520	25 mL 100 mL	
320-46392 328-46393	10% Nafion® Dispersion Solution DE1021	25 mL 100 mL	
323-46402 321-46403	10% Nafion® Dispersion Solution DE1020	25 mL 100 mL	
320-46412 328-46413	20% Nafion® Dispersion Solution DE2021	25 mL 100 mL	
327-46422 325-46423	20% Nafion® Dispersion Solution DE2020	25 mL 100 mL	

Zirconia, Sulfated ($\text{SO}_4 / \text{ZrO}_2$) Zirconia Tungstate ($\text{WO}_3 / \text{ZrO}_2$)

Zirconia, Sulfated and Zirconia Tungstate are solid yet have higher acid strength than sulfate and exhibit excellent catalytic efficiency in various acid catalysis reactions such as esterification, acylation¹⁾, isomerization, ether synthesis, alkylation, disproportionation, polymerization and degradation.

These products can be used at relatively high reaction temperature. Since these can be treated as solids, the corrosivity for apparatuses is low and the acid treatment is rarely required. Therefore these are superior catalysts from the perspective of environmental protection. Wako provides 2 types of products, powder and pellet, for use at various reaction conditions

[Example reactions]

Esterification




- (1) $\text{R}^1 = \text{methyl}$, $\text{R}^2 = \text{methyl, ethyl, propyl, butyl}$
- (2) $\text{R}^1 = \text{acetyl}$, $\text{R}^2 = \text{ethyl}$
- (3) $\text{R}^1 = \text{salicyl}$, $\text{R}^2 = \text{methyl}$
- (4) $\text{R}^1 = \text{phthal, terephthal}$, $\text{R}^2 = \text{octyl, 2-ethylhexyl}$
- (5) $\text{R}^1 = \text{methacryl}$, $\text{R}^2 = \text{methyl}^{*1}$
- (6) $\text{R}^1 = \text{heptyl}$, $\text{R}^2 = \text{methyl}$

*1


Temp.	130 °C
Methanol/Methacrylic Acid	3 mol/mol
LHSV	1 h ⁻¹
Methacrylic Acid Conv.	99
Methyl Methacrylate Select.	>99

Acylation



- | Catalyst | Yield (%) |
|--------------------------------------|-----------|
| SO_4/ZrO_2 | 31.7 |
| WO_3/ZrO_2 | 10.5 |
| $\text{SiO}_2\text{-Al}_2\text{O}_3$ | <0.1 |
- Reaction condition: *p*-chlorobenzoyl chloride 20 mmol, Chlorobenzene 200 mmol, Catalyst 2.0 g, Temp. 135°C, Time 3 h

Synthesis of ether by dehydration




Activity at just after start of reaction (GC) DME/(DME+MeOH): 50%
Activity at 20 hours after start of reaction (GC) DME/(DME+MeOH): 40%

Others





Isomerization



Catalyst	Conversion (wt %)	$2,2\text{-DMB}/\Sigma C_6$ (wt %)
Pt-SO ₄ /ZrO ₂ (Pellet)	89.8	27.5
Pt-SO ₄ /ZrO ₂ (Powder)	84.5	15.6

Reaction condition: Temp. 200°C, Press. 1.08 Mpa, LHSV 1.5 h⁻¹, H₂/Hexane 5 mol/mol

Alkylation



Just before use, this product should be left to dry for 1 hour at an air temperature of 300 ~ 500°C to obtain sufficient expression activity.

References

- 1) Matsuzawa, K. : *Prepr. Am. Chem. Soc. Div. Pet. Chem.*, **42**(4), 734 (1997).

Wako Catalog No	Description	Grade	Package Size
269-01471	Zirconia, Sulfated	Wako 1 st Grade	5 g
267-01472			25 g
268-01762	Zirconia, Sulfated, Pellet	for Organic Synthesis	25 g
260-01761			100 g
267-01771	Zirconia Tungstate	Wako 1 st Grade	5 g
265-01772			25 g
262-01782	Zirconia Tungstate, Pellet	for Organic Synthesis	25 g
264-01781			100 g




H. Biocatalysts

Catalysts for Asymmetric Synthesis


Amano Lipase / Acylase

Utilization of the biocatalysts that are enzymes or microorganisms is a notable technique from the green chemical standpoint of "preventing waste". It has been recognized as being useful for obtaining optically active substances without difficulty. Wako offers hydrolases such as lipase and acylase for producing optically active substance. Three kinds of reusable immobilized enzymes are also in our lineup.

[Reaction 1]



[Reaction 2]



Wako Catalog No.	Description	Source	Activity	Optimum pH	Optimum Temperature	Package Size
321-58331	Lipase AS Amano	<i>Aspergillus niger</i>	not less than 12,000units/g	6.5	45°C	10 g 50 g
327-58333						
328-58341	Lipase M Amano 10	<i>Mucor javanicus</i>	not less than 10,000units/g	7.0	35°C	10 g 50 g
324-58343						
325-58351	Lipase F-AP15	<i>Rhizopus oryzae</i>	not less than 150,000units/g	7.0	40°C	10 g 50 g
321-58353						
322-58361	Lipase G Amano 50	<i>Penicillium camembertii</i>	not less than 50,000units/g	5.0	40°C	10 g 50 g
328-58363						
329-58371	Lipase AYS Amano	<i>Candida rugosa</i>	not less than 30,000units/g	7.0	45°C	10 g 50 g
325-58373						
326-58381	Lipase PS Amano	<i>Burkholderia cepacia</i>	not less than 30,000units/g	7.0~8.0	45°C	10 g 50 g
322-58383						
323-58391	Lipase AK Amano	<i>Pseudomonas fluorescence</i>	not less than 25,000units/g	8.0	60°C	10 g 50 g
329-58393						
326-58401	Lipase PS-C Amano I, Immobilized on Ceramic	<i>Burkholderia cepacia</i>	not less than 1,000units/g	7.0~8.0	45°C	5 g 25 g
324-58402						
323-58411	Lipase PS-C Amano II, Immobilized on Ceramic	<i>Burkholderia cepacia</i>	not less than 600units/g	7.0~8.0	45°C	5 g 25g
321-58412						
320-58421	Lipase PS-D Amano I, Immobilized on Diatomaceous Earth	<i>Burkholderia cepacia</i>	not less than 500units/g	7.0~8.0	45°C	5 g 25 g
328-58422						
327-58431	Acylase Amano	<i>Aspergillus melleus</i>	not less than 30,000units/g	8.0	50°C	10 g 50 g
323-58433						
329-61061	D-Aminoacylase Amano	<i>Escherichia coli</i>	not less than 10.1Munits/g	8.0	45°C	10 Munits 50 Munits
325-61063						

I. Anthracene

Organic zeolite is formed from anthracenebis- resorcinol with metal alkoxides such as La(OPr-*i*)₃ or Zr(OBu-*t*)₄ based on hydrogen-bonded network. Pd-Nanocage is self-assembled from (en)Pd(NO₃)₂ (as an adhesive) and triazine derivatives (as a molecular panel). These unique materials are usefully effective catalysts for the organic reactions in aqueous medium. For example, the Diels-Alder reaction, Wacker oxidation, Michael reaction and nitroaldol reaction are accelerated in water by using these catalysts. These catalytic performances are caused by the function of metal and hydrophobic space containing their inner.

Anthracene-9,10-bis (5-resorcinol)

Cat. #018-18851 100 mg


#014-18853 1 g

RT, Solid (pale yellow, powder)

MW : 394.42 (C₂₆H₁₈O₄)

CAS : 153715-08-3

Assay (HPLC) : 98.0 + %



..... *WAKO PRODUCT UPDATE *

J. Ionic Liquid

Ionic liquid is quaternary ammonium salt with low melting point and high boiling point. With its high stability at from high to low temperatures.

[Features]

1. Nonvolatile liquid
2. Having low-viscosity although the liquid is ionic
3. Heat resistant and with wide liquid temperature range
4. Easy to collect and recycle because of its non-miscibility with organic solvent


1-Ethyl-3-methylimidazolium Trifluoromethanesulfonate

Cat. #059-07111 10 g


RT, Liquid (pale brown, liquid)

MW : 260.23 (C₇H₁₁F₃N₂O₃S)

CAS:145022-44-2



[Reaction 1] Aldol Reaction




A. Chiral Crown Ethers


Chiral crown ether, which forms complexes with free metals in water solution, has turned out to be an effective asymmetric ligand for catalytic asymmetric synthesis reaction.

[References]


S. Nagayama, S. and S. Kobayashi : *J. Am. Chem. Soc.*, **122**, 11531 (2000).

[Reaction]**(R)-2,2'-Binaphthyl-14-crwon-4**


Cat. #023-14681 1 g
RT, Solid (pale brown, powder)
CAS : 128778-82-5

**(R)-2,2'-Binaphthyl-17-crwon-5**


Cat. #023-14701 1 g
RT, Solid (pale yellow, powder)
CAS : 99630-51-0

**(R)-2,2'-Binaphthyl-20-crwon-6**

Cat. #027-14721 1 g
RT, Solid (colorless, viscous mass)
CAS : 75684-69-4

**(S)-2,2'-Binaphthyl-14-crwon-4**


Cat. #020-14691 1 g
RT, Solid (pale yellow, powder)
CAS : 55442-00-7



⟨ MW : 400.47 ($\text{C}_{26}\text{H}_{24}\text{O}_4$) ⟩

(S)-2,2'-Binaphthyl-17-crwon-5


Cat. #020-14711 1 g
RT, Solid (nearly pale yellow, powder)
CAS : 55442-01-8



⟨ MW : 444.52 ($\text{C}_{28}\text{H}_{28}\text{O}_5$) ⟩

(S)-2,2'-Binaphthyl-20-crwon-6

Cat. #024-14731 1 g
RT, Solid (colorless, viscous mass)
CAS : 41024-92-4





⟨ MW : 488.57 ($\text{C}_{30}\text{H}_{32}\text{O}_6$) ⟩


WAKO PRODUCT UPDATE

B. Ligands for Catalytic Asymmetric Synthesis: BINAP


BINAP is one of the most frequently used chiral phosphine compounds, which is used for asymmetric synthesis using various kinds of metal catalysts.




Wako Catalog No.	Product	Package Size	MW	CAS	Assay	Appearance	Storage
029-14301 025-14303 027-14302	(R)-(+)-2,2'-Bis(diphenylphosphino)-1,1'-binaphthyl [(R)-(+)-BINAP]	1 g 5 g 25 g	622.67 ($\text{C}_{44}\text{H}_{32}\text{P}_2$)	76189-55-4	97.0+% (HPLC)	Solid	RT
026-14311 022-14313 024-14312	(S)-(-)-2,2'-Bis(diphenylphosphino)-1,1'-binaphthyl [(S)-(-)-BINAP]	1 g 5 g 25 g	622.67 ($\text{C}_{44}\text{H}_{32}\text{P}_2$)	76189-56-5		Solid	RT

C. Ligands for Asymmetric Two-Center Catalysis : linked-BINOL




(R,R)-Ga-Li-linked-BINOL



The molecule models of Ga-Li-linked BINOL by X-ray crystal structure analysis

**[Features]**

Ligands for direct catalytic asymmetric aldol reaction with unmodified ketones and for the catalytic asymmetric Michael reaction

[Reaction 1] Asymmetric aldol reaction^{1) 2) 4)}**[Reaction 2] Asymmetric Michael reaction^{3) 4)}****[Reaction 3] Asymmetric Michael reaction^{5) 6)}****[References]**

- 1) N. Yoshikawa, N. Kumagai, S. Matsunaga, G. Moll, T. Ohshima, T. Suzuki, M. Shibasaki: *J. Am. Chem. Soc.*, **123**, 2466 (2001).
- 2) N. Kumagai, S. Matsunaga, N. Yoshikawa, T. Ohshima, M. Shibasaki: *Org. Lett.*, **3**, 1539 (2001).
- 3) N. Kumagai, S. Matsunaga, M. Shibasaki: *Org. Lett.*, **3**, 4251 (2001).
- 4) M. Shibasaki, M. Kanai, K. Funabashi: *Chem. Commun.*, 1989 (2002).
- 5) Y. S. Kim, S. Matsunaga, J. Das, A. Sekine, T. Ohshima, M. Shibasaki: *J. Am. Chem. Soc.*, **122**, 6506 (2000).
- 6) S. Matsunaga, T. Ohshima, M. Shibasaki: *Terathedron Lett.*, **41**, 8473 (2000).


Wako Catalog No.	Product	CAS No.	MW	Pkg. Size	Appearance	Storage
155-02421	3,3"-[Oxybis(methylene)] bis-(1<i>R</i>,1<i>"R</i>)-1,1'-bi-2-naphthol [(<i>R,R</i>)-linked-BINOL], 90.0+ % (HPLC)	265116-85-6	614.68 (C ₄₂ H ₃₀ O ₅)	200 mg	Solid	RT
152-02431	3,3"-[Oxybis(methylene)] bis-(1<i>S</i>,1<i>"S</i>)-1,1'-bi-2-naphthol [(<i>S,S</i>)-linked-BINOL], 95.0+ % (HPLC)	336800-79-4	614.68 (C ₄₂ H ₃₀ O ₅)	200 mg		

A. Sulfur Compounds with Less Odor

These sulfur compounds with less odor are made by replacing one of alkyl chains bound to a sulfur atom by a dodecyl radical. Each has less volatility and smells less disagreeable.


Recently, sulfur compounds with less odor have been investigated^{1) 2)}, and Sulfide ① and Sulfoxide ② are each applicable to alcohol oxidation reaction, typified by Corey-Kim oxidation and Swern oxidation. Especially in Corey-Kim oxidation, the reaction progresses even in solvents that are easy to treat. Moreover, Sulfide ① can be applied to dealkylation of ethers and esters. Sulfonium salt ③ can also be used as a methylation agent for micelle formation^{4) 5)} in addition to synthesis of oxirane³⁾.

Thiol ④ and ⑤ that introduced the trimethylsilyl (TMS) group, which readily induces functional group conversion, to the benzene ring will be released soon⁶⁾.




[Example reactions]

Corey-Kim Oxidation




Solvent	Yield (%)
CH_2C_2	96
Toluene	98
AcOEt	97
THF	99
CH_3CN	93
Acetone	100


Swern Oxidation



Dealkylation



Oxirane synthesis



[References]

- 1) K. Nishide, S. Ohsugi, H. Shiraki, H. Tamakita, M. Node: *Org. Lett.*, **3**, 3121 (2001).
- 2) M. Node, K. Kumar, K. Nishide, S. Ohsugi, T. Miyamoto: *Tetrahedron Lett.*, **42**, 9207 (2001).
- 3) Y. Yano, T. Okonogi, M. Sunaga, W. Tagaki: *J. Chem. Soc., Chem. Commun.*, 527 (1973)
- 4) K. Yamauchi, Y. Hisanaga, M. Kinoshita: *J. Am. Chem. Soc.*, **105**, 538 (1983),
- 5) K. Yamauchi, Y. Hisanaga, M. Kinoshita: *J. Chem. Soc. Perkin Trans.*, **1**, 1941 (1983).
- 6) K. Nishide, T. Miyamoto, K. Kumar, S. Ohsugi, M. Node: *Tetrahedron Lett.*, **43**, 8569 (2002).

Wako Catalog No	Description	Package Size
040-28581	Dodecyl Methyl Sulfide	10g
047-28591	Dodecyl Methyl Sulfoxide	10g
040-28601	Dodecyldimethylsulfonium Iodide	10g
209-15961	p -(Trimethylsilyl)benzenethiol	1g
		5g
206-15971	p -(Trimethylsilyl)phenylmethanethiol	10g

B. Condensing agent


Triazine-based Condensing Reagent

DMT-MM

[4-(4,6-Dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium Chloride n-Hydrate]

The most common technique to form amide and ester bonds is dehydrocondensation accompanied by activation of carboxylic acids. Carbodiimides (WSC, DCC) are widely used as a condensing reagent from reactive and economical standpoints. There are, however, problems that need to be addressed, which are induction of allergy and necessity of removal of byproducts. Wako introduces a new triazine-based condensing reagent, DMT-MM. DMT-MM can be used under the same reaction condition as for widely used condensing reagents, WSC and DCC. It also can be used in water, protic solvents, hydrated organic solvents and water-organic system. The solvents are not required to dehydrate.

[Structure]



4-(4,6-Dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium Chloride n-Hydrate

[CAS No. 3945-69-5]

[Reaction example]

Amidation



Run	Cndensing agent	Solvent	Amid (3) ^a	Ester ^a (%)	Amide / Ester
1	DMT-MM	MeOH	98	1.0 (R=Me)	98
2	DMT-MM	EtOH	99	0.7 (R=Et)	141
3	DMT-MM	<i>i</i> -PrOH	96	0.7 (R= <i>i</i> -Pr)	137
4	DCC	MeOH	27 ^b	7.0 (R=Me)	4
5	EDC	MeOH	53	16.0 (R=Me)	3.3

a : Isolated yield

b : *N*-Acylurea was isolated in 21%

Wako Catalog No	Description	Package Size
329-53751	4-(4,6-Dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium Chloride n-Hydrate	5g
327-53752		25g

[References]

M. Kunishima, C. Kawachi, K. Hioki, K. Terao, S. Tani : *Tetrahedron*, **57**, 1551 (2001).


C. Tosyloxylation Reaction

Poly [p-(hydroxy) (tosyloxy) iodostyrene]

Cat. #163-20461 500 mg

-20°C, Solid

Molecular Formula : $-\left[\text{C}_{15}\text{H}_{15}\text{IO}_4\text{S}\right]_n-$




[Features]

1. For oxidative α -tosyloxylation of alcohols, α -tosyloxylation of ketones and oxidative α -tosyloxylation of alkynes
2. Readily isolated from the reaction mixture by filtration

[References]

- 1) H. Togo and K. Sakuratani: *Synlett*, 1966 (2002)
- 2) S. Abe, K. Sakuratani and H. Togo: *Synlett*, 22 (2001)
- 3) S. Abe, K. Sakuratani and H. Togo: *J. Org. Chem.*, **66**, 6174 (2001)

[Reaction]



Oxidative α -tosyloxylation of alcohols

WAKO PRODUCT UPDATE

D. CFC Alternatives HCFC-225

HCFC-225 is a highly safe fluorinated solvent with no flash point. Since it also has excellent features such as low surface tension, low viscosity and selective solubility, it is widely used not only as a detergent, but also as a functional solvent in various fields of industry. HCFC-225 is highly appreciated internationally as the alternative to CFC-113 and is listed as an acceptable substance in the SNAP (Significant New Alternative Policy) Program under the US Clean Air Act.

[Features]

1. Physical properties is similar to CFC-113
2. Cleaning performance is similar to or better than CFC-113
3. Existing cleaning equipment can be used with no modification
4. Non-flammable
5. Compatible with most plastics and metals
6. Low boiling point and short drying time
7. Excellent permeability due to low surface tension
8. High chemical and heat stability
9. Easy solution management
10. Recyclable by distillation


Description	HCFC-225	
Chemical formula	$\text{CF}_3\text{CF}_2\text{CHCl}_2 / \text{CClF}_2\text{CF}_2\text{CHClF}$	
Boiling point	°C	54
Freezing point	°C	-131
Viscosity	$\text{mPa} \cdot \text{s}$	0.59
Surface tension	mN/m	16.2
Flash point	°C	nil

Reference values H₂O (25 °C) : Viscosity 1.0 mPa · s, Surface tension 72.0 mN/m

Wako Catalog No	Description	Package Size
321-39232	HCFC-225 (mixture of $\text{CF}_3\text{CF}_2\text{CHCl}_2$ and $\text{CClF}_2\text{CF}_2\text{CHClF}$)	25g
325-39235		500g

E. Dehydrated Solvents

Solvent for organic synthesis, containing its water content to its minimum. Can be used as a solvent for various organic synthesis reactions where water should be avoided. Each 100mL-and 500 mL-package has a special cap through which a syringe can be inserted. As shown in the figure, the solvent is collected while blowing nitrogen gas into a package.



Product Name	Water content	Package Size			Class	UN	Storage
		100 mL	500 mL	3 L			
Acetone, Dehydrated (99.5+ %)	max. 50 ppm	010-15533	016-15535	014-15531	3	1090	Protect from light
Acetonitrile, Dehydrated (99.0+ %)	max. 50 ppm	017-15543	013-15545	011-15541	3	1648	RT
Benzene, Dehydrated (99.5+ %)	max. 30 ppm	022-12853	028-12855	026-12851	3	1114	RT
1-Butanol, Dehydrated (99.0+ %)	max. 50 ppm	-	020-13035	028-13031	3	1120	RT
2-Butanone, Dehydrated (99.0+ %)	max. 50 ppm	-	027-13045	025-13041	3	1193	Protect from light
Butyl Acetate, Dehydrated (99.0+ %)	max. 50 ppm	027-13263	023-13265	-	3	1123	RT
Chloroform, Dehydrated, containing 0.3 ~ 1.0 % of Ethanol (99.0+ %)	max. 30 ppm	035-16283	031-16285	039-16281	6.1	1888	RT
Chloroform, Dehydrated, containing 150 ppm of Amylene (99.0+ %)	max. 30 ppm	032-16813	038-16815	036-16811	6.1	1888	RT
Cyclohexane, Dehydrated (99.5+ %)	max. 30 ppm	-	036-16595	034-16591	3	1145	RT
Dichloromethane, Dehydrated (99.0+ %)	max. 30 ppm	048-25503	044-25505	042-25501	6.1	1593	Protect from light
Diethyl Ether, Dehydrated (99.5+ %)	max. 50 ppm	-	041-25495	-	3	1155	Protect from light
N,N-Dimethylacetamide, Dehydrated (98.0+ %)	max. 50 ppm	-	042-25285	040-25281	-	-	RT
N,N-Dimethylformamide, Dehydrated (99.5+ %)	max. 50 ppm	041-25473	047-25475	045-25471	3	2265	Protect from light
Dimethyl Sulfoxide, Dehydrated (99.0+ %)	max. 50 ppm	046-26023	042-26025	-	-	-	RT
1,4-Dioxane, Dehydrated (99.0+ %)	max. 50 ppm	-	044-25485	042-25481	3	1165	RT
Ethanol, Dehydrated (99.5+ %)	max. 50 ppm	055-06133	051-06135	059-06131	3	1170	RT
Ethyl Acetate, Dehydrated (99.5+ %)	max. 50 ppm	050-06183	056-06185	054-06181	3	1173	RT
Ethylene Glycol, Dehydrated (99.5+ %)	max. 50 ppm	053-06313	059-06315	-	-	-	RT
Heptane, Dehydrated (99.0+ %)	max. 30 ppm	089-07273	085-07275	-	3	1206	RT
Hexane, Dehydrated (96.0+ %)	max. 30 ppm	089-07033	085-07035	083-07031	3	1208	RT
Methanol, Dehydrated (99.8+ %)	max. 50 ppm	136-12383	132-12385	130-12381	3 & 6.1	1230	RT
4-Methyl-2-pentanone, Dehydratede	max. 50 ppm	131-12713	137-12715	-	3	1245	Protect from light
1-Methyl-2-pyrrolidone, Dehydrated (97.0+ %)	max. 50 ppm	138-12723	134-12725	-	-	-	RT
1-Propanol, Dehydrated (99.5+ %)	max. 50 ppm	-	166-18305	164-18301	3	1274	RT
2-Propanol, Dehydrated, containing Iso-propanol (99.5+ %)	max. 50 ppm	165-17993	161-17995	169-17991	3	1219	RT
Pyridine, Dehydrated (99.5+ %)	max. 50 ppm	161-18453	167-18455	165-18451	3	1282	RT
Tetrahydrofuran, Dehydrated, containing 0.03% BHT (99.5+ %)	max. 50 ppm	206-13433	202-13435	200-13431	3	2056	Protect from light
Tetrahydrofuran, Dehydrated, containing no stabilizer (99.5+ %)	max. 50 ppm	207-13963	203-13965	201-13961	3	2056	Protect from light
Toluene, Dehydrated (99.5+ %)	max. 30 ppm	203-13443	209-13445	207-13441	3	1294	RT
Xylene, Dehydrated (80+ %) (mixture of o-, m- and p-Xylene)	max. 30 ppm	-	242-00685	240-00681	3	1307	RT

ALPHABETICAL INDEX

	page	Description	
A	27	Acetone, Dehydrated	Solvent
	1, 9	Acetonitrile Dichloropalladium (II)-Triphenylphosphine PE fibres	Polymer-supported Catalyst
	27	Acetonitrile, Dehydrated	Solvent
	20	Acylase, Amano	Biocatalyst
	21	Anthracene-9,10-bis(5-resorcinol)	Anthracene
	20	D-Aminoacylase Amano	Biocatalyst
	12	linear-Alkylbenzenesulfonic Acid	Reaction in Water
B	22	(R)-(+)-2,2'-Bis (diphenylphosphino)-1,1'-binaphthyl	Asymmetric Ligand
	22	(R)-(-)-BINAP, (S)-(-)-BINAP	Asymmetric Ligand
	22	(R)-2,2'-Binaphthyl-14-crown-4, (S)-2,2'-Binaphthyl-14-crown-4	Asymmetric Ligand
	22	(R)-2,2'-Binaphthyl-17-crown-5, (S)-2,2'-Binaphthyl-17-crown-5	Asymmetric Ligand
	22	(R)-2,2'-Binaphthyl-20-crown-6, (S)-2,2'-Binaphthyl-20-crown-6	Asymmetric Ligand
	13	(R,R)-3,5-Bistrifluoromethylphenyl-NAS Bromide	Chiral Phase-Transfer Catalyst
	23	(R,R)-linked-BINOL, (S,S)-linked-BINOL	Asymmetric Ligand
	22	(S)-(-)-2,2'-Bis (diphenylphosphino)-1,1'-binaphthyl	Asymmetric Ligand
	27	1-Butanol, Dehydrated	Solvent
	27	2-Butanone, Dehydrated	Solvent
	27	Benzene, Dehydrated	Solvent
	1, 9	cis-Bis(acetonitrile)dichloroplatinum (II)-Triphenylphosphine PE fibres	Polymer-supported Catalyst
	1, 16	Bis(butylcyclopentadienyl)hafnium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(butylcyclopentadienyl)titanium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(t-butylcyclopentadienyl)-titanium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(t-butylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(n-butylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(cyclopentadienyl)dimethyl-zirconium(IV)	Metallocene/Catalyst
	1, 16	Bis(cyclopentadienyl)hafnium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(cyclopentadienyl)titanium Dichloride	Metallocene/Catalyst
	2, 16	Bis(cyclopentadienyl)tungsten(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(cyclopentadienyl)zirconium(IV) Chloride Hydride	Metallocene/Catalyst
	2, 16	Bis(dodecylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(ethylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(indenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(methylcyclopentadienyl)-titanium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(methylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 9	Bis(η -norbornadiene)rhodium (I) Tetrafluoroborate-Triphenylphosphine PE fibres	Polymer-supported Catalyst
	2, 16	Bis(octylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	1, 16	Bis(pentamethylcyclopentadienyl)-chromium(II)	Metallocene/Catalyst
	2, 16	Bis(propylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Bis(trimethylsilylcyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	27	Butyl Acetate, Dehydrated	Solvent
C	2, 16	(Cyclopentadienyl)(indenyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	27	Chloroform, Dehydrated, 150 ppm Amylene added	Solvent
	27	Chloroform, Dehydrated, containing 0.3 ~ 1.0% Ethanol	Solvent
	28	Chromato Sheet	Analytical Chemistry
	1, 16	Cobaltocene	Metallocene/Catalyst
	27	Cyclohexane, Dehydrated	Solvent
D	12	DBSA	Reaction in Water
	1, 4	Dibromobis(triphenylphosphine)nickel (II), Supported PS Resin	Polymer-supported Catalyst
	1, 5	Di- μ -chlorobis [$(\eta$ -allyl) palladium (II)], Supported PEG-PS Resin	Polymer-supported Catalyst
	2, 9	Di- μ -chlorobis (η -norbonadiene)dirhodium (I)-Triphenylphosphine PE fibres	Polymer-supported Catalyst
	1, 4	Dichlorobis(triphenylphosphine)cobalt(II), Supported PS Resin	Polymer-supported Catalyst
	27	Dichloromethane, Dehydrated	Solvent
	27	Diethyl Ether, Dehydrated	Solvent
	27	1,4-Dioxane, Dehydrated	Solvent
	1, 9	Dihydrogen Hexachloroplatinate (IV) n-Hydrate-Pyridine PE fibres	Polymer-supported Catalyst
	25	4-(4,6-Dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium Chloride n-Hydrate	Condensing agent in Water
	27	N,N-Dimethylacetamide, Dehydrated	Solvent
	27	N,N-Dimethylformamide, Dehydrated	Solvent
	2, 16	Diphenylmethylidene(cyclopenta-diienyl)(9-fluorenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	27	Dimethyl Sulfoxide, Dehydrated	Solvent
	2, 16	Dimethylsilylenebis(cyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 17	rac-Dimethylsilylenebis(1-indenyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 17	meso-Dimethylsilylenebis(2-methyl-1-indenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 17	rac-Dimethylsilylenebis(2-methyl-1-indenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 17	rac-Dimethylsilylenebis(4,5,6,7-tetrahydro-1-indenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	1, 9	Dipotassium Dioxotetrahydroxoosmate (VIII)-Triethylamine PE fibres	Polymer-supported Catalyst
	28	met-DIRUX	Catalyst
	25	DMT-MM	Condensing agent in Water
	24	Dodecyl Methyl Sulfide	Alternative to dimethylsulfide


	page	Description	
D	24	Dodecyl Methyl Sulfoxide	Alternative to dimethylsulfoxide
	12	Dodecylbenzenesulfonic acid	Reaction in Water
	24	Dodecyldimethylsulfonium Iodide	Alternative to DMSO
E	27	Ethanol, Dehydrated	Solvent
	27	Ethyl Acetate, Dehydrated	Solvent
	2, 17	rac-Ethylenebis(1-indenyl)-dimethylzirconium(IV)	Metallocene/Catalyst
	2, 16	Ethylenebis(1-indenyl)titanium(IV) Dichloride	Metallocene/Catalyst
	2, 16	meso-Ethylenebis(1-indenyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 17	rac-1,2-Ethylenebis(2-methyl-1-indenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 17	rac-Ethylenebis(4,5,6,7-tetrahydro-1-indenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
F	27	Ethylene Glycol, Dehydrated	Solvent
	21	1-Ethyl-3-methylimidazolium Trifluoromethanesulfonate	Ionic Liquid
F	1, 9	FibreCat®	Polymer-supported Catalyst
H	26	HCFC-225 (mixture of CF ₃ CF ₂ CHCl ₂ and CClF ₂ CF ₂ CHClF)	CFC-Alternatives/Solvent
	27	Heptane, Dehydrated	Solvent
	27	Hexane, Dehydrated	Solvent
I	1, 15	ip-FOXAP	Catalyst with an Asymmetric Ligand
	2, 16	Isopropylidenebis(cyclopentadienyl)-zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Isopropylidene(cyclopentadienyl)-(9-fluorenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
	2, 16	Isopropylidene(3-methylcyclopenta-dienyl)(9-fluorenyl)zirconium(IV) Dichloride	Metallocene/Catalyst
L	23	(R,R)-linked-BINOL, (S,S)-linked-BINOL	Asymmetric Ligand
	20	Lipase, Amano	Biocatalyst
M	13	Maruoka Catalyst	Chiral Phase-Transfer Catalyst
	2, 8	met-DIRUX	Catalyst
	2, 8	Methanethiolate-Bridged Diruthenium Complex	Catalyst
	27	Methanol, Dehydrated	Solvent
N	27	4-Methyl-2-pentanone, Dehydrated	Solvent
	27	1-Methyl-2-pyrrolidone, Dehydrated	Solvent
	18	Nafion® Dispersion Solution	Solid Super-Acid Catalyst
O	18	Nafion® NR-50	Solid Super-Acid Catalyst
	23	3,3''-[Oxybis(methylene)]bis-(1 <i>R</i> ,1'' <i>R</i>)-1,1'-bi-2-naphthol, 3,3''-[Oxybis(methylene)]bis-(1 <i>S</i> ,1'' <i>S</i>)-1,1'-bi-2-naphthol	Asymmetric Ligand
P	1, 3	Osmium (VIII) Oxide, Microencapsulated	Polymer-supported Catalyst
	1, 9	Osmium (VIII) Oxide-Pyridine PE fibres	Polymer-supported Catalyst
P	1, 9	Palladium (II) Acetate-Dicyclohexylphenylphosphine PE fibres	Polymer-supported Catalyst
	1, 9	Palladium (II) Acetate-Triphenylphosphine PE fibres	Polymer-supported Catalyst
	1, 5	Palladium (II)-Hydrotalcite	Polymer-supported Catalyst
	1, 6	Palladium-Activated Carbon Ethylenediamine Complex (Pd: 5%)	Stable Pd/C Catalyst
	1, 11, 6	Palladium-Nanocage	Catalyst acted in Water
	1, 4	PEG-PS resin-supported phosphine-Palladium complex	Polymer-supported Catalyst
	1, 4	PEP-Pd	Polymer-supported Catalyst
	26	Poly[p-(hydroxy)(tosyloxy)iodostyrene]	Polymer-supported reagent
	27	1-Propanol, Dehydrated	Solvent
	27	2-Propanol, Dehydrated	Solvent
R	2, 7	Pyridine, Dehydrated	Solvent
	2, 7	RuHAP	Polymer-supported Catalyst
	2, 7	Ruthenium(III)-Hydroxyapatite	Polymer-supported Catalyst
S	2, 16	Ruthenocene	Metallocene/Catalyst
	1, 11	Scandium (III) Trifluoromethanesulfonate	Reaction in Water
	1, 12	Scandium Tris (dodecyl sulfate) Trihydrate	Reaction in Water
	2, 4	Scandium Trifluoromethanesulfonate, Microencapsulated	Polymer-supported Catalyst
	1, 9	Sodium Ruthenate (VII)-Triethylamine PE fibres	Polymer-supported Catalyst
T	2, 12	STDS	Reaction in Water
	14	TaDiAS-[(4 <i>R</i> ,5 <i>R</i>)-2-t-butyl-2-methyl-N,N,N'N'-tetrakis(4-methoxybenzyl)] Diiodide	Chiral Phase-Transfer Catalyst
	14	TaDiAS-[(4 <i>S</i> ,5 <i>S</i>)-2-t-butyl-2-methyl-N,N,N'N'-tetrakis(4-methoxybenzyl)] Diiodide	Chiral Phase-Transfer Catalyst
	14	TaDiAS-[(4 <i>R</i> ,5 <i>R</i>)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate)	Chiral Phase-Transfer Catalyst
	14	TaDiAS-[(4 <i>S</i> ,5 <i>S</i>)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate)	Chiral Phase-Transfer Catalyst
	27	Tetrahydrofuran, Dehydrated, containing 0.03% BHT	Solvent
	27	Tetrahydrofuran, Dehydrated, containing no stabilizer	Solvent
	27	Tetrahydrofuran, Dehydrated, containing no stabilizer	Solvent
	1, 4	Tetrakis(triphenylphosphine)palladium(0), Supported PS Resin	Polymer-supported Catalyst
	27	Toluene, Dehydrated	Dehydrated Solvent
X	13	(<i>R,R</i>)-3,4,5-Trifluorophenyl-NAS Bromide	Chiral Phase-Transfer Catalyst
	24	p-(Trimethylsilyl)benzenethiol	Less Odor Sulfur Compound
	24	p-(Trimethylsilyl)phenylmethanethiol	Less Odor Sulfur Compound
	27	Xylene, Dehydrated	Solvent
Z	2, 19	Zirconia, Sulfated	Solid Super-Acid Catalyst
	2, 19	Zirconia, Tungstate	Solid Super-Acid Catalyst
	2, 16	Zirconocene Dichloride	Metallocene/Catalyst

Chromato Sheet

Cat. #036-17151 25 sheets

RT


Chromato sheet is an environment-friendly product, making full use of advantages of "paper." It is light, easy to use and disposable with no damage to environment. Unlike a conventional thin-layer chromatoplate (TLC plate), this is made from paper fiber which is coated with silica gel and fluorescent F254. This product provides usability of paper and separability of silica gel (adsorption/partition mode).



Left : Chromato Sheet
Right : Silicagel 70F₂₅₄ Plate-wako
Sample : Wakogel B-Tester (#231-00051)
Developing solvent : Chloroform

[Specification]

Size : 20 x 20 cm
Weight : Approximately 7 g/sheet
Thickness : 0.3 mm
Silica gel : Wakogel C-500HG
containing F254 fluorescent
Indicator.



Detection at 254 nm
Sample : Brucine, oxypropyltheophylline
and caffeine
Developing solvent :
Chloroform(9) + Methanol(1)


[Features]

1. High resolution & reproducibility
2. Clipping, Wiring and Filing as a paper
3. Meet no-detached silica gel powder
4. Applicable to blotting
5. Applicable to fluorometric detection

[Limitation]

Inapplicable to use color-producing reagents containing strong acids and to treat carbonization by heating at a high temperature

**Chromato Sheet,
which is a TLC paper
can be treated
as a paper!**



- All products are sold for laboratory use only. They are not for use in humans.
- Please visit our online catalog to search for other products from Wako ; <http://search.wako-chem.com>
- This brochure may contain products that cannot be exported to your country due to regulations.
- Bulk quote requests for some products are welcomed. Please contact us.

04315IBK

Wako Pure Chemical Industries, Ltd.

1-2, Doshomachi 3-Chome
Chuo-Ku, Osaka 540-8605, Japan
Telephone :+81-6-6203-3741
Facsimile :+81-6-6201-5964
Online Cat. : <http://search.wako-chem.com>

Wako Chemicals USA, Inc.

1600 Bellwood Road
Richmond, VA 23237, U.S.A.
Toll-Free (U.S. only) :+1-877-714-1920
Telephone :+1-804-714-1920
Facsimile :+1-804-271-7791
<http://www.wakousa.com>

Wako Chemicals GmbH

Nissanstraße 2, D-41468
Neuss, Germany
Telephone :+49-(0)2131-311-0
Facsimile :+49-(0)2131-311100
E-mail : biochem@wako-chemicals.de