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**NH IMMUNOCHROMATO O111 <<Instruction Manual>>**

\* Read this manual before you start using this test kit

**[Intended use]**

NH Immunochromato O111 is intended for use in food testing laboratories for qualitative detection of *E. coli* O111 from various food products. Not intended for diagnostic use. This kit is intended for use by personnel familiar with the appropriate aseptic techniques for the isolation and identification of *E. coli* O111.


**[Introduction]**

Enterohemorrhagic *Escherichia coli* (EHEC) strains can cause bloody diarrhea or hemorrhagic colitis in humans. Hemorrhagic colitis can progress to potentially fatal hemolytic uremic syndrome. EHEC is transmitted to humans through consumption of contaminated foods and water. Although the most prevalent serogroup of EHEC is O157, non-O157 EHEC, such as such as O26, O111, and O103, has been increasingly detected all over the world.

NH Immunochromato O111 is an immunochromatographic test kit for simple and rapid testing of *E. coli* O111. The kit does not require skilled personnel or well-equipped laboratory.

**[Kit contents]**

- A: Test strip 2 tests × 10 packs (20 tests)  
 B: Instruction manual 1 sheet  
 C: Plastic bag with zipper 1 piece

**[Principle of assay]****1) Illustration of test strip**

- a. Sample application zone  
 b. Reagent zone  
 c. Nitrocellulose membrane  
 d. Absorbent pad  
 e. Test zone  
 f. Control zone

**Note 1:** Do not scratch or touch with bare fingers on the sample application zone (a) or nitrocellulose membrane (c).

**2) Principle of assay**

When a sample solution is applied to the sample application zone of the test strip, it is wicked through reagent zone. *E. coli* O111 (1) in the sample binds to the colloidal gold-labeled anti-*E. coli* O111 antibody (2) in the reagent zone. The pathogen-antibody complex flows through the nitrocellulose membrane which contains a zone of antibody specific to *E. coli* O111. The immune complex is captured and concentrated in the zone (test zone) (3), eventually forms a reddish purple line. The membrane also contains a control zone where the gold-labeled anti-*E. coli* O111 antibody is captured. The control line forms in the control zone (4) regardless of the presence of *E. coli* O111 in the sample solution. The control line ensures the test is working correctly.

**[Detection limit]**

The test using 7 strains of *E. coli* O111 revealed that enriched sample including  $1 \times 10^4$ – $10^6$  CFU/mL of *E. coli* O111 will yield a positive result.

**Note 2:** The lower detection limit of the kit may vary depending on the strain and components of the test solution.

**[Precautions]**

- E. coli* O111 has been shown to have low infective dose. Users of NH Immunochromato O111 should be familiar with the appropriate aseptic techniques for the isolation and identification of *E. coli* O111. Extreme care should be taken in handling samples, enriched cultures and used test strips.
- Do not use expired test strips. The expiration date is indicated on the aluminum pouches and the box.
- Good microbiological laboratory practices should be used to prevent contamination of samples, lab workers and the environment.


**[Additionally required materials and instruments]**

- Enrichment broth (e.g. mEC with novobiocin, mEC)
- Stomacher and stomacher bags (preferably with a filter)
- Autoclave
- Incubator (42±1°C)
- Balance (capable of weighing 25 g)
- Disposable plastic transfer pipettes and/or appropriate micro pipettes and disposable tips
- Timer
- Disposable polypropylene tubes for sterilization of samples (optional)

**[Test procedures]**

- Weigh 25-g sample into the stomacher bag and add 225-mL enrichment broth, e.g. mEC with novobiocin, mEC.
- Homogenize with a stomacher for 1 minute.
- Incubate the homogenized sample at 42±1°C for 22±2 hours.
- Remove the stomacher bag from the incubator and mix the culture gently using side-to-side motion.
- Transfer approximate 5-mL enriched sample into an appropriate test tube with a sterilized pipette (Note 3). Cover with a loose-fitting cap.
- Sterilize enriched sample by autoclaving (121°C, 20 minutes) or boiling (15 minutes). The sterilized sample solutions should be applied to test strips after cooling to room temperature (Note 4).
- Allow the test strips warm to room temperature (20–25°C) before opening the aluminum pouch. Open the pouch immediately before use (Notes 5, 6).

\*Enriched sample should be applied to a test strip by either of following procedures (8a or 8b).

- Place a test strip on flat surface and with a new pipette tip apply 100 µL of the enriched sample, drop-wise to the sample application zone. Careful not to overload the sample application zone; overloading it may cause spillage.
- Transfer 150 µL of the enriched sample into a test tube with a new pipette tip. Put a test strip in the test tube immersing the sample application zone in the sample solution (Note 7).
- Observe the test result 15 minutes after applying the sample solution to the test strip.


**Note 3:** Keep the remaining enriched samples unsterilized until the end of the test. Confirmatory tests, such as isolation of *E. coli* O111, might be required after testing with the kit. Unsterilized enriched samples can be used in confirmatory tests.

**Note 4:** Sterilization is not mandatory for the test performance but recommended to reduce the potential risk of infection and contamination.

**Note 5:** Do not open the aluminum pouch before it comes to the room temperature to prevent moisture absorption by test strips. Moist test strips can give incorrect test results. Unused test strips should be kept at 2–8°C in the plastic bag (provided with the kit) along with a desiccant sheet.

**Note 6:** Do not scratch or touch with bare fingers on the sample application zone or nitrocellulose membrane. Hold the absorbent pad when handling test strips.

**Note 7:** Keep the surface of sample solution below the arrows printed on the test strip.

**[Interpretation of the test results]**

- Positive results: Reddish purple lines appear in the test zone and the control zone in 15 minutes after application of sample solutions.
- Negative results: A reddish purple line appears in the control zone.
- Invalid results: No reddish purple line appears in the control zone.

**Note 8:** Presumptive positive results should be confirmed according to authorized methods <sup>1)</sup>. The enriched samples prepared for testing with the kit can be used in confirmatory tests.

**[Disposal]**

Decontaminate used test strips, test tubes, pipettes, and enriched culture by autoclave, etc. in strict compliance with local, state and federal disposal regulations.

**[Storage and expiration date]**

- Storage: NH immunochromato O111 should be stored at 2–8°C under protection from the light. NH immunochromato O111 must not be frozen.
- Expiration date: 12 months from the date of manufacture.

**[Warranty]**

NH Foods Ltd. makes no warranty, whether expressed or implied. If any materials are defective, Distributor or NH Foods Ltd. will provide a replacement of the product. There is no warranty of merchantability or fitness for a particular purpose. Each user is expected to verify the appropriateness of each operative procedure and the application to individual food materials. NH Foods Ltd. shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

**[References]**

- MHLW Notice No. 1120 Article 3, from the Inspection and Safety Division "Detection Methods for Enterohemorrhagic *E. coli* O26, O103, O111, O121, O145 and O157 from foods"

**[Technical assistance]**

For technical assistance, please contact distributor or manufacturer.

● **Distributor**  
(Labeled on the box)

● **Manufacturer**  
 R&D Center, NH Foods Ltd.  
 3-3 Midorigahara, Tsukuba, Ibaraki 300-2646, Japan  
 Phone: +81-29-847-7825 Fax: +81-29-848-1256  
 E-mail: nh.kit@nipponham.co.jp  
 http://www.rdc.nipponham.co.jp