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References

1 Coronaviridae Study Group of the International Committee on Taxonomy of Viruses. The

Species Severe Acute Respiratory Syndrome-Related Coronavirus: Classifying 2019-nCoV and Naming It SARS-CoV-2. *Nat Microbiol.* Mar 2 2020 (<https://doi.org/10.1038/s41564-020-0695-z>).

2 Lai C et al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *Int J Microb Agents.* Feb 20 2020 (<https://doi.org/10.1016/j.ijantimicag.2020.105924>).

3 Rabenau HF et al. SARS-coronavirus (SARS-CoV) and the safety of a solvent/detergent (S/D) treated immunoglobulin preparation. *Biologicals.* 33(2):95-99 (2005).

4 Darnell MER et al. Evaluation of Inactivation Methods for Severe Acute Respiratory Syndrome Coronavirus in Noncellular Blood Products. *Transfusion.* 46(10), 1770-7 (2006).

5 Howett, MK et al. A Broad-Spectrum Microbicide with Virucidal Activity against Sexually-Transmitted Viruses. *Antimicrob Agents Chem.* 43(2):314-321 (1999).

6 Piret, J et al. Sodium Lauryl Sulfate, a Microbicide Effective Against Enveloped and Nonenveloped Viruses. *Curr Drug Targets.*3(1):17-30 (2002).

7 Rabenau HF et al. Stability and inactivation of SARS coronavirus. *Med Microbiol Immunol.* 194:1-6 (2005).

8 Kariwa H et al. Inactivation of SARS coronavirus by means of povidone-iodine, physical conditions, and chemical reagents. *Jpn J Vet Res.* 52(3):105-112 (2004).

9 Paget A et al. Evaluation of SARS-Coronavirus Decontamination Procedures. *Applied Biosafety.* 12(2):100-8 (2007).

10 International Air Transport Association Dangerous Goods Regulations, 60th Edition. (2019).

TMDNA Genotek is a trademark and [®]ORAgene, OMNIgene and ORAcollect are registered trademarks of DNA Genotek Inc.

[®]OraSure is a registered trademark of OraSure Technologies, Inc.

Some OraSure Technologies, Inc and DNA Genotek products may not be available in all geographic regions, contact your sales representative for details.