TECHNICAL DATA SHEET

SALMONELLA ENRICHMENT

ENRICHMENT FOR *SALMONELLA*

1 INTENDED USE

Salmonella Enrichment is a special formulation of Buffered Peptone Water that has been created and controlled for optimal detection of *Salmonella* in food products and feed.

Salmonella Enrichment with Tween[®] 80 is used as enrichment medium for *Salmonella* analysis of products whose fat content exceeds 20%.

The *Salmonella* Enrichment line complies with NF EN ISO 6579-1 standard, Microbiology of the food chain - Horizontal method for the detection, enumeration and serotyping of *Salmonella* - Part 1 : detection of *Salmonella* spp., NF EN ISO 6887-1 Microbiology of food - General rules for the preparation of the initial suspension and decimal dilutions, Parts 2, 3, 4, & 5. *Salmonella* Enrichment can be used as Buffered Peptone Water in required methods but the inverse does not apply.

Salmonella Enrichment has been specially formulated for the validated methods IRIS *Salmonella*[®] and SESAME *Salmonella* TEST[®].

2 PRINCIPLES

The peptide composition and osmotic balance of *Salmonella* Enrichment medium have been optimized to allow an exceptional resuscitation level of *Salmonella* strains.

3 TYPICAL COMPOSITION

The composition can be adjusted in order to obtain optimal performance.

The *Salmonella* Enrichment formulation conforms to that of Buffered Peptone Water.

For 1 liter of *Salmonella* Enrichment :

- Peptone	10,0 g
- Sodium chloride	5,0 g
- Disodium phosphate, anhydrous	3,56 g
- Monopotassium phosphate	1,5 g

pH if the ready-to-use media at 25 °C : 7,0 ± 0,2.

For 1 liter of *Salmonella* Enrichment + Tween[®] 80 :

- Peptone	10,0 g
- Sodium chloride	5,0 g
- Disodium phosphate, anhydrous	3,56 g
- Monopotassium phosphate	1,5 g
- Tween [®] 80	10,0 g

pH of the ready-to-use media at 25 °C : 7,0 ± 0,2.

4 PREPARATION

- Dissolve 20,0 g of dehydrated media (BK194) in 1 liter of distilled or demineralized water.
- Stir slowly until complete dissolution.
- Dispense into tubes or vials.
- Sterilize in an autoclave at 121 °C for 15 minutes.
- Cool to room temperature.

✓ **Reconstitution :**
20,0 g/L

✓ **Sterilization :**
15 min at 121 °C

5 INSTRUCTIONS FOR USE

- Introduce aseptically **25 g** of the sample to be tested into **225 mL** ready-to-use **Salmonella Enrichment** in order to achieve a 1:10 dilution.
- or
- Introduce aseptically **X g** of the sample to be tested into **9 X mL** ready-to-use **Salmonella Enrichment** in order to respect the 1:10 dilution ratio of 1 part sample + 9 parts of diluent mL.
 - Mix well.
 - Incubate at temperatures and for the periods required by the analytical protocol chosen.

6 QUALITY CONTROL

Dehydrated media : cream-white powder, free-flowing and homogeneous.

Prepared media : amber solution, limpid, may present a slight precipitate after prolonged storage.

Typical culture response (NF EN ISO 11133) :

Microorganisms		Growth
⁽¹⁾ <i>Salmonella</i> Typhimurium	WDCM 00031	Positive, score 2
⁽¹⁾ <i>Salmonella</i> Enteritidis	WDCM 00030	Positive, score 2
⁽¹⁾ <i>Escherichia coli</i>	WDCM 00012	Positive, score 2
⁽²⁾ <i>Listeria monocytogenes</i> 4b	WDCM 00021	± 30 % colonies / T ₀
⁽²⁾ <i>Listeria monocytogenes</i> ½a	WDCM 00109	± 30 % colonies / T ₀
⁽³⁾ <i>Escherichia coli</i>	WDCM 00012	± 30 % colonies / T ₀
⁽³⁾ <i>Staphylococcus aureus</i>	WDCM 00034	± 30 % colonies / T ₀

⁽¹⁾ After 18 hours of incubation at 37 °C (inoculum ≤ 10² microorganisms)

⁽²⁾ After 60 minutes of incubation at 20 °C

⁽³⁾ After 45-60 minutes of incubation at 20-25 °C

7 STORAGE / SHELF LIFE

Dehydrated media : 2-30 °C.

Ready-to-use media in vials or flexible bags : 2-25 °C.

The expiration dates are indicated on the labels.

Prepared media in vials or tubes (*) : 180 days at 2-25 °C.

(*) Benchmark value determined under standard preparation conditions, following manufacturer's instructions.

8 PACKAGING

- **Salmonella Enrichment :**

Dehydrated media :

500 g bottle BK194HA

5 kg drum BK194GC

Ready-to-use media :

10 x 225 mL vials BM13608

3 x 3 L flexible bags BM13708

2 x 5 L flexible bags BM14408

- *Salmonella* Enrichment + Tween® 80 :

Ready-to-use media :

3 x 3 L flexible bags	BM16308
2 x 5 L flexible bags	BM19808

9 BIBLIOGRAPHY

Perry, D.F., and Quiring, C.. 1997. Fundamental aspects of enzyme/chromogenic substrate interactions in agar media formulations for esterase and glycosidase detection in *Salmonella*. In *Salmonella* and Salmonellosis-Proceedings. Ploufragan-France, 63-70.

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IRIS *Salmonella*®, **COMPASS®** and **SESAME *Salmonella* TEST®** are registered trademarks of SOLABIA S.A.S