

## TECHNICAL DATA SHEET

# LISTERIA ENRICHMENT BROTH (MODIFIED UVM I)

### SELECTIVE ENRICHMENT FOR *LISTERIA*

## 1 INTENDED USE

The enrichment broths UVM (University of Vermont) involve a 2 step process, leading to a higher degree of isolation of *Listeria monocytogenes* in meat products (including poultry).

UVM 1 enrichment broth allows for primary enrichment of *Listeria*.

## 2 HISTORY

UVM media were described by Donnelly and Baigent in 1986, based on the Dominguez-Rodriguez *et al.* (1984) enrichment medium for *Listeria*. The difference in acriflavin concentrations in UVM I and UVM II Broths, as well as the two-step enrichment procedure, lead to a very satisfactory recovery of *Listeria monocytogenes*.

In the modified formula currently used, the nalidixic acid content has been reduced from 40 to 20 mg/liter.

## 3 PRINCIPLES

Peptones, yeast extract and meat extract supply the essential nutrients for the growth of *Listeria*.

The sodium chloride concentration increases selectivity.

Phosphates buffer the pH of the medium.

Esculin is hydrolyzed to glucose and esculetin by *Listeria*.

Nalidixic acid blocks the DNA replication of bacteria sensitive to this antibacterial agent.

Accompanying Gram-positive microflora are inhibited by acriflavin.

## 4 TYPICAL COMPOSITION

The composition can be adjusted in order to obtain optimal performance.

For 1 liter of media :

- Proteose peptone .....	5,0 g
- Tryptone .....	5,0 g
- Yeast extract .....	5,0 g
- Bacteriological meat extract .....	5,0 g
- Sodium chloride .....	20,0 g
- Disodium phosphate, anhydrous .....	9,6* g
- Monopotassium phosphate .....	1,35 g
- Esculin .....	1,0 g
- Nadixic acid .....	20,0 mg
- Acriflavin (chlorhydrate) .....	12,0 mg

pH du of the ready-to-use media at °C : 7,4 ± 0,2.

\* 9.6 g of disodium phosphate anhydrous is equal to 12 g dibasic sodium phosphate.

## 5 PREPARATION

- Dissolve 52 g of dehydrated media (BK113) in 1 liter of distilled or demineralized water.
- Stir slowly until complete dissolution.
- Dispense in vials at 225 mL per vial.
- Sterilize in an autoclave at 115°C for 15 minutes.
- Cool to room temperature.

✓ **Reconstitution :**  
52 g/L

✓ **Sterilization :**  
15 min at 115 °C

## 6 INSTRUCTIONS FOR USE

- Aseptically add 25 g of the product to analyze to a tared flask containing 225 mL of media
- Mix well.
- Incubate at 30 ± 1 °C for 24 hours.

✓ **Inoculation :**  
25 g in 225 mL

✓ **Incubation :**  
24 h at 30 °C

## 7 RESULTS

Perform a secondary enrichment in UVM II broth or isolate on a selective media like **COMPASS® Listeria Agar** according to the reference in vigor.

## 8 QUALITY CONTROL

**Dehydrated media :** beige powder, free-flowing and homogeneous.

**Prepared media :** amber yellow solution, with a bluish reflection, limpid.

Typical culture response after 24 hours of incubation at 30 °C, followed by subculture on **COMPASS® Listeria Agar**.

Microorganisms		Growth
<i>Listeria monocytogenes</i> 4b + <i>Enterococcus faecalis</i> + <i>Escherichia coli</i>	WDCM 00021 WDCM 00087 WDCM 00013	> 10 characteristic colonies
<i>Listeria monocytogenes</i> ½ a + <i>Enterococcus faecalis</i> + <i>Escherichia coli</i>	WDCM 00109 WDCM 00087 WDCM 00013	> 10 characteristic colonies
<i>Enterococcus faecalis</i> <i>Escherichia coli</i>	WDCM 00087 WDCM 00013	< 100 colonies Inhibited, score 0

## 9 STORAGE / SHELF LIFE

**Dehydrated media :** 2-30 °C.

The expiration date is indicated on the label.

**Prepared media in vials (\*) :** 30 days at 2-8 °C, shielded from light.

(\*) Benchmark value determined under standard preparation conditions, following manufacturer's instructions.

## 10 PACKAGING

**Dehydrated media :**

500 g bottle ..... BK113HA

## 11 BIBLIOGRAPHY

Dominguez-Rodriguez, L., Suarez-Fernandez, G., Fernandez-Garayzabal, J.F., and Rodriguez-Ferri, E.. 1984. New methodology for the isolation of *Listeria* microorganisms from heavily contaminated environments. Applied and Environmental Microbiology, **47** : 1188-1190.

Donnelly, C.W., and Baigent, G.J.. 1986. Method for flow cytometric detection of *Listeria monocytogenes* in milk. Applied and Environmental Microbiology, **52** : 689-695.

Lee, W.H., and MacClain, D.. 1986. Improved *Listeria monocytogenes* selective agar. Applied and Environmental Microbiology, **52** : 1215-1217.

MacClain, D., and Lee, W.H. 1988. Development of USDA-FSIS method for isolation of *Listeria monocytogenes* from raw meat and poultry. Journal of Association of Official Analytical Chemists, **71** : 660-664.

## **12 ADDITIONAL INFORMATION**

---

The information provided on the labels take precedence over the formulations or instructions described in this document and are susceptible to modification at any time, without warning.

**COMPASS**<sup>®</sup> is a registered trademark of SOLABIA S.A.S

Document code : UVM1 MODIFIED BROTH\_ENv9.  
Creation date : 05-2003  
Updated : 02-2018  
Origin of revision : Modification of typical composition