FUJIFILM

Wako

Code No. 057-09731

# Exosomes, from Mesenchymal Stem Cells, Purified

This product is a highly purified exosome preparation obtained from the culture supernatant of human bone marrow-derived mesenchymal stem cells (MSC) using the PS affinity method. The purified exosomes have been confirmed to be positive not only for the exosome markers CD9, CD63, and CD81 but also for CD59 and CD73, which are known markers of MSC-derived exosomes. These can be used as positive control exosomes for the above markers in qualitative analyses such as Western blotting, flow cytometry, and ELISA.*1
Furthermore, when measuring exosomes using anti-CD63 antibodies with the PS Capture™ Exosome ELISA Kit series or the CD63-Capture Human Exosome ELISA Kit (Streptavidin HRP), this product can be used as a calibrator to generate a standard curve. This allows conversion of sample absorbance measurements to concentrations and helps to correct measurement variability between assays.
This product is for research use only.

*1 Expression of CD59 and CD73 has been confirmed only by Western blotting.

## 【Storage】

Store at -80℃ *1
It is recommended to aliquot and store the product in appropriate volumes.
Example : Dispense the product into 10 tubes of 5 μL each, store frozen at -80℃, and thaw at room temperature before use.

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10 μg/mL*2

*2 The concentration of this product is a protein concentration adjusted by CD63 signal values to minimize lot-to-lot variation when used as a positive control and calibrator.

## 【Precaution for use】

This product consists of exosomes derived from MSCs obtained from human bone marrow. The MSCs used have been confirmed negative for human immunodeficiency virus (HIV), hepatitis C virus (HCV), hepatitis B virus (HBV), and Treponema pallidum. However, tests for other pathogens have not been performed. Therefore, please handle this product as a specimen with potential infection risk.

## 【Procedure】

**〈Instructions for Use as Positive Control Exosomes for Human CD9, CD63, CD81, CD59, and CD73〉**

Dilute this product to the desired concentration and use it as a positive control in qualitative analyses such as Western blotting*3, flow cytometry, and ELISA.

*3 For Western blotting, the use of non-reducing samples is recommended. Please also refer to the following for the recommended antibodies.

■ Anti-CD9/CD63/CD81 Antibodies

We recommend using our company's antibodies. The available lineup can be found at:
https://labchem-wako.fujifilm.com/us/category/01115.html

■ Anti-CD59/CD73 Antibodies

The antibodies confirmed to have detection performance in our company are as follows:
・Anti-CD59 Antibody [MEM-43] (Abcam #ab9182)
・CD73 Monoclonal Antibody (1D7) (Thermo Fisher Scientific #MA5-15537)

**〈Instructions for Use as a Calibrator in PS Capture™ Exosome ELISA Kits and CD63-Capture Human Exosome ELISA Kit〉**

This product can be used as a calibrator when measuring exosomes with anti-CD63 antibody as the detection antibody in the following ELISA kits:

ELISA Kits :
① 297-79201 PS Capture™ Exosome ELISA Kit (Anti-Mouse IgG POD)
② 298-80601 PS Capture™ Exosome ELISA Kit (Streptavidin HRP)
③ 290-83601 CD63-Capture Human Exosome ELISA Kit (Streptavidin HRP)

By generating a standard curve with this product, sample absorbance readings can be converted to concentrations, allowing correction for measurement variability across assays. Prepare the standard curve by diluting this product to concentrations that yield ELISA absorbance at 450 nm (subtracting the absorbance at 620–650 nm as the reference wavelength) within the range of 0.1 to 3.0 O.D., using the dilution buffer specified by the kit.

Example of Use :
To prepare a standard curve for the PS Capture™ Exosome ELISA Kit (Streptavidin HRP), perform a two-fold serial dilution series over seven steps, starting from a concentration of 80 ng/mL :

(1) Mix 4 μL of this product with 496 μL of the Reaction Buffer supplied with the kit to prepare a 125-fold dilution (80 ng/mL).
(2) Prepare the six-point two-fold serial dilution series from 80 ng/mL by mixing with Reaction Buffer at the proportions shown in the table below.

| Final Conc. (ng/mL) | Starting Conc. (ng/mL) | Volume to dilute (μL) | Reaction Buffer (μL) |
|---|---|---|---|
| 80.0 | 10,000 | 4 | 496 |
| 40.0 | 80.0 | 250 | 250 |
| 20.0 | 40.0 | 250 | 250 |
| 10.0 | 20.0 | 250 | 250 |
| 5.0 | 10.0 | 250 | 250 |
| 2.5 | 5.0 | 250 | 250 |
| 1.25 | 2.5 | 250 | 250 |
| 0.00 (Blank) | – | – | 250 |

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