


SHIBAYAGI'S ALBUMIN ELISA KITS & TIA REAGENTS FOR PHARMACEUTICAL SAFETY & MONITORING

RESEARCH REAGENTS

Model Organism Specific Albumin Reagents

Drug safety plays a pivotal role during the collection, detection, assessment, monitoring and prevention of adverse effects in the manufacturing process of pharmaceutical products, involving the identification of hazards and minimizing risk of any harm that may come to patients.

FUJIFILM Wako Chemicals Europe GmbH provides a unique lineup of ELISA kits for evaluation of urine, serum and plasma albumin for rat and mouse model organism. For high throughput analysis TIA reagents for the application on clinical chemistry autoanalyser are available.



ELISA = Enzyme-Linked Immunosorbent Assay; TIA = Turbidimetric Immuno Assay

Features

- Reacts with rat or mouse albumin specifically
- Small sample volume (max. 7 μ L) (TIA reagents)
- Rapid assay (10 min) with automatic analyzers (TIA reagents)
- Easy ready-to-use reagents (TIA reagents)
- Wide assay range and excellent reproducibility
- Little influence by hemolysis

Ordering Information

Wako Code	Product name	Application	Assay range	Total reaction time
638-04309	LBIS Mouse Albumin ELISA Kit	Mouse Serum, Plasma, Urine	50-1,000 ng/mL	2 h 20 min
635-04319	LBIS Rat Albumin ELISA Kit	Rat Serum, Plasma, Urine	50-1,000 ng/mL	2 h 20 min
635-07099	LBIS Bovine Albumin ELISA Kit	Media	0.78 - 50 ng/mL	2 h 30 min
For application on clinical chemistry autoanalysers				
638-25309	LBIS Rat Urinary Albumin Assay Kit (S-type) (TIA)	Rat Serum, Plasma, Urine	6.17-500 µg/mL	10 min
632-25569	LBIS Mouse Urinary Albumin Assay Kit (S-type) (TIA)	Mouse Serum, Plasma, Urine	6.17-500 µg/mL	10 min

Related Products

Wako Code	Product name	Application	Assay range	Total reaction time
383-14491	High Sensitivity Human L-FABP ELISA Kit	Human Urine	0.3-60 ng/mL	3 h
386-14501	Rat L-FABP ELISA Kit	Rat Urine	1.56-200 ng/mL	3 h
383-14511	Primate L-FABP ELISA Kit	Primate Urine	2.6-525 ng/mL	3 h

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako; <https://labchem-wako.fujifilm.com> / This leaflet may contain products that can not be exported to your country due to regulations. / Bulk quote requests are welcome. Please contact us.

Manufacturer

FUJIFILM Wako Shibayagi Corporation

1062-1, Ishihara, Shibukawa, Gunma Pref., 377-0007, Japan

Telephone: +81-279-25-0279

Fax: +81-279-23-0313

E-mail: wksb-info@fujifilm.com

Distribution & Inquiries

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstr. 12 • 41468 Neuss, Germany

Telephone: +49-2131-311-271

Fax: +49-2131-311-110

E-mail: labchem_wkeu@fujifilm.com



<https://labchem-wako.fujifilm.com>

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome Chuo-Ku, Osaka 540-8605, Japan

Telephone: +81-6-6203-3741

Fax: +81-6-6203-1999

E-mail: ffwk-cservice@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road Richmond • VA 23237, USA

Telephone: +1-804-714-1920

Fax: +1-804-271-7791

E-mail: wkuslabchem@fujifilm.com