2024/05/17 Distributed by FUJIFILM Wako Pure Chemical Corporation Manufactured by FUJIFILM Wako Pure Chemical Corporation

SPECIFICATION

Distilled Water

for Amino Acid Sequence Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless, clear liquid
Absorbance (200~220nm)	max.0.02
Absorbance (230~400nm)	max.0.01
Chloride	to pass test(Cl : max. about0.1ppm)
Ammonium (NH4)	max.0.03ppm
Substances reducing permanganate	to pass test