

## SPECIFICATION

## Ethyl Acetate

for SHIMADZU Protein Peptide Sequencer

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (260nm)	max.0.200
Absorbance (270nm)	max.0.050
Absorbance (280nm)	max.0.015
Absorbance (300~400nm)	max.0.010
Water	max.0.02%
Aldehydes (as CH <sub>3</sub> CHO)	max.1ppm
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.1ppm
Suitability for amino acids sequence analysis	to pass test
Assay (capillary GC)	min.99.0%