

(90 × 190mm Size)

FUJIFILM

Wako

Code No. 013-22831 (50 μ L)

Anti Human Atg7, Rabbit

抗ヒト Atg7, ウサギ

Rabbit Anti Human Atg7 antiserum was raised against the synthetic peptides corresponding to the amino acid residues 556-571 (VVAPGDSTRDRTLDDQQ) of human Atg7. The synthetic peptides were chemically coupled with KLH and immunized in rabbits. The antiserum reacts with human, rat and mouse Atg7. For Western blotting, the antiserum can be diluted 1,000-fold in 1% BSA dissolved in 20 mmol/L Tris-HCl (pH 7.5)-0.15 mol/L NaCl-0.1% Tween 20-0.1% NaN_3 (TTBS- NaN_3). The diluted antibody solution can be used in Western blotting for at least 6 months. Inclusion of 0.1% Tween 20 is recommended to obtain lower background when the antibody reactivity is detected in chemiluminescence.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

Preparation : Antiserum contains no preservatives and stabilizers.

Specificity : Specific for human, rat and mouse.

Working Concentration : Western blotting 1 : 1,000

Diluent : 1% BSA dissolved in 20 mmol/L Tris-HCl (pH 7.5)-0.15 mol/L NaCl-0.1% Tween 20-0.1% NaN_3 (TTBS- NaN_3)

Storage : Store at -20°C .
After opening aliquot contents and freeze at -20°C .

Package : 50 μ L

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation
1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH
Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

1802KA1