

FUJIFILM

Wako

Code No. 015-27771 (20 μ L)
011-27773 (100 μ L)

Anti CD81, Monoclonal Antibody (17B1)

The 17B1 monoclonal antibody specifically binds to CD81, which is 26 kDa cell surface protein and is known as TAPA-1 (Target of the Antiproliferative Antibody 1) and Tetraspanin-28 (Tspan-28). CD81 interacts directly with immunoglobulin superfamily member 8 (IGSF8, CD316) and CD36 and forms a signal transduction complex with CD19, CD21 and Leu-13 (CD225) on the surface of B cell. In addition, on T cells CD81 associates with CD4 and CD8 and provides a costimulatory signal with CD3.

This product is intended for laboratory use only.

[Concentration]

Indicated on product label.

[Formulation]

1 \times TBS aqueous solution with 50% glycerol, 0.05w/v% sodium azide, pH 7.4

[Clone No.]

17B1

[Cross Reactivity]

Human/Bovine CD81

Note : This antibody does not recognize mouse and rat CD81.

[Subclass]

IgG2b

[Application]

Western blotting ^{※1}	1 : 1,000 – 1 : 8,000 (under non-reducing condition)
ELISA	1 : 2,000 – 1 : 16,000
Flow Cytometry	1 : 250 – 1 : 1,000
Immunoprecipitation	5 – 50 μ g/assay

Please optimize the most appropriate concentration for your analysis and samples.

※ 1 Use this antibody under non-reducing condition.

[Storage]

Store at -20°C.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>