

FUJIFILM**Wako**Code No. 017-27591 (10 μ L)
013-27593 (50 μ L)**Anti Human Iba1, Monoclonal
Antibody (NCNP27)****抗ヒトIba1, モノクローナル抗体
(NCNP27)****Background**

Iba1 is a protein highly expressed in microglia and macrophages with a molecular weight of about 16.7kDa¹⁾. The protein is a commonly known microglial marker in the nervous system. This item is a mouse monoclonal antibody that reacts with human Iba1. For research use only. Not for use in diagnostic procedure.

Description

[Purification] Purified from the hybridoma supernatant by Protein A-affinity chromatography

[Reactivity] Reacts with human Iba1

[Antigen] Synthetic peptide corresponding to the C-terminus of Iba1

[Clone No.] NCNP27

[Species cross reactivity] Human
(Other species have not been tested)

[Host] Mouse

[Isotype] IgG2b

[Concentration] Indicated on the label

[Formulation] 1 \times TBS (pH 7.4) aqueous solution with 50% glycerol containing 0.05% sodium azide.

Applications

Immunohistochemistry (paraffin section, DAB) 1 : 50~200

Optimal concentration should be determined by each laboratory for each application.

Storage

Store at -20 $^{\circ}$ C.
Avoid repeated freeze and thaw.

Package

Code No.	Packaging
017-27591	10 μ L
013-27593	50 μ L

- 1/2 -

Reference

- 1) Imai, Y., Ibata, I., Ito, D., Ohsawa, K. and Kohsaka, S. :
Biochem. Biophys. Res. Commun., **224**, 855 (1996).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : +81-6-6203-3741
Facsimile : +81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : +49-2131-311-0
Facsimile : +49-2131-311100
<http://www.wako-chemicals.de>

1904KA1

- 2/2 -