

CultureSure® Albumin, from Bovine Serum(BSA)

This is a bovine serum albumin (BSA) product for cell culture. It has been tested for mycoplasma test, bacterial endotoxins, and viable counts test, so it can be used safely for cell culture.

Features

- Has passed mycoplasma test, bacterial endotoxins test, and viable count test.
- Globulin-, fatty acids-, and protease-free

Product Information

Cat. No.	Product Name	Specifications	Package Size
034-25462	CultureSure® Albumin, from Bovine Serum(BSA)	For cell culture	25 g

Physical Properties

Requirement	Specifications
Appearance	White - pale yellowish brown or slightly greenish yellow, crystals - powder
Solubility in water	to pass test
Loss on drying over silica gel in vacuum	under investigation
pH (100g/L, 25°C)	6.5-7.5
Bacterial endotoxins	max. 3.0 EU/mg
Globulin	to pass test
Mycoplasma test	to pass test
Free fatty acids	to pass test
Protease	to pass test
Assay (electrophoresis)	min. 98%
Viable count test	max. 1,000 CFU/g

Bovine serum albumin products list

We have a product line of bovine serum albumin for selection by application, specification, and manufacturing method.

Purification method	Product Name	Cat. No.	Package Size	Oceania	Protease-free	Fatty acid-free	Globulin-free	Low-globulin	Bacterial endotoxins test	Low-sodium	RIA	pH5.2
Heat treatment/Solvent treatment	CultureSure® Albumin, from Bovine Serum (BSA)	034-25462	25 g		○	○	○		○			
	Albumin, from Bovine Serum (BSA), Protease Free	014-15151	5 g		○							
		010-15153	10 g									
		018-15154	50 g									
	Albumin, from Bovine Serum (BSA), Globulin Free	019-15101	10 g				○					
		015-15103	50 g									
		013-15104	100 g									
	Albumin, from Bovine Serum (BSA), Globulin Free-HG (High Grade)	016-15111	10 g				○					
		012-15113	50 g									
		010-15114	100 g									
	Albumin, from Bovine Serum (BSA), Fatty Acid Free	017-15146	5 g									
		017-15141	10 g			○						
		013-15143	50 g									
		011-15144	100 g									
	Albumin, from Bovine Serum, Protease Free, Low IgG	019-28391	10 g		○			○				
		015-28393	50 g									
	Albumin, from Bovine Serum, Fatty Acid/Protease Free, Low IgG	012-28381	10 g		○	○		○				
		018-28383	50 g									
	Albumin, from Bovine Serum (BSA), Low Salt	013-15121	10 g								○	
		019-15123	50 g									
		017-15124	100 g									
		015-15125	500 g									
	Albumin, from Bovine Serum (BSA), RIA Grade	010-15131	10 g								○	
		016-15133	50 g									
		014-15134	100 g									
	Albumin, from Bovine Serum (BSA), pH5.2(Fraction V)	011-21271	1 g									○
017-21273		10 g										
019-21272		25 g										
015-21274		100 g										
Albumin, from Bovine Serum (BSA), Crystallized	016-15091	1 g										
	016-15096	5 g										
	012-15093	10 g										
	010-15094	50 g										
Chromato graphy	Albumin, from Bovine Serum (BSA), pH7.0, New Zealand Origin	012-23381	5 g	○								
		010-23382	25 g									

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site
<https://labchem-wako.fujifilm.com>



FUJIFILM Wako Pure Chemical Corporation
 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
 Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
 ffwk-cservice@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation
 1600 Bellwood Road, Richmond, VA 23237, U.S.A.
 Toll-Free (U.S. only): +1 877 714 1920
 Tel: +1 804 271 7677 Fax: +1 804 271 7791
 wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH
 Fuggerstr 12, 41468 Neuss, Germany
 Tel: +49 2131 311 0 Fax: +49 2131 311 100
 labchem_wkeu@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited
 Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road,
 Tsuen Wan, N.T., Hong Kong
 Tel: +852-2799-9019 Fax: +852-2799-9808
 wkhk.info@fujifilm.com

FUJIFILM Wako (Guangzhou) Trading Corporation
 Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road,
 Guangzhou, 510095, China
 Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai)
 Tel: +86-10-6413-6388(Beijing)
 wkgz.info@fujifilm.com