

Catalyst

Flow Reaction Catalyst Pd/C, Beads


Features

- Pd/C with a particle size of about 200 μm
- Easy to use for flow reaction using heterogeneous catalyst
- Low flammability and safe
- Bulk supply possible
- No binder used




Applications


<Reduction of olefin>




<Reduction of Ar-NO₂>



<Deprotection of O-Bn>



Comparison of back pressure




Code No.	Product Name	Grade	Package Size
169-28861	Pd/C, Beads (Pd 5%)	for Organic Synthesis	5g
167-28862			25g

<Similar products>

Code No.	Product Name	Grade	Package Size
162-28851	Pt/C, Beads (Pt 5%)	for Organic Synthesis	5g
160-28852			25g

<Related products>

Hydrogenation of aromatic rings under mild conditions (50-80°C, 1 atm)



Code No.	Product Name	Grade	Package Size
186-03451	Rh-Pt/(DMPSi-Al ₂ O ₃)	for Organic Synthesis	1g
182-03453			5g

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site
<https://labchem-wako.fujifilm.com>



FUJIFILM Wako Pure Chemical Corporation
1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
ffwk-cserservice@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation
1600 Bellwood Road, Richmond, VA 23237, U.S.A.
Toll-Free (U.S. only): +1 877 714 1920
Tel: +1 804 271 7677 Fax: +1 804 271 7791
wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited
Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong
Tel: +852-2799-9019 Fax: +852-2799-9808
wkhk.info@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH
Fuggerstr 12, 41468 Neuss, Germany
Tel: +49 2131 311 0 Fax: +49 2131 311 100
labchem_wkeu@fujifilm.com

FUJIFILM Wako (Guangzhou) Trading Corporation
Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road, Guangzhou, 510095, China
Tel: +86-20-8732-6381 (Guangzhou) Tel: +86-21-6288-4751 (Shanghai)
Tel: +86-10-6413-6388 (Beijing)
wkgz.info@fujifilm.com