

SPECIFICATION

Methanol

for the Japanese Pharmacopoeia General Tests (for LC)

REQUIREMENT	SPECIFICATION
Description (JP)	Colorless clear liquid
Absorbance (210nm) (JP)	not more than 0.70
Absorbance (220nm) (JP)	not more than 0.30
Absorbance (230nm) (JP)	not more than 0.15
Absorbance (240nm) (JP)	not more than 0.07
Absorbance (254nm) (JP)	not more than 0.02
* Appearance	Colorless clear liquid
* Density (20°C)	0.789~0.792g/mL
* Refractive index 20°C	1.327~1.330
* Absorbance (210nm)	max.0.60
* Absorbance (220nm)	max.0.30
* Absorbance (230nm)	max.0.15
* Absorbance (240nm)	max.0.06
* Absorbance (254nm)	max.0.015
* Absorbance (260~400nm)	max.0.01
* Water	max.0.05%
* Residue after evaporation	max.5ppm
* Acidity (as HCOOH)	max.0.001%
* Peroxides (as H ₂ O ₂)	max.5ppm
* Fluorescence test	to pass test
* Assay (capillary GC)	min.99.7%
Description (USP-NF)	Clear, colorless liquid
Assay (USP-NF)	minimum 99.8%
Substances darkened by sulfuric acid (USP-NF)	passes test
Substances reducing permanganate (USP-NF)	passes test
Solubility in water (USP-NF)	passes test
Color (APHA) (USP-NF)	maximum 10
Water (H ₂ O) (USP-NF)	maximum 0.1%
Residue after evaporation (USP-NF)	maximum 0.001%

Carbonyl compounds (USP-NF)	maximum 0.001% each of acetone, formaldehyde, and acetaldehyde
Titration acid (USP-NF)	maximum 0.0003meq/g
Titration base (USP-NF)	maximum 0.0002meq/g
Absorbance (280~400nm) (USP-NF)	maximum 0.01
Absorbance (260nm) (USP-NF)	maximum 0.04
Absorbance (240nm) (USP-NF)	maximum 0.10
Absorbance (230nm) (USP-NF)	maximum 0.20
Absorbance (220nm) (USP-NF)	maximum 0.40
Absorbance (210nm) (USP-NF)	maximum 0.80
Absorbance (205nm) (USP-NF)	maximum 1.00
Gradient elution (USP-NF)	passes test
Description (Ph.Eur.)	Clear, colourless liquid
Sp. Gr.(20/20°C) (Ph.Eur.)	0.791~0.793
Absorbance (210nm) (Ph.Eur.)	maximum 0.70
Absorbance (220nm) (Ph.Eur.)	maximum 0.30
Absorbance (230nm) (Ph.Eur.)	maximum 0.13
Absorbance (250nm) (Ph.Eur.)	maximum 0.02
Absorbance (260nm) (Ph.Eur.)	maximum 0.01
Content (Ph.Eur.)	minimum 99.8%
Absorbance (225nm) (Ph.Eur.)	maximum 0.17

This product is suited for Methanol for liquid chromatography of general tests on The Japanese Pharmacopoeia. (*:Additional test performed by Wako)USP-NF Reagent Methanol, SpectrophotometricPh. Eur. Reagent Methanol, Methanol R1 and Methanol R2