




## QuickGene DNA Tissue Kit S (DT-S)

### Composition/Information on ingredients

Name of substance	Classification acc. to GHS	Pictograms
Proteinase K EDT-01	Resp. Sens. 1A / H334 Skin Sens. 1A / H317	
Lysis Buffer LDT-01	Acute Tox. 4 / H302 Skin Irrit. 2 / H315 Eye Irrit. 2 / H319	
Tissue lysis Buffer MDT-01	Eye Irrit. 2 / H319	
Wash Buffer WDT-02		
Elution Buffer CDT-01		

## Proteinase K EDT-01

Version number: 2.0  
Replaces version of: 2015-07-20

Revision: 2017-08-10  
First version: 20.07.2017

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

<b>Trade name</b>	<u>Proteinase K EDT-01</u>
<b>Product number</b>	EDT-01
<b>Registration number (REACH)</b>	not relevant (mixture)
<b>CAS number</b>	not relevant (mixture)

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

<b>Relevant identified uses</b>	Laboratory chemicals
---------------------------------	----------------------

#### 1.3 Details of the supplier of the safety data sheet

KURABO INDUSTRIES LTD.	Telephone: ++81-72-820-3079
Bio-Medical department	Telefax: ++81-72-820-3095
Neyagawa Techno Center, 14-5, Shimokida-cho	
Neyagawa, Osaka 572-0823 Japan	

#### Additional information

Supplier of the product					
Country	Name	Postal code/city	Telephone	Telefax	Website
Germany	Wako Chemicals GmbH	D-41468 Neuss	+49 (0) 2131 - 311-0	+49 (0) 2131 - 311-100	

**e-mail (competent person)** sdb@csb-online.de

Please do not use this e-mail address to ask for the latest safety data sheet. For this purpose contact Wako Chemicals GmbH.

#### 1.4 Emergency telephone number

As above or next toxicological information centre.

# Proteinase K EDT-01

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

#### Classification according to Regulation (EC) No 1272/2008 (CLP)

Classification acc. to GHS				
Section	Hazard class	Category	Hazard class and category	Hazard statement
3.4R	respiratory sensitisation	1A	Resp. Sens. 1A	H334
3.4S	skin sensitisation	1A	Skin Sens. 1A	H317

for full text of abbreviations: see SECTION 16

### 2.2 Label elements

#### Labelling according to Regulation (EC) No 1272/2008 (CLP)

**Signal word** danger

#### Pictograms

**GHS08**



#### Hazard statements

**H317** May cause an allergic skin reaction.

**H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled.

#### Precautionary statements

**P280** Wear protective gloves/protective clothing/eye protection/face protection.

**P303+P361+P353** IF ON SKIN (or hair): take off immediately all contaminated clothing. Rinse skin with water/shower.

**P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing.

**P308+P311** IF exposed or concerned: Call a POISON CENTER/doctor.

**Hazardous ingredients for labelling** Proteinase, Tritirachium album serine

### 2.3 Other hazards

There is no additional information.

#### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

# Proteinase K EDT-01


## SECTION 3: Composition/information on ingredients

### 3.1 Substances

not relevant (mixture)

### 3.2 Mixtures

#### Description of the mixture

Hazardous ingredients acc. to GHS					
Name of substance	Identifier	Wt%	Classification acc. to GHS	Pictograms	M-Factors
1,2,3-propanetriol	CAS No 56-81-5  EC No 200-289-5	40 – 60			
Proteinase K	CAS No 39450-01-6  EC No 254-457-8  Index No 647-014-00-9	1 – < 2.5	Skin Irrit. 2 / H315 Eye Irrit. 2 / H319 Resp. Sens. 1A / H334 Skin Sens. 1A / H317 STOT SE 3 / H335		

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General notes

Take off immediately all contaminated clothing.

In all cases of doubt, or when symptoms persist, seek medical advice.

#### Following inhalation

Provide fresh air.

If breathing is irregular or stopped, immediately seek medical assistance and start first aid actions.

#### Following skin contact

Wash with plenty of soap and water.

After contact with skin, take off immediately all contaminated clothing, and wash immediately with plenty of water.

#### Following eye contact

Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

# Proteinase K EDT-01

---

## **Following ingestion**

Rinse mouth. Do not induce vomiting.  
Get medical advice/attention if you feel unwell.

## **Notes for the doctor**

none

## **4.2 Most important symptoms and effects, both acute and delayed**

These information are not available.

## **4.3 Indication of any immediate medical attention and special treatment needed**

none

## **SECTION 5: Firefighting measures**

### **5.1 Extinguishing media**

#### **Suitable extinguishing media**

water spray, alcohol resistant foam, fire extinguishing powder, carbon dioxide (CO<sub>2</sub>)

#### **Unsuitable extinguishing media**

water jet

### **5.2 Special hazards arising from the substance or mixture**

Hazardous decomposition products: Section 10.

#### **Hazardous combustion products**

carbon monoxide (CO), carbon dioxide (CO<sub>2</sub>), hydrogen chloride (HCl)

### **5.3 Advice for firefighters**

In case of fire and/or explosion do not breathe fumes.  
Co-ordinate firefighting measures to the fire surroundings.  
Do not allow firefighting water to enter drains or water courses.  
Collect contaminated firefighting water separately.  
Fight fire with normal precautions from a reasonable distance.

#### **Special protective equipment for firefighters**

use suitable breathing apparatus

# Proteinase K EDT-01

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Remove persons to safety.

Ventilate affected area.

Wearing of suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing.

#### For emergency responders

Wear breathing apparatus if exposed to vapours/dust/spray/gases.

### 6.2 Environmental precautions

Keep away from drains, surface and ground water.

Retain contaminated washing water and dispose of it.

### 6.3 Methods and material for containment and cleaning up

#### Advices on how to clean up a spill

Collect spillage.

Absorbent material (e.g. sand, diatomaceous earth, acid binder, universal binder, sawdust, etc.).

#### Appropriate containment techniques

Use of adsorbent materials.

#### Other information relating to spills and releases

Place in appropriate containers for disposal.

Ventilate affected area.

### 6.4 Reference to other sections

Hazardous combustion products: see section 5.

Personal protective equipment: see section 8.

Incompatible materials: see section 10.

Disposal considerations: see section 13.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

#### Measures to prevent fire as well as aerosol and dust generation

Use local and general ventilation.

Keep away from sources of ignition - No smoking.

#### Specific notes/details

None.

# Proteinase K EDT-01

## Measures to protect the environment

Avoid release to the environment.

## Advice on general occupational hygiene

Do not eat, drink and smoke in work areas.

Wash hands after use.

Preventive skin protection (barrier creams/ointments) is recommended.

Remove contaminated clothing and protective equipment before entering eating areas.

## 7.2 Conditions for safe storage, including any incompatibilities

### Flammability hazards

None.

### Incompatible substances or mixtures

Incompatible materials: see section 10.

### Protect against external exposure, such as

frost

### Consideration of other advice

Keep away from food, drink and animal feedingstuffs.

### Ventilation requirements

Provision of sufficient ventilation.

### Packaging compatibilities

Keep only in original container.

## 7.3 Specific end use(s)

No information available.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

Occupational exposure limit values (Workplace Exposure Limits)							
Country	Name of agent	CAS No	Notation	Identifier	TWA [mg/m <sup>3</sup> ]	STEL [mg/m <sup>3</sup> ]	Source
GB	glycerol	56-81-5	mist	WEL	10		EH40/2005

#### Notation

mist as mists

STEL short-term exposure limit: a limit value above which exposure should not occur and which is related to a 15-minute period unless otherwise specified

TWA time-weighted average (long-term exposure limit): measured or calculated in relation to a reference period of 8 hours time-weighted average

# Proteinase K EDT-01

Relevant DNELs of components of the mixture						
Name of substance	CAS No	Endpoint	Threshold level	Protection goal, route of exposure	Used in	Exposure time
1,2,3-propanetriol	56-81-5	DNEL	229 mg/kg	human, oral	consumer (private households)	chronic - systemic effects
1,2,3-propanetriol	56-81-5	DNEL	56 mg/m <sup>3</sup>	human, inhalatory	worker (industry)	chronic - local effects
1,2,3-propanetriol	56-81-5	DNEL	33 mg/m <sup>3</sup>	human, inhalatory	consumer (private households)	chronic - local effects

Relevant PNECs of components of the mixture				
Name of substance	CAS No	Endpoint	Threshold level	Environmental compartment
1,2,3-propanetriol	56-81-5	PNEC	0.885 mg/l	freshwater
1,2,3-propanetriol	56-81-5	PNEC	0.0885 mg/l	marine water
1,2,3-propanetriol	56-81-5	PNEC	1,000 mg/l	sewage treatment plant (STP)
1,2,3-propanetriol	56-81-5	PNEC	3.3 mg/kg	freshwater sediment
1,2,3-propanetriol	56-81-5	PNEC	0.33 mg/kg	marine sediment
1,2,3-propanetriol	56-81-5	PNEC	0.141 mg/kg	soil
1,2,3-propanetriol	56-81-5	PNEC	8.85 mg/l	water

## 8.2 Exposure controls

### Appropriate engineering controls

General ventilation.

### Individual protection measures (personal protective equipment)

#### Eye/face protection

Wear eye/face protection.

#### Hand protection

Wear suitable gloves.

Chemical protection gloves are suitable, which are tested according to EN 374.

Check leak-tightness/impermeability prior to use.

In the case of wanting to use the gloves again, clean them before taking off and air them well.

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.



# Proteinase K EDT-01

## Respiratory protection

In case of inadequate ventilation wear respiratory protection.

## Environmental exposure controls

Use appropriate container to avoid environmental contamination.  
Keep away from drains, surface and ground water.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

#### Appearance

Physical state	liquid
Form	fluid
Colour	colourless
Odour	odourless
Odour threshold	these information are not available

#### Other safety parameters

pH (value)	7.5
Melting point/freezing point	these information are not available
Initial boiling point and boiling range	(unknown)
Flash point	these information are not available
Evaporation rate	these information are not available
Flammability (solid, gas)	not relevant (fluid)

#### Explosive limits

Lower explosion limit (LEL)	these information are not available
Upper explosion limit (UEL)	these information are not available
Vapour pressure	these information are not available
Density	1.1 g/cm <sup>3</sup>
Vapour density	these information are not available
Relative density	these information are not available

#### Solubility(ies)

Water solubility	miscible in any proportion
------------------	----------------------------

# Proteinase K EDT-01

## Partition coefficient

n-octanol/water (log KOW)	these information are not available
Auto-ignition temperature	these information are not available
Relative self-ignition temperature for solids	not relevant (Fluid)
Decomposition temperature	these information are not available

## Viscosity

Kinematic viscosity	these information are not available
Dynamic viscosity	these information are not available
Explosive properties	not explosive
Oxidising properties	shall not be classified as oxidising

## 9.2 Other information

None

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material is not reactive under normal ambient conditions.

### 10.2 Chemical stability

The material is stable under normal ambient and anticipated storage and handling conditions of temperature and pressure.

### 10.3 Possibility of hazardous reactions

No known hazardous reactions.

### 10.4 Conditions to avoid

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

### 10.5 Incompatible materials

oxidisers

### 10.6 Hazardous decomposition products

Reasonably anticipated hazardous decomposition products produced as a result of use, storage, spill and heating are not known.

Hazardous combustion products: see section 5.

# Proteinase K EDT-01

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### Classification procedure

If not otherwise specified the classification is based on:  
Ingredients of the mixture (additivity formula).

#### Classification according to GHS (1272/2008/EC, CLP)

#### Acute toxicity

Acute toxicity of components of the mixture					
Name of substance	CAS No	Exposure route	Endpoint	Value	Species
1,2,3-propanetriol	56-81-5	oral	LD50	27,200 mg/kg	rat, female
1,2,3-propanetriol	56-81-5	oral	LD50	23,000 mg/kg	mouse, male
1,2,3-propanetriol	56-81-5	oral	LD50	≥10,000 mg/kg	guinea pig
1,2,3-propanetriol	56-81-5	dermal	LD50	56,750 mg/kg	guinea pig
1,2,3-propanetriol	56-81-5	inhalation: dust/mist	LC0	>11 mg/l/1h	rat, male
1,2,3-propanetriol	56-81-5	inhalation: dust/mist	LC50	>2.75 mg/l/4h	rat, male

#### Skin corrosion/irritation

Shall not be classified as corrosive/irritant to skin.

#### Serious eye damage/eye irritation

Shall not be classified as seriously damaging to the eye or eye irritant.

#### Respiratory or skin sensitisation

##### Skin sensitisation

May cause an allergic skin reaction.

##### Respiratory sensitisation

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

#### Germ cell mutagenicity

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

#### Carcinogenicity

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

# Proteinase K EDT-01

## Reproductive toxicity

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Specific target organ toxicity - single exposure

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Specific target organ toxicity - repeated exposure

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Aspiration hazard

Shall not be classified as presenting an aspiration hazard.

## SECTION 12: Ecological information

### 12.1 Toxicity

#### Aquatic toxicity (acute)

Test data are not available for the complete mixture.

#### Aquatic toxicity (acute) of components of the mixture

Aquatic toxicity (acute) of components of the mixture					
Name of substance	CAS No	Endpoint	Value	Species	Exposure time
1,2,3-propanetriol	56-81-5	LC50	54,000 mg/l	rainbow trout (Oncorhynchus mykiss)	96 h
1,2,3-propanetriol	56-81-5	LC50	1,955 mg/l	daphnia magna	48 h

#### Aquatic toxicity (chronic)

Test data are not available for the complete mixture.

### 12.2 Persistence and degradability

#### Degradability of components of the mixture

Degradability of components of the mixture					
Name of substance	CAS No	Process	Degradation rate	Time	Source
1,2,3-propanetriol	56-81-5	DOC removal	94 %	1 d	ECHA

# Proteinase K EDT-01

## Biodegradation

The relevant substances of the mixture are readily biodegradable.

## Persistence

Data are not available.

## 12.3 Bioaccumulative potential

Data are not available.

### Bioaccumulative potential of components of the mixture

Bioaccumulative potential of components of the mixture			
Name of substance	CAS No	BCF	Log KOW
1,2,3-propanetriol	56-81-5		-1.75 (pH value: 7.4, 25 °C)

## 12.4 Mobility in soil

Data are not available.

## 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

## 12.6 Other adverse effects

Data are not available.

### Endocrine disrupting potential

None of the ingredients are listed.

### Remarks

Water hazard class - WHC (Wassergefährdungsklasse): 1 (Slightly hazardous to water)

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

This material and its container must be disposed of as hazardous waste.

#### Sewage disposal-relevant information

Do not empty into drains.

#### Waste treatment of containers/packagings

Handle contaminated packages in the same way as the substance itself.

### Remarks

Please consider the relevant national or regional provisions.

# Proteinase K EDT-01

## SECTION 14: Transport information

- 14.1 UN number** not subject to transport regulations
- 14.2 UN proper shipping name** -
- 14.3 Transport hazard class(es)**
- Class** -
- 14.4 Packing group** -
- 14.5 Environmental hazards** non-environmentally hazardous acc. to the dangerous goods regulations
- 14.6 Special precautions for user**  
There is no additional information.
- 14.7 Transport in bulk according to Annex II of MARPOL and the IBC Code**  
The cargo is not intended to be carried in bulk.
- 14.8 Information for each of the UN Model Regulations**
- Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN)**  
Not subject to ADR, RID and ADN.
- International Maritime Dangerous Goods Code (IMDG)**  
Not subject to IMDG.
- International Civil Aviation Organization (ICAO-IATA/DGR)**  
Not subject to ICAO-IATA.

## SECTION 15: Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

**Relevant provisions of the European Union (EU)**

**Restrictions according to REACH, Annex XVII**

none of the ingredients are listed

<b>Dangerous substances with restrictions (REACH, Annex XVII)</b>				
<b>Name of substance</b>	<b>Name acc. to inventory</b>	<b>Type of registration</b>	<b>Restriction</b>	<b>No</b>
Proteinase K EDT-01	this product meets the criteria for classification in accordance with Regulation No 1272/2008/EC	1907/2006/EC annex XVII	R3	3

**Legend**

- R3 1. Shall not be used in:  
- ornamental articles intended to produce light or colour effects by means of different phases, for example in ornamental lamps and ashtrays,

# Proteinase K EDT-01

## Legend

- tricks and jokes,
  - games for one or more participants, or any article intended to be used as such, even with ornamental aspects,
2. Articles not complying with paragraph 1 shall not be placed on the market.
  3. Shall not be placed on the market if they contain a colouring agent, unless required for fiscal reasons, or perfume, or both, if they:
    - can be used as fuel in decorative oil lamps for supply to the general public, and,
    - present an aspiration hazard and are labelled with R65 or H304,
  4. Decorative oil lamps for supply to the general public shall not be placed on the market unless they conform to the European Standard on Decorative oil lamps (EN 14059) adopted by the European Committee for Standardisation (CEN).
  5. Without prejudice to the implementation of other Community provisions relating to the classification, packaging and labelling of dangerous substances and mixtures, suppliers shall ensure, before the placing on the market, that the following requirements are met:
    - (a) lamp oils, labelled with R65 or H304, intended for supply to the general public are visibly, legibly and indelibly marked as follows: 'Keep lamps filled with this liquid out of the reach of children'; and, by 1 December 2010, 'Just a sip of lamp oil - or even sucking the wick of lamps - may lead to life-threatening lung damage';
    - (b) grill lighter fluids, labelled with R65 or H304, intended for supply to the general public are legibly and indelibly marked by 1 December 2010 as follows: 'Just a sip of grill lighter may lead to life threatening lung damage';
    - (c) lamp oils and grill lighters, labelled with R65 or H304, intended for supply to the general public are packaged in black opaque containers not exceeding 1 litre by 1 December 2010.
  6. No later than 1 June 2014, the Commission shall request the European Chemicals Agency to prepare a dossier, in accordance with Article 69 of the present Regulation with a view to ban, if appropriate, grill lighter fluids and fuel for decorative lamps, labelled R65 or H304, intended for supply to the general public.
  7. Natural or legal persons placing on the market for the first time lamp oils and grill lighter fluids, labelled with R65 or H304, shall by 1 December 2011, and annually thereafter, provide data on alternatives to lamp oils and grill lighter fluids labelled R65 or H304 to the competent authority in the Member State concerned. Member States shall make those data available to the Commission.

## List of substances subject to authorisation (REACH, Annex XIV)

none of the ingredients are listed

## Seveso Directive

2012/18/EU (Seveso III)			
No	Dangerous substance/hazard categories	Qualifying quantity (tonnes) for the application of lower and upper-tier requirements	Notes
	not assigned		

## Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) - Annex II

none of the ingredients are listed

# Proteinase K EDT-01

## Regulation 166/2006/EC concerning the establishment of a European Pollutant Release and Transfer Register (PRTR)

none of the ingredients are listed

## Directive 2000/60/EC establishing a framework for Community action in the field of water policy (WFD)

none of the ingredients are listed

## Regulation 98/2013/EU on the marketing and use of explosives precursors

none of the ingredients are listed

## SECTION 16: Other information

### Abbreviations and acronyms

Abbreviations and acronyms	
Abbr.	Descriptions of used abbreviations
ADN	Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	Accord européen relatif au transport international des marchandises dangereuses par route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BCF	Bioconcentration factor
CAS	Chemical Abstracts Service (service that maintains the most comprehensive list of chemical substances)
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DGR	Dangerous Goods Regulations (see IATA/DGR)
DNEL	Derived No-Effect Level
EC No	The EC Inventory (EINECS, ELINCS and the NLP-list) is the source for the seven-digit EC number, an identifier of substances commercially available within the EU (European Union)
EH40/2005	EH40/2005 Workplace exposure limits ( <a href="http://www.nationalarchives.gov.uk/doc/open-government-licence/">http://www.nationalarchives.gov.uk/doc/open-government-licence/</a> )
EINECS	European Inventory of Existing Commercial Chemical Substances
ELINCS	European List of Notified Chemical Substances
Eye Dam.	Seriously damaging to the eye
Eye Irrit.	Irritant to the eye
GHS	"Globally Harmonized System of Classification and Labelling of Chemicals" developed by the United Nations
IATA	International Air Transport Association
IATA/DGR	Dangerous Goods Regulations (DGR) for the air transport (IATA)
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods Code



# Proteinase K EDT-01

<b>Abbreviations and acronyms</b>	
<b>Abbr.</b>	<b>Descriptions of used abbreviations</b>
index No	The Index number is the identification code given to the substance in Part 3 of Annex VI to Regulation (EC) No 1272/2008
log KOW	n-Octanol/water
MARPOL	International Convention for the Prevention of Pollution from Ships (abbr. of "Marine Pollutant")
NLP	No-Longer Polymer
PBT	Persistent, Bioaccumulative and Toxic
PNEC	Predicted No-Effect Concentration
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
Resp. Sens.	Respiratory sensitisation
RID	Règlement concernant le transport International ferroviaire des marchandises Dangereuses (Regulations concerning the International carriage of Dangerous goods by Rail)
Skin Corr.	Corrosive to skin
Skin Irrit.	Irritant to skin
Skin Sens.	Skin sensitisation
STEL	Short-term exposure limit
STOT SE	Specific target organ toxicity - single exposure
TWA	Time-weighted average
vPvB	Very Persistent and very Bioaccumulative
WEL	Workplace exposure limit

## Key literature references and sources for data

Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures.

Regulation (EC) No. 1907/2006 (REACH), amended by 2015/830/EU.

Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN).

International Maritime Dangerous Goods Code (IMDG).

Dangerous Goods Regulations (DGR) for the air transport (IATA).

## Classification procedure

Physical and chemical properties.

Health hazards.

Environmental hazards.

The method for classification of the mixture is based on ingredients of the mixture (additivity formula).

# Proteinase K EDT-01

---

## List of relevant phrases (code and full text as stated in chapter 2 and 3)

List of relevant phrases (code and full text as stated in chapter 2 and 3)	
Code	Text
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.

### Responsible for the safety data sheet

C.S.B. GmbH  
Düsseldorfer Str. 113  
47809 Krefeld

Telephone: +49 (0) 2151 - 652086 - 0  
Telefax: +49 (0) 2151 - 652086 - 9  
e-Mail: [info@csb-online.de](mailto:info@csb-online.de)  
Website: [www.csb-online.de](http://www.csb-online.de)

### Disclaimer

This information is based upon the present state of our knowledge.  
This SDS has been compiled and is solely intended for this product.



# Lysis Buffer LDT-01

---

for full text of abbreviations: see SECTION 16

## 2.2 Label elements

### Labelling according to Regulation (EC) No 1272/2008 (CLP)

**Signal word** warning

### Pictograms

**GHS07**



### Hazard statements

**H302** Harmful if swallowed.  
**H315** Causes skin irritation.  
**H319** Causes serious eye irritation.

### Precautionary statements

**P280** Wear protective gloves/protective clothing/eye protection/face protection.  
**P301+P330+P331** IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
**P303+P361+P353** IF ON SKIN (or hair): take off immediately all contaminated clothing. Rinse skin with water/shower.  
**P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
**P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
**P311** Call a POISON CENTER/doctor.

### Supplemental hazard information

**EUH208** Contains 2,4,7,9-tetramethyldec-5-yne-4,7-diol. May produce an allergic reaction.

**Hazardous ingredients for labelling** Guanadine hydrochloride

## 2.3 Other hazards

There is no additional information.

### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

# Lysis Buffer LDT-01



## SECTION 3: Composition/information on ingredients

### 3.1 Substances

not relevant (mixture)

### 3.2 Mixtures

#### Description of the mixture

Hazardous ingredients acc. to GHS				
Name of substance	Identifier	Wt%	Classification acc. to GHS	Pictograms
Guanadine hydrochloride	CAS No 50-01-1  EC No 200-002-3  Index No 607-148-00-0	50 – < 75	Acute Tox. 4 / H302 Acute Tox. 4 / H332 Skin Irrit. 2 / H315 Eye Irrit. 2 / H319	
2,4,7,9-Tetramethyldec-5-in-4,7-diol	CAS No 126-86-3  EC No 204-809-1	< 1	Acute Tox. 4 / H302 Eye Dam. 1 / H318 Skin Sens. 1B / H317 Aquatic Chronic 3 / H412	

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General notes

Take off immediately all contaminated clothing.

In all cases of doubt, or when symptoms persist, seek medical advice.

#### Following inhalation

Provide fresh air.

If breathing is irregular or stopped, immediately seek medical assistance and start first aid actions.

#### Following skin contact

Wash with plenty of soap and water.

#### Following eye contact

Rinse cautiously with water for several minutes.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Remove contact lenses, if present and easy to do. Continue rinsing.

# Lysis Buffer LDT-01

---

## **Following ingestion**

Rinse mouth. Do not induce vomiting.  
Get medical advice/attention if you feel unwell.

## **Notes for the doctor**

none

## **4.2 Most important symptoms and effects, both acute and delayed**

These information are not available.

## **4.3 Indication of any immediate medical attention and special treatment needed**

none

## **SECTION 5: Firefighting measures**

### **5.1 Extinguishing media**

#### **Suitable extinguishing media**

water spray, alcohol resistant foam, fire extinguishing powder, carbon dioxide (CO<sub>2</sub>)

#### **Unsuitable extinguishing media**

water jet

### **5.2 Special hazards arising from the substance or mixture**

Hazardous decomposition products: Section 10.

### **5.3 Advice for firefighters**

In case of fire and/or explosion do not breathe fumes.  
Co-ordinate firefighting measures to the fire surroundings.  
Do not allow firefighting water to enter drains or water courses.  
Collect contaminated firefighting water separately.  
Fight fire with normal precautions from a reasonable distance.

#### **Special protective equipment for firefighters**

use suitable breathing apparatus

## **SECTION 6: Accidental release measures**

### **6.1 Personal precautions, protective equipment and emergency procedures**

#### **For non-emergency personnel**

Remove persons to safety.  
Ventilate affected area.  
Wearing of suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing.

# Lysis Buffer LDT-01

---

## **For emergency responders**

Wear breathing apparatus if exposed to vapours/dust/spray/gases.

## **6.2 Environmental precautions**

Keep away from drains, surface and ground water.

Retain contaminated washing water and dispose of it.

## **6.3 Methods and material for containment and cleaning up**

### **Advices on how to clean up a spill**

Collect spillage.

Absorbent material (e.g. sand, diatomaceous earth, acid binder, universal binder, sawdust, etc.).

### **Appropriate containment techniques**

Use of adsorbent materials.

### **Other information relating to spills and releases**

Place in appropriate containers for disposal.

Ventilate affected area.

## **6.4 Reference to other sections**

Personal protective equipment: see section 8.

Incompatible materials: see section 10.

Disposal considerations: see section 13.

## **SECTION 7: Handling and storage**

## **7.1 Precautions for safe handling**

### **Measures to prevent fire as well as aerosol and dust generation**

Use local and general ventilation.

Keep away from sources of ignition - No smoking.

### **Specific notes/details**

None.

### **Measures to protect the environment**

Avoid release to the environment.

### **Advice on general occupational hygiene**

Do not eat, drink and smoke in work areas.

Wash hands after use.

Preventive skin protection (barrier creams/ointments) is recommended.

Remove contaminated clothing and protective equipment before entering eating areas.

## **7.2 Conditions for safe storage, including any incompatibilities**

### **Flammability hazards**

None.

# Lysis Buffer LDT-01

## Incompatible substances or mixtures

Incompatible materials: see section 10.

## Protect against external exposure, such as

frost

## Consideration of other advice

Keep away from food, drink and animal feedingstuffs.

## Ventilation requirements

Provision of sufficient ventilation.

## Specific designs for storage rooms or vessels

### Storage temperature

15 – 28 °C

### Packaging compatibilities

Keep only in original container.

## 7.3 Specific end use(s)

No information available.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

Relevant DNELs of components of the mixture						
Name of substance	CAS No	Endpoint	Threshold level	Protection goal, route of exposure	Used in	Exposure time
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	DNEL	1.76 mg/m <sup>3</sup>	human, inhalatory	worker (industry)	chronic - systemic effects
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	DNEL	0.5 mg/kg bw/day	human, dermal	worker (industry)	chronic - systemic effects
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	DNEL	0.43 mg/m <sup>3</sup>	human, inhalatory	consumer (private households)	chronic - systemic effects
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	DNEL	0.25 mg/kg bw/day	human, dermal	consumer (private households)	chronic - systemic effects
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	DNEL	0.25 mg/kg bw/day	human, oral	consumer (private households)	chronic - systemic effects

Relevant PNECs of components of the mixture				
Name of substance	CAS No	Endpoint	Threshold level	Environmental compartment
Guanadine hydrochloride	50-01-1	PNEC	0.29 mg/l	freshwater
Guanadine hydrochloride	50-01-1	PNEC	0.029 mg/l	marine water



## Lysis Buffer LDT-01

<b>Relevant PNECs of components of the mixture</b>				
<b>Name of substance</b>	<b>CAS No</b>	<b>Endpoint</b>	<b>Threshold level</b>	<b>Environmental compartment</b>
Guanadine hydrochloride	50-01-1	PNEC	7,125 mg/l	sewage treatment plant (STP)
Guanadine hydrochloride	50-01-1	PNEC	1.08 mg/kg	freshwater sediment
Guanadine hydrochloride	50-01-1	PNEC	0.108 mg/kg	marine sediment
Guanadine hydrochloride	50-01-1	PNEC	4.35 mg/kg	soil
Guanadine hydrochloride	50-01-1	PNEC	0.335 mg/l	water
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.032 mg/cm <sup>3</sup>	marine sediment
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.004 mg/cm <sup>3</sup>	marine water
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.028 mg/cm <sup>3</sup>	soil
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.04 mg/cm <sup>3</sup>	freshwater
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.32 mg/cm <sup>3</sup>	freshwater sediment
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.4 mg/cm <sup>3</sup>	water
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	7 mg/cm <sup>3</sup>	sewage treatment plant (STP)
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.04 mg/l	freshwater
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.004 mg/l	marine water
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.4 mg/l	water
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	7 mg/l	sewage treatment plant (STP)
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.32 mg/kg	freshwater sediment
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.032 mg/kg	marine sediment
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	PNEC	0.028 mg/kg	soil

# Lysis Buffer LDT-01

## 8.2 Exposure controls

### Appropriate engineering controls

General ventilation.

### Individual protection measures (personal protective equipment)

#### Eye/face protection

Wear eye/face protection.

#### Hand protection

Material	Material thickness	Breakthrough times of the glove material
no information available		

Wear suitable gloves.

Chemical protection gloves are suitable, which are tested according to EN 374.

Check leak-tightness/impermeability prior to use.

In the case of wanting to use the gloves again, clean them before taking off and air them well.

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

#### Environmental exposure controls

Use appropriate container to avoid environmental contamination.

Keep away from drains, surface and ground water.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

#### Appearance

Physical state	liquid
Form	fluid
Colour	according to product specification
Odour	characteristic
Odour threshold	these information are not available

#### Other safety parameters

pH (value)	these information are not available
Melting point/freezing point	these information are not available
Initial boiling point and boiling range	(unknown)
Flash point	321 °C
Evaporation rate	these information are not available

## Lysis Buffer LDT-01

---

Flammability (solid, gas)	not relevant (fluid)
<b>Explosive limits</b>	
Lower explosion limit (LEL)	these information are not available
Upper explosion limit (UEL)	these information are not available
Vapour pressure	these information are not available
Density	these information are not available
Vapour density	these information are not available
Relative density	these information are not available
<b>Solubility(ies)</b>	
Water solubility	miscible in any proportion
<b>Partition coefficient</b>	
n-octanol/water (log KOW)	these information are not available
Auto-ignition temperature	these information are not available
Relative self-ignition temperature for solids	not relevant (Fluid)
Decomposition temperature	these information are not available
<b>Viscosity</b>	
Kinematic viscosity	these information are not available
Dynamic viscosity	these information are not available
Explosive properties	not explosive
Oxidising properties	shall not be classified as oxidising
<b>9.2 Other information</b>	
None	

# Lysis Buffer LDT-01

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material is not reactive under normal ambient conditions.

### 10.2 Chemical stability

The material is stable under normal ambient and anticipated storage and handling conditions of temperature and pressure.

### 10.3 Possibility of hazardous reactions

No known hazardous reactions.

### 10.4 Conditions to avoid

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

### 10.5 Incompatible materials

There is no additional information.

### 10.6 Hazardous decomposition products

Reasonably anticipated hazardous decomposition products produced as a result of use, storage, spill and heating are not known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### Classification procedure

If not otherwise specified the classification is based on:  
Ingredients of the mixture (additivity formula).

#### Classification according to GHS (1272/2008/EC, CLP)

#### Acute toxicity

Harmful if swallowed.

Acute toxicity estimate (ATE) of components of the mixture			
Name of substance	CAS No	Exposure route	ATE
Guanadine hydrochloride	50-01-1	oral	556.5 mg/kg
Guanadine hydrochloride	50-01-1	inhalation: dust/mist	3.181 mg/l/4h
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	oral	500 mg/kg

Acute toxicity of components of the mixture					
Name of substance	CAS No	Exposure route	Endpoint	Value	Species
Guanadine hydrochloride	50-01-1	oral	LD50	556.5 mg/kg	rat
Guanadine hydrochloride	50-01-1	inhalation: dust/mist	LC50	3.181 mg/l/4h	rat
Guanadine hydrochloride	50-01-1	dermal	LD50	>2,000 mg/kg	rabbit

## Lysis Buffer LDT-01

<b>Acute toxicity of components of the mixture</b>					
<b>Name of substance</b>	<b>CAS No</b>	<b>Exposure route</b>	<b>Endpoint</b>	<b>Value</b>	<b>Species</b>
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	oral	LD50	>500 mg/kg	rat
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	dermal	LD50	>2,000 mg/kg	rat

### **Skin corrosion/irritation**

Causes skin irritation.

### **Serious eye damage/eye irritation**

Causes serious eye irritation.

### **Respiratory or skin sensitisation**

Contains 2,4,7,9-tetramethyldec-5-yne-4,7-diol. May produce an allergic reaction.

### **Skin sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Respiratory sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Germ cell mutagenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Carcinogenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Reproductive toxicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Specific target organ toxicity - single exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Specific target organ toxicity - repeated exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

# Lysis Buffer LDT-01

## Aspiration hazard

Shall not be classified as presenting an aspiration hazard.

## SECTION 12: Ecological information

### 12.1 Toxicity

#### Aquatic toxicity (acute)

Test data are not available for the complete mixture.

#### Aquatic toxicity (acute) of components of the mixture

Aquatic toxicity (acute) of components of the mixture					
Name of substance	CAS No	Endpoint	Value	Species	Exposure time
Guanadine hydrochloride	50-01-1	LC50	1,758 mg/l	fish	48 h
Guanadine hydrochloride	50-01-1	EC50	11.8 mg/l	algae	72 h
Guanadine hydrochloride	50-01-1	ErC50	33.5 mg/l	algae	72 h
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	EC50	91 mg/l	daphnia magna	48 h
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	LC50	36 mg/l	fathead minnow (pimephales promelas)	96 h

#### Aquatic toxicity (chronic)

Test data are not available for the complete mixture.

#### Aquatic toxicity (chronic) of components of the mixture

Aquatic toxicity (chronic) of components of the mixture					
Name of substance	CAS No	Endpoint	Value	Species	Exposure time
Guanadine hydrochloride	50-01-1	LC50	1,872 mg/l	fish	24 h
Guanadine hydrochloride	50-01-1	growth (EbCx) 10%	7,125 mg/l	microorganisms	18 h

### 12.2 Persistence and degradability

# Lysis Buffer LDT-01

## Degradability of components of the mixture

Degradability of components of the mixture						
Name of substance	CAS No	Process	Degradation rate	Time	Method	Source
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3	carbon dioxide generation	5 %	29 d	OECD Guideline 301	ECHA
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3		12 %	60 d	ISO DIS 9439	ECHA
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3		8 %	60 d	ISO DIS 9439	ECHA

### Biodegradation

Data are not available.

### Persistence

Data are not available.

## 12.3 Bioaccumulative potential

Data are not available.

### Bioaccumulative potential of components of the mixture

Bioaccumulative potential of components of the mixture			
Name of substance	CAS No	BCF	Log KOW
Guanadine hydrochloride	50-01-1		<-1.7 (pH value: 7.4, 20 °C)
2,4,7,9-Tetramethyldec-5-in-4,7-diol	126-86-3		2.64

## 12.4 Mobility in soil

Data are not available.

## 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

## 12.6 Other adverse effects

Data are not available.

# Lysis Buffer LDT-01

## Endocrine disrupting potential

None of the ingredients are listed.

## Remarks

Water hazard class - WHC (Wassergefährdungsklasse): 1 (Slightly hazardous to water)

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

This material and its container must be disposed of as hazardous waste.

### Sewage disposal-relevant information

Do not empty into drains.

### Waste treatment of containers/packagings

Handle contaminated packages in the same way as the substance itself.

## Remarks

Please consider the relevant national or regional provisions.

## SECTION 14: Transport information

**14.1 UN number** not subject to transport regulations

**14.2 UN proper shipping name** -

**14.3 Transport hazard class(es)**  
**Class** -

**14.4 Packing group** -

**14.5 Environmental hazards** -

**14.6 Special precautions for user**  
There is no additional information.

**14.7 Transport in bulk according to Annex II of MARPOL and the IBC Code**  
The cargo is not intended to be carried in bulk.

### 14.8 Information for each of the UN Model Regulations

**Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN)**  
Not subject to ADR, RID and ADN.

**International Maritime Dangerous Goods Code (IMDG)**  
Not subject to IMDG.



# Lysis Buffer LDT-01

## International Civil Aviation Organization (ICAO-IATA/DGR)

Not subject to ICAO-IATA.

### SECTION 15: Regulatory information

#### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

##### Relevant provisions of the European Union (EU)

##### Restrictions according to REACH, Annex XVII

none of the ingredients are listed

Dangerous substances with restrictions (REACH, Annex XVII)				
Name of substance	Name acc. to inventory	Type of registration	Restriction	No
Lysis Buffer LDT-01	this product meets the criteria for classification in accordance with Regulation No 1272/2008/EC	1907/2006/EC annex XVII	R3	3

##### Legend

- R3
1. Shall not be used in:
    - ornamental articles intended to produce light or colour effects by means of different phases, for example in ornamental lamps and ashtrays,
    - tricks and jokes,
    - games for one or more participants, or any article intended to be used as such, even with ornamental aspects,
  2. Articles not complying with paragraph 1 shall not be placed on the market.
  3. Shall not be placed on the market if they contain a colouring agent, unless required for fiscal reasons, or perfume, or both, if they:
    - can be used as fuel in decorative oil lamps for supply to the general public, and,
    - present an aspiration hazard and are labelled with R65 or H304,
  4. Decorative oil lamps for supply to the general public shall not be placed on the market unless they conform to the European Standard on Decorative oil lamps (EN 14059) adopted by the European Committee for Standardisation (CEN).
  5. Without prejudice to the implementation of other Community provisions relating to the classification, packaging and labelling of dangerous substances and mixtures, suppliers shall ensure, before the placing on the market, that the following requirements are met:
    - (a) lamp oils, labelled with R65 or H304, intended for supply to the general public are visibly, legibly and indelibly marked as follows: 'Keep lamps filled with this liquid out of the reach of children'; and, by 1 December 2010, 'Just a sip of lamp oil - or even sucking the wick of lamps - may lead to life-threatening lung damage';
    - (b) grill lighter fluids, labelled with R65 or H304, intended for supply to the general public are legibly and indelibly marked by 1 December 2010 as follows: 'Just a sip of grill lighter may lead to life threatening lung damage';
    - (c) lamp oils and grill lighters, labelled with R65 or H304, intended for supply to the general public are packaged in black opaque containers not exceeding 1 litre by 1 December 2010.
  6. No later than 1 June 2014, the Commission shall request the European Chemicals Agency to prepare a dossier, in accordance with Article 69 of the present Regulation with a view to ban, if appropriate, grill lighter fluids and fuel for decorative lamps, labelled R65 or H304, intended for supply to the general public.
  7. Natural or legal persons placing on the market for the first time lamp oils and grill lighter fluids, labelled with R65 or H304, shall by 1 December 2011, and annually thereafter, provide data on alternatives to lamp oils and grill lighter fluids labelled R65 or H304 to the competent authority in the Member State concerned. Member States shall make those data available to the Commission.

# Lysis Buffer LDT-01

## List of substances subject to authorisation (REACH, Annex XIV)

none of the ingredients are listed

## Seveso Directive

2012/18/EU (Seveso III)			
No	Dangerous substance/hazard categories	Qualifying quantity (tonnes) for the application of lower and upper-tier requirements	Notes
	not assigned		

## Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) - Annex II

none of the ingredients are listed

## Regulation 166/2006/EC concerning the establishment of a European Pollutant Release and Transfer Register (PRTR)

none of the ingredients are listed

## Directive 2000/60/EC establishing a framework for Community action in the field of water policy (WFD)

none of the ingredients are listed

## Regulation 98/2013/EU on the marketing and use of explosives precursors

none of the ingredients are listed

## SECTION 16: Other information

### Indication of changes (revised safety data sheet)

Indication of changes: Section 7, 8, 15

### Abbreviations and acronyms

Abbreviations and acronyms	
Abbr.	Descriptions of used abbreviations
Acute Tox.	Acute toxicity
ADN	Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	Accord européen relatif au transport international des marchandises dangereuses par route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
Aquatic Chronic	Hazardous to the aquatic environment - chronic hazard
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CAS	Chemical Abstracts Service (service that maintains the most comprehensive list of chemical substances)

## Lysis Buffer LDT-01

<b>Abbreviations and acronyms</b>	
<b>Abbr.</b>	<b>Descriptions of used abbreviations</b>
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DGR	Dangerous Goods Regulations (see IATA/DGR)
DNEL	Derived No-Effect Level
EC No	The EC Inventory (EINECS, ELINCS and the NLP-list) is the source for the seven-digit EC number, an identifier of substances commercially available within the EU (European Union)
EINECS	European Inventory of Existing Commercial Chemical Substances
ELINCS	European List of Notified Chemical Substances
Eye Dam.	Seriously damaging to the eye
Eye Irrit.	Irritant to the eye
GHS	"Globally Harmonized System of Classification and Labelling of Chemicals" developed by the United Nations
IATA	International Air Transport Association
IATA/DGR	Dangerous Goods Regulations (DGR) for the air transport (IATA)
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods Code
index No	The Index number is the identification code given to the substance in Part 3 of Annex VI to Regulation (EC) No 1272/2008
log KOW	n-Octanol/water
MARPOL	International Convention for the Prevention of Pollution from Ships (abbr. of "Marine Pollutant")
NLP	No-Longer Polymer
PBT	Persistent, Bioaccumulative and Toxic
PNEC	Predicted No-Effect Concentration
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Règlement concernant le transport International ferroviaire des marchandises Dangereuses (Regulations concerning the International carriage of Dangerous goods by Rail)
Skin Corr.	Corrosive to skin
Skin Irrit.	Irritant to skin
Skin Sens.	Skin sensitisation
vPvB	Very Persistent and very Bioaccumulative

# Lysis Buffer LDT-01

## Key literature references and sources for data

Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures.  
Regulation (EC) No. 1907/2006 (REACH), amended by 2015/830/EU.  
Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN).  
International Maritime Dangerous Goods Code (IMDG).  
Dangerous Goods Regulations (DGR) for the air transport (IATA).

## Classification procedure

Physical and chemical properties.  
Health hazards.  
Environmental hazards.  
The method for classification of the mixture is based on ingredients of the mixture (additivity formula).

## List of relevant phrases (code and full text as stated in chapter 2 and 3)

List of relevant phrases (code and full text as stated in chapter 2 and 3)	
Code	Text
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.
H412	Harmful to aquatic life with long lasting effects.

## Responsible for the safety data sheet

C.S.B. GmbH  
Düsseldorfer Str. 113  
47809 Krefeld

Telephone: +49 (0) 2151 - 652086 - 0  
Telefax: +49 (0) 2151 - 652086 - 9  
e-Mail: [info@csb-online.de](mailto:info@csb-online.de)  
Website: [www.csb-online.de](http://www.csb-online.de)

## Disclaimer

This information is based upon the present state of our knowledge.  
This SDS has been compiled and is solely intended for this product.

## Tissue lysis Buffer MDT-01

Version number: 2.1  
Replaces version of: 2015-10-06

Revision: 2017-08-09  
First version: 06.10.2015

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

<b>Trade name</b>	<u>Tissue lysis Buffer MDT-01</u>
<b>Product number</b>	MDT-01
<b>Registration number (REACH)</b>	not relevant (mixture)
<b>CAS number</b>	not relevant (mixture)

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** Chemicals for various applications

#### 1.3 Details of the supplier of the safety data sheet

KURABO INDUSTRIES LTD. Telephone: ++81-72-820-3079  
Bio-Medical department Telefax: ++81-72-820-3095  
Neyagawa Techno Center, 14-5, Shimokida-cho  
Neyagawa, Osaka 572-0823 Japan

**e-mail (competent person)** sdb@csb-online.de

Please do not use this e-mail adress to ask for the latest safety data sheet. For this purpose contact

#### 1.4 Emergency telephone number

As above or next toxicological information centre.

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 (CLP)

Classification acc. to GHS				
Section	Hazard class	Category	Hazard class and category	Hazard state-ment
3.3	serious eye damage/eye irritation	2	Eye Irrit. 2	H319

for full text of abbreviations: see SECTION 16

# Tissue lysis Buffer MDT-01

The most important adverse physicochemical, human health and environmental effects

Spillage and fire water can cause pollution of watercourses.

## 2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008 (CLP)

Signal word warning

Pictograms

GHS07



Hazard statements

H319 Causes serious eye irritation.

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

## 2.3 Other hazards

There is no additional information.

Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.



## SECTION 3: Composition/information on ingredients

### 3.1 Substances


not relevant (mixture)

### 3.2 Mixtures

Description of the mixture

Hazardous ingredients acc. to GHS					
Name of substance	Identifier	Wt%	Classification acc. to GHS	Pictograms	M-Factors
edetetic acid	CAS No 60-00-4  EC No 200-449-4  Index No 607-429-00-8	1 - <5	Acute Tox. 4 / H332 Eye Irrit. 2 / H319 STOT RE 2 / H373	 	

# Tissue lysis Buffer MDT-01

Hazardous ingredients acc. to GHS					
Name of substance	Identifier	Wt%	Classification acc. to GHS	Pictograms	M-Factors
sodium dodecyl sulphate	CAS No 151-21-3  EC No 205-788-1	1 - < 5	Acute Tox. 4 / H302 Skin Irrit. 2 / H315 Eye Dam. 1 / H318 Aquatic Chronic 3 / H412		

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General notes

Take off immediately all contaminated clothing.  
In all cases of doubt, or when symptoms persist, seek medical advice.

#### Following inhalation

Provide fresh air.  
If breathing is irregular or stopped, immediately seek medical assistance and start first aid actions.

#### Following skin contact

Wash with plenty of soap and water.

#### Following eye contact

Rinse cautiously with water for several minutes.  
In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.  
Remove contact lenses, if present and easy to do. Continue rinsing.

#### Following ingestion

Rinse mouth. Do not induce vomiting.  
Get medical advice/attention if you feel unwell.

#### Notes for the doctor

none

### 4.2 Most important symptoms and effects, both acute and delayed

These information are not available.

### 4.3 Indication of any immediate medical attention and special treatment needed

none

# Tissue lysis Buffer MDT-01

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### Suitable extinguishing media

water spray, alcohol resistant foam, fire extinguishing powder, carbon dioxide (CO<sub>2</sub>)

#### Unsuitable extinguishing media

water jet

### 5.2 Special hazards arising from the substance or mixture

Hazardous decomposition products: Section 10.

### 5.3 Advice for firefighters

In case of fire and/or explosion do not breathe fumes.

Co-ordinate firefighting measures to the fire surroundings.

Do not allow firefighting water to enter drains or water courses.

Collect contaminated firefighting water separately.

Fight fire with normal precautions from a reasonable distance.

#### Special protective equipment for firefighters

use suitable breathing apparatus

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Remove persons to safety.

Ventilate affected area.

Wearing of suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing.

#### For emergency responders

Wear breathing apparatus if exposed to vapours/dust/spray/gases.

### 6.2 Environmental precautions

Keep away from drains, surface and ground water.

Retain contaminated washing water and dispose of it.

If substance has entered a water course or sewer, inform the responsible authority.

### 6.3 Methods and material for containment and cleaning up

#### Advices on how to clean up a spill

Collect spillage.

Absorbent material (e.g. sand, diatomaceous earth, acid binder, universal binder, sawdust, etc.).



# Tissue lysis Buffer MDT-01

---

## **Appropriate containment techniques**

Use of adsorbent materials.

## **Other information relating to spills and releases**

Place in appropriate containers for disposal.

Ventilate affected area.

## **6.4 Reference to other sections**

Personal protective equipment: see section 8.

Incompatible materials: see section 10.

Disposal considerations: see section 13.

## **SECTION 7: Handling and storage**

### **7.1 Precautions for safe handling**

#### **Measures to prevent fire as well as aerosol and dust generation**

Use local and general ventilation.

Keep away from sources of ignition - No smoking.

#### **Specific notes/details**

None.

#### **Measures to protect the environment**

Avoid release to the environment.

#### **Advice on general occupational hygiene**

Do not eat, drink and smoke in work areas.

Wash hands after use.

Preventive skin protection (barrier creams/ointments) is recommended.

Remove contaminated clothing and protective equipment before entering eating areas.

### **7.2 Conditions for safe storage, including any incompatibilities**

#### **Flammability hazards**

None.

#### **Incompatible substances or mixtures**

Incompatible materials: see section 10.

#### **Protect against external exposure, such as**

frost

#### **Consideration of other advice**

Keep away from food, drink and animal feedingstuffs.

#### **Ventilation requirements**

Provision of sufficient ventilation.

# Tissue lysis Buffer MDT-01

## Specific designs for storage rooms or vessels

**Storage temperature** 15 – 28 °C

## Packaging compatibilities

Keep only in original container.

## 7.3 Specific end use(s)

No information available.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

Relevant DNELs of components of the mixture						
Name of substance	CAS No	End-point	Threshold level	Protection goal, route of exposure	Used in	Exposure time
edetic acid	60-00-4	DNEL	1.5 mg/m <sup>3</sup>	human, inhalatory	worker (industry)	chronic - local effects
edetic acid	60-00-4	DNEL	0.6 mg/m <sup>3</sup>	human, inhalatory	consumer (private households)	chronic - local effects
edetic acid	60-00-4	DNEL	25 mg/kg bw/day	human, oral	consumer (private households)	chronic - systemic effects
sodium dodecyl sulphate	151-21-3	DNEL	85 mg/m <sup>3</sup>	human, inhalatory	consumer (private households)	chronic - systemic effects
sodium dodecyl sulphate	151-21-3	DNEL	24 mg/kg	human, oral	consumer (private households)	chronic - systemic effects
sodium dodecyl sulphate	151-21-3	DNEL	4,060 mg/kg	human, dermal	worker (industry)	chronic - systemic effects
sodium dodecyl sulphate	151-21-3	DNEL	285 mg/m <sup>3</sup>	human, inhalatory	worker (industry)	chronic - systemic effects
sodium dodecyl sulphate	151-21-3	DNEL	2,440 mg/kg	human, dermal	consumer (private households)	chronic - systemic effects

Relevant PNECs of components of the mixture				
Name of substance	CAS No	Endpoint	Threshold level	Environmental compartment
edetic acid	60-00-4	PNEC	2.2 mg/l	freshwater
edetic acid	60-00-4	PNEC	0.22 mg/l	marine water
edetic acid	60-00-4	PNEC	1.2 mg/l	water

# Tissue lysis Buffer MDT-01

Relevant PNECs of components of the mixture				
Name of substance	CAS No	Endpoint	Threshold level	Environmental compartment
edetiac acid	60-00-4	PNEC	43 mg/l	sewage treatment plant (STP)
edetiac acid	60-00-4	PNEC	0.72 mg/kg	soil

## 8.2 Exposure controls

### Appropriate engineering controls

General ventilation.

### Individual protection measures (personal protective equipment)

#### Eye/face protection

Wear eye/face protection.

#### Hand protection

Material
these information are not available

Wear suitable gloves.

Chemical protection gloves are suitable, which are tested according to EN 374.

Check leak-tightness/impermeability prior to use.

In the case of wanting to use the gloves again, clean them before taking off and air them well.

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

#### Environmental exposure controls

Use appropriate container to avoid environmental contamination.

Keep away from drains, surface and ground water.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

#### Appearance

Physical state	liquid
Form	fluid
Colour	colourless

## Tissue lysis Buffer MDT-01

---

Odour	odourless
Odour threshold	these information are not available
<b>Other safety parameters</b>	
pH (value)	8.5
Melting point/freezing point	these information are not available
Initial boiling point and boiling range	100 °C
Flash point	170 °C
Evaporation rate	these information are not available
Flammability (solid, gas)	not relevant (fluid)
<b>Explosive limits</b>	
Lower explosion limit (LEL)	these information are not available
Upper explosion limit (UEL)	these information are not available
Vapour pressure	these information are not available
Density	0.9835 g/cm <sup>3</sup>
Vapour density	these information are not available
Relative density	these information are not available
<b>Solubility(ies)</b>	
Water solubility	miscible in any proportion
<b>Partition coefficient</b>	
n-octanol/water (log KOW)	these information are not available
Auto-ignition temperature	these information are not available
Relative self-ignition temperature for solids	not relevant (Fluid)
Decomposition temperature	these information are not available
<b>Viscosity</b>	
Kinematic viscosity	these information are not available
Dynamic viscosity	these information are not available
Explosive properties	not explosive
Oxidising properties	shall not be classified as oxidising

# Tissue lysis Buffer MDT-01

## 9.2 Other information

None

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material is not reactive under normal ambient conditions.

### 10.2 Chemical stability

The material is stable under normal ambient and anticipated storage and handling conditions of temperature and pressure.

### 10.3 Possibility of hazardous reactions

No known hazardous reactions.

### 10.4 Conditions to avoid

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

### 10.5 Incompatible materials

There is no additional information.

### 10.6 Hazardous decomposition products

Reasonably anticipated hazardous decomposition products produced as a result of use, storage, spill and heating are not known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### Classification procedure

If not otherwise specified the classification is based on:  
Ingredients of the mixture (additivity formula).

#### Classification according to GHS (1272/2008/EC, CLP)

#### Acute toxicity

Acute toxicity of components of the mixture					
Name of substance	CAS No	Exposure route	Endpoint	Value	Species
edetic acid	60-00-4	oral	LD50	4,500 mg/kg	rat
sodium dodecyl sulphate	151-21-3	oral	LD50	1,288 mg/kg	rat
sodium dodecyl sulphate	151-21-3	dermal	LD50	>2,000 mg/kg	rabbit

#### Skin corrosion/irritation

Shall not be classified as corrosive/irritant to skin.

# Tissue lysis Buffer MDT-01

## Serious eye damage/eye irritation

Causes serious eye irritation.

## Respiratory or skin sensitisation

### Skin sensitisation

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### Respiratory sensitisation

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Germ cell mutagenicity

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Carcinogenicity

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Reproductive toxicity

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Specific target organ toxicity - single exposure

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Specific target organ toxicity - repeated exposure

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## Aspiration hazard

Shall not be classified as presenting an aspiration hazard.

## SECTION 12: Ecological information

### 12.1 Toxicity

#### Aquatic toxicity (acute)

Test data are not available for the complete mixture.

Aquatic toxicity (acute) of components of the mixture					
Name of substance	CAS No	Endpoint	Value	Species	Exposure time
edetic acid	60-00-4	LC50	41 mg/l	bluegill ( <i>Lepomis macrochirus</i> )	96 h
edetic acid	60-00-4	EC50	625 mg/l	daphnia magna	24 h
edetic acid	60-00-4	EC50	610 mg/l	daphnia magna	24 h
edetic acid	60-00-4	EC50	140 mg/l	daphnia magna	48 h

## Tissue lysis Buffer MDT-01

Aquatic toxicity (acute) of components of the mixture					
Name of substance	CAS No	Endpoint	Value	Species	Exposure time
sodium dodecyl sulphate	151-21-3	EC50	5.55 mg/l	daphnia	48 h
sodium dodecyl sulphate	151-21-3	LC50	>120 mg/l	algae (Desmod-esmus subspicatus)	72 h

### Aquatic toxicity (chronic)

Test data are not available for the complete mixture.

Aquatic toxicity (chronic) of components of the mixture					
Name of substance	CAS No	Endpoint	Value	Species	Exposure time
edetic acid	60-00-4	NOEC	25.7 mg/l	zebra fish (danio rerio)	35 d
edetic acid	60-00-4	LOEC	50 mg/l	daphnia magna	21 d

## 12.2 Persistence and degradability

### Biodegradation

Data are not available.

### Persistence

Data are not available.

## 12.3 Bioaccumulative potential

Data are not available.

Bioaccumulative potential of components of the mixture			
Name of substance	CAS No	BCF	Log KOW
edetic acid	60-00-4	1.8	-3.86

## 12.4 Mobility in soil

Data are not available.

## 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

# Tissue lysis Buffer MDT-01

## 12.6 Other adverse effects

Data are not available.

### Endocrine disrupting potential

None of the ingredients are listed.

### Remarks

Water hazard class - WHC (Wassergefährdungsklasse): 2 (Hazardous to water)

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

This material and its container must be disposed of as hazardous waste.

#### Sewage disposal-relevant information

Do not empty into drains.

#### Waste treatment of containers/packagings

Handle contaminated packages in the same way as the substance itself.

### Remarks

Please consider the relevant national or regional provisions.

## SECTION 14: Transport information

14.1	<b>UN number</b>	not subject to transport regulations
14.2	<b>UN proper shipping name</b>	-
14.3	<b>Transport hazard class(es)</b>	
	<b>Class</b>	-
14.4	<b>Packing group</b>	-
14.5	<b>Environmental hazards</b>	non-environmentally hazardous acc. to the dangerous goods regulations
14.6	<b>Special precautions for user</b>	
		There is no additional information.
14.7	<b>Transport in bulk according to Annex II of MARPOL and the IBC Code</b>	
		The cargo is not intended to be carried in bulk.
14.8	<b><u>Information for each of the UN Model Regulations</u></b>	
	<b>Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN)</b>	
		Not subject to ADR, RID and ADN.



# Tissue lysis Buffer MDT-01

## International Maritime Dangerous Goods Code (IMDG)

Not subject to IMDG.

## International Civil Aviation Organization (ICAO-IATA/DGR)

Not subject to ICAO-IATA.

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Relevant provisions of the European Union (EU)

#### Restrictions according to REACH, Annex XVII

none of the ingredients are listed

Dangerous substances with restrictions (REACH, Annex XVII)					
Name of substance	Name acc. to inventory	CAS No	Type of registration	Conditions of restriction	No
Tissue lysis Buffer MDT-01	this product meets the criteria for classification in accordance with Regulation No 1272/2008/EC		1907/2006/EC annex XVII	R3	3

#### Legend

- R3
1. Shall not be used in:
    - ornamental articles intended to produce light or colour effects by means of different phases, for example in ornamental lamps and ashtrays,
    - tricks and jokes,
    - games for one or more participants, or any article intended to be used as such, even with ornamental aspects,
  2. Articles not complying with paragraph 1 shall not be placed on the market.
  3. Shall not be placed on the market if they contain a colouring agent, unless required for fiscal reasons, or perfume, or both, if they:
    - can be used as fuel in decorative oil lamps for supply to the general public, and,
    - present an aspiration hazard and are labelled with R65 or H304,
  4. Decorative oil lamps for supply to the general public shall not be placed on the market unless they conform to the European Standard on Decorative oil lamps (EN 14059) adopted by the European Committee for Standardisation (CEN).
  5. Without prejudice to the implementation of other Community provisions relating to the classification, packaging and labelling of dangerous substances and mixtures, suppliers shall ensure, before the placing on the market, that the following requirements are met:
    - (a) lamp oils, labelled with R65 or H304, intended for supply to the general public are visibly, legibly and indelibly marked as follows: 'Keep lamps filled with this liquid out of the reach of children'; and, by 1 December 2010, 'Just a sip of lamp oil - or even sucking the wick of lamps - may lead to life-threatening lung damage';
    - (b) grill lighter fluids, labelled with R65 or H304, intended for supply to the general public are legibly and indelibly marked by 1 December 2010 as follows: 'Just a sip of grill lighter may lead to life threatening lung damage';
    - (c) lamp oils and grill lighters, labelled with R65 or H304, intended for supply to the general public are packaged in black opaque containers not exceeding 1 litre by 1 December 2010.
  6. No later than 1 June 2014, the Commission shall request the European Chemicals Agency to prepare a dossier, in accordance with Article 69 of the present Regulation with a view to ban, if appropriate, grill lighter fluids and fuel for decorative lamps, labelled R65 or H304, intended for supply to the general public.
  7. Natural or legal persons placing on the market for the first time lamp oils and grill lighter fluids, labelled with

# Tissue lysis Buffer MDT-01

## Legend

R65 or H304, shall by 1 December 2011, and annually thereafter, provide data on alternatives to lamp oils and grill lighter fluids labelled R65 or H304 to the competent authority in the Member State concerned. Member States shall make those data available to the Commission.

## List of substances subject to authorisation (REACH, Annex XIV)

none of the ingredients are listed

## Seveso Directive

2012/18/EU (Seveso III)			
No	Dangerous substance/hazard categories	Qualifying quantity (tonnes) for the application of lower and upper-tier requirements	Notes
	not assigned		

## Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) - Annex II

none of the ingredients are listed

## Regulation 166/2006/EC concerning the establishment of a European Pollutant Release and Transfer Register (PRTR)

none of the ingredients are listed

## Directive 2000/60/EC establishing a framework for Community action in the field of water policy (WFD)

none of the ingredients are listed

## Regulation 98/2013/EU on the marketing and use of explosives precursors

none of the ingredients are listed

## SECTION 16: Other information

### Indication of changes (revised safety data sheet)

Indication of changes: Section 7, 8, 15

### Abbreviations and acronyms

Abbreviations and acronyms	
Abbr.	Descriptions of used abbreviations
Acute Tox.	Acute toxicity
ADN	Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	Accord européen relatif au transport international des marchandises dangereuses par route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
Aquatic Chronic	Hazardous to the aquatic environment - chronic hazard

## Tissue lysis Buffer MDT-01

<b>Abbreviations and acronyms</b>	
<b>Abbr.</b>	<b>Descriptions of used abbreviations</b>
BCF	Bioconcentration factor
CAS	Chemical Abstracts Service (service that maintains the most comprehensive list of chemical substances)
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DGR	Dangerous Goods Regulations (see IATA/DGR)
DNEL	Derived No-Effect Level
EC No	The EC Inventory (EINECS, ELINCS and the NLP-list) is the source for the seven-digit EC number, an identifier of substances commercially available within the EU (European Union)
EINECS	European Inventory of Existing Commercial Chemical Substances
ELINCS	European List of Notified Chemical Substances
Eye Dam.	Seriously damaging to the eye
Eye Irrit.	Irritant to the eye
GHS	"Globally Harmonized System of Classification and Labelling of Chemicals" developed by the United Nations
IATA	International Air Transport Association
IATA/DGR	Dangerous Goods Regulations (DGR) for the air transport (IATA)
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods Code
index No	The Index number is the identification code given to the substance in Part 3 of Annex VI to Regulation (EC) No 1272/2008
log KOW	n-Octanol/water
MARPOL	International Convention for the Prevention of Pollution from Ships (abbr. of "Marine Pollutant")
NLP	No-Longer Polymer
PBT	Persistent, Bioaccumulative and Toxic
PNEC	Predicted No-Effect Concentration
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Règlement concernant le transport International ferroviaire des marchandises Dangereuses (Regulations concerning the International carriage of Dangerous goods by Rail)
Skin Corr.	Corrosive to skin
Skin Irrit.	Irritant to skin
STOT RE	Specific target organ toxicity - repeated exposure
vPvB	Very Persistent and very Bioaccumulative

# Tissue lysis Buffer MDT-01

## Key literature references and sources for data

Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures.  
Regulation (EC) No. 1907/2006 (REACH), amended by 2015/830/EU.  
Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN).  
International Maritime Dangerous Goods Code (IMDG).  
Dangerous Goods Regulations (DGR) for the air transport (IATA).

## Classification procedure

Physical and chemical properties.  
Health hazards.  
Environmental hazards.  
The method for classification of the mixture is based on ingredients of the mixture (additivity formula).

## List of relevant phrases (code and full text as stated in chapter 2 and 3)

List of relevant phrases (code and full text as stated in chapter 2 and 3)	
Code	Text
H302	Harmful if swallowed.
H315	Causes skin irritation.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.
H373	May cause damage to organs through prolonged or repeated exposure.
H412	Harmful to aquatic life with long lasting effects.

## Responsible for the safety data sheet

C.S.B. GmbH  
Düsseldorfer Str. 113  
47809 Krefeld

Telephone: +49 (0) 2151 - 652086 - 0  
Telefax: +49 (0) 2151 - 652086 - 9  
e-Mail: [info@csb-online.de](mailto:info@csb-online.de)  
Website: [www.csb-online.de](http://www.csb-online.de)

## Disclaimer

This information is based upon the present state of our knowledge.  
This SDS has been compiled and is solely intended for this product.

## Wash Buffer WDT-02

Version number: 2.0  
Replaces version of: 2015-07-21

Revision: 2017-08-18  
First version: 21.07.2015

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

<b>Trade name</b>	<u>Wash Buffer WDT-02</u>
<b>Product number</b>	WDT-02
<b>Registration number (REACH)</b>	not relevant (mixture)
<b>CAS number</b>	not relevant (mixture)

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** Laboratory chemicals

#### 1.3 Details of the supplier of the safety data sheet

KURABO INDUSTRIES LTD. Telephone: ++81-72-820-3079  
Bio-Medical department Telefax: ++81-72-820-3095  
Neyagawa Techno Center, 14-5, Shimokida-cho  
Neyagawa, Osaka 572-0823 Japan

**e-mail (competent person)** sdb@csb-online.de

Please do not use this e-mail adress to ask for the latest safety data sheet. For this purpose contact

#### 1.4 Emergency telephone number

As above or next toxicological information centre.

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

##### Classification according to Regulation (EC) No 1272/2008 (CLP)

This mixture does not meet the criteria for classification in accordance with Regulation No 1272/2008/EC.

#### 2.2 Label elements

##### Labelling according to Regulation (EC) No 1272/2008 (CLP)

not required

## Wash Buffer WDT-02

---

### 2.3 Other hazards

There is no additional information.

#### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

## SECTION 3: Composition/information on ingredients

### 3.1 Substances

not relevant (mixture)

### 3.2 Mixtures

#### Hazardous ingredients acc. to EU regulation

None

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General notes

In all cases of doubt, or when symptoms persist, seek medical advice.

#### Following inhalation

Provide fresh air.

#### Following skin contact

Wash with plenty of soap and water.

#### Following eye contact

Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

#### Following ingestion

Rinse mouth. Do not induce vomiting.

Get medical advice/attention if you feel unwell.

#### Notes for the doctor

none

### 4.2 Most important symptoms and effects, both acute and delayed

These information are not available.

### 4.3 Indication of any immediate medical attention and special treatment needed

none

# Wash Buffer WDT-02

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### Suitable extinguishing media

water spray, alcohol resistant foam, fire extinguishing powder, carbon dioxide (CO<sub>2</sub>)

#### Unsuitable extinguishing media

water jet

### 5.2 Special hazards arising from the substance or mixture

Hazardous decomposition products: Section 10.

### 5.3 Advice for firefighters

In case of fire and/or explosion do not breathe fumes.

Co-ordinate firefighting measures to the fire surroundings.

Do not allow firefighting water to enter drains or water courses.

Collect contaminated firefighting water separately.

Fight fire with normal precautions from a reasonable distance.

#### Special protective equipment for firefighters

use suitable breathing apparatus

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Ventilate affected area.

Wearing of suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing.

#### For emergency responders

Wear breathing apparatus if exposed to vapours/dust/spray/gases.

### 6.2 Environmental precautions

Keep away from drains, surface and ground water.

Retain contaminated washing water and dispose of it.

### 6.3 Methods and material for containment and cleaning up

#### Advices on how to clean up a spill

Collect spillage.

Absorbent material (e.g. sand, diatomaceous earth, acid binder, universal binder, sawdust, etc.).

#### Appropriate containment techniques

Use of adsorbent materials.

# Wash Buffer WDT-02

---

## Other information relating to spills and releases

Place in appropriate containers for disposal.  
Ventilate affected area.

## 6.4 Reference to other sections

Personal protective equipment: see section 8.  
Incompatible materials: see section 10.  
Disposal considerations: see section 13.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

#### Measures to prevent fire as well as aerosol and dust generation

Use local and general ventilation.

#### Specific notes/details

None.

#### Measures to protect the environment

Avoid release to the environment.

#### Advice on general occupational hygiene

Do not eat, drink and smoke in work areas.  
Wash hands after use.  
Preventive skin protection (barrier creams/ointments) is recommended.  
Remove contaminated clothing and protective equipment before entering eating areas.

### 7.2 Conditions for safe storage, including any incompatibilities

#### Flammability hazards

None.

#### Incompatible substances or mixtures

Incompatible materials: see section 10.

#### Protect against external exposure, such as

frost

#### Consideration of other advice

Keep away from food, drink and animal feedingstuffs.

#### Ventilation requirements

Provision of sufficient ventilation.

#### Specific designs for storage rooms or vessels

**Storage temperature** 15 – 28 °C

#### Packaging compatibilities

Keep only in original container.



# Wash Buffer WDT-02

## 7.3 Specific end use(s)

No information available.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

No data available.

### 8.2 Exposure controls

#### Appropriate engineering controls

General ventilation.

#### Individual protection measures (personal protective equipment)

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

Material
no information available

Wear suitable gloves.

Chemical protection gloves are suitable, which are tested according to EN 374.

Check leak-tightness/impermeability prior to use.

In the case of wanting to use the gloves again, clean them before taking off and air them well.

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

##### Environmental exposure controls

Use appropriate container to avoid environmental contamination.

Keep away from drains, surface and ground water.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

#### Appearance

Physical state	liquid
Form	fluid
Colour	colourless

## Wash Buffer WDT-02

---

Odour	recognizable
Odour threshold	these information are not available
<b>Other safety parameters</b>	
pH (value)	these information are not available
Melting point/freezing point	these information are not available
Initial boiling point and boiling range	(unknown)
Flash point	not applicable
Evaporation rate	these information are not available
Flammability (solid, gas)	not relevant (fluid)
<b>Explosive limits</b>	
Lower explosion limit (LEL)	these information are not available
Upper explosion limit (UEL)	these information are not available
Vapour pressure	these information are not available
Density	these information are not available
Vapour density	these information are not available
Relative density	these information are not available
<b>Solubility(ies)</b>	
Water solubility	miscible in any proportion
<b>Partition coefficient</b>	
n-octanol/water (log KOW)	these information are not available
Auto-ignition temperature	these information are not available
Relative self-ignition temperature for solids	not relevant (Fluid)
Decomposition temperature	these information are not available
<b>Viscosity</b>	
Kinematic viscosity	these information are not available
Dynamic viscosity	these information are not available
Explosive properties	not explosive
Oxidising properties	shall not be classified as oxidising

## Wash Buffer WDT-02

---

### 9.2 Other information

None

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material is not reactive under normal ambient conditions.

### 10.2 Chemical stability

The material is stable under normal ambient and anticipated storage and handling conditions of temperature and pressure.

### 10.3 Possibility of hazardous reactions

No known hazardous reactions.

### 10.4 Conditions to avoid

There are no specific conditions known which have to be avoided.

### 10.5 Incompatible materials

There is no additional information.

### 10.6 Hazardous decomposition products

Reasonably anticipated hazardous decomposition products produced as a result of use, storage, spill and heating are not known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### Classification procedure

If not otherwise specified the classification is based on:  
Ingredients of the mixture (additivity formula).

#### Classification according to GHS (1272/2008/EC, CLP)

This mixture does not meet the criteria for classification in accordance with Regulation No 1272/2008/EC.

#### Acute toxicity

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

#### Skin corrosion/irritation

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Serious eye damage/eye irritation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Respiratory or skin sensitisation**

#### **Skin sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

#### **Respiratory sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Germ cell mutagenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Carcinogenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Reproductive toxicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Specific target organ toxicity - single exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Specific target organ toxicity - repeated exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Aspiration hazard**

Shall not be classified as presenting an aspiration hazard.

## SECTION 12: Ecological information

### **12.1 Toxicity**

#### **Aquatic toxicity (acute)**

Test data are not available for the complete mixture.

#### **Aquatic toxicity (chronic)**

Test data are not available for the complete mixture.

### **12.2 Persistence and degradability**

#### **Biodegradation**

The relevant substances of the mixture are readily biodegradable.

## Wash Buffer WDT-02

### Persistence

Data are not available.

### 12.3 Bioaccumulative potential

Data are not available.

### 12.4 Mobility in soil

Data are not available.

### 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### 12.6 Other adverse effects

Data are not available.

### Endocrine disrupting potential

None of the ingredients are listed.

### Remarks

Water hazard class - WHC (Wassergefährdungsklasse): 1 (Slightly hazardous to water)

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local/regional/national/international regulations.

### Sewage disposal-relevant information

Do not empty into drains.

### Waste treatment of containers/packagings

Handle contaminated packages in the same way as the substance itself.

### Remarks

Please consider the relevant national or regional provisions.

## SECTION 14: Transport information

14.1	UN number	not subject to transport regulations
14.2	UN proper shipping name	-
14.3	Transport hazard class(es)	
	Class	-
14.4	Packing group	-
14.5	Environmental hazards	-

## Wash Buffer WDT-02

### 14.6 Special precautions for user

There is no additional information.

### 14.7 Transport in bulk according to Annex II of MARPOL and the IBC Code

The cargo is not intended to be carried in bulk.

### 14.8 Information for each of the UN Model Regulations

#### Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN)

Not subject to ADR, RID and ADN.

#### International Maritime Dangerous Goods Code (IMDG)

Not subject to IMDG.

#### International Civil Aviation Organization (ICAO-IATA/DGR)

Not subject to ICAO-IATA.

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Relevant provisions of the European Union (EU)

#### Restrictions according to REACH, Annex XVII

none of the ingredients are listed

#### List of substances subject to authorisation (REACH, Annex XIV)

none of the ingredients are listed

#### Seveso Directive

2012/18/EU (Seveso III)			
No	Dangerous substance/hazard categories	Qualifying quantity (tonnes) for the application of lower and upper-tier requirements	Notes
	not assigned		

#### Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) - Annex II

none of the ingredients are listed

#### Regulation 166/2006/EC concerning the establishment of a European Pollutant Release and Transfer Register (PRTR)

none of the ingredients are listed

## Wash Buffer WDT-02

### Directive 2000/60/EC establishing a framework for Community action in the field of water policy (WFD)

none of the ingredients are listed

### Regulation 98/2013/EU on the marketing and use of explosives precursors

none of the ingredients are listed

## SECTION 16: Other information

### Indication of changes (revised safety data sheet)

Indication of changes: Section 7, 8, 15

### Abbreviations and acronyms

Abbreviations and acronyms	
Abbr.	Descriptions of used abbreviations
ADN	Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	Accord européen relatif au transport international des marchandises dangereuses par route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
CAS	Chemical Abstracts Service (service that maintains the most comprehensive list of chemical substances)
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DGR	Dangerous Goods Regulations (see IATA/DGR)
GHS	"Globally Harmonized System of Classification and Labelling of Chemicals" developed by the United Nations
IATA	International Air Transport Association
IATA/DGR	Dangerous Goods Regulations (DGR) for the air transport (IATA)
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods Code
MARPOL	International Convention for the Prevention of Pollution from Ships (abbr. of "Marine Pollutant")
PBT	Persistent, Bioaccumulative and Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Règlement concernant le transport International ferroviaire des marchandises Dangereuses (Regulations concerning the International carriage of Dangerous goods by Rail)
vPvB	Very Persistent and very Bioaccumulative

## Wash Buffer WDT-02

---

### Key literature references and sources for data

Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures.  
Regulation (EC) No. 1907/2006 (REACH), amended by 2015/830/EU.  
Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN).  
International Maritime Dangerous Goods Code (IMDG).  
Dangerous Goods Regulations (DGR) for the air transport (IATA).

### Classification procedure

Physical and chemical properties.  
Health hazards.  
Environmental hazards.  
The method for classification of the mixture is based on ingredients of the mixture (additivity formula).

### Responsible for the safety data sheet

C.S.B. GmbH	Telephone: +49 (0) 2151 - 652086 - 0
Düsseldorfer Str. 113	Telefax: +49 (0) 2151 - 652086 - 9
47809 Krefeld	e-Mail: <a href="mailto:info@csb-online.de">info@csb-online.de</a>
	Website: <a href="http://www.csb-online.de">www.csb-online.de</a>

### Disclaimer

This information is based upon the present state of our knowledge.  
This SDS has been compiled and is solely intended for this product.



## Elution Buffer CDT-01

Version number: 2.0  
Replaces version of: 2015-10-06

Revision: 2017-08-09  
First version: 06.10.2015

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

<b>Trade name</b>	<u>Elution Buffer CDT-01</u>
<b>Product number</b>	CDT-017
<b>Registration number (REACH)</b>	not relevant (mixture)
<b>CAS number</b>	not relevant (mixture)

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** Chemicals for various applications

#### 1.3 Details of the supplier of the safety data sheet

KURABO INDUSTRIES LTD. Telephone: ++81-72-820-3079  
Bio-Medical department Telefax: ++81-72-820-3095  
Neyagawa Techno Center, 14-5, Shimokida-cho  
Neyagawa, Osaka 572-0823 Japan

**e-mail (competent person)** sdb@csb-online.de

Please do not use this e-mail adress to ask for the latest safety data sheet. For this purpose contact

#### 1.4 Emergency telephone number

As above or next toxicological information centre.

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

##### Classification according to Regulation (EC) No 1272/2008 (CLP)

This mixture does not meet the criteria for classification in accordance with Regulation No 1272/2008/EC.

#### 2.2 Label elements

##### Labelling according to Regulation (EC) No 1272/2008 (CLP)

not required

# Elution Buffer CDT-01

---

## 2.3 Other hazards

There is no additional information.

### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

## SECTION 3: Composition/information on ingredients

### 3.1 Substances

not relevant (mixture)

### 3.2 Mixtures

#### Hazardous ingredients acc. to EU regulation

None

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General notes

In all cases of doubt, or when symptoms persist, seek medical advice.

#### Following inhalation

Provide fresh air.

#### Following skin contact

Wash with plenty of soap and water.

#### Following eye contact

Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

#### Following ingestion

Rinse mouth. Do not induce vomiting.

Get medical advice/attention if you feel unwell.

#### Notes for the doctor

none

### 4.2 Most important symptoms and effects, both acute and delayed

These information are not available.

### 4.3 Indication of any immediate medical attention and special treatment needed

none

# Elution Buffer CDT-01

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### Suitable extinguishing media

water spray, alcohol resistant foam, fire extinguishing powder, carbon dioxide (CO<sub>2</sub>)

#### Unsuitable extinguishing media

water jet

### 5.2 Special hazards arising from the substance or mixture

Hazardous decomposition products: Section 10.

### 5.3 Advice for firefighters

In case of fire and/or explosion do not breathe fumes.

Co-ordinate firefighting measures to the fire surroundings.

Collect contaminated firefighting water separately.

Fight fire with normal precautions from a reasonable distance.

#### Special protective equipment for firefighters

use suitable breathing apparatus

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Ventilate affected area.

Wearing of suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing.

#### For emergency responders

Wear breathing apparatus if exposed to vapours/dust/spray/gases.

### 6.2 Environmental precautions

Keep away from drains, surface and ground water.

Retain contaminated washing water and dispose of it.

### 6.3 Methods and material for containment and cleaning up

#### Advices on how to clean up a spill

Collect spillage.

Absorbent material (e.g. sand, diatomaceous earth, acid binder, universal binder, sawdust, etc.).

#### Appropriate containment techniques

Use of adsorbent materials.

# Elution Buffer CDT-01

---

## Other information relating to spills and releases

Place in appropriate containers for disposal.  
Ventilate affected area.

## 6.4 Reference to other sections

Personal protective equipment: see section 8.  
Incompatible materials: see section 10.  
Disposal considerations: see section 13.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

#### Measures to prevent fire as well as aerosol and dust generation

Use local and general ventilation.

#### Specific notes/details

None.

#### Measures to protect the environment

Avoid release to the environment.

#### Advice on general occupational hygiene

Do not eat, drink and smoke in work areas.  
Wash hands after use.  
Preventive skin protection (barrier creams/ointments) is recommended.  
Remove contaminated clothing and protective equipment before entering eating areas.

### 7.2 Conditions for safe storage, including any incompatibilities

#### Flammability hazards

None.

#### Incompatible substances or mixtures

Incompatible materials: see section 10.

#### Protect against external exposure, such as

frost

#### Consideration of other advice

Keep away from food, drink and animal feedingstuffs.

#### Ventilation requirements

Provision of sufficient ventilation.

#### Specific designs for storage rooms or vessels

**Storage temperature** 15 – 28 °C

#### Packaging compatibilities

Keep only in original container.

# Elution Buffer CDT-01

## 7.3 Specific end use(s)

No information available.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

No data available.

### 8.2 Exposure controls

#### Appropriate engineering controls

General ventilation.

#### Individual protection measures (personal protective equipment)

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

Material
these information are not available

Wear suitable gloves.

Chemical protection gloves are suitable, which are tested according to EN 374.

Check leak-tightness/impermeability prior to use.

In the case of wanting to use the gloves again, clean them before taking off and air them well.

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

##### Environmental exposure controls

Use appropriate container to avoid environmental contamination.

Keep away from drains, surface and ground water.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

#### Appearance

Physical state	liquid
Form	fluid
Colour	according to product specification

## Elution Buffer CDT-01

---

Odour	characteristic
Odour threshold	these information are not available
<b>Other safety parameters</b>	
pH (value)	these information are not available
Melting point/freezing point	these information are not available
Initial boiling point and boiling range	(unknown)
Flash point	not applicable
Evaporation rate	these information are not available
Flammability (solid, gas)	not relevant (fluid)
<b>Explosive limits</b>	
Lower explosion limit (LEL)	these information are not available
Upper explosion limit (UEL)	these information are not available
Vapour pressure	these information are not available
Density	these information are not available
Vapour density	these information are not available
Relative density	these information are not available
<b>Solubility(ies)</b>	
Water solubility	not miscible in any proportion
<b>Partition coefficient</b>	
n-octanol/water (log KOW)	these information are not available
Auto-ignition temperature	these information are not available
Relative self-ignition temperature for solids	not relevant (Fluid)
Decomposition temperature	these information are not available
<b>Viscosity</b>	
Kinematic viscosity	these information are not available
Dynamic viscosity	these information are not available
Explosive properties	not explosive
Oxidising properties	shall not be classified as oxidising

# Elution Buffer CDT-01

---

## 9.2 Other information

None

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material is not reactive under normal ambient conditions.

### 10.2 Chemical stability

The material is stable under normal ambient and anticipated storage and handling conditions of temperature and pressure.

### 10.3 Possibility of hazardous reactions

No known hazardous reactions.

### 10.4 Conditions to avoid

There are no specific conditions known which have to be avoided.

### 10.5 Incompatible materials

There is no additional information.

### 10.6 Hazardous decomposition products

Reasonably anticipated hazardous decomposition products produced as a result of use, storage, spill and heating are not known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### Classification procedure

If not otherwise specified the classification is based on:  
Ingredients of the mixture (additivity formula).

#### Classification according to GHS (1272/2008/EC, CLP)

This mixture does not meet the criteria for classification in accordance with Regulation No 1272/2008/EC.

#### Acute toxicity

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

#### Skin corrosion/irritation

Classification could not be established because:  
Data are lacking, inconclusive, or conclusive but not sufficient for classification.

# Elution Buffer CDT-01

---

## **Serious eye damage/eye irritation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Respiratory or skin sensitisation**

### **Skin sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

### **Respiratory sensitisation**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Germ cell mutagenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Carcinogenicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Reproductive toxicity**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Specific target organ toxicity - single exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Specific target organ toxicity - repeated exposure**

Classification could not be established because:

Data are lacking, inconclusive, or conclusive but not sufficient for classification.

## **Aspiration hazard**

Shall not be classified as presenting an aspiration hazard.

## **SECTION 12: Ecological information**

### **12.1 Toxicity**

#### **Aquatic toxicity (acute)**

Test data are not available for the complete mixture.

#### **Aquatic toxicity (chronic)**

Test data are not available for the complete mixture.

### **12.2 Persistence and degradability**

#### **Biodegradation**

The relevant substances of the mixture are readily biodegradable.



# Elution Buffer CDT-01

## Persistence

Data are not available.

### 12.3 Bioaccumulative potential

Data are not available.

### 12.4 Mobility in soil

Data are not available.

### 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### 12.6 Other adverse effects

Data are not available.

#### Endocrine disrupting potential

None of the ingredients are listed.

#### Remarks

Water hazard class - WHC (Wassergefährdungsklasse): nwg (Non-hazardous to water)

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local/regional/national/international regulations.

#### Sewage disposal-relevant information

Do not empty into drains.

#### Waste treatment of containers/packagings

Handle contaminated packages in the same way as the substance itself.

#### Remarks

Please consider the relevant national or regional provisions.

## SECTION 14: Transport information

14.1	UN number	not subject to transport regulations
14.2	UN proper shipping name	-
14.3	Transport hazard class(es)	
	Class	-
14.4	Packing group	-
14.5	Environmental hazards	non-environmentally hazardous acc. to the dangerous goods regulations

# Elution Buffer CDT-01

## 14.6 Special precautions for user

There is no additional information.

## 14.7 Transport in bulk according to Annex II of MARPOL and the IBC Code

The cargo is not intended to be carried in bulk.

## 14.8 Information for each of the UN Model Regulations

### Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN)

Not subject to ADR, RID and ADN.

### International Maritime Dangerous Goods Code (IMDG)

Not subject to IMDG.

### International Civil Aviation Organization (ICAO-IATA/DGR)

Not subject to ICAO-IATA.

## SECTION 15: Regulatory information

## 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

### Relevant provisions of the European Union (EU)

#### Restrictions according to REACH, Annex XVII

none of the ingredients are listed

#### List of substances subject to authorisation (REACH, Annex XIV)

none of the ingredients are listed

#### Seveso Directive

2012/18/EU (Seveso III)			
No	Dangerous substance/hazard categories	Qualifying quantity (tonnes) for the application of lower and upper-tier requirements	Notes
	not assigned		

#### Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) - Annex II

none of the ingredients are listed

#### Regulation 166/2006/EC concerning the establishment of a European Pollutant Release and Transfer Register (PRTR)

none of the ingredients are listed

# Elution Buffer CDT-01

## Directive 2000/60/EC establishing a framework for Community action in the field of water policy (WFD)

none of the ingredients are listed

## Regulation 98/2013/EU on the marketing and use of explosives precursors

none of the ingredients are listed

### SECTION 16: Other information

#### Indication of changes (revised safety data sheet)

Indication of changes: Section 7, 8, 15

#### Abbreviations and acronyms

Abbreviations and acronyms	
Abbr.	Descriptions of used abbreviations
ADN	Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	Accord européen relatif au transport international des marchandises dangereuses par route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
CAS	Chemical Abstracts Service (service that maintains the most comprehensive list of chemical substances)
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures
DGR	Dangerous Goods Regulations (see IATA/DGR)
GHS	"Globally Harmonized System of Classification and Labelling of Chemicals" developed by the United Nations
IATA	International Air Transport Association
IATA/DGR	Dangerous Goods Regulations (DGR) for the air transport (IATA)
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods Code
MARPOL	International Convention for the Prevention of Pollution from Ships (abbr. of "Marine Pollutant")
PBT	Persistent, Bioaccumulative and Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Règlement concernant le transport International ferroviaire des marchandises Dangereuses (Regulations concerning the International carriage of Dangerous goods by Rail)
vPvB	Very Persistent and very Bioaccumulative

# Elution Buffer CDT-01

---

## Key literature references and sources for data

Regulation (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures.  
Regulation (EC) No. 1907/2006 (REACH), amended by 2015/830/EU.  
Transport of dangerous goods by road, rail and inland waterway (ADR/RID/ADN).  
International Maritime Dangerous Goods Code (IMDG).  
Dangerous Goods Regulations (DGR) for the air transport (IATA).

## Classification procedure

Physical and chemical properties.  
Health hazards.  
Environmental hazards.  
The method for classification of the mixture is based on ingredients of the mixture (additivity formula).

## Responsible for the safety data sheet

C.S.B. GmbH	Telephone: +49 (0) 2151 - 652086 - 0
Düsseldorfer Str. 113	Telefax: +49 (0) 2151 - 652086 - 9
47809 Krefeld	e-Mail: <a href="mailto:info@csb-online.de">info@csb-online.de</a>
	Website: <a href="http://www.csb-online.de">www.csb-online.de</a>

## Disclaimer

This information is based upon the present state of our knowledge.  
This SDS has been compiled and is solely intended for this product.