

**FUJIFILM**

**Wako**

Code No. 018-22141 (100  $\mu$ L)

**Anti SQSTM1/A170/p62, Rabbit**  
**抗SQSTM1/A170/p62, ウサギ**

Anti SQSTM1/A170/p62 antibody is raised in rabbits using recombinant murine SQSTM1/A170/p62 PEST domain (T7 tag + amino acid 254-333 + His tag) as the immunogen.

**Preparation** : The antiserum is diluted with an equal amount of PBS and absorbed with SQSTM1/A170/p62 knockout mouse tissue proteins.

**Formulation** : Diluted antiserum containing 0.1 % Sodium azide

**Specificity** : Specific for mouse and rat SQSTM1/A170/p62. Slightly reactive with human SQSTM1/A170/p62.

**Working dilution** :

Western blot	1 : 200 ~ 1 : 1,000
Immunohistochemistry	1 : 400 ~ 1 : 1,000
Immunofluorescence	1 : 400 ~ 1 : 1,000

**[Storage]** Store at  $-20^{\circ}\text{C}$ .

**[Package]** 100  $\mu$ L

**[Reference]**

1) Ishii, T., Itoh, K., Takahashi, S., Sato, H., Yanagawa, T., Katoh, Y., Bannai, S. and Yamamoto, M. : *J. Biol. Chem.*, **275**, 16023 (2000).

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1802KA1