



FUJIFILM Cellular Dynamics, Inc.  
525 Science Drive  
Madison, WI 53711 USA

Toll Free in US (877) 320-6688 / (608) 310-5100 **T**  
(608) 310-5101 **F**  
fcdi-support@fujifilm.com **E**  
www.fujifilmcdi.com **W**

## **Biosafety Documentation:** *iCell® Microglia*

Catalog Number(s): C1110  
Donor ID Number: 01279

### **Cell Source and Biosafety Level Classification**

iCell® Microglia are human cells differentiated from a master bank of stably induced pluripotent stem (iPS) cells. FUJIFILM Cellular Dynamics, Inc. (FCDI), classifies these cells as Biosafety Level 1 (BSL1) based on the United States Centers for Disease Control and Prevention publication: *Biosafety in Microbiological and Biomedical Laboratories*. We recommend handling iCell Microglia according to the biosafety guidelines applicable in your region.

### **Reprogramming**

The iPS cell lines were generated from human peripheral blood through ectopic expression of reprogramming factors (i.e. Oct4, Sox2, Nanog, Lin28, Klf4, L-Myc, SV40LT) by episomal transfection. Following reprogramming, no episomal plasmids were detected by PCR in the iPS cell line.

### **Engineering**

No engineering was performed on the reprogrammed iPS cell line.

### **Infectious Disease Testing**

The iPS cell line is negative for HBV, HCV, HIV-1, HIV-2, HTLV-1, and HTLV-2.

### **Reference(s)**

FCDI Technical Support: fcdi-support@fujifilm.com; US and Canada toll-free +1 (877) 320-6688 or +1 (608) 310-5100.