

SAFETY DATA SHEET

Issuing Date: 15-Jan-2019

Version 1

MyCell® Disease and Diversity Products

1. PRODUCT AND COMPANY IDENTIFICATION

Product Name	MyCell® Disease and Diversity Products
Product Code(s)	C1081, C1085, C1119, DDP-CMC-0.5X,DDP-CMC-1X; C1120, C1121, DDP-NRC-1X, DDP-NRC-0.5X; C1112, C1113
Product Use	For research use only.
<u>Manufactured by</u> FUJIFILM Cellular Dynamics, Inc. 525 Science Drive Madison, WI 53711	
SDSs are available at the following website(s):	https://fujifilmcdi.com
Company Phone Number	U.S.A: 877-320-6688 U.S.A: 608-310-5100 Singapore: 6862-2116
Emergency Telephone	Transport-CHEMTREC Inside NA: 800-424-9300 Transport-CHEMTREC Outside NA: 703-527-3887 Transport-CANUTEC Inside Canada: 613-996-66666 Medical Emergency (24 hour): 877-935-7387
E-mail	fcdi-support@fujifilm.com

2. HAZARDS IDENTIFICATION

Classification

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS)

GHS Label elements, including precautionary statements

Not classified

Precautionary Statements

Prevention

Not applicable

Response Not applicable

Storage Not applicable

Disposal Not applicable

Hazards not otherwise classified (HNOC) Not classified

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No	Weight-%
Dimethyl Sulfoxide	67-68-5	7-13%

*The exact percentage (concentration) of composition has been withheld as a trade secret

4. FIRST AID MEASURES

First aid measures for different exposure routes

General advice	If symptoms persist, call a physician.
Eye contact	Rinse thoroughly with plenty of water, also under the eyelids. Get medical attention if irritation persists.
Skin contact	Wash off immediately with plenty of water. Get medical attention if irritation develops and persists.
Inhalation	Move to fresh air. If symptoms persist, call a physician.
Ingestion	If swallowed, do not induce vomiting - seek medical advice.

Most important symptoms/effects, acute and delayed

None known.

Indication of immediate medical attention and special treatment needed, if necessary

Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Dry chemical, CO₂, alcohol-resistant foam or water spray.

Unsuitable Extinguishing Media

None known.

Specific hazards arising from the chemical

None known.

Hazardous Combustion Products Carbon oxides. Sulfur oxides.

Explosion Data

Sensitivity to Mechanical Impact none

Sensitivity to Static Discharge none

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Use personal protective equipment. Avoid contact with skin, eyes or clothing.

Environmental precautions

Prevent product from entering drains.

Methods and materials for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Use personal protective equipment. Soak up with inert absorbent material. Pick up and transfer to properly labeled containers. Clean contaminated surface thoroughly. After cleaning, flush away traces with water.

7. HANDLING AND STORAGE

Precautions for safe handling

Wear personal protective equipment. Avoid contact with skin, eyes or clothing. Avoid breathing vapors or mists. Ensure adequate ventilation.

Conditions for safe storage, including any incompatibilities

Keep in properly labeled containers.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH	AIHA - Workplace Environmental Exposure Levels (WEELs) - TWAs
Dimethyl Sulfoxide				250 ppm TWA

Exposure controls

Engineering Measures Ventilation systems

Individual protection measures, such as personal protective equipment

Eye/Face Protection	Tightly fitting safety goggles.
Skin and body protection	Wear protective gloves/clothing.
Respiratory protection	If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.
General Hygiene Considerations	When using do not eat, drink or smoke. Take off contaminated clothing and wash before reuse. Regular cleaning of equipment, work area and clothing is recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance Odor Threshold pH Specific Gravity Flash point Decomposition temperature Melting point / melting range Flammability Limit in Air	clear Not available Not available Not available Not available 32 °F / 0 °C Not available	Odor Physical State @20°C Molecular Weight Autoignition temperature Boiling point / boiling range Freezing Point	Slight Liquid Not available Not available Not available 24.8 °F / -4 °C
Oxidizing Properties Solubility Evaporation rate Vapor density VOC (Ib/gal) Dynamic viscosity	Not available Soluble in water Not available Not available None Not available	Explosive Properties Partition coefficient Vapor Pressure Density VOC (g/I)	Not available Not available Not available Not available None

10. STABILITY AND REACTIVITY

Reactivity

Stable under recommended storage conditions.

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

None under normal processing.

Conditions to Avoid

Excessive heat. Freezing.

Incompatible Materials

Strong oxidizing agents. Strong acids. Strong bases. Acid chlorides.

Hazardous Decomposition Products

None known.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information

Acute toxicity	
Inhalation	Non-irritating during normal use.
Eyes	Non-irritating during normal use.
Skin	Causes mild skin irritation.
Ingestion	Ingestion may cause stomach discomfort.

Component Information

Chemical Name	Oral LD50	Dermal LD50	LC50 (lethal concentration)
Dimethyl Sulfoxide	= 14500 mg/kg (Rat) = 28300	= 40 g/kg (Rat)	> 5.33 mg/L (Rat)4 h
	mg/kg (Rat)		

Information on toxicological effects

No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Irritation Corrosivity Sensitization Mutagenic Effects Reproductive Toxicity Carcinogenicity	Mild skin irritant. No information available. No information available. No information available. No information available. None known.	
STOT - single exposure	No information available.	
STOT - repeated exposure	No information available.	
Target Organ Effects	Skin.	
Aspiration hazard	No information available.	
Numerical measures of toxicity - Product Information		
ATEmix (oral)	>5000 mg/kg mg/l	

ATE: Acute toxicity estimate

12. ECOLOGICAL INFORMATION

Ecotoxicity

The environmental impact of this product has not been fully investigated.

Chemical Name	Algae toxicity	Toxicity to fish	Toxicity to microorganisms	Toxicity to daphnia and other aquatic invertebrates
Dimethyl Sulfoxide		Pimephales promelas: 34000 mg/L at 96 h		

Persistence and degradability

No information available.

Bioaccumulation

No information available.

Chemical Name	Octanol Water Partition Coefficient (log pow)
Dimethyl Sulfoxide	-2.03

Mobility

No information available.

Other adverse effects

No information available.

13. DISPOSAL CONSIDERATIONS

Waste Disposal Methods

Dispose of in accordance with local regulations.

Contaminated packaging

Do not re-use empty containers.

14. TRANSPORT INFORMATION		
DOT	Not regulated	
TDG	Not regulated.	
<u>MEX</u>	Not regulated	
ICAO	Not regulated	
IATA	Not regulated	
IMDG	Not regulated	
ADR/RID	Not regulated	
ADN	Not regulated	

15. REGULATORY INFORMATION

International Inventories

TSCA	Yes
DSL/NDSL	Yes
PICCS	Yes
EINECS/ELINCS	Yes
ENCS	Yes
IECSC	Yes
KECL	Yes

AICS

Yes

*Yes - All component(s) of this product are included or are exempt from listing on the inventory.

*No - Indicates the component(s) of this product are either not listed or have not been determined to be listed on the inventory.

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory
DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List
EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
ENCS - Japan Existing and New Chemical Substances
IECSC - China Inventory of Existing Chemical Substances
KECL - Korean Existing and Evaluated Chemical Substances
PICCS - Philippines Inventory of Chemicals and Chemical Substances
Australian Inventory of Chemical Substances

U.S. Federal Regulations

TSCA Sections 4, 5 and 12(b)

This product does not contain any chemicals regulated by TSCA Sections 4, 5 or 12(b).

<u>SARA 313</u>

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

SARA 311/312 Hazard Categories

Classification is shown in section 2 of this SDS

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals

U.S. State Right-to-Know Regulations

Chemical Name	Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Dimethyl Sulfoxide		Х			

International Regulations

Canada - NDSL

This product does not contain any NDSL chemicals.

<u> Mexico - Grade</u>

No information available

Mexico - Carcinogen Status and Exposure Limits

No information available

Other Regulations

16. OTHER INFORMATION							
NFPA	Health Hazard 1	Flammability 0	Instability 0	Physical and chemical hazards -			
HMIS	Health Hazard 1	Flammability 0	Physical Hazard 0	Personal protection B			
Prepared By	FUJIFILM	FUJIFILM Environment, Health and Safety, phone: (01) (800) 473-3854					
Revision Date	15-Jan-20	15-Jan-2019					
Revision Note	No information available						
Disclaimer	knowledg given is o transport specifica may not l	The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.					

end