SAFETY DATA SHEET

According to Regulation 1907/2006/EC

PAPAIC DIGEST OF LIVER

Created: Last update: April 20, 2007 October 5, 2012 (v.2)

1. PRODUCT AND COMPANY IDENTIFICATION

Trade name Reference Papaic digest of liver A1911

Recommended use

Fermentation - Cell culture - Diagnostics

Supplier

SOLABIA 29 rue Delizy 93698 - Pantin Cedex FRANCE www.solabia.com

 Tel.
 00 33 1 48 10 19 40

 Fax.
 00 33 1 48 91 18 77

 Mail
 info.fds@solabia.fr

00.33.1.45.42.59.59

Emergency telephone (ORFILA)

| ר | | | IDENIT | | TION |
|----------|----------------|-----|--------------|-------|------|
| Ζ. Γ | 1 A / A | RDS | IDENT | IFICA | |
| | .,, | | | | |

| Human health hazards | Not classified as hazardous | |
|--------------------------|--|--|
| Eye hazards: | May cause slight transient irritation | |
| Skin hazards: | Not harmful | |
| Ingestion hazards: | Not harmful | |
| Inhalation hazards: | May be harmful if inhaled (irritation of mucous membranes and upper respiratory tract) | |
| Environment hazards | Not classified as hazardous | |
| Physico-chemical hazards | None | |

3. INFORMATION ON INGREDIENTS

Composition

See technical data sheet

This product does not contain, from supplier knowledge, any dangerous substances in quantity justifying their inscription in the present section, according to European Union regulation.

4. FIRST AID MEASURES

| Inhalation | Remove to fresh air. |
|--------------|---|
| Ingestion | Wash out mouth with water. |
| Skin contact | Wash skin with water |
| Eye contact | Rinse opened eyes with plenty of water. |





5. FIRE FIGHTING MEASURES

| Suitable extinguishing media | Water spray, carbon dioxide, dry chemical powder or appropriate foam |
|---|--|
| Extinguishing media which must not be used for safety reasons | Not concerned |
| Special hazards | None |

6. ACCIDENTAL RELEASE MEASURES

| Personal protection/precautions | Wear protective clothing |
|---------------------------------|--|
| Environmental precautions | Avoid raising dust. Aerate the workstation and wash the dissipating area after having completely recuperated the product |
| Methods for cleaning up | Sweep up and place in an appropriate sack before elimination |

7. HANDLING AND STORAGE

Handling

Do not breathe dust

Storage

| Storage conditions: | Store at room temperature |
|--------------------------------|---------------------------|
| Recommended package materials: | Plastic |

Specific uses

Not listed

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure limit values

Not concerned

Occupational exposure controls

Respiratory protection:Wear respiratory maskHand protection:Wear chemical resistant glovesEye protection:Wear safety glasses or chemical gogglesSkin protection:Wear protective clothingHygiene measures:Wash thoroughly after handling

Environmental exposure controls Not available

9. PHYSICAL AND CHEMICAL PROPERTIES

General information

Physical state at 20°C: Powder Color: Light brown

Properties

Density: 0.8



| 10. STABILITY AND REACTIVITY | | |
|---|---|--|
| Stability | Stable product | |
| Conditions to avoid | None | |
| Materials to avoid | None | |
| Hazardous decomposition products | Not available | |
| 44 - | | |
| 11. TOXICOLOGICAL INFORMATI | ION | |
| Acute toxicity | No available data | |
| 12. ECOLOGICAL INFORMATION | | |
| Ecotoxicity | Biodegradable product | |
| | | |
| 13. DISPOSAL CONSIDERATIONS | | |
| Disposal of residues | Comply with national and Community regulations in force | |
| Disposal of soiled packaging | Comply with national and Community regulations in force | |
| | | |
| 14. TRANSPORT INFORMATION | | |
| Not subject to the transport' regulations | | |

15. REGULATORY INFORMATIONClassification and labelling according to European directives' into force:Not classified as hazardousSymbolNot concernedRisk phrases RNot concernedSafety phrases SNot concerned

16. OTHER INFORMATION

This sheet completes the instruction notice but does not replace it. The data are based on the current state of our knowledge. They are given in good faith. User attention is directed, among other points, to eventual risks incurred when the product is used in ways contrary to its original purpose.

This document does not in any event exempt the user from knowing and applying all norms and texts regulating its activity. The user alone assumes sole responsibility for the precautions linked to the use of the product.

Reason for the revision:

Compliance with REACH Regulation n° 1907/2006/EC