

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
 Product name : MKTTn BROTH
 Product code : BM078
 Product group : Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Industrial use, Professional use
 Use of the substance/mixture : Microbiological control

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

BIOKAR DIAGNOSTICS - SOLABIA GROUP
 Rue des Quarante Mines
 60002 Beauvais Cedex - FRANCE
 T 0033 3 44 14 33 33 - F 0033 3 44 14 33 34
info.fds@solabia.fr - www.biokar-diagnostics.fr

1.4. Emergency telephone number

No additional information available

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

2.2. Label elements

Safety data sheet available on request.

2.3. Other hazards

Adverse physicochemical, human health and environmental effects : May be slightly irritating to eyes. Repeated or prolonged contact may cause skin irritation. May be harmful if swallowed.

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Calcium carbonate	(CAS No) 471-34-1 (EC no) 207-439-9	3,9	Not classified
Bile salts	(CAS No) - (EC no) -	< 1	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319
Iode	(CAS No) 7553-56-2 (EC no) 231-442-4 (EC index no) 053-001-00-3	0,5	Acute Tox. 4 (Inhalation), H332 Acute Tox. 4 (Dermal), H312 Aquatic Acute 1, H400
Novobiocine	(CAS No) 1476-53-5 (EC no) 216-023-6	0,04 - 0,07	Eye Irrit. 2, H319 Skin Sens. 1, H317
Brilliant green	(CAS No) 633-03-4 (EC no) 211-190-1	0,001	Acute Tox. 4 (Oral), H302 Eye Irrit. 2, H319

Full text of H-phrases: see section 16

MKTTn BROTH

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) - Annex II

SECTION 4: First aid measures

4.1. Description of first aid measures

- First-aid measures after inhalation** : Not applicable.
- First-aid measures after skin contact** : Wash with plenty of soap and water.
- First-aid measures after eye contact** : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persist.
- First-aid measures after ingestion** : Rinse mouth. Do NOT induce vomiting. Call a POISON CENTER or doctor/physician if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Not combustible.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : None.

5.3. Advice for firefighters

Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment : Wear personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

6.2. Environmental precautions

Avoid release to the environment. Prevent from entering sewers, basements and workpits, or any place where its accumulation can be dangerous.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up liquid spill into absorbent material. Contain leaking substance in suitable containers.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Wear personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

Hygiene measures : Always wash hands after handling the product.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a dark area. Keep container tightly closed.

Storage temperature : 2 - 8 °C

Packaging materials : Plastic.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Calcium carbonate (471-34-1)

France	Local name	Calcium (carbonate de)
France	VME (mg/m ³)	10 mg/m ³

MKTTn BROTH

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) - Annex II

8.2. Exposure controls

Appropriate engineering controls	: Ensure good ventilation of the work station.
Materials for protective clothing	: Wear suitable protective clothing
Hand protection	: Wear suitable gloves
Eye protection	: Safety glasses with side guards should be worn to prevent injury from airborne particles and/or other eye contact with this product
Skin and body protection	: Wear suitable protective clothing
Respiratory protection	: Not applicable

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: suspension. Presence of precipitate.
Colour	: bluish.
Odour	: No data available
Odour threshold	: No data available
pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

No additional information available

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Moisture.

10.5. Incompatible materials

Acids. Strong oxidizing agents. Magnesium. Aluminium.

10.6. Hazardous decomposition products

No additional information available

MKTTn BROTH

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) - Annex II

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity	: Not classified No data available
Skin corrosion/irritation	: Not classified No data available
Serious eye damage/irritation	: Not classified No data available
Respiratory or skin sensitisation	: Not classified No data available
Germ cell mutagenicity	: Not classified No data available
Carcinogenicity	: Not classified No data available
Reproductive toxicity	: Not classified No data available
Specific target organ toxicity (single exposure)	: Not classified No data available
Specific target organ toxicity (repeated exposure)	: Not classified No data available
Aspiration hazard	: Not classified No data available

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : No data available.

12.2. Persistence and degradability

Persistence and degradability No data available.

12.3. Bioaccumulative potential

Bioaccumulative potential Not established.

12.4. Mobility in soil

Ecology - soil No data available.

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Other adverse effects : No data available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Incineration, disposal or recycling at specific offsite provider.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

Not regulated for transport

14.2. UN proper shipping name

Proper Shipping Name (ADR)	: Not applicable
Proper Shipping Name (IMDG)	: Not applicable
Proper Shipping Name (IATA)	: Not applicable
Proper Shipping Name (ADN)	: Not applicable

MKTTn BROTH

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) - Annex II

Proper Shipping Name (RID) : Not applicable

14.3. Transport hazard class(es)

ADR

Transport hazard class(es) (ADR) : Not applicable

IMDG

Transport hazard class(es) (IMDG) : Not applicable

IATA

Transport hazard class(es) (IATA) : Not applicable

ADN

Transport hazard class(es) (ADN) : Not applicable

RID

Transport hazard class(es) (RID) : Not applicable

14.4. Packing group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

14.5. Environmental hazards

Dangerous for the environment : No

Marine pollutant : No

Other information : No supplementary information available

14.6. Special precautions for user

- Overland transport

No data available

- Transport by sea

No data available

- Air transport

No data available

- Inland waterway transport

Carriage prohibited (ADN) : No

Not subject to ADN : No

- Rail transport

Carriage prohibited (RID) : No

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No additional information available

MKTTn BROTH

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) - Annex II

SECTION 16: Other information

Full text of H- and EUH-statements:

Acute Tox. 4 (Dermal)	Acute toxicity (dermal), Category 4
Acute Tox. 4 (Inhalation)	Acute toxicity (inhal.), Category 4
Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Aquatic Acute 1	Hazardous to the aquatic environment — Acute Hazard, Category 1
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
Skin Sens. 1	Sensitisation — Skin, category 1
H302	Harmful if swallowed
H312	Harmful in contact with skin
H315	Causes skin irritation
H317	May cause an allergic skin reaction
H319	Causes serious eye irritation
H332	Harmful if inhaled
H400	Very toxic to aquatic life

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product