

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form	: Mixture
Product name	: IRIS SALMONELLA® SUPPLEMENT (TABLETS) 10 gr
Product code	: BS093
Product group	: Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category	: Industrial use, Professional use
Use of the substance/mixture	: Microbiological control

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

BIOKAR DIAGNOSTICS - SOLABIA GROUP
 Rue des Quarante Mines
 ZAC de Ther - Allonne
 P.O. Box 10245
 60002 Beauvais Cedex - FRANCE
 T 0033 3 44 14 33 33 - F 0033 3 44 14 33 34
info.fds@solabia.fr - www.biokar-diagnostics.fr

1.4. Emergency telephone number

No additional information available

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Resp. Sens. 1	H334
Skin Sens. 1	H317

Full text of H statements : see section 16

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS08

Signal word (CLP) :

Danger

Hazard statements (CLP) :

H317 - May cause an allergic skin reaction.
 H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precautionary statements (CLP) :

P304+P341 - IF INHALED: if breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
 P302+P352 - IF ON SKIN: Wash with plenty of water/...
 P312 - Call a POISON CENTER/doctor if you feel unwell.

EUH-statements :

EUH208 - Contains . May produce an allergic reaction.

2.3. Other hazards

Adverse physicochemical, human health and environmental effects : May cause an allergic skin reaction. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause moderate irritation to the eyes. Repeated or prolonged contact may cause skin irritation. May be harmful if swallowed.

IRIS SALMONELLA® SUPPLEMENT (TABLETS) 10 gr

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Sulfadiazine sodium	(CAS-No.) 547-32-0 (EC-No.) 208-919-0	1 - 20	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319 Resp. Sens. 1, H334 Skin Sens. 1, H317 STOT SE 3, H335
Novobiocine	(CAS-No.) 1476-53-5 (EC-No.) 216-023-6	0,1 - 1	Eye Irrit. 2, H319 Skin Sens. 1, H317

Full text of H-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: Call a poison center or a doctor if you feel unwell.
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a poison center or a doctor.
First-aid measures after skin contact	: Wash with plenty of water/.... Take off contaminated clothing. If skin irritation or rash occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persists.
First-aid measures after ingestion	: Rinse mouth. Do NOT induce vomiting. Call a POISON CENTER/doctor if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after inhalation	: May cause allergy or asthma symptoms or breathing difficulties if inhaled.
Symptoms/effects after skin contact	: May cause an allergic skin reaction.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Water spray. Dry powder. Foam.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : Toxic fumes may be released.

5.3. Advice for firefighters

Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment : Wear personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

Emergency procedures : Ventilate spillage area. Avoid contact with skin and eyes.

Measures in case of dust release : Avoid raising dust. Avoid breathing dust.

6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

6.2. Environmental precautions

Avoid release to the environment. Prevent from entering sewers, basements and workpits, or any place where its accumulation can be dangerous.

IRIS SALMONELLA® SUPPLEMENT (TABLETS) 10 gr

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

6.3. Methods and material for containment and cleaning up

- For containment** : Collect spillage.
- Methods for cleaning up** : Mechanically recover the product. Take up mechanically (sweeping, shovelling) and collect in suitable container for disposal. Minimise generation of dust.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

- Precautions for safe handling** : Ensure good ventilation of the work station. Do not breathe dust. Avoid contact with skin and eyes. Wear personal protective equipment. For further information refer to section 8: "Exposure controls/personal protection".
- Hygiene measures** : Always wash hands after handling the product. Wash contaminated clothing before reuse.

7.2. Conditions for safe storage, including any incompatibilities

- Storage conditions** : Keep container tightly closed.
- Storage temperature** : 2 - 8 °C
- Packaging materials** : aluminium. PVC.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

- Appropriate engineering controls** : Ensure good ventilation of the work station.
- Materials for protective clothing** : Wear suitable protective clothing
- Hand protection** : Wear suitable gloves
- Eye protection** : Safety glasses with side guards should be worn to prevent injury from airborne particles and/or other eye contact with this product
- Skin and body protection** : Wear suitable protective clothing
- Respiratory protection** : Suitable respiratory equipment
- Environmental exposure controls** : Avoid release to the environment.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

- Physical state** : Solid
- Appearance** : tablet.
- Colour** : light blue.
- Odour** : No data available
- Odour threshold** : No data available
- pH** : No data available
- Relative evaporation rate (butylacetate=1)** : No data available
- Melting point** : No data available
- Freezing point** : Not applicable
- Boiling point** : No data available
- Flash point** : Not applicable
- Auto-ignition temperature** : Not applicable
- Decomposition temperature** : No data available
- Flammability (solid, gas)** : Non flammable.
- Vapour pressure** : No data available
- Relative vapour density at 20 °C** : No data available
- Relative density** : Not applicable
- Solubility** : Soluble in water.
- Log Pow** : No data available

IRIS SALMONELLA® SUPPLEMENT (TABLETS) 10 gr

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

Viscosity, kinematic	: Not applicable
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: Not applicable

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

Strong oxidizing agents.

10.6. Hazardous decomposition products

Nitrogen oxides. carbon oxides. sulphur oxides. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity	: Not classified
	No data available

Sulfadiazine sodium (547-32-0)

RTECS	WP2000000
-------	-----------

Novobiocine (1476-53-5)

LD50 oral rat	3500 mg/kg
---------------	------------

RTECS	RD5425000
-------	-----------

Skin corrosion/irritation	: Not classified
	No data available
Serious eye damage/irritation	: Not classified
	No data available
Respiratory or skin sensitisation	: May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction.
	No data available
Germ cell mutagenicity	: Not classified
	No data available
Carcinogenicity	: Not classified
	No data available
Reproductive toxicity	: Not classified
	No data available
STOT-single exposure	: Not classified
	No data available
STOT-repeated exposure	: Not classified
	No data available
Aspiration hazard	: Not classified
	No data available

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general	: No data available.
--------------------------	----------------------

IRIS SALMONELLA® SUPPLEMENT (TABLETS) 10 gr

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

12.2. Persistence and degradability

Persistence and degradability No data available.

12.3. Bioaccumulative potential

Bioaccumulative potential Not established.

12.4. Mobility in soil

Ecology - soil No data available.

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Other adverse effects : No data available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Incineration, disposal or recycling at specific offsite provider.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

Not regulated for transport

14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not applicable

Proper Shipping Name (IMDG) : Not applicable

Proper Shipping Name (IATA) : Not applicable

Proper Shipping Name (ADN) : Not applicable

Proper Shipping Name (RID) : Not applicable

14.3. Transport hazard class(es)

ADR

Transport hazard class(es) (ADR) : Not applicable

IMDG

Transport hazard class(es) (IMDG) : Not applicable

IATA

Transport hazard class(es) (IATA) : Not applicable

ADN

Transport hazard class(es) (ADN) : Not applicable

RID

Transport hazard class(es) (RID) : Not applicable

14.4. Packing group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

14.5. Environmental hazards

Dangerous for the environment : No

Marine pollutant : No

Other information : No supplementary information available

IRIS SALMONELLA® SUPPLEMENT (TABLETS) 10 gr

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

14.6. Special precautions for user

- Overland transport

No data available

- Transport by sea

No data available

- Air transport

No data available

- Inland waterway transport

Carriage prohibited (ADN) : No

Not subject to ADN : No

- Rail transport

Carriage prohibited (RID) : No

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes:

	Date of issue	Modified	
	Revision date	Modified	
	Supersedes	Modified	
3	Composition/information on ingredients	Modified	

Full text of H- and EUH-statements:

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Resp. Sens. 1	Respiratory sensitisation, Category 1
Skin Irrit. 2	Skin corrosion/irritation, Category 2
Skin Sens. 1	Skin sensitisation, Category 1
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
EUH208	Contains . May produce an allergic reaction.

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.