2024/05/03 Distributed by FUJIFILM Wako Pure Chemical Corporation Manufactured by FUJIFILM Wako Pure Chemical Corporation

SPECIFICATION

1-Chlorobutane

for Amino Acid Sequence Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (230nm)	max.0.25
Absorbance (240nm)	max.0.05
Absorbance (254nm)	max.0.02
Absorbance (260 \sim 400nm)	max.0.01
Peroxides	to pass test
Tollens reactive substances	to pass test
Assay (capillary GC)	min.99.0%