

SPECIFICATION

25% Trifluoroacetic Acid Solution (for Gradient Elution)

for SHIMADZU Protein Peptide Sequencer

| REQUIREMENT | SPECIFICATION |
|---|------------------------|
| Appearance | Colorless clear liquid |
| Absorbance (269nm) | max.0.040 |
| Absorbance (280nm) | max.0.020 |
| Absorbance (300nm) | max.0.010 |
| Absorbance (350nm) | max.0.005 |
| Suitability for amino acids sequence analysis | to pass test |
| Assay | 25.0~25.5% |