

2017年9月作成

品名	D-PBS(-)	10×D-PBS(-)	PBS(-)	10×PBS(-)
Code No.	045-29795 043-29791 049-29793	048-29805	166-23555 164-23551	163-25265
<b>Components</b>	(mg/L)	(mg/L)	(mg/L)	(mg/L)
<b><i>Inorganic salts:</i></b>				
KCl	200.00	2,000.00	-	-
NaCl	8,000.00	80,000.00	9,000.00	90,000.00
KH <sub>2</sub> PO <sub>4</sub>	200.00	2,000.00	144.00	1,440.00
Na <sub>2</sub> HPO <sub>4</sub> (anhyd.)	1,150.00	11,500.00	421.00	4,211.13

品名	HBSS(-) with Phenol Red	HBSS(-) without Phenol Red	10×HBSS(-) without Phenol Red	HBSS(+) with Phenol Red	HBSS(+) without Phenol Red
Code No.	084-08345	085-09355	082-09865	082-09365	084-08965 082-08961
<b>Components</b>	(mg/L)	(mg/L)	(mg/L)	(mg/L)	(mg/L)
<b><i>Inorganic salts:</i></b>					
KCl	400.00	400.00	4,000.00	400.00	400.00
NaCl	8,000.00	8,000.00	80,000.00	8,000.00	8,000.00
NaHCO <sub>3</sub>	350.00	350.00	-	350.00	350.00
KH <sub>2</sub> PO <sub>4</sub>	60.00	60.00	600.00	60.00	60.00
Na <sub>2</sub> HPO <sub>4</sub> (anhyd.)	47.90	47.90	479.00	47.80	47.80
CaCl <sub>2</sub> (anhyd.)	-	-	-	140.00	140.00
MgCl <sub>2</sub> 6H <sub>2</sub> O	-	-	-	100.00	100.00
MgSO <sub>4</sub> 7H <sub>2</sub> O	-	-	-	100.00	100.00
<b><i>Other components:</i></b>					
D-Glucose	1,000.00	1,000.00	10,000.00	1,000.00	1,000.00
Phenol Red	10.00	-	-	10.00	-