

2025年7月 改訂

品名	Penicillin-Streptomycin Solution (×50)	Penicillin-Streptomycin Solution (×100)	Penicillin-Streptomycin-Amphotericin B Suspension (×100)	Penicillin-Streptomycin-L-Glutamine Solution (×100)
Code No.	164-25251	168-23191	161-23181	161-23201
Components	(mg/L)	(mg/L)	(mg/L)	(mg/L)
Sodium Penicillin G	※5,000,000.00	※10,000,000.00	※10,000,000.00	※10,000,000.00
Streptomycin Sulfate	5,000.00	10,000.00	10,000.00	-
Streptomycin	-	-	-	10,000.00
Amphotericin B	-	-	25.00	-
Sodium Deoxycholate	-	-	205.00	-
L-Glutamine	-	-	-	29,200.00
NaCl	8,500.00	8,500.00	8,500.00	8,500.00

※ : units/L

品名	Gentamicin Sulfate Solution
Code No.	078-06061
Components	(mg/L)
Gentamicin Sulfate	50,000.00

品名	Kanamycin Sulfate Solution
Code No.	117-00961
Components	(mg/L)
Kanamycin Sulfate	50,000.00
NaCl	8,500.00

品名	Amphotericin B Suspension
Code No.	019-23891
Components	(mg/L)
Amphotericin B	250.00
Sodium Deoxycholate	205.00

品名	1mg/ml Mitomycin C Solution
Code No.	133-15931
Components	(mg/L)
Mitomycin C	1,000.00

10v/v%エタノール、90v/v%エチレングリコールで調製されています。

品名	50mg/mL G-418 Sulfate Solution
Code No.	074-06801 070-06803
Components	(mg/L)
G-418 2H ₂ SO ₄	69,752.29

品名	10mg/mL Blasticidin S Hydrochloride Solution
Code No.	026-18711 022-18713
Components	(mg/L)
Blasticidin S Hydrochloride	10,000.00
HEPES	4,766.20