

90×190mmサイズ

**FUJIFILM**

**Wako**

Code No. 015-22411 (50  $\mu$ L)

**Anti Ago1, Monoclonal Antibody (2A7)**  
**抗Ago1, モノクローナル抗体(2A7)**

Monoclonal Anti Ago1 (Argonaute 1) antibody is produced by hybridoma clone 2A7, established by fusing mouse myeloma cells P3U1 and lymphocytes of BALB/c mice immunized with synthetic partial peptides of Ago1 protein which was coupled with bovine thyroglobulin (BTG).

This antibody was purified from mouse ascites by the Protein A affinity chromatography.

This product is for laboratory research use only ; use in any other application is the responsibility of the user.

**[Concentration]**

Indicated on the label.

**[Formulation]**

1 × TBS (pH 7.4) aqueous solution with 10% glycerol containing 0.05% sodium azide.

**[Clone No.]**

2A7

**[Cross Reactivity]**

Human and mouse Ago1

**[Subclass]**

IgG<sub>2a</sub>  $\kappa$

**[Working Dilution]**

Immunoprecipitation 5~10  $\mu$ g/20  $\mu$ L 10% Protein G beads slurry

**[Storage]**

Keep at 4 °C in the dark.

Avoid repeated freeze and thaw.

**[Note]**

Even if some aggregates can be seen, the quality of this antibody solution was confirmed in normal use.

**[Package]**

50  $\mu$ L

- 1/2 -

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-311-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1802KA1

- 2/2 -