

**FUJIFILM**

**Wako**

Code No. 018-23241 (50  $\mu$ L)

**Anti Human Ago3, Monoclonal Antibody (1C12)**

**抗ヒトAgo3, モノクローナル抗体 (1C12)**

Monoclonal Anti Human Ago3 (Argonaute 3) antibody is produced by hybridoma clone 1C12, established by fusing mouse myeloma cells P3U1 and lymphocytes of BALB/c mice immunized with synthetic partial peptides of human Ago3 protein which was coupled with bovine thyroglobulin (BTG).

This antibody was purified from mouse ascites by using affinity chromatography.

This product is for laboratory research use only ; use in any other application is the responsibility of the user.

**[Concentration]**

Indicated on the label.

**[Formulation]**

1  $\times$  TBS (pH 7.4) aqueous solution with 10% glycerol containing 0.05% sodium azide.

**[Clone No.]**

1C12

**[Subclass]**

IgG<sub>1</sub>

**[Specificity]**

Specific for human Ago3

**[Working Dilution]**

Immunoprecipitation 5~10  $\mu$ g/20  $\mu$ L Protein G beads slurry

**[Storage]**

Keep at 2~10  $^{\circ}$ C in the dark.

Avoid repeated freeze and thaw.

**[Package]**

50  $\mu$ L

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-311-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1802KA1