

90 × 190mm サイズ

FUJIFILM

Wako

< for Research Use Only >

Code No. 011-25431 (1 mg)

017-25433 (1 mg × 5)

ACG Lectin, recombinant, Solution
ACG レクチン, 組換え体, 溶液

ACG (*Agrocybe cylindracea* galectin) lectin is isolated from *Agrocybe cylindracea*, the mushroom. ACG is categorized into galectin family and recognizes not only α 2-3 linked sialic acid containing glycoconjugates but also various β -galactosides.

Source : *E. coli* expressed *Agrocybe cylindracea* galectin lectin (ACG).

Appearance : PBS (pH 7.4).

Protein concentration : Indicated on the label.

Carbohydrate specificity : α 2-3Sia

Manufacture : Affinity purification by lactose agarose.

[Expiration date] Indicated on the label.

[Storage] Store at -20°C . After thawing, keep at $2-10^{\circ}\text{C}$.
Avoid the freeze-thaw cycles.

[Package]

Code No.	Packaging
011-25431	1 mg
017-25433	1 mg × 5 vials

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-3111100
<http://www.wako-chemicals.de>

1801KA1