

90 × 190mm サイズ

FUJIFILM

Wako

< for Research Use Only >

Code No. 033-23771 (1 mg)

039-23773 (1 mg × 5)

CGL2 Lectin, recombinant, Solution
CGL2レクチン, 組換え体, 溶液

CGL (*Coprinopsis cinerea* lectin) is isolated from *Coprinopsis cinerea*, the mushroom. CGL2 is categorized into galectin and recognizes galactose.

Source : *E. coli* expressed *Coprinopsis cinerea* lectin (CGL2).

Appearance : PBS (pH 7.4).

Protein concentration : Indicated on the label.

Carbohydrate specificity : β Gal, GalNAc α 1-3Gal (Blood Group A), Gal α 1-3Gal (Blood Group B)

Manufacture : Affinity purification by lactose agarose.

[**Expiration date**] Indicated on the label.

[**Storage**] Store at -20°C . After thawing, keep at $2-10^{\circ}\text{C}$.
Avoid the freeze-thaw cycles.

[**Package**]

Code No.	Packaging
033-23771	1 mg
039-23773	1 mg × 5 vials

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

1802KA1