

90 × 190mm サイズ

**FUJIFILM**

**Wako**

< for Research Use Only >

Code No. 136-18221 (1 mg)  
132-18223 (1 mg × 5)

**Malectin Lectin, recombinant, Solution**  
**マレクチンレクチン, 組換え体, 溶液**

Malectin is an endoplasmic reticulum resident lectin involving in glycoprotein quality control and known that it is highly conserved in animals. It is known that specially bind to Glc<sub>2</sub>-high-mannose *N*-glycan.

**Source** : *E. coli* expressed human malectin.

**Appearance** : PBS (pH 7.4).

**Protein concentration** : Indicated on the label.

**Carbohydrate specificity** : Glc<sub>2</sub>-*N*-glycan

**Manufacture** : Affinity purification by maltose agarose.

[**Expiration date**] Indicated on the label.

[**Storage**] Store at -20 °C. After thawing, keep at 2-10 °C.  
Avoid the freeze-thaw cycles.

[**Package**]

Code No.	Packaging
136-18221	1 mg
132-18223	1 mg × 5 vials

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggenstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-311-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1803KA1