

90 × 190mm サイズ

FUJIFILM

Wako

< for Research Use Only >

Code No. 072-06341 (1 mg)
078-06343 (1 mg × 5)

Gal1-S Lectin, recombinant, Solution
Gal1-Sレクチン, 組換え体, 溶液
(Galectin-1)

Galectins are a family of proteins which typically bind β -galactose-containing glycoconjugates and conserve the primary structural homology in their carbohydrate recognition domains (CRDs). The galectin family is classified into three groups according to their structure: the prototypical-type, chimeric-type and tandem-repeat-type.

This galectin-1 is the prototypical-type. Galectin-1 shows lectin activity only in the reduced form. The oxidized form of it exhibits no lectin activities, and shows several biological activities.

Source : *E. coli* expressed human galectin1-S (Gal1-S).

Appearance : PBS (pH 7.4).

Protein concentration : Indicated on the label.

Carbohydrate specificity : Gal, branched LacNAc

Manufacture : Affinity purification by lactose agarose.

[**Expiration date**] Indicated on the label.

[**Storage**] Store at -20°C . After thawing, keep at $2-10^{\circ}\text{C}$.
Avoid the freeze-thaw cycles.

[**Package**]

Code No.	Packaging
072-06341	1 mg
078-06343	1 mg × 5 vials

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>