

90×190mmサイズ

**FUJIFILM**

**Wako**

< for Research Use Only >

Code No. 168-26751 (0.5 mg)

**PPL Lectin, recombinant, Solution**

**PPLレクチン, 組換え体, 溶液**

---

PPL is isolated from *Pleurocybella porrigens*, the mushroom. It is known as the galactose-binding lectins and binds to blood group A, B and O.

**Source :** *E. coli* expressed *Pleurocybella porrigens* lectin (PPL).

**Appearance :** PBS (pH 7.4).

**Protein concentration :** Indicated on the label.

**Carbohydrate specificity :**  $\alpha/\beta$  GalNAc

**Manufacture :** Affinity purification by lactose agarose.

[**Expiration date**] Indicated on the label.

[**Storage**] Store at  $-20^{\circ}\text{C}$ . After thawing, keep at  $2-10^{\circ}\text{C}$ .  
Avoid the freeze-thaw cycles.

[**Package**] 0.5 mg

---

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-311-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1803KA1