

90 × 210mm サイズ

**FUJIFILM**

**Wako**

< for Research Use Only >

Code No. 010-26861 (100  $\mu$ g)

**Anti TERT, Monoclonal Antibody (TMab-6)**  
**抗 TERT, モノクローナル抗体 (TMab-6)**

TERT (Teromerase Reverse Transcriptase) is one of the subunit of teromerase, which consists of TERT, TERC, dyskerin and more subunits.

The mutations in the promoter region of *TERT* gene have been found in various cancers including diffuse gliomas, and have been demonstrated the importance of the TERT promoter mutations in the classification of diffuse gliomas.

Anti TERT, monoclonal antibody (TMab-6) is established against the synthetic peptide of hTERT. This antibody is useful for immunohistochemistry.

This product is for laboratory use only ; use in any such application is the responsibility of the users.

**Description**

[Reactivity] Reacts with human TERT  
[Species cross reactivity] Human  
[Clone No.] TMab-6  
[Host] Mouse  
[Isotype] IgM  $\cdot$   $\kappa$   
[Concentration] Indicated on the label.  
[Formulation] PBS (pH 7.4) containing sodium azide

**Application**

Western Blot 1  $\mu$ g/mL ~  
Immunohistochemistry 1  $\mu$ g/mL ~  
ELISA 1  $\mu$ g/mL ~

\*Please optimize the most appropriate concentration for your analysis conditions and samples.

**Storage**

Store at -20  $^{\circ}$ C.  
Avoid repeated freeze and thaw.

**Package**

100  $\mu$ g

**Reference**

1) Masui, k. *et al.*, in press.

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1801KA1