

FUJIFILM**Wako**Code No. 013-27951 (20 μ L)
019-27953 (100 μ L)**Anti CD9, Rat Monoclonal Antibody (30B),
Biotin Conjugated**
**抗 CD9, ラットモノクローナル抗体 (30B),
ビオチン結合**

The 30B monoclonal antibody specifically binds to CD9. CD9 is a cell surface glycoprotein which belongs to the tetraspanin superfamily and is also known to form a complex with integrins and other transmembrane modulate cell adhesion and migration and also trigger platelet activation and aggregation.

This product is antibody conjugated biotin to Anti CD9, Monoclonal antibody (30B).

This product is intended for laboratory use only.

[Formulation]

1 \times TBS aqueous solution with 50% glycerol, 0.05 w/v% sodium azide, pH 7.2.

[Clone No.]

30B

[Host]

Rat

[Cross Reactivity]

Human CD9

Note : It shows a little reactivity with rat CD9.

[Subclass]

IgG2b

[Application]Western blotting^{※1} 1 : 1,000 - 1 : 4,000 (under non-reducing condition)

ELISA 1 : 2,000 - 1 : 8,000

Please optimize the most appropriate concentration for your analysis and samples.

※1 Use this antibody under non-reducing condition.

[Storage]

Store at -20°C.

[Package]

Code No.	Packaging
013-27951	20 μ L
019-27953	100 μ L

FUJIFILM Wako Pure Chemical Corporation1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>**FUJIFILM Wako Chemicals U.S.A. Corporation**1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>**FUJIFILM Wako Chemicals Europe GmbH**Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-3111100
<http://www.wako-chemicals.de>