

**FUJIFILM****Wako**Code No. 010-26883 (50  $\mu$ L)**Anti Amyloid  $\beta$ , Monoclonal Antibody (BNT77)**  
**抗アミロイド $\beta$ , モノクローナル抗体 (BNT77)****[Background]**

Amyloid  $\beta$  (A $\beta$ ) is a peptide which consists of about 40 amino acids. A $\beta$  is considered the cause of Alzheimer's disease, because it is a major component of senile plaques formed in brains of Alzheimer's disease patients. This item is a mouse monoclonal antibody that reacts with A $\beta$ <sup>1)</sup>. The antibody was developed by Takeda Pharmaceutical Company Limited.

**[Description]**

**[Reactivity]** Reacts with Amyloid  $\beta$   
**[Antigen]** 11-28 a.a of human Amyloid  $\beta$   
**[Species cross reactivity]** Human, mouse and rat  
**[Clone No.]** BNT77  
**[Host]** Mouse  
**[Isotype]** IgA $\cdot\kappa$   
**[Purification]** Protein A-affinity purified  
**[Concentration]** Indicated on the label  
**[Formulation]** PBS with 0.05% Sodium Azide

**[Applications]**

Immunoprecipitation<sup>2) 3)</sup> 1 : 500  
Immunohistochemistry 1 : 500-10,000  
ELISA<sup>1) 2) 3)</sup> 1 : 250-1,000

Optimal concentration should be determined by each laboratory for each application.

**[Storage]**

Store at -20°C.  
Avoid repeated freeze and thaw.

**[Package]**

Code No.	Packaging
010-26883	50 $\mu$ L

**[References]**

- 1) Asami-Odaka, A. *et. al.* : *Biochemistry*, **34** (32), 10272 (1995).
- 2) Tomita, T. *et. al.* : *J. Neurosci.*, **19** (24), 10627 (1999).
- 3) Sudoh, S. *et. al.* : *J. Neurochem.*, **71** (4), 1535 (1998).

**[Acknowledgement]**

The Purchaser shall use the Product at its research facility solely for its internal research. A separate license by Takeda Pharmaceutical Company Limited (1-1, Doshomachi 4-chome, Chuo-ku, Osaka 540-8645, Japan) is required for Commercial Use of the Product. Information about a license for the Commercial Use of the Product may be obtained from :

FUJIFILM Wako Pure Chemical Corporation. Laboratory Chemicals Division, 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>