

FUJIFILM**Wako**Code No. 010-27221 (10 μ L)
016-27223 (50 μ L)**Anti Fibrin, Monoclonal Antibody (FIB1-11)**
抗フィブリン, モノクローナル抗体(FIB1-11)**Background**

Fibrin is a major component of blood clots produced in response to bleeding. It is formed from a soluble protein in plasma, fibrinogen. When bleeding occurs due to tissue damage, fibrinogen is converted to Fibrin at the wound site by activity of thrombin which is the clotting enzyme. This product is an antibody that specifically recognizes Fibrin¹⁾⁻³⁾.

Description

[Reactivity]	Reacts with Fibrin
[Antigen]	Human Fibrin N-terminal peptide
[Species cross reactivity]	Human, rabbit and baboon
[Clone No.]	FIB1-11
[Host]	Mouse
[Isotype]	IgG1
[Concentration]	Indicated on the label
[Formulation]	PBS

Applications

ELISA 1 : 100-10,000

Optimal concentration should be determined by each laboratory for each application.

Storage

Store at -20 °C.
Avoid repeated freeze and thaw.

Package

Code No.	Packaging
010-27221	10 μ L
016-27223	50 μ L

References

- 1) Kurokawa, T., *et al.* : *Thromb. Haemost.*, **66**, 6, 684 (1991).
- 2) Kurokawa, T., *et al.* : *Thromb. Res. Suppl.*, **10**, 83 (1990).
- 3) Tada, H., *et al.* : *J. Biotechnol.*, **33**, 2, 157 (1994).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA3