

FUJIFILM

Wako

Code No. 010-27721 (100 μ L)

Anti PA tag, Rat Monoclonal Antibody, Fluorescein Conjugated

PA tag is consisted of 12 amino acids (GVAMPGAEDDVV) and used for purification of PA tag fusion recombinant proteins.

Anti PA tag, Rat Monoclonal Antibody, Fluorescein Conjugated is the antibody derived from hybridoma (Clone No. NZ-1) and recognizes PA tag. This product is conjugated to fluorescein.

【Concentration】

Indicated on the label

【Formulation】

1 $\times$ PBS aqueous solution containing 0.05% Sodium Azide

【Storage】

Store in the dark at 2-10°C

【Clone No.】

NZ-1

【Subclass】

IgG_{2a}

【Host】

Rat

【Specificity】

Specific for GVAMPGAEDDVV peptide.

【Conjugate】

Fluorescein

【Spectral Property】

Excitation : 494 nm, Emission : 521 nm

【Application】

Flow cytometry : 1:10-1:1,000

【Intended use】

For research use only

【Package】

Code No.	Packaging
010-27721	100 μ L

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2304KA2